

---

# PyMeasure Documentation

*Release 0.9.0*

**PyMeasure Developers**

**Nov 28, 2021**




# LEARNING PYMEASURE

<b>1</b>	<b>Introduction</b>	<b>3</b>
1.1	Instrument ready . . . . .	3
1.2	Graphical displays . . . . .	3
<b>2</b>	<b>Quick start</b>	<b>5</b>
2.1	Setting up Python . . . . .	5
2.2	Installing PyMeasure . . . . .	5
<b>3</b>	<b>Tutorials</b>	<b>7</b>
3.1	Connecting to an instrument . . . . .	7
3.2	Making a measurement . . . . .	9
3.3	Using a graphical interface . . . . .	16
<b>4</b>	<b>pymeasure.adapters</b>	<b>31</b>
4.1	Adapter base class . . . . .	31
4.2	Fake adapter . . . . .	32
4.3	Serial adapter . . . . .	33
4.4	Prologix adapter . . . . .	34
4.5	VISA adapter . . . . .	36
4.6	VXI-11 adapter . . . . .	38
4.7	Telnet adapter . . . . .	39
<b>5</b>	<b>pymeasure.experiment</b>	<b>41</b>
5.1	Experiment class . . . . .	41
5.2	Listener class . . . . .	42
5.3	Procedure class . . . . .	43
5.4	Parameter classes . . . . .	44
5.5	Worker class . . . . .	46
5.6	Results class . . . . .	47
<b>6</b>	<b>pymeasure.display</b>	<b>49</b>
6.1	Browser classes . . . . .	49
6.2	Curves classes . . . . .	49
6.3	Inputs classes . . . . .	50
6.4	Listeners classes . . . . .	52
6.5	Log classes . . . . .	52
6.6	Manager classes . . . . .	52
6.7	Plotter class . . . . .	53
6.8	Qt classes . . . . .	54
6.9	Thread classes . . . . .	54
6.10	Widget classes . . . . .	54

6.11	Windows classes . . . . .	57
<b>7</b>	<b>pymeasure.instruments</b>	<b>61</b>
7.1	Instrument classes . . . . .	61
7.2	Validator functions . . . . .	64
7.3	Comedi data acquisition . . . . .	65
7.4	Resource Manager . . . . .	65
7.5	Advantest . . . . .	66
7.6	Agilent . . . . .	66
7.7	Ametek . . . . .	103
7.8	AMI . . . . .	105
7.9	Anaheim Automation . . . . .	107
7.10	Anapico . . . . .	109
7.11	Anritsu . . . . .	110
7.12	Attocube . . . . .	112
7.13	BK Precision . . . . .	115
7.14	Danfysik . . . . .	116
7.15	Delta Elektronika . . . . .	119
7.16	Edwards . . . . .	120
7.17	Fluke . . . . .	120
7.18	F.W. Bell . . . . .	121
7.19	Hewlett Packard . . . . .	122
7.20	Keithley . . . . .	126
7.21	Keysight . . . . .	158
7.22	Lake Shore Cryogenics . . . . .	164
7.23	Newport . . . . .	167
7.24	National Instruments . . . . .	168
7.25	Oxford Instruments . . . . .	180
7.26	Parker . . . . .	182
7.27	Pendulum . . . . .	184
7.28	Razorbill . . . . .	185
7.29	Rohde & Schwarz . . . . .	186
7.30	Signal Recovery . . . . .	200
7.31	Stanford Research Systems . . . . .	203
7.32	Tektronix . . . . .	213
7.33	Thorlabs . . . . .	214
7.34	Yokogawa . . . . .	215
<b>8</b>	<b>Contributing</b>	<b>219</b>
8.1	Using the development version . . . . .	219
8.2	Working on a new feature . . . . .	220
8.3	Making a pull-request . . . . .	220
8.4	Unit testing . . . . .	220
<b>9</b>	<b>Reporting an error</b>	<b>223</b>
<b>10</b>	<b>Adding instruments</b>	<b>225</b>
10.1	File structure . . . . .	225
10.2	Instrument file . . . . .	226
10.3	Writing properties . . . . .	227
10.4	Advanced properties . . . . .	228
<b>11</b>	<b>Coding Standards</b>	<b>235</b>
11.1	Python style guides . . . . .	235
11.2	Documentation . . . . .	235

11.3 Usage of getter and setter functions . . . . .	235
<b>12 Authors</b>	<b>237</b>
<b>13 License</b>	<b>239</b>
<b>Python Module Index</b>	<b>241</b>
<b>Index</b>	<b>243</b>





# PyMeasure

PyMeasure makes scientific measurements easy to set up and run. The package contains a repository of instrument classes and a system for running experiment procedures, which provides graphical interfaces for graphing live data and managing queues of experiments. Both parts of the package are independent, and when combined provide all the necessary requirements for advanced measurements with only limited coding.

Installing Python and PyMeasure are demonstrated in the [Quick Start guide](#). From there, checkout the existing [instruments that are available for use](#).

PyMeasure is currently under active development, so please report any issues you experience on our [Issues page](#).

The main documentation for the site is organized into a couple sections:

- [Learning PyMeasure](#)
- [API References](#)
- [About PyMeasure](#)

Information about development is also available:

- [Getting involved](#)





## INTRODUCTION

PyMeasure uses an object-oriented approach for communicating with scientific instruments, which provides an intuitive interface where the low-level SCPI and GPIB commands are hidden from normal use. Users can focus on solving the measurement problems at hand, instead of re-inventing how to communicate with instruments.

Instruments with VISA (GPIB, Serial, etc) are supported through the [PyVISA package](#) under the hood. [Prologix GPIB](#) adapters are also supported. Communication protocols can be swapped, so that instrument classes can be used with all supported protocols interchangeably.

Before using PyMeasure, you may find it helpful to be acquainted with [basic Python programming for the sciences](#) and understand the concept of objects.

### 1.1 Instrument ready

The package includes a number of *instruments already defined*. Their definitions are organized based on the manufacturer name of the instrument. For example the class that defines the *Keithley 2400 SourceMeter* can be imported by calling:


```
from pymeasure.instruments.keithley import Keithley2400
```

The [Tutorials](#) section will go into more detail on *connecting to an instrument*. If you don't find the instrument you are looking for, but are interested in contributing, see the documentation on *adding an instrument*.

### 1.2 Graphical displays

Graphical user interfaces (GUIs) can be easily generated to manage execution of measurement procedures with PyMeasure. This includes live plotting for data, and a queue system for managing large numbers of experiments.

These features are explored in the *Using a graphical interface* tutorial.



## QUICK START

This section provides instructions for getting up and running quickly with PyMeasure.

### 2.1 Setting up Python

The easiest way to install the necessary Python environment for PyMeasure is through the [Anaconda distribution](#), which includes 720 scientific packages. The advantage of using this approach over just relying on the `pip` installer is that it Anaconda correctly installs the required Qt libraries.

Download and install the appropriate Python version of [Anaconda](#) for your operating system.

### 2.2 Installing PyMeasure

#### 2.2.1 Install with conda

If you have the [Anaconda distribution](#) you can use the conda package mangager to easily install PyMeasure and all required dependencies.

Open a terminal and type the following commands (on Windows look for the *Anaconda Prompt* in the Start Menu):

```
conda config --add channels conda-forge
conda install pymeasure
```

This will install PyMeasure and all the required dependencies.

#### 2.2.2 Install with pip

PyMeasure can also be installed with `pip`.

```
pip install pymeasure
```

Depending on your operating system, using this method may require additional work to install the required dependencies, which include the Qt libraries.

### 2.2.3 Checking the version

Now that you have Python and PyMeasure installed, open up a “Jupyter Notebook” to test which version you have installed. Execute the following code into a notebook cell.

```
import pymeasure
pymeasure.__version__
```

You should see the version of PyMeasure printed out. At this point you have PyMeasure installed, and you are ready to start using it! Are you ready to *connect to an instrument*?

## TUTORIALS

The following sections provide instructions for getting started with PyMeasure.

### 3.1 Connecting to an instrument

After following the *Quick Start* section, you now have a working installation of PyMeasure. This section describes connecting to an instrument, using a Keithley 2400 SourceMeter as an example. To follow the tutorial, open a command prompt, IPython terminal, or Jupyter notebook.

First import the instrument of interest.

```
from pymeasure.instruments.keithley import Keithley2400
```

Then construct an object by passing the GPIB address. For this example we connect to the instrument over GPIB (using VISA) with an address of 4. See the *adapters* section below for more details.

```
sourcemeter = Keithley2400("GPIB::4")
```

For instruments with standard SCPI commands, an `id` property will return the results of a `*IDN?` SCPI command, identifying the instrument.

```
sourcemeter.id
```

This is equivalent to manually calling the SCPI command.

```
sourcemeter.ask("*IDN?")
```

Here the `ask` method writes the SCPI command, reads the result, and returns that result. This is further equivalent to calling the methods below.

```
sourcemeter.write("*IDN?")  
sourcemeter.read()
```

This example illustrates that the top-level methods like `id` are really composed of many lower-level methods. Both can be called depending on the operation that is desired. PyMeasure hides the complexity of these lower-level operations, so you can focus on the bigger picture.

### 3.1.1 Using adapters

PyMeasure supports a number of adapters, which are responsible for communicating with the underlying hardware. In the example above, we passed the string “GPIB::4” when constructing the instrument. By default this constructs a `VISAAdapter` class to connect to the instrument using VISA. Instead of passing a string, we could equally pass an adapter object.

```
from pymeasure.adapters import VISAAdapter

adapter = VISAAdapter("GPIB::4")
sourcemeter = Keithely2400(adapter)
```

To instead use a Prologix GPIB device connected on `/dev/ttyUSB0` (proper permissions are needed in Linux, see [PrologixAdapter](#)), the adapter is constructed in a similar way. Unlike the VISA adapter which is specific to each instrument, the Prologix adapter can be shared by many instruments. Therefore, they are addressed separately based on the GPIB address number when passing the adapter into the instrument construction.

```
from pymeasure.adapters import PrologixAdapter

adapter = PrologixAdapter('/dev/ttyUSB0')
sourcemeter = Keithley2400(adapter.gpib(4))
```

For instruments using serial communication that have particular settings that need to be matched, a custom *Adapter* sub-class can be made. For example, the LakeShore 425 Gaussmeter connects via USB, but uses particular serial communication settings. Therefore, a *LakeShoreUSBAdapter* class enables these requirements in the background.

```
from pymeasure.instruments.lakeshore import LakeShore425

gaussmeter = LakeShore425('/dev/lakeshore425')
```

Behind the scenes the `/dev/lakeshore425` port is passed to the *LakeShoreUSBAdapter*.

Some equipment may require the vxi-11 protocol for communication. An example would be a Agilent E5810B ethernet to GPIB bridge. To use this type equipment the `python-vxi11` library has to be installed which is part of the extras package requirements.

```
from pymeasure.adapters import VXI11Adapter
from pymeasure.instruments import Instrument

adapter = VXI11Adapter("TCPIP::192.168.0.100::inst0::INSTR")
instr = Instrument(adapter, "my_instrument")
```

The above examples illustrate different methods for communicating with instruments, using adapters to keep instrument code independent from the communication protocols. Next we present the methods for setting up measurements.

## 3.2 Making a measurement

This tutorial will walk you through using PyMeasure to acquire a current-voltage (IV) characteristic using a Keithley 2400. Even if you don't have access to this instrument, this tutorial will explain the method for making measurements with PyMeasure. First we describe using a simple script to make the measurement. From there, we show how *Procedure* objects greatly simplify the workflow, which leads to making the measurement with a graphical interface.

### 3.2.1 Using scripts

Scripts are a quick way to get up and running with a measurement in PyMeasure. For our IV characteristic measurement, we perform the following steps:

- 1) Import the necessary packages
- 2) Set the input parameters to define the measurement
- 3) Connect to the Keithley 2400
- 4) Set up the instrument for the IV characteristic
- 5) Allocate arrays to store the resulting measurements
- 6) Loop through the current points, measure the voltage, and record
- 7) Save the final data to a CSV file
- 8) Shutdown the instrument

These steps are expressed in code as follows.

```
# Import necessary packages
from pymeasure.instruments.keithley import Keithley2400
import numpy as np
import pandas as pd
from time import sleep

# Set the input parameters
data_points = 50
averages = 50
max_current = 0.01
min_current = -max_current

# Connect and configure the instrument
sourcemeter = Keithley2400("GPIB::4")
sourcemeter.reset()
sourcemeter.use_front_terminals()
sourcemeter.measure_voltage()
sourcemeter.config_current_source()
sleep(0.1) # wait here to give the instrument time to react
sourcemeter.set_buffer(averages)

# Allocate arrays to store the measurement results
currents = np.linspace(min_current, max_current, num=data_points)
voltages = np.zeros_like(currents)
voltage_stds = np.zeros_like(currents)
```

(continues on next page)

(continued from previous page)

```

# Loop through each current point, measure and record the voltage
for i in range(data_points):
    sourcemeter.current = currents[i]
    sourcemeter.reset_buffer()
    sleep(0.1)
    sourcemeter.start_buffer()
    sourcemeter.wait_for_buffer()

    # Record the average and standard deviation
    voltages[i] = sourcemeter.means
    voltage_stds[i] = sourcemeter.standard_devs

# Save the data columns in a CSV file
data = pd.DataFrame({
    'Current (A)': currents,
    'Voltage (V)': voltages,
    'Voltage Std (V)': voltage_stds,
})
data.to_csv('example.csv')

sourcemeter.shutdown()

```

Running this example script will execute the measurement and save the data to a CSV file. While this may be sufficient for very basic measurements, this example illustrates a number of issues that PyMeasure solves. The issues with the script example include:

- The progress of the measurement is not transparent
- Input parameters are not associated with the data that is saved
- Data is not plotted during the execution (nor at all in this case)
- Data is only saved upon successful completion, which is otherwise lost
- Canceling a running measurement causes the system to end in an undetermined state
- Exceptions also end the system in an undetermined state

The *Procedure* class allows us to solve all of these issues. The next section introduces the *Procedure* class and shows how to modify our script example to take advantage of these features.

### 3.2.2 Using Procedures

The Procedure object bundles the sequence of steps in an experiment with the parameters required for its successful execution. This simple structure comes with huge benefits, since a number of convenient tools for making the measurement use this common interface.

Let's start with a simple example of a procedure which loops over a certain number of iterations. We make the SimpleProcedure object as a sub-class of Procedure, since SimpleProcedure *is a* Procedure.

```

from time import sleep
from pymeasure.experiment import Procedure
from pymeasure.experiment import IntegerParameter

class SimpleProcedure(Procedure):

```

(continues on next page)



(continued from previous page)

```

# a Parameter that defines the number of loop iterations
iterations = IntegerParameter('Loop Iterations')

# a list defining the order and appearance of columns in our data file
DATA_COLUMNS = ['Iteration']

def execute(self):
    """ Loops over each iteration and emits the current iteration,
    before waiting for 0.01 sec, and then checking if the procedure
    should stop """
    for i in range(self.iterations):
        self.emit('results', {'Iteration': i})
        sleep(0.01)
        if self.should_stop():
            break

```

At the top of the SimpleProcedure class we define the required Parameters. In this case, `iterations` is a `IntegerParameter` that defines the number of loops to perform. Inside our Procedure class we reference the value in the iterations Parameter by the class variable where the Parameter is stored (`self.iterations`). PyMeasure swaps out the Parameters with their values behind the scene, which makes accessing the values of parameters very convenient.

We define the data columns that will be recorded in a list stored in `DATA_COLUMNS`. This sets the order by which columns are stored in the file. In this example, we will store the Iteration number for each loop iteration.

The `execute` methods defines the main body of the procedure. Our example method consists of a loop over the number of iterations, in which we emit the data to be recorded (the Iteration number). The data is broadcast to any number of listeners by using the `emit` method, which takes a topic as the first argument. Data with the 'results' topic and the proper data columns will be recorded to a file. The sleep function in our example provides two very useful features. The first is to delay the execution of the next lines of code by the time argument in units of seconds. The seconds is that during this delay time, the CPU is free to perform other code. Successful measurements often require the intelligent use of sleep to deal with instrument delays and ensure that the CPU is not hogged by a single script. After our delay, we check to see if the Procedure should stop by calling `self.should_stop()`. By checking this flag, the Procedure will react to a user canceling the procedure execution.

This covers the basic requirements of a Procedure object. Now let's construct our SimpleProcedure object with 100 iterations.

```

procedure = SimpleProcedure()
procedure.iterations = 100

```

Next we will show how to run the procedure.

## Running Procedures

A Procedure is run by a Worker object. The Worker executes the Procedure in a separate Python thread, which allows other code to execute in parallel to the procedure (e.g. a graphical user interface). In addition to performing the measurement, the Worker spawns a Recorder object, which listens for the ‘results’ topic in data emitted by the Procedure, and writes those lines to a data file. The Results object provides a convenient abstraction to keep track of where the data should be stored, the data in an accessible form, and the Procedure that pertains to those results.

We first construct a Results object for our Procedure.

```
from pymeasure.experiment import Results

data_filename = 'example.csv'
results = Results(procedure, data_filename)
```

Constructing the Results object for our Procedure creates the file using the data\_filename, and stores the Parameters for the Procedure. This allows the Procedure and Results objects to be reconstructed later simply by loading the file using Results.load(data\_filename). The Parameters in the file are easily readable.

We now construct a Worker with the Results object, since it contains our Procedure.

```
from pymeasure.experiment import Worker

worker = Worker(results)
```

The Worker publishes data and other run-time information through specific queues, but can also publish this information over the local network on a specific TCP port (using the optional port argument. Using TCP communication allows great flexibility for sharing information with Listener objects. Queues are used as the standard communication method because they preserve the data order, which is of critical importance to storing data accurately and reacting to the measurement status in order.

Now we are ready to start the worker.

```
worker.start()
```

This method starts the worker in a separate Python thread, which allows us to perform other tasks while it is running. When writing a script that should block (wait for the Worker to finish), we need to join the Worker back into the main thread.

```
worker.join(timeout=3600) # wait at most 1 hr (3600 sec)
```

Let’s put all the pieces together. Our SimpleProcedure can be run in a script by the following.

```
from time import sleep
from pymeasure.experiment import Procedure, Results, Worker
from pymeasure.experiment import IntegerParameter

class SimpleProcedure(Procedure):

    # a Parameter that defines the number of loop iterations
    iterations = IntegerParameter('Loop Iterations')

    # a list defining the order and appearance of columns in our data file
    DATA_COLUMNS = ['Iteration']

    def execute(self):
```

(continues on next page)

(continued from previous page)

```

""" Loops over each iteration and emits the current iteration,
before waiting for 0.01 sec, and then checking if the procedure
should stop
"""

for i in range(self.iterations):
    self.emit('results', {'Iteration': i})
    sleep(0.01)
    if self.should_stop():
        break

if __name__ == "__main__":
    procedure = SimpleProcedure()
    procedure.iterations = 100

    data_filename = 'example.csv'
    results = Results(procedure, data_filename)

    worker = Worker(results)
    worker.start()

    worker.join(timeout=3600) # wait at most 1 hr (3600 sec)

```

Here we have included an if statement to only run the script if the `__name__` is `__main__`. This precaution allows us to import the `SimpleProcedure` object without running the execution.

## Using Logs

Logs keep track of important details in the execution of a procedure. We describe the use of the Python logging module with PyMeasure, which makes it easy to document the execution of a procedure and provides useful insight when diagnosing issues or bugs.

Let's extend our `SimpleProcedure` with logging.

```

import logging
log = logging.getLogger(__name__)
log.addHandler(logging.NullHandler())

from time import sleep
from pymeasure.log import console_log
from pymeasure.experiment import Procedure, Results, Worker
from pymeasure.experiment import IntegerParameter

class SimpleProcedure(Procedure):

    iterations = IntegerParameter('Loop Iterations')

    DATA_COLUMNS = ['Iteration']

    def execute(self):
        log.info("Starting the loop of %d iterations" % self.iterations)
        for i in range(self.iterations):
            data = {'Iteration': i}

```

(continues on next page)

(continued from previous page)

```

        self.emit('results', data)
        log.debug("Emitting results: %s" % data)
        sleep(0.01)
        if self.should_stop():
            log.warning("Caught the stop flag in the procedure")
            break

if __name__ == "__main__":
    console_log(log)

    log.info("Constructing a SimpleProcedure")
    procedure = SimpleProcedure()
    procedure.iterations = 100

    data_filename = 'example.csv'
    log.info("Constructing the Results with a data file: %s" % data_filename)
    results = Results(procedure, data_filename)

    log.info("Constructing the Worker")
    worker = Worker(results)
    worker.start()
    log.info("Started the Worker")

    log.info("Joining with the worker in at most 1 hr")
    worker.join(timeout=3600) # wait at most 1 hr (3600 sec)
    log.info("Finished the measurement")

```

First, we have imported the Python logging module and grabbed the logger using the `__name__` argument. This gives us logging information specific to the current file. Conversely, we could use the `'` argument to get all logs, including those of `pymeasure`. We use the `console_log` function to conveniently output the log to the console. Further details on how to use the logger are addressed in the Python logging documentation.

## Modifying our script

Now that you have a background on how to use the different features of the `Procedure` class, and how they are run, we will revisit our IV characteristic measurement using Procedures. Below we present the modified version of our example script, now as a `IVProcedure` class.

```

# Import necessary packages
from pymeasure.instruments.keithley import Keithley2400
from pymeasure.experiment import Procedure
from pymeasure.experiment import IntegerParameter, FloatParameter
from time import sleep

class IVProcedure(Procedure):

    data_points = IntegerParameter('Data points', default=50)
    averages = IntegerParameter('Averages', default=50)
    max_current = FloatParameter('Maximum Current', units='A', default=0.01)
    min_current = FloatParameter('Minimum Current', units='A', default=-0.01)

```

(continues on next page)

(continued from previous page)

```

DATA_COLUMNS = ['Current (A)', 'Voltage (V)', 'Voltage Std (V)']

def startup(self):
    log.info("Connecting and configuring the instrument")
    self.sourcemeter = Keithley2400("GPIB::4")
    self.sourcemeter.reset()
    self.sourcemeter.use_front_terminals()
    self.sourcemeter.measure_voltage()
    self.sourcemeter.config_current_source()
    sleep(0.1) # wait here to give the instrument time to react
    self.sourcemeter.set_buffer(averages)

def execute(self):
    currents = np.linspace(
        self.min_current,
        self.max_current,
        num=self.data_points
    )

    # Loop through each current point, measure and record the voltage
    for current in currents:
        log.info("Setting the current to %g A" % current)
        self.sourcemeter.current = current
        self.sourcemeter.reset_buffer()
        sleep(0.1)
        self.sourcemeter.start_buffer()
        log.info("Waiting for the buffer to fill with measurements")
        self.sourcemeter.wait_for_buffer()

        self.emit('results', {
            'Current (A)': current,
            'Voltage (V)': self.sourcemeter.means,
            'Voltage Std (V)': self.sourcemeter.standard_devs
        })
        sleep(0.01)
        if self.should_stop():
            log.info("User aborted the procedure")
            break

def shutdown(self):
    self.sourcemeter.shutdown()
    log.info("Finished measuring")

if __name__ == "__main__":
    console_log(log)

    log.info("Constructing an IVProcedure")
    procedure = IVProcedure()
    procedure.data_points = 100
    procedure.averages = 50
    procedure.max_current = -0.01
    procedure.min_current = 0.01

```

(continues on next page)

(continued from previous page)

```

data_filename = 'example.csv'
log.info("Constructing the Results with a data file: %s" % data_filename)
results = Results(procedure, data_filename)

log.info("Constructing the Worker")
worker = Worker(results)
worker.start()
log.info("Started the Worker")

log.info("Joining with the worker in at most 1 hr")
worker.join(timeout=3600) # wait at most 1 hr (3600 sec)
log.info("Finished the measurement")

```

At this point, you are familiar with how to construct a Procedure sub-class. The next section shows how to put these procedures to work in a graphical environment, where will have live-plotting of the data and the ability to easily queue up a number of experiments in sequence. All of these features come from using the Procedure object.

## 3.3 Using a graphical interface

In the previous tutorial we measured the IV characteristic of a sample to show how we can set up a simple experiment in PyMeasure. The real power of PyMeasure comes when we also use the graphical tools that are included to turn our simple example into a full-fledged user interface.

### 3.3.1 Using the Plotter

While it lacks the nice features of the ManagedWindow, the Plotter object is the simplest way of getting live-plotting. The Plotter takes a Results object and plots the data at a regular interval, grabbing the latest data each time from the file.

**Warning:** The example in this section is known to raise issues when executed: a *QApplication was not created in the main thread / nextEventMatchingMask should only be called from the Main Thread* warning is raised. While the example works without issues on some operating systems and python configurations, users are advised not to rely on the plotter while this issue is unresolved. Users can hence skip this example and continue with the *Using the ManagedWindow* section.

Let's extend our SimpleProcedure with a RandomProcedure, which generates random numbers during our loop. This example does not include instruments to provide a simpler example.

```

import logging
log = logging.getLogger(__name__)
log.addHandler(logging.NullHandler())

import random
from time import sleep
from pymeasure.log import console_log
from pymeasure.display import Plotter
from pymeasure.experiment import Procedure, Results, Worker

```

(continues on next page)

(continued from previous page)

```

from pymeasure.experiment import IntegerParameter, FloatParameter, Parameter

class RandomProcedure(Procedure):

    iterations = IntegerParameter('Loop Iterations')
    delay = FloatParameter('Delay Time', units='s', default=0.2)
    seed = Parameter('Random Seed', default='12345')

    DATA_COLUMNS = ['Iteration', 'Random Number']

    def startup(self):
        log.info("Setting the seed of the random number generator")
        random.seed(self.seed)

    def execute(self):
        log.info("Starting the loop of %d iterations" % self.iterations)
        for i in range(self.iterations):
            data = {
                'Iteration': i,
                'Random Number': random.random()
            }
            self.emit('results', data)
            log.debug("Emitting results: %s" % data)
            sleep(self.delay)
            if self.should_stop():
                log.warning("Caught the stop flag in the procedure")
                break

if __name__ == "__main__":
    console_log(log)

    log.info("Constructing a RandomProcedure")
    procedure = RandomProcedure()
    procedure.iterations = 100

    data_filename = 'random.csv'
    log.info("Constructing the Results with a data file: %s" % data_filename)
    results = Results(procedure, data_filename)

    log.info("Constructing the Plotter")
    plotter = Plotter(results)
    plotter.start()
    log.info("Started the Plotter")

    log.info("Constructing the Worker")
    worker = Worker(results)
    worker.start()
    log.info("Started the Worker")

    log.info("Joining with the worker in at most 1 hr")
    worker.join(timeout=3600) # wait at most 1 hr (3600 sec)

```

(continues on next page)


(continued from previous page)

```
log.info("Finished the measurement")
```

The important addition is the construction of the Plotter from the Results object.

```
plotter = Plotter(results)
plotter.start()
```

The Plotter is started in a different process so that it can be run on a separate CPU for higher performance. The Plotter launches a Qt graphical interface using pyqtgraph which allows the Results data to be viewed based on the columns in the data.



### 3.3.2 Using the ManagedWindow

The ManagedWindow is the most convenient tool for running measurements with your Procedure. This has the major advantage of accepting the input parameters graphically. From the parameters, a graphical form is automatically generated that allows the inputs to be typed in. With this feature, measurements can be started dynamically, instead of defined in a script.

Another major feature of the ManagedWindow is its support for running measurements in a sequential queue. This allows you to set up a number of measurements with different input parameters, and watch them unfold on the live-plot. This is especially useful for long running measurements. The ManagedWindow achieves this through the Manager



object, which coordinates which Procedure the Worker should run and keeps track of its status as the Worker progresses.

Below we adapt our previous example to use a ManagedWindow.

```
import logging
log = logging.getLogger(__name__)
log.addHandler(logging.NullHandler())

import sys
import tempfile
import random
from time import sleep
from pymeasure.log import console_log
from pymeasure.display.Qt import QtGui
from pymeasure.display.windows import ManagedWindow
from pymeasure.experiment import Procedure, Results
from pymeasure.experiment import IntegerParameter, FloatParameter, Parameter

class RandomProcedure(Procedure):

    iterations = IntegerParameter('Loop Iterations')
    delay = FloatParameter('Delay Time', units='s', default=0.2)
    seed = Parameter('Random Seed', default='12345')

    DATA_COLUMNS = ['Iteration', 'Random Number']

    def startup(self):
        log.info("Setting the seed of the random number generator")
        random.seed(self.seed)

    def execute(self):
        log.info("Starting the loop of %d iterations" % self.iterations)
        for i in range(self.iterations):
            data = {
                'Iteration': i,
                'Random Number': random.random()
            }
            self.emit('results', data)
            log.debug("Emitting results: %s" % data)
            sleep(self.delay)
            if self.should_stop():
                log.warning("Caught the stop flag in the procedure")
                break

class MainWindow(ManagedWindow):

    def __init__(self):
        super(MainWindow, self).__init__(
            procedure_class=RandomProcedure,
            inputs=['iterations', 'delay', 'seed'],
            displays=['iterations', 'delay', 'seed'],
            x_axis='Iteration',
            y_axis='Random Number'
```

(continues on next page)

(continued from previous page)

```

    )
    self.setWindowTitle('GUI Example')

    def queue(self):
        filename = tempfile.mktemp()


        procedure = self.make_procedure()
        results = Results(procedure, filename)
        experiment = self.new_experiment(results)

        self.manager.queue(experiment)

if __name__ == "__main__":
    app = QtGui.QApplication(sys.argv)
    window = MainWindow()
    window.show()
    sys.exit(app.exec_())


```

This results in the following graphical display.




In the code, the `MainWindow` class is a sub-class of the `ManagedWindow` class. We override the constructor to provide information about the procedure class and its options. The `inputs` are a list of `Parameters` class-variable names, which the display will generate graphical fields for. When the list of inputs is long, a boolean key-word argument `inputs_in_scrollarea` is provided that adds a scrollbar to the input area. The `displays` is a list similar to the `inputs` list, which instead defines the parameters to display in the browser window. This browser keeps track of the experiments being run in the sequential queue.


The `queue` method establishes how the `Procedure` object is constructed. We use the `self.make_procedure` method to create a `Procedure` based on the graphical input fields. Here we are free to modify the procedure before putting it on the queue. In this context, the `Manager` uses an `Experiment` object to keep track of the `Procedure`, `Results`, and its associated graphical representations in the browser and live-graph. This is then given to the `Manager` to queue the experiment.



By default the `Manager` starts a measurement when its procedure is queued. The abort button can be pressed to stop an experiment. In the `Procedure`, the `self.should_stop` call will catch the abort event and halt the measurement. It is important to check this value, or the `Procedure` will not be responsive to the abort event.



If you abort a measurement, the resume button must be pressed to continue the next measurement. This allows you to adjust anything, which is presumably why the abort was needed.



Now that you have learned about the `ManagedWindow`, you have all of the basics to get up and running quickly with a measurement and produce an easy to use graphical interface with PyMeasure.

**Note:** For performance reasons, the default linewidth of all the graphs has been set to 1. If performance is not an issue, the linewidth can be changed to 2 (or any other value) for better visibility by using the `linewidth` keyword-argument in the `Plotter` or the `ManagedWindow`. Whenever a linewidth of 2 is preferred and a better performance is required, it is possible to enable using OpenGL in the import section of the file:

```
import pyqtgraph as pg
pg.setConfigOption("useOpenGL", True)
```

### 3.3.3 Customising the plot options

For both the `PlotterWindow` and `ManagedWindow`, plotting is provided by the `pyqtgraph` library. This library allows you to change various plot options, as you might expect: axis ranges (by default auto-ranging), logarithmic and semilogarithmic axes, downsampling, grid display, FFT display, etc. There are two main ways you can do this:

1. You can right click on the plot to manually change any available options. This is also a good way of getting an overview of what options are available in `pyqtgraph`. Option changes will, of course, not persist across a restart of your program.
2. You can programmatically set these options using `pyqtgraph`'s `PlotItem` API, so that the window will open with these display options already set, as further explained below.

For `Plotter`, you can make a sub-class that overrides the `setup_plot()` method. This method will be called when the `Plotter` constructs the window. As an example

```
class LogPlotter(Plotter):
    def setup_plot(self, plot):
        # use logarithmic x-axis (e.g. for frequency sweeps)
        plot.setLogMode(x=True)
```

For `ManagedWindow`, the mechanism to customize plots is much more flexible by using specialization via inheritance. Indeed `ManagedWindowBase` is the base class for `ManagedWindow` and `ManagedImageWindow` which are subclasses ready to use for GUI.

### 3.3.4 Defining your own ManagedWindow's widgets

The parameter `widget_list` in `ManagedWindowBase` constructor allow to introduce user's defined widget in the GUI results display area. The user's widget should inherit from `TabWidget` and could reimplement any of the methods that needs customization. In order to get familiar with the mechanism, users can check the following widgets already provided:

- `LogWidget`
- `PlotWidget`
- `ImageWidget`

### 3.3.5 Using the sequencer

As an extension to the way of graphically inputting parameters and executing multiple measurements using the `ManagedWindow`, `SequencerWidget` is provided which allows users to queue a series of measurements with varying one, or more, of the parameters. This sequencer thereby provides a convenient way to scan through the parameter space of the measurement procedure.

To activate the sequencer, two additional keyword arguments are added to `ManagedWindow`, namely `sequencer` and `sequencer_inputs`. `sequencer` accepts a boolean stating whether or not the sequencer has to be included into the window and `sequencer_inputs` accepts either `None` or a list of the parameter names are to be scanned over. If no list of parameters is given, the parameters displayed in the manager queue are used.

In order to be able to use the sequencer, the `ManagedWindow` class is required to have a `queue` method which takes a keyword (or better keyword-only for safety reasons) argument `procedure`, where a `procedure` instance can be passed. The sequencer will use this method to queue the parameter scan.

In order to implement the sequencer into the previous example, only the `MainWindow` has to be modified slightly (where modified lines are marked):

```

class MainWindow(ManagedWindow):

    def __init__(self):
        super(MainWindow, self).__init__(
            procedure_class=TestProcedure,
            inputs=['iterations', 'delay', 'seed'],
            displays=['iterations', 'delay', 'seed'],
            x_axis='Iteration',
            y_axis='Random Number',
            sequencer=True,                                # Added line
            sequencer_inputs=['iterations', 'delay', 'seed'], # Added line
            sequence_file="gui_sequencer_example_sequence.txt", # Added line, optional
        )
        self.setWindowTitle('GUI Example')

    def queue(self, procedure=None):                        # Modified line
        filename = tempfile.mktemp()

        if procedure is None:                              # Added line
            procedure = self.make_procedure()              # Indented

        results = Results(procedure, filename)
        experiment = self.new_experiment(results)

        self.manager.queue(experiment)

```

This adds the sequencer underneath the the input panel.

Input Parameters

Loop Iterations:

100

Delay Time:

0.2 s

Random Seed:

12345

Queue

Abort

Sequencer

Level	Parameter	Sequence
0	Delay Time	arange(0.25, 1, 0.25)
1	Random Seed	[1, 4, 8]
2	Loop Iterations	exp(linspace(1, 5, 3))
1	Random Seed	arange(10, 100, 10)

Add root item

Add item

Remove item

Load sequence

Queue sequence

The widget contains a tree-view where you can build the sequence. It has three columns: `level` (indicated how deep an item is nested), `parameter` (a drop-down menu to select which parameter is being sequenced by that item), and `sequence` (the text-box where you can define the sequence). While the two former columns are rather straightforward, filling in the later requires some explanation.

In order to maintain flexibility, the sequence is defined in a text-box, allowing the user to enter any list-generating single-line piece of code. To assist in this, a number of functions is supported, either from the main python library (namely `range`, `sorted`, and `list`) or the numpy library. The supported numpy functions (prepending `numpy.` or any abbreviation is not required) are: `arange`, `linspace`, `arccos`, `arcsin`, `arctan`, `arctan2`, `ceil`, `cos`, `cosh`, `degrees`, `e`, `exp`, `fabs`, `floor`, `fmod`, `frexp`, `hypot`, `ldexp`, `log`, `log10`, `modf`, `pi`, `power`, `radians`, `sin`, `sinh`, `sqrt`, `tan`, and `tanh`.

As an example, `arange(0, 10, 1)` generates a list increasing with steps of 1, while using `exp(arange(0, 10, 1))` generates an exponentially increasing list. This way complex sequences can be entered easily.

The sequences can be extended and shortened using the buttons `Add root item`, `Add item`, and `Remove item`. The later two either add a item as a child of the currently selected item or remove the selected item, respectively. To queue the entered sequence the button `Queue sequence` can be used. If an error occurs in evaluating the sequence text-boxes, this is mentioned in the logger, and nothing is queued.

Finally, it is possible to write a simple text file to quickly load a pre-defined sequence with the `Load sequence` button, such that the user does not need to write the sequence again each time. In the sequence file each line adds one item



to the sequence tree, starting with a number of dashes (-) to indicate the level of the item (starting with 1 dash for top level), followed by the name of the parameter and the sequence string, both as a python string between parentheses. An example of such a sequence file is given below, resulting in the sequence shown in the figure above.

```
- "Delay Time", "arange(0.25, 1, 0.25)"
-- "Random Seed", "[1, 4, 8]"
--- "Loop Iterations", "exp(linspace(1, 5, 3))"
-- "Random Seed", "arange(10, 100, 10)"
```

This file can also be automatically loaded at the start of the program by adding the key-word argument `sequence_file="filename.txt"` to the `super(MainWindow, self).__init__` call, as was done in the example.

### 3.3.6 Using the directory input

It is possible to add a directory input in order to choose where the experiment's result will be saved. This option is activated by passing a boolean key-word argument `directory_input` during the [ManagedWindow](#) init. The value of the directory can be retrieved using the property `directory`.

Only the `MainWindow` needs to be modified in order to use this option (modified lines are marked).

```
class MainWindow(ManagedWindow):


    def __init__(self):
        super(MainWindow, self).__init__(
            procedure_class=TestProcedure,
            inputs=['iterations', 'delay', 'seed'],
            displays=['iterations', 'delay', 'seed'],
            x_axis='Iteration',
            y_axis='Random Number',
            directory_input=True,           # Added line
        )
        self.setWindowTitle('GUI Example')

    def queue(self):
        directory = self.directory
        filename = unique_filename(directory)  # Modified line

        results = Results(procedure, filename)
        experiment = self.new_experiment(results)

        self.manager.queue(experiment)
```

This adds the input line above the Queue and Abort buttons.



A completer is implemented allowing to quickly select an existing folder, and a button on the right side of the input widget opens a browse dialog.

### 3.3.7 Using the estimator widget

In order to provide estimates of the measurement procedure, an *EstimatorWidget* is provided that allows the user to define and calculate estimates. The widget is automatically activated when the `get_estimates` method is added in the Procedure.

The quickest and most simple implementation of the `get_estimates` function simply returns the estimated duration of the measurement in seconds (as an `int` or a `float`). As an example, in the example provided in the *Using the ManagedWindow* section, the Procedure is changed to:

```
class RandomProcedure(Procedure):

    # ...

    def get_estimates(self, sequence_length=None, sequence=None):

        return self.iterations * self.delay
```

This will add the estimator widget at the dock on the left. The duration and finishing-time of a single measurement is always displayed in this case. Depending on whether the *SequencerWidget* is also used, the length, duration and finishing-time of the full sequence is also shown.

For maximum flexibility (e.g. for showing multiple and other types of estimates, such as the duration, filesize, finishing-time, etc.) it is also possible that the `get_estimates` returns a list of tuples. Each of these tuple consists of two strings: the first is the name (label) of the estimate, the second is the estimate itself.

As an example, in the example provided in the *Using the ManagedWindow* section, the Procedure is changed to:

```
class RandomProcedure(Procedure):

    # ...

    def get_estimates(self, sequence_length=None, sequence=None):
```

(continues on next page)

(continued from previous page)

```


duration = self.iterations * self.delay

estimates = [
    ("Duration", "%d s" % int(duration)),
    ("Number of lines", "%d" % int(self.iterations)),
    ("Sequence length", str(sequence_length)),
    ('Measurement finished at', str(datetime.now() +
→timedelta(seconds=duration))),
]

return estimates

```

This will add the estimator widget at the dock on the left.



Note that after the initialisation of the widget both the label of the estimate as of course the estimate itself can be modified, but the amount of estimates is fixed.

The keyword arguments are not required in the implementation of the function, but are passed if asked for (i.e. `def get_estimates(self)` does also works). Keyword arguments that are accepted are `sequence`, which contains the full sequence of the sequencer (if present), and `sequence_length`, which gives the length of the sequence as integer (if present). If the sequencer is not present or the sequence cannot be parsed, both `sequence` and `sequence_length` will contain `None`.

The estimates are automatically updated every 2 seconds. Changing this update interval is possible using the “Update continuously”-checkbox, which can be toggled between three states: off (i.e. no updating), auto-update every two seconds (default) or auto-update every 100 milliseconds. Manually updating the estimates (useful whenever continuous updating is turned off) is also possible using the “update”-button.

### 3.3.8 Flexible hiding of inputs

There can be situations when it may be relevant to turn on or off a number of inputs (e.g. when a part of the measurement script is skipped upon turning of a single `BooleanParameter`). For these cases, it is possible to assign a `Parameter` to a controlling `Parameter`, which will hide or show the `Input` of the `Parameter` depending on the value of the `Parameter`. This is done with the `group_by` key-word argument.

```

toggle = BooleanParameter("toggle", default=True)
param = FloatParameter('some parameter', group_by='toggle')

```

When both the `toggle` and `param` are visible in the `InputsWidget` (via `inputs=['iterations', 'delay', 'seed']` as demonstrated above) one can control whether the input-field of `param` is visible by checking and unchecking the checkbox of `toggle`. By default, the group will be visible if the value of the `group_by` `Parameter` is `True`

(which is only relevant for a `BooleanParameter`), but it is possible to specify other value as conditions using the `group_condition` keyword argument.

```
iterations = IntegerParameter('Loop Iterations', default=100)
param = FloatParameter('some parameter', group_by='iterations', group_condition=99)
```

Here the input of `param` is only visible if `iterations` has a value of 99. This works with any type of `Parameter` as `group_by` parameter.

To allow for even more flexibility, it is also possible to pass a (lambda)function as a condition:

```
iterations = IntegerParameter('Loop Iterations', default=100)
param = FloatParameter('some parameter', group_by='iterations', group_condition=lambda
    ↪ v: 50 < v < 100)
```

Now the input of `param` is only shown if the value of `iterations` is between 51 and 99.

Using the `hide_groups` keyword-argument of the `ManagedWindow` you can choose between hiding the groups (`hide_groups = True`) and disabling / graying-out the groups (`hide_groups = False`).

Finally, it is also possible to provide multiple parameters to the `group_by` argument, in which case the input will only be visible if all of the conditions are true. Multiple parameters for grouping can either be passed as a dict of string: condition pairs, or as a list of strings, in which case the `group_condition` can be either a single condition or a list of conditions:

```
iterations = IntegerParameter('Loop Iterations', default=100)
toggle = BooleanParameter('A checkbox')
param_A = FloatParameter('some parameter', group_by=['iterations', 'toggle'], group_
    ↪ condition=[lambda v: 50 < v < 100, True])
param_B = FloatParameter('some parameter', group_by={'iterations': lambda v: 50 < v <
    ↪ 100, 'toggle': True})
```

Note that in this example, `param_A` and `param_B` are identically grouped: they're only visible if `iterations` is between 51 and 99 and if the `toggle` checkbox is checked (i.e. `True`).

## PYMEASURE.ADAPTERS

The adapter classes allow the instruments to be independent of the communication method used.

Adapters for specific instruments should be grouped in an `adapters.py` file in the corresponding manufacturer's folder of `pymeasure.instruments`. For example, the adapter for communicating with LakeShore instruments over USB, [\*LakeShoreUSBAdapter\*](#), is found in `pymeasure.instruments.lakeshore.adapters`.

### 4.1 Adapter base class

**class** `pymeasure.adapters.Adapter`(*preprocess\_reply=None, \*\*kwargs*)

Base class for Adapter child classes, which adapt between the Instrument object and the connection, to allow flexible use of different connection techniques.

This class should only be inherited from.

#### Parameters

- **preprocess\_reply** – optional callable used to preprocess strings received from the instrument. The callable returns the processed string.
- **kwargs** – all other keyword arguments are ignored.

**ask**(*command*)

Writes the command to the instrument and returns the resulting ASCII response

**Parameters** **command** – SCPI command string to be sent to the instrument

**Returns** String ASCII response of the instrument

**binary\_values**(*command, header\_bytes=0, dtype=<class 'numpy.float32'>*)

Returns a numpy array from a query for binary data

#### Parameters

- **command** – SCPI command to be sent to the instrument
- **header\_bytes** – Integer number of bytes to ignore in header
- **dtype** – The NumPy data type to format the values with

**Returns** NumPy array of values

**read**()

Reads until the buffer is empty and returns the resulting ASCII response

**Returns** String ASCII response of the instrument.

**values**(*command, separator=', ', cast=<class 'float'>, preprocess\_reply=None*)

Writes a command to the instrument and returns a list of formatted values from the result

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **separator** – A separator character to split the string into a list
- **cast** – A type to cast the result
- **preprocess\_reply** – optional callable used to preprocess values received from the instrument. The callable returns the processed string. If not specified, the Adapter default is used if available, otherwise no preprocessing is done.

**Returns** A list of the desired type, or strings where the casting fails

**write**(*command*)

Writes a command to the instrument

**Parameters** **command** – SCPI command string to be sent to the instrument

## 4.2 Fake adapter

**class** `pymeasure.adapters.FakeAdapter`(*preprocess\_reply=None, \*\*kwargs*)

Bases: `pymeasure.adapters.adapter.Adapter`

Provides a fake adapter for debugging purposes, which bounces back the command so that arbitrary values testing is possible.

```
a = FakeAdapter()
assert a.read() == ""
a.write("5")
assert a.read() == "5"
assert a.read() == ""
assert a.ask("10") == "10"
assert a.values("10") == [10]
```

**ask**(*command*)

Writes the command to the instrument and returns the resulting ASCII response

**Parameters** **command** – SCPI command string to be sent to the instrument

**Returns** String ASCII response of the instrument

**binary\_values**(*command, header\_bytes=0, dtype=<class 'numpy.float32'>*)

Returns a numpy array from a query for binary data

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **header\_bytes** – Integer number of bytes to ignore in header
- **dtype** – The NumPy data type to format the values with

**Returns** NumPy array of values

**read**()

Returns the last commands given after the last read call.

**values**(*command, separator=', ', cast=<class 'float'>, preprocess\_reply=None*)

Writes a command to the instrument and returns a list of formatted values from the result

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **separator** – A separator character to split the string into a list
- **cast** – A type to cast the result
- **preprocess\_reply** – optional callable used to preprocess values received from the instrument. The callable returns the processed string. If not specified, the Adapter default is used if available, otherwise no preprocessing is done.

**Returns** A list of the desired type, or strings where the casting fails

**write**(*command*)

Writes the command to a buffer, so that it can be read back.

## 4.3 Serial adapter

**class** `pymeasure.adapters.SerialAdapter`(*port*, *preprocess\_reply=None*, *\*\*kwargs*)

Bases: `pymeasure.adapters.adapter.Adapter`

Adapter class for using the Python Serial package to allow serial communication to instrument

**Parameters**

- **port** – Serial port
- **preprocess\_reply** – optional callable used to preprocess strings received from the instrument. The callable returns the processed string.
- **kwargs** – Any valid key-word argument for `serial.Serial`

**\_format\_binary\_values**(*values*, *datatype='f'*, *is\_big\_endian=False*, *header\_fmt='ieee'*)

Format values in binary format, used internally in `write_binary_values()`.

**Parameters**

- **values** – data to be written to the device.
- **datatype** – the format string for a single element. See `struct` module.
- **is\_big\_endian** – boolean indicating endianness.
- **header\_fmt** – Format of the header prefixing the data (“ieee”, “hp”, “empty”).

**Returns** binary string.

**Return type** bytes

**ask**(*command*)

Writes the command to the instrument and returns the resulting ASCII response

**Parameters** **command** – SCPI command string to be sent to the instrument

**Returns** String ASCII response of the instrument

**binary\_values**(*command*, *header\_bytes=0*, *dtype=<class 'numpy.float32'>*)

Returns a numpy array from a query for binary data

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **header\_bytes** – Integer number of bytes to ignore in header
- **dtype** – The NumPy data type to format the values with

**Returns** NumPy array of values

**read()**

Reads until the buffer is empty and returns the resulting ASCII response

**Returns** String ASCII response of the instrument.

**values**(*command*, *separator*=' ', *cast*=<class 'float'>, *preprocess\_reply*=None)

Writes a command to the instrument and returns a list of formatted values from the result

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **separator** – A separator character to split the string into a list
- **cast** – A type to cast the result
- **preprocess\_reply** – optional callable used to preprocess values received from the instrument. The callable returns the processed string. If not specified, the Adapter default is used if available, otherwise no preprocessing is done.

**Returns** A list of the desired type, or strings where the casting fails

**write**(*command*)

Writes a command to the instrument

**Parameters** **command** – SCPI command string to be sent to the instrument

**write\_binary\_values**(*command*, *values*, *\*\*kwargs*)

Write binary data to the instrument, e.g. waveform for signal generators

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **values** – iterable representing the binary values
- **kwargs** – Key-word arguments to pass onto `_format_binary_values()`

**Returns** number of bytes written

## 4.4 Prologix adapter

**class** `pymeasure.adapters.PrologixAdapter`(*port*, *address*=None, *rw\_delay*=None, *serial\_timeout*=0.5, *preprocess\_reply*=None, *\*\*kwargs*)

Bases: `pymeasure.adapters.serial.SerialAdapter`

Encapsulates the additional commands necessary to communicate over a Prologix GPIB-USB Adapter, using the `SerialAdapter`.

Each `PrologixAdapter` is constructed based on a serial port or connection and the GPIB address to be communicated to. Serial connection sharing is achieved by using the `gpiplib()` method to spawn new `PrologixAdapters` for different GPIB addresses.

**Parameters**

- **port** – The Serial port name or a `serial.Serial` object
- **address** – Integer GPIB address of the desired instrument
- **rw\_delay** – An optional delay to set between a write and read call for slow to respond instruments.



- **preprocess\_reply** – optional callable used to preprocess strings received from the instrument. The callable returns the processed string.
- **kwargs** – Key-word arguments if constructing a new serial object

**Variables** **address** – Integer GPIB address of the desired instrument

To allow user access to the Prologix adapter in Linux, create the file: `/etc/udev/rules.d/51-prologix.rules`, with contents:

```
SUBSYSTEMS=="usb",ATTRS{idVendor}=="0403",ATTRS{idProduct}=="6001",MODE="0666"
```

Then reload the udev rules with:

```
sudo udevadm control --reload-rules
sudo udevadm trigger
```

**\_format\_binary\_values**(*values*, *datatype='f'*, *is\_big\_endian=False*, *header\_fmt='ieee'*)

Format values in binary format, used internally in [write\\_binary\\_values\(\)](#).

#### Parameters

- **values** – data to be written to the device.
- **datatype** – the format string for a single element. See struct module.
- **is\_big\_endian** – boolean indicating endianness.
- **header\_fmt** – Format of the header prefixing the data (“ieee”, “hp”, “empty”).

**Returns** binary string.

**Return type** bytes

**ask**(*command*)

Ask the Prologix controller, include a forced delay for some instruments.

**Parameters** **command** – SCPI command string to be sent to instrument

**binary\_values**(*command*, *header\_bytes=0*, *dtype=<class 'numpy.float32'>*)

Returns a numpy array from a query for binary data

#### Parameters

- **command** – SCPI command to be sent to the instrument
- **header\_bytes** – Integer number of bytes to ignore in header
- **dtype** – The NumPy data type to format the values with

**Returns** NumPy array of values

**gpiib**(*address*, *rw\_delay=None*)

Returns and PrologixAdapter object that references the GPIB address specified, while sharing the Serial connection with other calls of this function

#### Parameters

- **address** – Integer GPIB address of the desired instrument
- **rw\_delay** – Set a custom Read/Write delay for the instrument

**Returns** PrologixAdapter for specific GPIB address

**read**()

Reads the response of the instrument until timeout

**Returns** String ASCII response of the instrument

**set\_defaults()**

Sets up the default behavior of the Prologix-GPIB adapter

**values**(*command*, *separator*=' ', *cast*=<class 'float'>, *preprocess\_reply*=None)

Writes a command to the instrument and returns a list of formatted values from the result

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **separator** – A separator character to split the string into a list
- **cast** – A type to cast the result
- **preprocess\_reply** – optional callable used to preprocess values received from the instrument. The callable returns the processed string. If not specified, the Adapter default is used if available, otherwise no preprocessing is done.

**Returns** A list of the desired type, or strings where the casting fails

**wait\_for\_srq**(*timeout*=25, *delay*=0.1)

Blocks until a SRQ, and leaves the bit high

**Parameters**

- **timeout** – Timeout duration in seconds
- **delay** – Time delay between checking SRQ in seconds

**write**(*command*)

Writes the command to the GPIB address stored in the [address](#)

**Parameters** **command** – SCPI command string to be sent to the instrument

**write\_binary\_values**(*command*, *values*, *\*\*kwargs*)

Write binary data to the instrument, e.g. waveform for signal generators.

values are encoded in a binary format according to IEEE 488.2 Definite Length Arbitrary Block Response Data block.

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **values** – iterable representing the binary values
- **kwargs** – Key-word arguments to pass onto [\\_format\\_binary\\_values\(\)](#)

**Returns** number of bytes written

## 4.5 VISA adapter

**class** `pymeasure.adapters.VISAAdapter`(*resource\_name*, *visa\_library*="", *preprocess\_reply*=None, *\*\*kwargs*)

Bases: [pymeasure.adapters.adapter.Adapter](#)

Adapter class for the VISA library using PyVISA to communicate with instruments.

**Parameters**

- **resource** – VISA resource name that identifies the address
- **visa\_library** – VisaLibrary Instance, path of the VISA library or VisaLibrary spec string (@py or @ni). if not given, the default for the platform will be used.

- **preprocess\_reply** – optional callable used to preprocess strings received from the instrument. The callable returns the processed string.
- **kwargs** – Any valid key-word arguments for constructing a PyVISA instrument

**ask**(*command*)

Writes the command to the instrument and returns the resulting ASCII response

**Parameters** **command** – SCPI command string to be sent to the instrument

**Returns** String ASCII response of the instrument

**ask\_values**(*command*, *\*\*kwargs*)

Writes a command to the instrument and returns a list of formatted values from the result. This leverages the *query\_ascii\_values* method in PyVISA.

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **kwargs** – Key-word arguments to pass onto *query\_ascii\_values*

**Returns** Formatted response of the instrument.

**binary\_values**(*command*, *header\_bytes=0*, *dtype=<class 'numpy.float32'>*)

Returns a numpy array from a query for binary data

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **header\_bytes** – Integer number of bytes to ignore in header
- **dtype** – The NumPy data type to format the values with

**Returns** NumPy array of values

**static has\_supported\_version()**

Returns True if the PyVISA version is greater than 1.8

**read()**

Reads until the buffer is empty and returns the resulting ASCII response

**Returns** String ASCII response of the instrument.

**read\_bytes**(*size*)

Reads specified number of bytes from the buffer and returns the resulting ASCII response

**Parameters** **size** – Number of bytes to read from the buffer

**Returns** String ASCII response of the instrument.

**values**(*command*, *separator=''*, *cast=<class 'float'>*, *preprocess\_reply=None*)

Writes a command to the instrument and returns a list of formatted values from the result

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **separator** – A separator character to split the string into a list
- **cast** – A type to cast the result
- **preprocess\_reply** – optional callable used to preprocess values received from the instrument. The callable returns the processed string. If not specified, the Adapter default is used if available, otherwise no preprocessing is done.

**Returns** A list of the desired type, or strings where the casting fails

**wait\_for\_srq**(*timeout=25, delay=0.1*)

Blocks until a SRQ, and leaves the bit high

**Parameters**

- **timeout** – Timeout duration in seconds
- **delay** – Time delay between checking SRQ in seconds

**write**(*command*)

Writes a command to the instrument

**Parameters** **command** – SCPI command string to be sent to the instrument

**write\_binary\_values**(*command, values, \*\*kwargs*)

Write binary data to the instrument, e.g. waveform for signal generators

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **values** – iterable representing the binary values
- **kwargs** – Key-word arguments to pass onto *write\_binary\_values*

**Returns** number of bytes written

## 4.6 VXI-11 adapter

**class** `pymeasure.adapters.VXI11Adapter`(*host, preprocess\_reply=None, \*\*kwargs*)

Bases: `pymeasure.adapters.adapter.Adapter`

**VXI11 Adapter class. Provides a adapter object that** wraps around the read, write and ask functionality of the vx11 library.

**Parameters**

- **host** – string containing the visa connection information.
- **preprocess\_reply** – optional callable used to preprocess strings received from the instrument. The callable returns the processed string.

**ask**(*command*)

Wrapper function for the ask command using the vx11 interface.

**Parameters** **command** – string with the command that will be transmitted to the instrument.

:returns string containing a response from the device.

**ask\_raw**(*command*)

Wrapper function for the ask\_raw command using the vx11 interface.

**Parameters** **command** – binary string with the command that will be transmitted to the instrument

:returns binary string containing the response from the device.

**binary\_values**(*command, header\_bytes=0, dtype=<class 'numpy.float32'>*)

Returns a numpy array from a query for binary data

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **header\_bytes** – Integer number of bytes to ignore in header

- **dtype** – The NumPy data type to format the values with

**Returns** NumPy array of values

**read()**

Wrapper function for the read command using the vx11 interface.

:return string containing a response from the device.

**read\_raw()**

Wrapper function for the read\_raw command using the vx11 interface.

:returns binary string containing the response from the device.

**values**(*command*, *separator*=' ', *cast*=<class 'float'>, *preprocess\_reply*=None)

Writes a command to the instrument and returns a list of formatted values from the result

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **separator** – A separator character to split the string into a list
- **cast** – A type to cast the result
- **preprocess\_reply** – optional callable used to preprocess values received from the instrument. The callable returns the processed string. If not specified, the Adapter default is used if available, otherwise no preprocessing is done.

**Returns** A list of the desired type, or strings where the casting fails

**write**(*command*)

Wrapper function for the write command using the vx11 interface.

**Parameters** **command** – string with command the that will be transmitted to the instrument.

**write\_raw**(*command*)

Wrapper function for the write\_raw command using the vx11 interface.

**Parameters** **command** – binary string with the command that will be transmitted to the instrument

## 4.7 Telnet adapter

**class** `pymeasure.adapters.TelnetAdapter`(*host*, *port*=0, *query\_delay*=0, *preprocess\_reply*=None, *\*\*kwargs*)

Bases: `pymeasure.adapters.adapter.Adapter`

Adapter class for using the Python telnetlib package to allow communication to instruments

**Parameters**

- **host** – host address of the instrument
- **port** – TCPIP port
- **query\_delay** – delay in seconds between write and read in the ask method
- **preprocess\_reply** – optional callable used to preprocess strings received from the instrument. The callable returns the processed string.
- **kwargs** – Valid keyword arguments for telnetlib.Telnet, currently this is only 'timeout'

**ask**(*command*)

Writes a command to the instrument and returns the resulting ASCII response

**Parameters** **command** – command string to be sent to the instrument

**Returns** String ASCII response of the instrument

**binary\_values**(*command*, *header\_bytes*=0, *dtype*=<class 'numpy.float32'>)

Returns a numpy array from a query for binary data

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **header\_bytes** – Integer number of bytes to ignore in header
- **dtype** – The NumPy data type to format the values with

**Returns** NumPy array of values

**read()**

Read something even with blocking the I/O. After something is received check again to obtain a full reply.

**Returns** String ASCII response of the instrument.

**values**(*command*, *separator*=' ', *cast*=<class 'float'>, *preprocess\_reply*=None)

Writes a command to the instrument and returns a list of formatted values from the result

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **separator** – A separator character to split the string into a list
- **cast** – A type to cast the result
- **preprocess\_reply** – optional callable used to preprocess values received from the instrument. The callable returns the processed string. If not specified, the Adapter default is used if available, otherwise no preprocessing is done.

**Returns** A list of the desired type, or strings where the casting fails

**write**(*command*)

Writes a command to the instrument

**Parameters** **command** – command string to be sent to the instrument

## PYMEASURE.EXPERIMENT

This section contains specific documentation on the classes and methods of the package.

### 5.1 Experiment class

The Experiment class is intended for use in the Jupyter notebook environment.

**class** `pymasure.experiment.experiment.Experiment`(*title, procedure, analyse=<function Experiment.<lambda>>>*)

Bases: object

Class which starts logging and creates/runs the results and worker processes.

```
procedure = Procedure()
experiment = Experiment(title, procedure)
experiment.start()
experiment.plot_live('x', 'y', style='.-')
```

**for** a multi-subplot graph:

```
import pylab as pl
ax1 = pl.subplot(121)
experiment.plot('x', 'y', ax=ax1)
ax2 = pl.subplot(122)
experiment.plot('x', 'z', ax=ax2)
experiment.plot_live()
```

**Variables value** – The value of the parameter

#### Parameters

- **title** – The experiment title
- **procedure** – The procedure object
- **analyse** – Post-analysis function, which takes a pandas dataframe as input and returns it with added (analysed) columns. The analysed results are accessible via `experiment.data`, as opposed to `experiment.results.data` for the 'raw' data.
- **\_data\_timeout** – Time limit for how long live plotting should wait for datapoints.

**clear\_plot()**

Clear the figures and plot lists.

**property data**

Data property which returns analysed data, if an analyse function is defined, otherwise returns the raw data.

**pcolor**(*xname, yname, zname, \*args, \*\*kwargs*)

Plot the results from the experiment.data pandas dataframe in a pcolor graph. Store the plots in a plots list attribute.

**plot**(*\*args, \*\*kwargs*)

Plot the results from the experiment.data pandas dataframe. Store the plots in a plots list attribute.

**plot\_live**(*\*args, \*\*kwargs*)

Live plotting loop for jupyter notebook, which automatically updates (an) in-line matplotlib graph(s). Will create a new plot as specified by input arguments, or will update (an) existing plot(s).

**start**()

Start the worker

**update\_line**(*ax, hl, xname, yname*)

Update a line in a matplotlib graph with new data.

**update\_pcolor**(*ax, xname, yname, zname*)

Update a pcolor graph with new data.

**update\_plot**()

Update the plots in the plots list with new data from the experiment.data pandas dataframe.

**wait\_for\_data**()

Wait for the data attribute to fill with datapoints.

**pymeasure.experiment.experiment.create\_filename**(*title*)

Create a new filename according to the style defined in the config file. If no config is specified, create a temporary file.

**pymeasure.experiment.experiment.get\_array**(*start, stop, step*)

Returns a numpy array from start to stop

**pymeasure.experiment.experiment.get\_array\_steps**(*start, stop, numsteps*)

Returns a numpy array from start to stop in numsteps

**pymeasure.experiment.experiment.get\_array\_zero**(*maxval, step*)

Returns a numpy array from 0 to maxval to -maxval to 0

## 5.2 Listener class

**class** **pymeasure.experiment.listeners.Listener**(*port, topic="", timeout=0.01*)

Bases: **pymeasure.thread.StoppableThread**

Base class for Threads that need to listen for messages on a ZMQ TCP port and can be stopped by a thread-safe method call

**message\_waiting**()

Check if we have a message, wait at most until timeout.

**receive**(*flags=0*)**class** **pymeasure.experiment.listeners.Monitor**(*results, queue*)

Bases: **pymeasure.log.QueueListener**

**class** **pymeasure.experiment.listeners.Recorder**(*results, queue, \*\*kwargs*)

Bases: **pymeasure.log.QueueListener**



Recorder loads the initial Results for a filepath and appends data by listening for it over a queue. The queue ensures that no data is lost between the Recorder and Worker.

#### **stop()**

Stop the listener.

This asks the thread to terminate, and then waits for it to do so. Note that if you don't call this before your application exits, there may be some records still left on the queue, which won't be processed.

## 5.3 Procedure class

**class** `pymeasure.experiment.procedure.Procedure(**kwargs)`

Provides the base class of a procedure to organize the experiment execution. Procedures should be run by Workers to ensure that asynchronous execution is properly managed.

```
procedure = Procedure()
results = Results(procedure, data_filename)
worker = Worker(results, port)
worker.start()
```

Inheriting classes should define the startup, execute, and shutdown methods as needed. The shutdown method is called even with a software exception or abort event during the execute method.

If keyword arguments are provided, they are added to the object as attributes.

#### **check\_parameters()**

Raises an exception if any parameter is missing before calling the associated function. Ensures that each value can be set and got, which should cast it into the right format. Used as a decorator `@check_parameters` on the startup method

#### **execute()**

Performs the commands needed for the measurement itself. During execution the shutdown method will always be run following this method. This includes when Exceptions are raised.

#### **gen\_measurement()**

Create MEASURE and DATA\_COLUMNS variables for get\_datapoint method.

#### **get\_estimates()**

Function that returns estimates that are to be displayed by the EstimatorWidget. Must be reimplemented by subclasses. Should return an int or float representing the duration in seconds, or a list with a tuple for each estimate. The tuple should consists of two strings: the first will be used as the label of the estimate, the second as the displayed estimate.

#### **parameter\_objects()**

Returns a dictionary of all the Parameter objects and grabs any current values that are not in the default definitions

#### **parameter\_values()**

Returns a dictionary of all the Parameter values and grabs any current values that are not in the default definitions

#### **parameters\_are\_set()**

Returns True if all parameters are set

#### **refresh\_parameters()**

Enforces that all the parameters are re-cast and updated in the meta dictionary

**set\_parameters**(*parameters*, *except\_missing=True*)

Sets a dictionary of parameters and raises an exception if additional parameters are present if *except\_missing* is True

**shutdown**()

Executes the commands necessary to shut down the instruments and leave them in a safe state. This method is always run at the end.

**startup**()

Executes the commands needed at the start-up of the measurement

**class** `pymeasure.experiment.procedure.UnknownProcedure`(*parameters*)

Handles the case when a *Procedure* object can not be imported during loading in the *Results* class

**startup**()

Executes the commands needed at the start-up of the measurement

## 5.4 Parameter classes

The parameter classes are used to define input variables for a *Procedure*. They each inherit from the *Parameter* base class.

**class** `pymeasure.experiment.parameters.BooleanParameter`(*name*, *default=None*, *ui\_class=None*,  
*group\_by=None*, *group\_condition=True*)

*Parameter* sub-class that uses the boolean type to store the value.

**Variables** *value* – The boolean value of the parameter

**Parameters**

- **name** – The parameter name
- **default** – The default boolean value
- **ui\_class** – A Qt class to use for the UI of this parameter

**class** `pymeasure.experiment.parameters.FloatParameter`(*name*, *units=None*, *minimum=- 1000000000.0*,  
*maximum=1000000000.0*, *decimals=15*,  
*\*\*kwargs*)

*Parameter* sub-class that uses the floating point type to store the value.

**Variables** *value* – The floating point value of the parameter

**Parameters**

- **name** – The parameter name
- **units** – The units of measure for the parameter
- **minimum** – The minimum allowed value (default: -1e9)
- **maximum** – The maximum allowed value (default: 1e9)
- **decimals** – The number of decimals considered (default: 15)
- **default** – The default floating point value
- **ui\_class** – A Qt class to use for the UI of this parameter

**class** `pymeasure.experiment.parameters.IntegerParameter`(*name*, *units=None*, *minimum=-*  
*1000000000.0*, *maximum=1000000000.0*,  
*\*\*kwargs*)

*Parameter* sub-class that uses the integer type to store the value.

**Variables** **value** – The integer value of the parameter

#### Parameters

- **name** – The parameter name
- **units** – The units of measure for the parameter
- **minimum** – The minimum allowed value (default: -1e9)
- **maximum** – The maximum allowed value (default: 1e9)
- **default** – The default integer value
- **ui\_class** – A Qt class to use for the UI of this parameter

**class** `pymeasure.experiment.parameters.ListParameter(name, choices=None, units=None, **kwargs)`  
*Parameter* sub-class that stores the value as a list. String representation of choices must be unique.

#### Parameters

- **name** – The parameter name
- **choices** – An explicit list of choices, which is disregarded if None
- **units** – The units of measure for the parameter
- **default** – The default value
- **ui\_class** – A Qt class to use for the UI of this parameter

#### property choices

Returns an immutable iterable of choices, or None if not set.

**class** `pymeasure.experiment.parameters.Measurable(name, fget=None, units=None, measure=True, default=None, **kwargs)`

Encapsulates the information for a measurable experiment parameter with information about the name, fget function and units if supplied. The value property is called when the procedure retrieves a datapoint and calls the fget function. If no fget function is specified, the value property will return the latest set value of the parameter (or default if never set).

**Variables** **value** – The value of the parameter

#### Parameters

- **name** – The parameter name
- **fget** – The parameter fget function (e.g. an instrument parameter)
- **default** – The default value

**class** `pymeasure.experiment.parameters.Parameter(name, default=None, ui_class=None, group_by=None, group_condition=True)`

Encapsulates the information for an experiment parameter with information about the name, and units if supplied.

**Variables** **value** – The value of the parameter

#### Parameters

- **name** – The parameter name
- **default** – The default value
- **ui\_class** – A Qt class to use for the UI of this parameter
- **group\_by** – Defines the Parameter(s) that controls the visibility of the associated input; can be a string containing the Parameter name, a list of strings with multiple Parameter names, or a dict containing {"Parameter name": condition} pairs.

- **group\_condition** – The condition for the group\_by Parameter that controls the visibility of this parameter, provided as a value or a (lambda)function. If the group\_by argument is provided as a list of strings, this argument can be either a single condition or a list of conditions. If the group\_by argument is provided as a dict this argument is ignored.

**is\_set()**

Returns True if the Parameter value is set

**class** pymeasure.experiment.parameters.**PhysicalParameter**(*name, uncertaintyType='absolute', \*\*kwargs*)

*VectorParameter* sub-class of 2 dimensions to store a value and its uncertainty.

**Variables** **value** – The value of the parameter as a list of 2 floating point numbers

**Parameters**

- **name** – The parameter name
- **uncertainty\_type** – Type of uncertainty, 'absolute', 'relative' or 'percentage'
- **units** – The units of measure for the parameter
- **default** – The default value
- **ui\_class** – A Qt class to use for the UI of this parameter

**class** pymeasure.experiment.parameters.**VectorParameter**(*name, length=3, units=None, \*\*kwargs*)

*Parameter* sub-class that stores the value in a vector format.

**Variables** **value** – The value of the parameter as a list of floating point numbers

**Parameters**

- **name** – The parameter name
- **length** – The integer dimensions of the vector
- **units** – The units of measure for the parameter
- **default** – The default value
- **ui\_class** – A Qt class to use for the UI of this parameter

## 5.5 Worker class

**class** pymeasure.experiment.workers.**Worker**(*results, log\_queue=None, log\_level=20, port=None*)

Bases: pymeasure.thread.StoppableThread

Worker runs the procedure and emits information about the procedure and its status over a ZMQ TCP port. In a child thread, a Recorder is run to write the results to

**emit**(*topic, record*)

Emits data of some topic over TCP

**handle\_abort()**

**handle\_error()**

**join**(*timeout=0*)

Joins the current thread and forces it to stop after the timeout if necessary

**Parameters** **timeout** – Timeout duration in seconds

**run()**

Method representing the thread's activity.

You may override this method in a subclass. The standard `run()` method invokes the callable object passed to the object's constructor as the target argument, if any, with sequential and keyword arguments taken from the `args` and `kwargs` arguments, respectively.

**shutdown()****update\_status**(*status*)

## 5.6 Results class

**class** `pymeasure.experiment.results.CSVFormatter`(*columns*, *delimiter*=',')

Formatter of data results

**format**(*record*)

Formats a record as csv.

**Parameters** **record** (*dict*) – record to format.**Returns** a string**class** `pymeasure.experiment.results.Results`(*procedure*, *data\_filename*)

The Results class provides a convenient interface to reading and writing data in connection with a [Procedure](#) object.

### Variables

- **COMMENT** – The character used to identify a comment (default: #)
- **DELIMITER** – The character used to delimit the data (default: ,)
- **LINE\_BREAK** – The character used for line breaks (default n)
- **CHUNK\_SIZE** – The length of the data chunk that is read

### Parameters

- **procedure** – Procedure object
- **data\_filename** – The data filename where the data is or should be stored

**format**(*data*)

Returns a formatted string containing the data to be written to a file

**header**()

Returns a text header to accompany a datafile so that the procedure can be reconstructed

**labels**()

Returns the columns labels as a string to be written to the file

**static load**(*data\_filename*, *procedure\_class*=None)

Returns a Results object with the associated Procedure object and data

**parse**(*line*)

Returns a dictionary containing the data from the line

**static parse\_header**(*header*, *procedure\_class*=None)

Returns a Procedure object with the parameters as defined in the header text.

**reload**()

Performs a full reloading of the file data, neglecting any changes in the comments

```
pymethods.experiment.results.replace_placeholders(string, procedure, date_format='%Y-%m-%d',  
                                                  time_format='%H:%M:%S')
```

Replace placeholders in string with values from procedure parameters.

Replaces the placeholders in the provided string with the values of the associated parameters, as provided by the procedure. This uses the standard python string.format syntax. Apart from the parameter in the procedure (which should be called by their full names) “date” and “time” are also added as optional placeholders.

#### Parameters

- **string** – The string in which the placeholders are to be replaced. Python string.format syntax is used, e.g. “{Parameter Name}” to insert a FloatParameter called “Parameter Name”, or “{Parameter Name:.2f}” to also specifically format the parameter.
- **procedure** – The procedure from which to get the parameter values.
- **date\_format** – A string to represent how the additional placeholder “date” will be formatted.
- **time\_format** – A string to represent how the additional placeholder “time” will be formatted.

```
pymethods.experiment.results.unique_filename(directory, prefix='DATA', suffix='', ext='csv',  
                                             dated_folder=False, index=True,  
                                             datetimeformat='%Y-%m-%d', procedure=None)
```

Returns a unique filename based on the directory and prefix

## PYMEASURE.DISPLAY

This section contains specific documentation on the classes and methods of the package.

### 6.1 Browser classes

**class** `pymeasure.display.browser.Browser(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtGui`.

Graphical list view of [Experiment](#) objects allowing the user to view the status of queued Experiments as well as loading and displaying data from previous runs.

In order that different Experiments be displayed within the same Browser, they must have entries in `DATA_COLUMNS` corresponding to the *measured\_quantities* of the Browser.

**add**(*experiment*)

Add a [Experiment](#) object to the Browser. This function checks to make sure that the Experiment measures the appropriate quantities to warrant its inclusion, and then adds a `BrowserItem` to the Browser, filling all relevant columns with Parameter data.

**class** `pymeasure.display.browser.BrowserItem(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtGui`.

Represent a row in the [Browser](#) tree widget

### 6.2 Curves classes

**class** `pymeasure.display.curves.BufferCurve(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph`.

Creates a curve based on a predefined buffer size and allows data to be added dynamically, in addition to supporting error bars

**append**(*x*, *y*, *xError=None*, *yError=None*)

Appends data to the curve with optional errors

**prepare**(*size*, *dtype=<class 'numpy.float32'>*)

Prepares the buffer based on its size, data type

**class** `pymeasure.display.curves.Crosshairs(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtCore`.

Attaches crosshairs to the a plot and provides a signal with the x and y graph coordinates

**mouseMoved**(*event=None*)

Updates the mouse position upon mouse movement

**update**()

Updates the mouse position based on the data in the plot. For dynamic plots, this is called each time the data changes to ensure the x and y values correspond to those on the display.

**class** `pymeasure.display.curves.ResultsCurve`(\*args: Any, \*\*kwargs: Any)

Bases: `pyqtgraph`.

Creates a curve loaded dynamically from a file through the Results object and supports error bars. The data can be forced to fully reload on each update, useful for cases when the data is changing across the full file instead of just appending.

**update\_data**()

Updates the data by polling the results

**class** `pymeasure.display.curves.ResultsImage`(\*args: Any, \*\*kwargs: Any)

Bases: `pyqtgraph`.

Creates an image loaded dynamically from a file through the Results object.

**find\_img\_index**(*x, y*)

Finds the integer image indices corresponding to the closest x and y points of the data given some x and y data.

**round\_up**(*x*)

Convenience function since numpy rounds to even

## 6.3 Inputs classes

**class** `pymeasure.display.inputs.BooleanInput`(\*args: Any, \*\*kwargs: Any)

Bases: `pymeasure.display.inputs.Input`, `pyqtgraph.Qt.QtGui`.

Checkbox for boolean values, connected to a `BooleanParameter`.

**set\_parameter**(*parameter*)

Connects a new parameter to the input box, and initializes the box value.

**Parameters** *parameter* – parameter to connect.

**class** `pymeasure.display.inputs.FloatInput`(\*args: Any, \*\*kwargs: Any)

Bases: `pymeasure.display.inputs.Input`, `pyqtgraph.Qt.QtGui`.

Spin input box for floating-point values, connected to a `FloatParameter`.

**See also:**

**Class** `ScientificInput` For inputs in scientific notation.

**set\_parameter**(*parameter*)

Connects a new parameter to the input box, and initializes the box value.

**Parameters** *parameter* – parameter to connect.

**class** `pymeasure.display.inputs.Input`(*parameter*, \*\*kwargs)

Bases: `object`

Mix-in class that connects a `Parameter` object to a GUI input box.

**Parameters** *parameter* – The parameter to connect to this input box.



**Attr parameter** Read-only property to access the associated parameter.

**property parameter**

The connected parameter object. Read-only property; see [set\\_parameter\(\)](#).

Note that reading this property will have the side-effect of updating its value from the GUI input box.

**set\_parameter(*parameter*)**

Connects a new parameter to the input box, and initializes the box value.

**Parameters** *parameter* – parameter to connect.

**update\_parameter()**

Update the parameter value with the Input GUI element's current value.

**class** `pymeasure.display.inputs.IntegerInput(*args: Any, **kwargs: Any)`

Bases: [pymeasure.display.inputs.Input](#), `pyqtgraph.Qt.QtGui`.

Spin input box for integer values, connected to a `IntegerParameter`.

**set\_parameter(*parameter*)**

Connects a new parameter to the input box, and initializes the box value.

**Parameters** *parameter* – parameter to connect.

**class** `pymeasure.display.inputs.ListInput(*args: Any, **kwargs: Any)`

Bases: [pymeasure.display.inputs.Input](#), `pyqtgraph.Qt.QtGui`.

Dropdown for list values, connected to a `ListParameter`.

**set\_parameter(*parameter*)**

Connects a new parameter to the input box, and initializes the box value.

**Parameters** *parameter* – parameter to connect.

**class** `pymeasure.display.inputs.ScientificInput(*args: Any, **kwargs: Any)`

Bases: [pymeasure.display.inputs.Input](#), `pyqtgraph.Qt.QtGui`.

Spinner input box for floating-point values, connected to a `FloatParameter`. This box will display and accept values in scientific notation when appropriate.

See also:

**Class** [FloatInput](#) For a non-scientific floating-point input box.

**set\_parameter(*parameter*)**

Connects a new parameter to the input box, and initializes the box value.

**Parameters** *parameter* – parameter to connect.

**class** `pymeasure.display.inputs.StringInput(*args: Any, **kwargs: Any)`

Bases: [pymeasure.display.inputs.Input](#), `pyqtgraph.Qt.QtGui`.

String input box connected to a `Parameter`. Parameter subclasses that are string-based may also use this input, but non-string parameters should use more specialised input classes.

## 6.4 Listeners classes

**class** `pymeasure.display.listeners.Monitor(*args: Any, **kwargs: Any)`  
Bases: `pyqtgraph.Qt.QtCore.`

Monitor listens for status and progress messages from a Worker through a queue to ensure no messages are losts

**class** `pymeasure.display.listeners.QListener(*args: Any, **kwargs: Any)`  
Bases: `pyqtgraph.Qt.QtCore.`

Base class for QThreads that need to listen for messages on a ZMQ TCP port and can be stopped by a thread- and process-safe method call

## 6.5 Log classes

**class** `pymeasure.display.log.LogHandler`  
Bases: `logging.Handler`

**class** `Emitter(*args: Any, **kwargs: Any)`  
Bases: `pyqtgraph.Qt.QtCore.`

**record**

**emit**(*record*)

Do whatever it takes to actually log the specified logging record.

This version is intended to be implemented by subclasses and so raises a `NotImplementedError`.

## 6.6 Manager classes

**class** `pymeasure.display.manager.Experiment(*args: Any, **kwargs: Any)`  
Bases: `pyqtgraph.Qt.QtCore.`

The Experiment class helps group the [Procedure](#), [Results](#), and their display functionality. Its function is only a convenient container.

### Parameters

- **results** – [Results](#) object
- **curve\_list** – [ResultsCurve](#) list. List of curves associated with an experiment. They could represent different views of the same experiment.
- **browser\_item** – [BrowserItem](#) object

**class** `pymeasure.display.manager.ExperimentQueue(*args: Any, **kwargs: Any)`  
Bases: `pyqtgraph.Qt.QtCore.`

Represents a Queue of Experiments and allows queries to be easily preformed

**has\_next**()

Returns True if another item is on the queue

**next**()

Returns the next experiment on the queue

**class** `pymeasure.display.manager.Manager(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtCore.`

Controls the execution of [Experiment](#) classes by implementing a queue system in which Experiments are added, removed, executed, or aborted. When instantiated, the Manager is linked to a [Browser](#) and a PyQtGraph *PlotItem* within the user interface, which are updated in accordance with the execution status of the Experiments.

**abort()**

Aborts the currently running Experiment, but raises an exception if there is no running experiment

**clear()**

Remove all Experiments

**is\_running()**

Returns True if a procedure is currently running

**load(*experiment*)**

Load a previously executed Experiment

**next()**

Initiates the start of the next experiment in the queue as long as no other experiments are currently running and there is a procedure in the queue.

**queue(*experiment*)**

Adds an experiment to the queue.

**remove(*experiment*)**

Removes an Experiment

**resume()**

Resume processing of the queue.

## 6.7 Plotter class

**class** `pymeasure.display.plotter.Plotter(results, refresh_time=0.1, linewidth=1)`

Bases: `pymeasure.thread.StoppableThread`

Plotter dynamically plots data from a file through the Results object and supports error bars.

**See also:**

**Tutorial** [Using the Plotter](#) A tutorial and example on using the Plotter and PlotterWindow.

**run()**

Method representing the thread's activity.

You may override this method in a subclass. The standard `run()` method invokes the callable object passed to the object's constructor as the target argument, if any, with sequential and keyword arguments taken from the `args` and `kwargs` arguments, respectively.

**setup\_plot(*plot*)**

This method does nothing by default, but can be overridden by the child class in order to set up custom options for the plot window, via its [PlotItem](#).

**Parameters** `plot` – This window's [PlotItem](#) instance.

## 6.8 Qt classes

All Qt imports should reference `pymeasure.display.Qt`, for consistant importing from either PySide or PyQt4.

`Qt.fromUi(**kwargs)`

Returns a Qt object constructed using `loadUiType` based on its arguments. All `QWidget` objects in the form class are set in the returned object for easy accessibility.

## 6.9 Thread classes

`class pymeasure.display.thread.StoppableQThread(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtCore.`

Base class for QThreads which require the ability to be stopped by a thread-safe method call

`join(timeout=0)`

Joins the current thread and forces it to stop after the timeout if necessary

**Parameters** `timeout` – Timeout duration in seconds

## 6.10 Widget classes

`class pymeasure.display.widgets.BrowserWidget(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtGui.`

`class pymeasure.display.widgets.DirectoryLineEdit(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtGui.`

Widget that allows to choose a directory path. A completer is implemented for quick completion. A browse button is available.

`class pymeasure.display.widgets.EstimatorThread(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtCore.`

`class pymeasure.display.widgets.EstimatorWidget(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtGui.`

Widget that allows to display up-front estimates of the measurement procedure.

This widget relies on a `get_estimates` method of the *Procedure* class. `get_estimates` is expected to return a list of tuples, where each tuple contains two strings: a label and the estimate.

If the *SequencerWidget* is also used, it is possible to ask for the current sequencer or its length by asking for two keyword arguments in the Implementation of the `get_estimates` function: *sequence* and *sequence\_length*, respectively.

`check_get_estimates_signature()`

Method that checks the signature of the `get_estimates` function. It checks which input arguments are allowed and, if the output is correct for the *EstimatorWidget*, stores the number of estimates.

`display_estimates(estimates)`

Method that updates the shown estimates for the given set of estimates.

**Parameters** `estimates` – The set of estimates to be shown in the form of a list of tuples of (2) strings

**get\_estimates()**  
Method that makes a procedure with the currently entered parameters and returns the estimates for these parameters.

**update\_estimates()**  
Method that gets and displays the estimates. Implemented for connecting to the 'update'-button.

**class** `pymeasure.display.widgets.ImageFrame(*args: Any, **kwargs: Any)`  
Bases: `pyqtgraph.Qt.QtGui`.  
Extends `PlotFrame` to plot also axis Z using colors

**ResultsClass**  
alias of `pymeasure.display.curves.ResultsImage`

**class** `pymeasure.display.widgets.ImageWidget(*args: Any, **kwargs: Any)`  
Bases: `pymeasure.display.widgets.TabWidget`, `pyqtgraph.Qt.QtGui`.  
Extends the `ImageFrame` to allow different columns of the data to be dynamically choosen

**load**(*curve*)  
Add curve to widget

**new\_curve**(*results*, *color*=`pyqtgraphintColor`, \*\*kwargs)  
Creates a new image

**remove**(*curve*)  
Remove curve from widget

**class** `pymeasure.display.widgets.InputsWidget(*args: Any, **kwargs: Any)`  
Bases: `pyqtgraph.Qt.QtGui`.

**get\_procedure()**  
Returns the current procedure

**class** `pymeasure.display.widgets.LogWidget(*args: Any, **kwargs: Any)`  
Bases: `pymeasure.display.widgets.TabWidget`, `pyqtgraph.Qt.QtGui`.  
Widget to display logging information in GUI  
It is recommended to include this widget in all subclasses of `ManagedWindowBase`

**class** `pymeasure.display.widgets.PlotFrame(*args: Any, **kwargs: Any)`  
Bases: `pyqtgraph.Qt.QtGui`.  
Combines a `PyQtGraph Plot` with `Crosshairs`. Refreshes the plot based on the `refresh_time`, and allows the axes to be changed on the fly, which updates the plotted data

**ResultsClass**  
alias of `pymeasure.display.curves.ResultsCurve`

**parse\_axis**(*axis*)  
Returns the units of an axis by searching the string

**class** `pymeasure.display.widgets.PlotWidget(*args: Any, **kwargs: Any)`  
Bases: `pymeasure.display.widgets.TabWidget`, `pyqtgraph.Qt.QtGui`.  
Extends the `PlotFrame` to allow different columns of the data to be dynamically choosen

**load**(*curve*)  
Add curve to widget

**new\_curve**(*results*, *color*=`pyqtgraphintColor`, \*\*kwargs)  
Create a new curve

**remove**(*curve*)

Remove curve from widget

**set\_color**(*curve, color*)

Change the color of the pen of the curve

**class** pymeasure.display.widgets.**ResultsDialog**(\*args: Any, \*\*kwargs: Any)

Bases: PyQtgraph.Qt.QtGui.

**exception** pymeasure.display.widgets.**SequenceEvaluationException**

Bases: Exception

Raised when the evaluation of a sequence string goes wrong.

**class** pymeasure.display.widgets.**SequencerWidget**(\*args: Any, \*\*kwargs: Any)

Bases: PyQtgraph.Qt.QtGui.

Widget that allows to generate a sequence of measurements with varying parameters. Moreover, one can write a simple text file to easily load a sequence.

Currently requires a queue function of the ManagedWindow to have a “procedure” argument.

**static eval\_string**(*string, name=None, depth=None*)

Evaluate the given string. The string is evaluated using a list of pre-defined functions that are deemed safe to use, to prevent the execution of malicious code. For this purpose, also any built-in functions or global variables are not available.

#### Parameters

- **string** – String to be interpreted.
- **name** – Name of the to-be-interpreted string, only used for error messages.
- **depth** – Depth of the to-be-interpreted string, only used for error messages.

**get\_sequence\_from\_tree**()

Generate a list of parameters from the sequence tree.

**load\_sequence**(\*, *fileName=None*)

Load a sequence from a .txt file.

**Parameters** **fileName** – Filename (string) of the to-be-loaded file.

**queue\_sequence**()

Obtain a list of parameters from the sequence tree, enter these into procedures, and queue these procedures.

**class** pymeasure.display.widgets.**TabWidget**(*name, \*args, \*\*kwargs*)

Bases: object

Utility class to define default implementation for some basic methods.

When defining a widget to be used in subclasses of ManagedWindowBase, users should inherit from this class and provide the specialized implementation of these method’s

**load**(*curve*)

Add curve to widget

**new\_curve**(*\*args, \*\*kwargs*)

Create a new curve

**remove**(*curve*)

Remove curve from widget

**set\_color**(*curve, color*)

Set color for widget

## 6.11 Windows classes

**class** `pymeasure.display.windows.ManagedImageWindow(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtGui`.

Display experiment output with an `ImageWidget` class.

### Parameters

- **procedure\_class** – procedure class describing the experiment (see [Procedure](#))
- **x\_axis** – the data-column for the x-axis of the plot, cannot be changed afterwards for the image-plot
- **y\_axis** – the data-column for the y-axis of the plot, cannot be changed afterwards for the image-plot
- **z\_axis** – the initial data-column for the z-axis of the plot, can be changed afterwards
- **\*\*kwargs** – optional keyword arguments that will be passed to [ManagedWindow](#)

**class** `pymeasure.display.windows.ManagedWindow(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtGui`.

Display experiment output with an `PlotWidget` class.

See also:

**Tutorial [Using the ManagedWindow](#)** A tutorial and example on the basic configuration and usage of `ManagedWindow`.

### Parameters

- **procedure\_class** – procedure class describing the experiment (see [Procedure](#))
- **x\_axis** – the initial data-column for the x-axis of the plot
- **y\_axis** – the initial data-column for the y-axis of the plot
- **linewidth** – linewidth for the displayed curves, default is 1
- **\*\*kwargs** – optional keyword arguments that will be passed to [ManagedWindowBase](#)

**class** `pymeasure.display.windows.ManagedWindowBase(*args: Any, **kwargs: Any)`

Bases: `pyqtgraph.Qt.QtGui`.

Base class for GUI experiment management .

The `ManagedWindowBase` provides an interface for inputting experiment parameters, running several experiments ([Procedure](#)), plotting result curves, and listing the experiments conducted during a session.

The `ManagedWindowBase` uses a `Manager` to control `Workers` in a `Queue`, and provides a simple interface. The `queue()` method must be overridden by the child class.

The `ManagedWindowBase` allow user to define a set of widget that display information about the experiment. The information displayed may include: plots, tabular view, logging information, etc.

This class is not intended to be used directly, but it should be subclassed to provide some appropriate widget list. Example of classes usable as element of widget list are:

- [LogWidget](#)
- [PlotWidget](#)
- [ImageWidget](#)

Of course, users can define its own widget making sure that inherits from *TabWidget*.

Examples of ready to use classes inherited from *ManagedWindowBase* are:

- *ManagedWindow*
- *ManagedImageWindow*

See also:

**Tutorial *Using the ManagedWindow*** A tutorial and example on the basic configuration and usage of *ManagedWindow*.

Parameters for `__init__` constructor.

#### Parameters

- **procedure\_class** – procedure class describing the experiment (see *Procedure*)
- **widget\_list** – list of widget to be displayed in the GUI
- **inputs** – list of *Parameter* instance variable names, which the display will generate graphical fields for
- **displays** – list of *Parameter* instance variable names displayed in the browser window
- **log\_channel** – logging.Logger instance to use for logging output
- **log\_level** – logging level
- **parent** – Parent widget or None
- **sequencer** – a boolean stating whether or not the sequencer has to be included into the window
- **sequencer\_inputs** – either None or a list of the parameter names to be scanned over. If no list of parameters is given, the parameters displayed in the manager queue are used
- **sequence\_file** – simple text file to quickly load a pre-defined sequence with the Load sequence button
- **inputs\_in\_scrollarea** – boolean that display or hide a scrollbar to the input area
- **directory\_input** – specify, if present, where the experiment’s result will be saved.
- **hide\_groups** – a boolean controlling whether parameter groups are hidden (True, default) or disabled/grayed-out (False) when the group conditions are not met.

#### **open\_file\_externally**(filename)

Method to open the datafile using an external editor or viewer. Uses the default application to open a datafile of this filetype, but can be overridden by the child class in order to open the file in another application of choice.

#### **queue**(procedure=None)

Abstract method, which must be overridden by the child class.

Implementations must call `self.manager.queue(experiment)` and pass an experiment (*Experiment*) object which contains the *Results* and *Procedure* to be run.

The optional *procedure* argument is not required for a basic implementation, but is required when the *SequencerWidget* is used.

For example:



```

def queue(self):
    filename = unique_filename('results', prefix="data") # from pymeasure.
    ↪experiment

    procedure = self.make_procedure() # Procedure class was passed at
    ↪construction
    results = Results(procedure, filename)
    experiment = self.new_experiment(results)

    self.manager.queue(experiment)

```

**set\_parameters(parameters)**

This method should be overwritten by the child class. The parameters argument is a dictionary of Parameter objects. The Parameters should overwrite the GUI values so that a user can click “Queue” to capture the same parameters.

**class** pymeasure.display.windows.**PlotterWindow**(\*args: Any, \*\*kwargs: Any)

Bases: pyqtgraph.Qt.QtGui.

A window for plotting experiment results. Should not be instantiated directly, but only via the [Plotter](#) class.

**See also:**

**Tutorial** [Using the Plotter](#) A tutorial and example code for using the Plotter and PlotterWindow.

**check\_stop()**

Checks if the Plotter should stop and exits the Qt main loop if so



## PYMEASURE.INSTRUMENTS

This section contains documentation on the instrument classes.

### 7.1 Instrument classes

**class** `pymeasure.instruments.Instrument`(*adapter, name, includeSCPI=True, \*\*kwargs*)

This provides the base class for all Instruments, which is independent of the particular Adapter used to connect for communication to the instrument. It provides basic SCPI commands by default, but can be toggled with `includeSCPI`.

**Parameters**

- **adapter** – An *Adapter* object
- **name** – A string name
- **includeSCPI** – A boolean, which toggles the inclusion of standard SCPI commands

**ask**(*command*)

Writes the command to the instrument through the adapter and returns the read response.

**Parameters** **command** – command string to be sent to the instrument

**check\_errors**()

Read all errors from the instrument.

**Returns** list of error entries

**clear**()

Clears the instrument status byte

**property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device's Output Queue when all pending selected device operations have been finished.

**static control**(*get\_command, set\_command, docs, validator=<function Instrument.<lambda>>, values=(), map\_values=False, get\_process=<function Instrument.<lambda>>, set\_process=<function Instrument.<lambda>>, check\_set\_errors=False, check\_get\_errors=False, \*\*kwargs*)

Returns a property for the class based on the supplied commands. This property may be set and read from the instrument.

**Parameters**

- **get\_command** – A string command that asks for the value

- **set\_command** – A string command that writes the value
- **docs** – A docstring that will be included in the documentation
- **validator** – A function that takes both a value and a group of valid values and returns a valid value, while it otherwise raises an exception
- **values** – A list, tuple, range, or dictionary of valid values, that can be used as to map values if **map\_values** is True.
- **map\_values** – A boolean flag that determines if the values should be interpreted as a map
- **get\_process** – A function that take a value and allows processing before value mapping, returning the processed value
- **set\_process** – A function that takes a value and allows processing before value mapping, returning the processed value
- **check\_set\_errors** – Toggles checking errors after setting
- **check\_get\_errors** – Toggles checking errors after getting

#### **property id**

Requests and returns the identification of the instrument.

**static measurement**(*get\_command*, *docs*, *values*=(), *map\_values*=None, *get\_process*=<function Instrument.<lambda>>, *command\_process*=<function Instrument.<lambda>>, *check\_get\_errors*=False, *\*\*kwargs*)

Returns a property for the class based on the supplied commands. This is a measurement quantity that may only be read from the instrument, not set.

#### **Parameters**

- **get\_command** – A string command that asks for the value
- **docs** – A docstring that will be included in the documentation
- **values** – A list, tuple, range, or dictionary of valid values, that can be used as to map values if **map\_values** is True.
- **map\_values** – A boolean flag that determines if the values should be interpreted as a map
- **get\_process** – A function that take a value and allows processing before value mapping, returning the processed value
- **command\_process** – A function that take a command and allows processing before executing the command, for both getting and setting
- **check\_get\_errors** – Toggles checking errors after getting

#### **property options**

Requests and returns the device options installed.

#### **read()**

Reads from the instrument through the adapter and returns the response.

#### **reset()**

Resets the instrument.

**static setting**(*set\_command*, *docs*, *validator*=<function Instrument.<lambda>>, *values*=(), *map\_values*=False, *set\_process*=<function Instrument.<lambda>>, *check\_set\_errors*=False, *\*\*kwargs*)

Returns a property for the class based on the supplied commands. This property may be set, but raises an exception when being read from the instrument.

**Parameters**

- **set\_command** – A string command that writes the value
- **docs** – A docstring that will be included in the documentation
- **validator** – A function that takes both a value and a group of valid values and returns a valid value, while it otherwise raises an exception
- **values** – A list, tuple, range, or dictionary of valid values, that can be used as to map values if **map\_values** is True.
- **map\_values** – A boolean flag that determines if the values should be interpreted as a map
- **set\_process** – A function that takes a value and allows processing before value mapping, returning the processed value
- **check\_set\_errors** – Toggles checking errors after setting

**shutdown()**

Brings the instrument to a safe and stable state

**property status**

Requests and returns the status byte and Master Summary Status bit.

**values(command, \*\*kwargs)**

Reads a set of values from the instrument through the adapter, passing on any key-word arguments.

**write(command)**

Writes the command to the instrument through the adapter.

**Parameters** **command** – command string to be sent to the instrument

**class** `pymeasure.instruments.Mock(wait=0.1, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Mock instrument for testing.

**get\_time()**

Get elapsed time

**get\_voltage()**

Get the voltage.

**get\_wave()**

Get wave.

**reset\_time()**

Reset the timer to 0 s.

**set\_output\_voltage(value)**

Set the voltage.

**set\_time(value)**

Wait for the timer to reach the specified time. If value = 0, reset.

**property time**

Get elapsed time

**property voltage**

Get the voltage.

**property wave**

Get wave.

## 7.2 Validator functions

Validators are used in conjunction with the [\*Instrument.control\*](#) function to allow properties with complex restrictions for valid values. They are described in more detail in the [\*Advanced properties\*](#) section.

`pymeasure.instruments.validators.discreteTruncate(number, discreteSet)`

Truncates the number to the closest element in the positive discrete set. Returns False if the number is larger than the maximum value or negative.

`pymeasure.instruments.validators.joined_validators(*validators)`

Join a list of validators together as a single. Expects a list of validator functions and values.

**Parameters** `validators` – an iterable of other validators

`pymeasure.instruments.validators.modular_range(value, values)`

Provides a validator function that returns the value if it is in the range. Otherwise it returns the value, modulo the max of the range.

**Parameters**

- **value** – a value to test
- **values** – A set of values that are valid

`pymeasure.instruments.validators.modular_range_bidirectional(value, values)`

Provides a validator function that returns the value if it is in the range. Otherwise it returns the value, modulo the max of the range. Allows negative values.

**Parameters**

- **value** – a value to test
- **values** – A set of values that are valid

`pymeasure.instruments.validators.strict_discrete_range(value, values, step)`

Provides a validator function that returns the value if its value is less than the maximum and greater than the minimum of the range and is a multiple of step. Otherwise it raises a `ValueError`.

**Parameters**

- **value** – A value to test
- **values** – A range of values (range, list, etc.)
- **step** – Minimum stepsize (resolution limit)

**Raises** `ValueError` if the value is out of the range

`pymeasure.instruments.validators.strict_discrete_set(value, values)`

Provides a validator function that returns the value if it is in the discrete set. Otherwise it raises a `ValueError`.

**Parameters**

- **value** – A value to test
- **values** – A set of values that are valid

**Raises** `ValueError` if the value is not in the set

`pymeasure.instruments.validators.strict_range(value, values)`

Provides a validator function that returns the value if its value is less than the maximum and greater than the minimum of the range. Otherwise it raises a `ValueError`.

**Parameters**

- **value** – A value to test

- **values** – A range of values (range, list, etc.)

**Raises** ValueError if the value is out of the range

`pymeasure.instruments.validators.truncated_discrete_set(value, values)`

Provides a validator function that returns the value if it is in the discrete set. Otherwise, it returns the smallest value that is larger than the value.

**Parameters**

- **value** – A value to test
- **values** – A set of values that are valid

`pymeasure.instruments.validators.truncated_range(value, values)`

Provides a validator function that returns the value if it is in the range. Otherwise it returns the closest range bound.

**Parameters**

- **value** – A value to test
- **values** – A set of values that are valid

## 7.3 Comedi data acquisition

The Comedi libraries provide a convenient method for interacting with data acquisition cards, but are restricted to Linux compatible operating systems.

`pymeasure.instruments.comedi.getAI(device, channel, range=None)`

Returns the analog input channel as specified for a given device

`pymeasure.instruments.comedi.getAO(device, channel, range=None)`

Returns the analog output channel as specified for a given device

`pymeasure.instruments.comedi.readAI(device, channel, range=None, count=1)`

Reads a single measurement (count==1) from the analog input channel of the device specified. Multiple readings can be preformed with count not equal to one, which are seperated by an arbitrary time

`pymeasure.instruments.comedi.writeAO(device, channel, voltage, range=None)`

Writes a single voltage to the analog output channel of the device specified

## 7.4 Resource Manager

The `list_resources` function provides an interface to check connected instruments interactively.

`pymeasure.instruments.list_resources()`

Prints the available resources, and returns a list of VISA resource names

```
resources = list_resources()
#prints (e.g.)
#0 : GPIB0::22::INSTR : Agilent Technologies,34410A,*****
#1 : GPIB0::26::INSTR : Keithley Instruments Inc., Model 2612, *****
dmm = Agilent34410(resources[0])
```

Instruments by manufacturer:

## 7.5 Advantest

This section contains specific documentation on the Advantest instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.5.1 Advantest R3767CG Vector Network Analyzer

**class** pymeasure.instruments.advantest.advantestR3767CG.**AdvantestR3767CG**(resourceName,  
\*\*kwargs)

Bases: *pymeasure.instruments.instrument.Instrument*

Represents the Advantest R3767CG VNA. Implements controls to change the analysis range and to retrieve the data for the trace.

**property** center\_frequency

Center Frequency in Hz

**property** id

Reads the instrument identification

**property** span\_frequency

Span Frequency in Hz

**property** start\_frequency

Starting frequency in Hz

**property** stop\_frequency

Stopping frequency in Hz

**property** trace\_1

Reads the Data array from trace 1 after formatting

## 7.6 Agilent

This section contains specific documentation on the Agilent instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.6.1 Agilent 8257D Signal Generator

**class** pymeasure.instruments.agilent.**Agilent8257D**(adapter, \*\*kwargs)

Bases: *pymeasure.instruments.instrument.Instrument*

Represents the Agilent 8257D Signal Generator and provides a high-level interface for interacting with the instrument.

```
generator = Agilent8257D("GPIB::1")

generator.power = 0           # Sets the output power to 0 dBm
generator.frequency = 5       # Sets the output frequency to 5 GHz
generator.enable()            # Enables the output
```

**property** amplitude\_depth

A floating point property that controls the amplitude modulation in percent, which can take values from 0 to 100 %.



**property amplitude\_source**

A string property that controls the source of the amplitude modulation signal, which can take the values: 'internal', 'internal 2', 'external', and 'external 2'.

**property center\_frequency**

A floating point property that represents the center frequency in Hz. This property can be set.

**config\_amplitude\_modulation**(*frequency=1000.0, depth=100.0, shape='sine'*)

Configures the amplitude modulation of the output signal.

**Parameters**

- **frequency** – A modulation frequency for the internal oscillator
- **depth** – A linear depth percentage
- **shape** – A string that describes the shape for the internal oscillator

**config\_low\_freq\_out**(*source='internal', amplitude=3*)

Configures the low-frequency output signal.

**Parameters**

- **source** – The source for the low-frequency output signal.
- **amplitude** – Amplitude of the low-frequency output

**config\_pulse\_modulation**(*frequency=1000.0, input='square'*)

Configures the pulse modulation of the output signal.

**Parameters**

- **frequency** – A pulse rate frequency in Hertz
- **input** – A string that describes the internal pulse input

**config\_step\_sweep**()

Configures a step sweep through frequency

**disable**()

Disables the output of the signal.

**disable\_amplitude\_modulation**()

Disables amplitude modulation of the output signal.

**disable\_low\_freq\_out**()

Disables low frequency output

**disable\_modulation**()

Disables the signal modulation.

**disable\_pulse\_modulation**()

Disables pulse modulation of the output signal.

**property dwell\_time**

A floating point property that represents the settling time in seconds at the current frequency or power setting. This property can be set.

**enable**()

Enables the output of the signal.

**enable\_amplitude\_modulation**()

Enables amplitude modulation of the output signal.

**enable\_low\_freq\_out**()

Enables low frequency output

**enable\_pulse\_modulation()**

Enables pulse modulation of the output signal.

**property frequency**

A floating point property that represents the output frequency in Hz. This property can be set.

**property has\_amplitude\_modulation**

Reads a boolean value that is True if the amplitude modulation is enabled.

**property has\_modulation**

Reads a boolean value that is True if the modulation is enabled.

**property has\_pulse\_modulation**

Reads a boolean value that is True if the pulse modulation is enabled.

**property internal\_frequency**

A floating point property that controls the frequency of the internal oscillator in Hertz, which can take values from 0.5 Hz to 1 MHz.

**property internal\_shape**

A string property that controls the shape of the internal oscillations, which can take the values: 'sine', 'triangle', 'square', 'ramp', 'noise', 'dual-sine', and 'swept-sine'.

**property is\_enabled**

Reads a boolean value that is True if the output is on.

**property low\_freq\_out\_amplitude**

A floating point property that controls the peak voltage (amplitude) of the low frequency output in volts, which can take values from 0-3.5V

**property low\_freq\_out\_source**

A string property which controls the source of the low frequency output, which can take the values 'internal [2]' for the internal source, or 'function [2]' for an internal function generator which can be configured.

**property power**

A floating point property that represents the output power in dBm. This property can be set.

**property pulse\_frequency**

A floating point property that controls the pulse rate frequency in Hertz, which can take values from 0.1 Hz to 10 MHz.

**property pulse\_input**

A string property that controls the internally generated modulation input for the pulse modulation, which can take the values: 'square', 'free-run', 'triggered', 'doublet', and 'gated'.

**property pulse\_source**

A string property that controls the source of the pulse modulation signal, which can take the values: 'internal', 'external', and 'scalar'.

**shutdown()**

Shuts down the instrument by disabling any modulation and the output signal.

**property start\_frequency**

A floating point property that represents the start frequency in Hz. This property can be set.

**property start\_power**

A floating point property that represents the start power in dBm. This property can be set.

**start\_step\_sweep()**

Starts a step sweep.

**property step\_points**

An integer number of points in a step sweep. This property can be set.

**property stop\_frequency**

A floating point property that represents the stop frequency in Hz. This property can be set.

**property stop\_power**

A floating point property that represents the stop power in dBm. This property can be set.

**stop\_step\_sweep()**

Stops a step sweep.

## 7.6.2 Agilent 8722ES Vector Network Analyzer

**class** `pymeasure.instruments.agilent.Agilent8722ES(resourceName, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Agilent8722ES Vector Network Analyzer and provides a high-level interface for taking scans of the scattering parameters.

**property averages**

An integer representing the number of averages to take. Note that averaging must be enabled for this to take effect. This property can be set.

**property averaging\_enabled**

A bool that indicates whether or not averaging is enabled. This property can be set.

**property data**

Returns the real and imaginary data from the last scan

**property data\_complex**

Returns the complex power from the last scan

**property data\_log\_magnitude**

Returns the absolute magnitude values in dB from the last scan

**property data\_magnitude**

Returns the absolute magnitude values from the last scan

**property data\_phase**

Returns the phase in degrees from the last scan

**disable\_averaging()**

Disables averaging

**enable\_averaging()**

Enables averaging

**property frequencies**

Returns a list of frequencies from the last scan

**is\_averaging()**

Returns True if averaging is enabled

**log\_magnitude(real, imaginary)**

Returns the magnitude in dB from a real and imaginary number or numpy arrays

**magnitude(real, imaginary)**

Returns the magnitude from a real and imaginary number or numpy arrays

**phase(real, imaginary)**

Returns the phase in degrees from a real and imaginary number or numpy arrays

**scan(averages=None, blocking=None, timeout=None, delay=None)**

Initiates a scan with the number of averages specified and blocks until the operation is complete.

**scan\_continuous()**

Initiates a continuous scan

**property scan\_points**

Gets the number of scan points

**scan\_single()**

Initiates a single scan

**set\_IF\_bandwidth(*bandwidth*)**

Sets the resolution bandwidth (IF bandwidth)

**set\_averaging(*averages*)**

Sets the number of averages and enables/disables averaging. Should be between 1 and 999

**set\_fixed\_frequency(*frequency*)**

Sets the scan to be of only one frequency in Hz

**property start\_frequency**

A floating point property that represents the start frequency in Hz. This property can be set.

**property stop\_frequency**

A floating point property that represents the stop frequency in Hz. This property can be set.

**property sweep\_time**

A floating point property that represents the sweep time in seconds. This property can be set.

### 7.6.3 Agilent E4408B Spectrum Analyzer

**class** `pymeasure.instruments.agilent.AgilentE4408B(resourceName, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the AgilentE4408B Spectrum Analyzer and provides a high-level interface for taking scans of high-frequency spectrums

**property center\_frequency**

A floating point property that represents the center frequency in Hz. This property can be set.

**property frequencies**

Returns a numpy array of frequencies in Hz that correspond to the current settings of the instrument.

**property frequency\_points**

An integer property that represents the number of frequency points in the sweep. This property can take values from 101 to 8192.

**property frequency\_step**

A floating point property that represents the frequency step in Hz. This property can be set.

**property start\_frequency**

A floating point property that represents the start frequency in Hz. This property can be set.

**property stop\_frequency**

A floating point property that represents the stop frequency in Hz. This property can be set.

**property sweep\_time**

A floating point property that represents the sweep time in seconds. This property can be set.

**trace(*number=1*)**

Returns a numpy array of the data for a particular trace based on the trace number (1, 2, or 3).

**trace\_df**(*number=1*)

Returns a pandas DataFrame containing the frequency and peak data for a particular trace, based on the trace number (1, 2, or 3).

## 7.6.4 Agilent E4980 LCR Meter

**class** `pymeasure.instruments.agilent.AgilentE4980`(*adapter*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents LCR meter E4980A/AL

**property** `ac_current`

AC current level, in Amps

**property** `ac_voltage`

AC voltage level, in Volts

**aperture**(*time=None*, *averages=1*)

Set and get aperture.

### Parameters

- **time** – integration time as string: SHORT, MED, LONG (case insensitive); if None, get values
- **averages** – number of averages, numeric

**freq\_sweep**(*freq\_list*, *return\_freq=False*)

Run frequency list sweep using sequential trigger.

### Parameters

- **freq\_list** – list of frequencies
- **return\_freq** – if True, returns the frequencies read from the instrument

Returns values as configured with `mode`

**property** `frequency`

AC frequency (range depending on model), in Hertz

**property** `impedance`

Measured data A and B, according to `mode`

**property** `mode`

Select quantities to be measured:

- CPD: Parallel capacitance [F] and dissipation factor [number]
- CPQ: Parallel capacitance [F] and quality factor [number]
- CPG: Parallel capacitance [F] and parallel conductance [S]
- CPRP: Parallel capacitance [F] and parallel resistance [Ohm]
- CSD: Series capacitance [F] and dissipation factor [number]
- CSQ: Series capacitance [F] and quality factor [number]
- CSRS: Series capacitance [F] and series resistance [Ohm]
- LPD: Parallel inductance [H] and dissipation factor [number]
- LPQ: Parallel inductance [H] and quality factor [number]

- LPG: Parallel inductance [H] and parallel conductance [S]
- LPRP: Parallel inductance [H] and parallel resistance [Ohm]
- LSD: Series inductance [H] and dissipation factor [number]
- LSQ: Series inductance [H] and quality factor [number]
- LSRS: Series inductance [H] and series resistance [Ohm]
- RX: Resistance [Ohm] and reactance [Ohm]
- ZTD: Impedance, magnitude [Ohm] and phase [deg]
- ZTR: Impedance, magnitude [Ohm] and phase [rad]
- GB: Conductance [S] and susceptance [S]
- YTD: Admittance, magnitude [Ohm] and phase [deg]
- YTR: Admittance magnitude [Ohm] and phase [rad]

**property trigger\_source**

Select trigger source; accept the values:

- HOLD: manual
- INT: internal
- BUS: external bus (GPIB/LAN/USB)
- EXT: external connector

## 7.6.5 Agilent 34410A Multimeter

**class** `pymeasure.instruments.agilent.Agilent34410A(adapter, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represent the HP/Agilent/Keysight 34410A and related multimeters.

Implemented measurements: voltage\_dc, voltage\_ac, current\_dc, current\_ac, resistance, resistance\_4w

**property current\_ac**

AC current, in Amps

**property current\_dc**

DC current, in Amps

**property resistance**

Resistance, in Ohms

**property resistance\_4w**

Four-wires (remote sensing) resistance, in Ohms

**property voltage\_ac**

AC voltage, in Volts

**property voltage\_dc**

DC voltage, in Volts

## 7.6.6 HP/Agilent/Keysight 34450A Digital Multimeter

**class** `pymeasure.instruments.agilent.Agilent34450A(adapter, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represent the HP/Agilent/Keysight 34450A and related multimeters.

```
dmm = Agilent34450A("USB0:...")
dmm.reset()
dmm.configure_voltage()
print(dmm.voltage)
dmm.shutdown()
```

**beep()**

Sounds a system beep.

**property capacitance**

Reads a capacitance measurement in Farads, based on the active mode.

**property capacitance\_auto\_range**

A boolean property that toggles auto ranging for capacitance.

**property capacitance\_range**

A property that controls the capacitance range in Farads, which can take values 1E-9, 10E-9, 100E-9, 1E-6, 10E-6, 100E-6, 1E-3, 10E-3, as well as “MIN”, “MAX”, or “DEF” (1E-6). Auto-range is disabled when this property is set.

**configure\_capacitance(capacitance\_range='AUTO')**

Configures the instrument to measure capacitance.

**Parameters** **capacitance\_range** – A capacitance in Farads to set the capacitance range, can be 1E-9, 10E-9, 100E-9, 1E-6, 10E-6, 100E-6, 1E-3, 10E-3, as well as “MIN”, “MAX”, “DEF” (1E-6), or “AUTO”.

**configure\_continuity()**

Configures the instrument to measure continuity.

**configure\_current(current\_range='AUTO', ac=False, resolution='DEF')**

Configures the instrument to measure current.

**Parameters**

- **current\_range** – A current in Amps to set the current range. DC values can be 100E-6, 1E-3, 10E-3, 100E-3, 1, 10, as well as “MIN”, “MAX”, “DEF” (100 mA), or “AUTO”. AC values can be 10E-3, 100E-3, 1, 10, as well as “MIN”, “MAX”, “DEF” (100 mA), or “AUTO”.
- **ac** – False for DC current, and True for AC current
- **resolution** – Desired resolution, can be 3.00E-5, 2.00E-5, 1.50E-6 (5 1/2 digits), as well as “MIN”, “MAX”, or “DEF” (1.50E-6).

**configure\_diode()**

Configures the instrument to measure diode voltage.

**configure\_frequency(measured\_from='voltage\_ac', measured\_from\_range='AUTO', aperture='DEF')**

Configures the instrument to measure frequency.

**Parameters**

- **measured\_from** – “voltage\_ac” or “current\_ac”

- **measured\_from\_range** – range of measured\_from. AC voltage can have ranges 100E-3, 1, 10, 100, 750, as well as “MIN”, “MAX”, “DEF” (10 V), or “AUTO”. AC current can have ranges 10E-3, 100E-3, 1, 10, as well as “MIN”, “MAX”, “DEF” (100 mA), or “AUTO”.
- **aperture** – Aperture time in Seconds, can be 100 ms, 1 s, as well as “MIN”, “MAX”, or “DEF” (1 s).

**configure\_resistance**(*resistance\_range='AUTO', wires=2, resolution='DEF'*)

Configures the instrument to measure resistance.

#### Parameters

- **resistance\_range** – A resistance in Ohms to set the resistance range, can be 100, 1E3, 10E3, 100E3, 1E6, 10E6, 100E6, as well as “MIN”, “MAX”, “DEF” (1E3), or “AUTO”.
- **wires** – Number of wires used for measurement, can be 2 or 4.
- **resolution** – Desired resolution, can be 3.00E-5, 2.00E-5, 1.50E-6 (5 1/2 digits), as well as “MIN”, “MAX”, or “DEF” (1.50E-6).

**configure\_temperature**()

Configures the instrument to measure temperature.

**configure\_voltage**(*voltage\_range='AUTO', ac=False, resolution='DEF'*)

Configures the instrument to measure voltage.

#### Parameters

- **voltage\_range** – A voltage in Volts to set the voltage range. DC values can be 100E-3, 1, 10, 100, 1000, as well as “MIN”, “MAX”, “DEF” (10 V), or “AUTO”. AC values can be 100E-3, 1, 10, 100, 750, as well as “MIN”, “MAX”, “DEF” (10 V), or “AUTO”.
- **ac** – False for DC voltage, True for AC voltage
- **resolution** – Desired resolution, can be 3.00E-5, 2.00E-5, 1.50E-6 (5 1/2 digits), as well as “MIN”, “MAX”, or “DEF” (1.50E-6).

**property continuity**

Reads a continuity measurement in Ohms, based on the active mode.

**property current**

Reads a DC current measurement in Amps, based on the active mode.

**property current\_ac**

Reads an AC current measurement in Amps, based on the active mode.

**property current\_ac\_auto\_range**

A boolean property that toggles auto ranging for AC current.

**property current\_ac\_range**

A property that controls the AC current range in Amps, which can take values 10E-3, 100E-3, 1, 10, as well as “MIN”, “MAX”, or “DEF” (100 mA). Auto-range is disabled when this property is set.

**property current\_ac\_resolution**

An property that controls the resolution in the AC current readings, which can take values 3.00E-5, 2.00E-5, 1.50E-6 (5 1/2 digits), as well as “MIN”, “MAX”, or “DEF” (1.50E-6).

**property current\_auto\_range**

A boolean property that toggles auto ranging for DC current.



**property current\_range**

A property that controls the DC current range in Amps, which can take values 100E-6, 1E-3, 10E-3, 100E-3, 1, 10, as well as “MIN”, “MAX”, or “DEF” (100 mA). Auto-range is disabled when this property is set.

**property current\_resolution**

A property that controls the resolution in the DC current readings, which can take values 3.00E-5, 2.00E-5, 1.50E-6 (5 1/2 digits), as well as “MIN”, “MAX”, and “DEF” (3.00E-5).

**property diode**

Reads a diode measurement in Volts, based on the active mode.

**property frequency**

Reads a frequency measurement in Hz, based on the active mode.

**property frequency\_aperture**

A property that controls the frequency aperture in seconds, which sets the integration period and measurement speed. Takes values 100 ms, 1 s, as well as “MIN”, “MAX”, or “DEF” (1 s).

**property frequency\_current\_auto\_range**

A boolean property that toggles auto ranging for AC current in frequency measurements.

**property frequency\_current\_range**

A property that controls the current range in Amps for frequency on AC current measurements, which can take values 10E-3, 100E-3, 1, 10, as well as “MIN”, “MAX”, or “DEF” (100 mA). Auto-range is disabled when this property is set.

**property frequency\_voltage\_auto\_range**

A boolean property that toggles auto ranging for AC voltage in frequency measurements.

**property frequency\_voltage\_range**

A property that controls the voltage range in Volts for frequency on AC voltage measurements, which can take values 100E-3, 1, 10, 100, 750, as well as “MIN”, “MAX”, or “DEF” (10 V). Auto-range is disabled when this property is set.

**property resistance**

Reads a resistance measurement in Ohms for 2-wire configuration, based on the active mode.

**property resistance\_4w**

Reads a resistance measurement in Ohms for 4-wire configuration, based on the active mode.

**property resistance\_4w\_auto\_range**

A boolean property that toggles auto ranging for 4-wire resistance.

**property resistance\_4w\_range**

A property that controls the 4-wire resistance range in Ohms, which can take values 100, 1E3, 10E3, 100E3, 1E6, 10E6, 100E6, as well as “MIN”, “MAX”, or “DEF” (1E3). Auto-range is disabled when this property is set.

**property resistance\_4w\_resolution**

A property that controls the resolution in the 4-wire resistance readings, which can take values 3.00E-5, 2.00E-5, 1.50E-6 (5 1/2 digits), as well as “MIN”, “MAX”, or “DEF” (1.50E-6).

**property resistance\_auto\_range**

A boolean property that toggles auto ranging for 2-wire resistance.

**property resistance\_range**

A property that controls the 2-wire resistance range in Ohms, which can take values 100, 1E3, 10E3, 100E3, 1E6, 10E6, 100E6, as well as “MIN”, “MAX”, or “DEF” (1E3). Auto-range is disabled when this property is set.

**property resistance\_resolution**

A property that controls the resolution in the 2-wire resistance readings, which can take values 3.00E-5, 2.00E-5, 1.50E-6 (5 1/2 digits), as well as “MIN”, “MAX”, or “DEF” (1.50E-6).

**property temperature**

Reads a temperature measurement in Celsius, based on the active mode.

**property voltage**

Reads a DC voltage measurement in Volts, based on the active mode.

**property voltage\_ac**

Reads an AC voltage measurement in Volts, based on the active mode.

**property voltage\_ac\_auto\_range**

A boolean property that toggles auto ranging for AC voltage.

**property voltage\_ac\_range**

A property that controls the AC voltage range in Volts, which can take values 100E-3, 1, 10, 100, 750, as well as “MIN”, “MAX”, or “DEF” (10 V). Auto-range is disabled when this property is set.

**property voltage\_ac\_resolution**

A property that controls the resolution in the AC voltage readings, which can take values 3.00E-5, 2.00E-5, 1.50E-6 (5 1/2 digits), as well as “MIN”, “MAX”, or “DEF” (1.50E-6).

**property voltage\_auto\_range**

A boolean property that toggles auto ranging for DC voltage.

**property voltage\_range**

A property that controls the DC voltage range in Volts, which can take values 100E-3, 1, 10, 100, 1000, as well as “MIN”, “MAX”, or “DEF” (10 V). Auto-range is disabled when this property is set.

**property voltage\_resolution**

A property that controls the resolution in the DC voltage readings, which can take values 3.00E-5, 2.00E-5, 1.50E-6 (5 1/2 digits), as well as “MIN”, “MAX”, or “DEF” (1.50E-6).

## 7.6.7 Agilent 4155/4156 Semiconductor Parameter Analyzer

**class** `pymeasure.instruments.agilent.agilent4156.Agilent4156`(*resourceName*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Agilent 4155/4156 Semiconductor Parameter Analyzer and provides a high-level interface for taking current-voltage (I-V) measurements.

```
from pymeasure.instruments.agilent import Agilent4156

# explicitly define r/w terminations; set sufficiently large timeout or None.
smu = Agilent4156("GPIB0::25", read_termination = '\n', write_termination = '\n',
    ↪ timeout=None)

# reset the instrument
smu.reset()

# define configuration file for instrument and load config
smu.configure("configuration_file.json")

# save data variables, some or all of which are defined in the json config file.
smu.save(['VC', 'IC', 'VB', 'IB'])
```

(continues on next page)

(continued from previous page)

```
# take measurements
status = smu.measure()

# measured data is a pandas dataframe and can be exported to csv.
data = smu.get_data(path='./t1.csv')
```

The JSON file is an ascii text configuration file that defines the settings of each channel on the instrument. The JSON file is used to configure the instrument using the convenience function `configure()` as shown in the example above. For example, the instrument setup for a bipolar transistor measurement is shown below.

```
{
  "SMU1": {
    "voltage_name" : "VC",
    "current_name" : "IC",
    "channel_function" : "VAR1",
    "channel_mode" : "V",
    "series_resistance" : "0OHM"
  },
  "SMU2": {
    "voltage_name" : "VB",
    "current_name" : "IB",
    "channel_function" : "VAR2",
    "channel_mode" : "I",
    "series_resistance" : "0OHM"
  },
  "SMU3": {
    "voltage_name" : "VE",
    "current_name" : "IE",
    "channel_function" : "CONS",
    "channel_mode" : "V",
    "constant_value" : 0,
    "compliance" : 0.1
  },
  "SMU4": {
    "voltage_name" : "VS",
    "current_name" : "IS",
    "channel_function" : "CONS",
    "channel_mode" : "V",
    "constant_value" : 0,
    "compliance" : 0.1
  },
  "VAR1": {
    "start" : 1,
    "stop" : 2,
    "step" : 0.1,
    "spacing" : "LINEAR",
    "compliance" : 0.1
  }
```

(continues on next page)

(continued from previous page)

```

    },
    "VAR2": {
        "start" : 0,
        "step" : 10e-6,
        "points" : 3,
        "compliance" : 2
    }
}

```

**property analyzer\_mode**

A string property that controls the instrument operating mode.

- Values: SWEEP, SAMPLING

```
smu.analyzer_mode = "SWEEP"
```

**configure(*config\_file*)**

Convenience function to configure the channel setup and sweep using a **JSON (JavaScript Object Notation)** configuration file.

**Parameters** *config\_file* – JSON file to configure instrument channels.

```
instr.configure('config.json')
```

**property data\_variables**

Gets a string list of data variables for which measured data is available. This looks for all the variables saved by the *save()* and *save\_var()* methods and returns it. This is useful for creation of dataframe headers.

**Returns** List

```
header = instr.data_variables
```

**property delay\_time**

A floating point property that measurement delay time in seconds, which can take the values from 0 to 65s in 0.1s steps.

```
instr.delay_time = 1 # delay time of 1-sec
```

**disable\_all()**

Disables all channels in the instrument.

```
instr.disable_all()
```

**get\_data(*path=None*)**

Gets the measurement data from the instrument after completion. If the measurement period is set to INF in the *measure()* method, then the measurement must be stopped using *stop()* before getting valid data.

**Parameters** *path* – Path for optional data export to CSV.

**Returns** Pandas Dataframe

```
df = instr.get_data(path='./datafolder/data1.csv')
```

**property hold\_time**

A floating point property that measurement hold time in seconds, which can take the values from 0 to 655s in 1s steps.

```
instr.hold_time = 2 # hold time of 2-secs.
```

**property integration\_time**

A string property that controls the integration time.

- Values: SHORT, MEDIUM, LONG

```
instr.integration_time = "MEDIUM"
```

**measure**(*period='INF', points=100*)

Performs a single measurement and waits for completion in sweep mode. In sampling mode, the measurement period and number of points can be specified.

**Parameters**

- **period** – Period of sampling measurement from 6E-6 to 1E11 seconds. Default setting is INF.
- **points** – Number of samples to be measured, from 1 to 10001. Default setting is 100.

**save**(*trace\_list*)

Save the voltage or current in the instrument display list

**Parameters** **trace\_list** – A list of channel variables whose measured data should be saved.

A maximum of 8 variables are allowed. If only one variable is being saved, a string can be specified.

```
instr.save(['IC', 'IB', 'VC', 'VB']) #for list of variables
instr.save('IC') #for single variable
```

**save\_var**(*trace\_list*)

Save the voltage or current in the instrument variable list. This is useful if one or two more variables need to be saved in addition to the 8 variables allowed by [save\(\)](#).

**Parameters** **trace\_list** – A list of channel variables whose measured data should be saved.

A maximum of 2 variables are allowed. If only one variable is being saved, a string can be specified.

```
instr.save_var(['VA', 'VB'])
```

**stop()**

Stops the ongoing measurement

```
instr.stop()
```

**class** `pymeasure.instruments.agilent.agilent4156.SMU`(*resourceName, channel, \*\*kwargs*)

Bases: [pymeasure.instruments.instrument.Instrument](#)

**property channel\_function**

A string property that controls the SMU<n> channel function.

- Values: VAR1, VAR2, VARD or CONS.

```
instr.smu1.channel_function = "VAR1"
```

**property channel\_mode**

A string property that controls the SMU<n> channel mode.

- Values: V, I or COMM

VPULSE AND IPULSE are not yet supported.

```
instr.smul.channel_mode = "V"
```

**property compliance**

This command sets the *constant* compliance value of SMU<n>. If the SMU channel is setup as a variable (VAR1, VAR2, VARD) then compliance limits are set by the variable definition.

- Value: Voltage in (-200V, 200V) and current in (-1A, 1A) based on [channel\\_mode\(\)](#).

```
instr.smul.compliance = 0.1
```

**property constant\_value**

This command sets the constant source value of SMU<n>. You use this command only if [channel\\_function\(\)](#) is CONS and also [channel\\_mode\(\)](#) should not be COMM.

**Parameters const\_value** – Voltage in (-200V, 200V) and current in (-1A, 1A). Voltage or current depends on if [channel\\_mode\(\)](#) is set to V or I.

```
instr.smul.constant_value = 1
```

**property current\_name**

Define the current name of the channel.

If input is greater than 6 characters long or starts with a number, the name is autocorrected and prepended with 'a'. Event is logged.

```
instr.smul.current_name = "Ibase"
```

**property disable**

This command deletes the settings of SMU<n>.

```
instr.smul.disable()
```

**property series\_resistance**

This command controls the series resistance of SMU<n>.

- Values: 00HM, 10KOHM, 100KOHM, or 1MOHM

```
instr.smul.series_resistance = "10KOHM"
```

**property voltage\_name**

Define the voltage name of the channel.

If input is greater than 6 characters long or starts with a number, the name is autocorrected and prepended with 'a'. Event is logged.

```
instr.smul.voltage_name = "Vbase"
```

**class** pymeasure.instruments.agilent.agilent4156.VAR1(*resourceName*, *\*\*kwargs*)

Bases: [pymeasure.instruments.agilent.agilent4156.VARX](#)

Class to handle all the specific definitions needed for VAR1. Most common methods are inherited from base class.

**property spacing**

This command selects the sweep type of VAR1.

- Values: LINEAR, LOG10, LOG25, LOG50.

**class** `pymeasure.instruments.agilent.agilent4156.VAR2(resource_name, **kwargs)`

Bases: `pymeasure.instruments.agilent.agilent4156.VARX`

Class to handle all the specific definitions needed for VAR2. Common methods are imported from base class.

**property points**

This command sets the number of sweep steps of VAR2. You use this command only if there is an SMU or VSU whose function (FCTN) is VAR2.

```
instr.var2.points = 10
```

**class** `pymeasure.instruments.agilent.agilent4156.VARD(resource_name, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Class to handle all the definitions needed for VARD. VARD is always defined in relation to VAR1.

**property compliance**

This command sets the sweep COMPLIANCE value of VARD.

```
instr.vard.compliance = 0.1
```

**property offset**

This command sets the OFFSET value of VARD. For each step of sweep, the output values of VAR1' are determined by the following equation:  $VARD = VAR1 \times RATIO + OFFSET$ . You use this command only if there is an SMU or VSU whose function is VARD.

```
instr.vard.offset = 1
```

**property ratio**

This command sets the RATIO of VAR1'. For each step of sweep, the output values of VAR1' are determined by the following equation:  $VAR1' = VAR1 \times RATIO + OFFSET$ . You use this command only if there is an SMU or VSU whose function (FCTN) is VAR1'.

```
instr.vard.ratio = 1
```

**class** `pymeasure.instruments.agilent.agilent4156.VARX(resource_name, var_name, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Base class to define sweep variable settings

**property compliance**

Sets the sweep COMPLIANCE value.

```
instr.var1.compliance = 0.1
```

**property start**

Sets the sweep START value.

```
instr.var1.start = 0
```

**property step**

Sets the sweep STEP value.

```
instr.var1.step = 0.1
```

**property stop**

Sets the sweep STOP value.

```
instr.var1.stop = 3
```

**class** `pymeasure.instruments.agilent.agilent4156.VMU(resourceName, channel, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

**property channel\_mode**

A string property that controls the VMU<n> channel mode.

- Values: V, DVOL

**property disable**

This command disables the settings of VMU<n>.

```
instr.vmu1.disable()
```

**property voltage\_name**

Define the voltage name of the VMU channel.

If input is greater than 6 characters long or starts with a number, the name is autocorrected and prepended with 'a'. Event is logged.

```
instr.vmu1.voltage_name = "Vanode"
```

**class** `pymeasure.instruments.agilent.agilent4156.VSU(resourceName, channel, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

**property channel\_function**

A string property that controls the VSU channel function.

- Value: VAR1, VAR2, VARD or CONS.

**property channel\_mode**

Get channel mode of VSU<n>.

**property constant\_value**

This command sets the constant source value of VSU<n>.

```
instr.vsu1.constant_value = 0
```

**property disable**

This command deletes the settings of VSU<n>.

```
instr.vsu1.disable()
```

**property voltage\_name**

Define the voltage name of the VSU channel

If input is greater than 6 characters long or starts with a number, the name is autocorrected and prepended with 'a'. Event is logged.

```
instr.vsu1.voltage_name = "Ve"
```



## 7.6.8 Agilent 33220A Arbitrary Waveform Generator

**class** `pymeasure.instruments.agilent.Agilent33220A`(*adapter*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Agilent 33220A Arbitrary Waveform Generator.

```
# Default channel for the Agilent 33220A
wfg = Agilent33220A("GPIB::10")

wfg.shape = "SINUSOID"           # Sets a sine waveform
wfg.frequency = 4.7e3             # Sets the frequency to 4.7 kHz
wfg.amplitude = 1                 # Set amplitude of 1 V
wfg.offset = 0                   # Set the amplitude to 0 V

wfg.burst_state = True           # Enable burst mode
wfg.burst_ncycles = 10           # A burst will consist of 10 cycles
wfg.burst_mode = "TRIGGERED"     # A burst will be applied on a trigger
wfg.trigger_source = "BUS"       # A burst will be triggered on TRG*

wfg.output = True                # Enable output of waveform generator
wfg.trigger()                    # Trigger a burst
wfg.wait_for_trigger()           # Wait until the triggering is finished
wfg.beep()                       # "beep"

print(wfg.check_errors())        # Get the error queue
```

### property `amplitude`

A floating point property that controls the voltage amplitude of the output waveform in V, from 10e-3 V to 10 V. Can be set.

### property `amplitude_unit`

A string property that controls the units of the amplitude. Valid values are Vpp (default), Vrms, and dBm. Can be set.

### `beep()`

Causes a system beep.

### property `beeper_state`

A boolean property that controls the state of the beeper. Can be set.

### property `burst_mode`

A string property that controls the burst mode. Valid values are: TRIG<GERED>, GAT<ED>. This setting can be set.

### property `burst_ncycles`

An integer property that sets the number of cycles to be output when a burst is triggered. Valid values are 1 to 50000. This can be set.

### property `burst_state`

A boolean property that controls whether the burst mode is on (True) or off (False). Can be set.

### property `frequency`

A floating point property that controls the frequency of the output waveform in Hz, from 1e-6 (1 uHz) to 20e+6 (20 MHz), depending on the specified function. Can be set.

### property `offset`

A floating point property that controls the voltage offset of the output waveform in V, from 0 V to 4.995 V,

depending on the set voltage amplitude (maximum offset =  $(10 - \text{voltage}) / 2$ ). Can be set.

**property output**

A boolean property that turns on (True) or off (False) the output of the function generator. Can be set.

**property pulse\_dutycycle**

A floating point property that controls the duty cycle of a pulse waveform function in percent. Can be set.

**property pulse\_hold**

A string property that controls if either the pulse width or the duty cycle is retained when changing the period or frequency of the waveform. Can be set to: WIDT<H> or DCYC<LE>.

**property pulse\_period**

A floating point property that controls the period of a pulse waveform function in seconds, ranging from 200 ns to 2000 s. Can be set and overwrites the frequency for *all* waveforms. If the period is shorter than the pulse width + the edge time, the edge time and pulse width will be adjusted accordingly.

**property pulse\_transition**

A floating point property that controls the the edge time in seconds for both the rising and falling edges. It is defined as the time between 0.1 and 0.9 of the threshold. Valid values are between 5 ns to 100 ns. The transition time has to be smaller than  $0.625 * \text{the pulse width}$ . Can be set.

**property pulse\_width**

A floating point property that controls the width of a pulse waveform function in seconds, ranging from 20 ns to 2000 s, within a set of restrictions depending on the period. Can be set.

**property ramp\_symmetry**

A floating point property that controls the symmetry percentage for the ramp waveform. Can be set.

**property remote\_local\_state**

A string property that controls the remote/local state of the function generator. Valid values are: LOC<AL>, REM<OTE>, RWL<OCK>. This setting can only be set.

**property shape**

A string property that controls the output waveform. Can be set to: SIN<USOID>, SQU<ARE>, RAMP, PULS<E>, NOIS<E>, DC, USER.

**property square\_dutycycle**

A floating point property that controls the duty cycle of a square waveform function in percent. Can be set.

**trigger()**

Send a trigger signal to the function generator.

**property trigger\_source**

A string property that controls the trigger source. Valid values are: IMM<EDIATE> (internal), EXT<ERNAL> (rear input), BUS (via trigger command). This setting can be set.

**property trigger\_state**

A boolean property that controls whether the output is triggered (True) or not (False). Can be set.

**property voltage\_high**

A floating point property that controls the upper voltage of the output waveform in V, from -4.990 V to 5 V (must be higher than low voltage). Can be set.

**property voltage\_low**

A floating point property that controls the lower voltage of the output waveform in V, from -5 V to 4.990 V (must be lower than high voltage). Can be set.

**wait\_for\_trigger**(*timeout=3600, should\_stop=<function Agilent33220A.<lambda>>*)

Wait until the triggering has finished or timeout is reached.

**Parameters**

- **timeout** – The maximum time the waiting is allowed to take. If timeout is exceeded, a `TimeoutError` is raised. If timeout is set to zero, no timeout will be used.
- **should\_stop** – Optional function (returning a bool) to allow the waiting to be stopped before its end.

### 7.6.9 Agilent 33500 Function/Arbitrary Waveform Generator Family

**class** `pymeasure.instruments.agilent.Agilent33500`(*adapter*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Agilent 33500 Function/Arbitrary Waveform Generator family. Individual devices are represented by subclasses.

```
generator = Agilent33500("GPIB::1")

generator.shape = 'SIN'           # Sets the output signal shape to sine
generator.frequency = 1e3         # Sets the output frequency to 1 kHz
generator.amplitude = 5           # Sets the output amplitude to 5 Vpp
generator.output = 'on'          # Enables the output

generator.shape = 'ARB'           # Set shape to arbitrary
generator.arb_srate = 1e6         # Set sample rate to 1MSa/s

generator.data_volatile_clear()   # Clear volatile internal memory
generator.data_arb(               # Send data points of arbitrary waveform
    'test',
    range(-10000, 10000, +20),    # In this case a simple ramp
    data_format='DAC'            # Data format is set to 'DAC'
)
generator.arb_file = 'test'       # Select the transmitted waveform 'test'
```

#### **property amplitude**

A floating point property that controls the voltage amplitude of the output waveform in V, from 10e-3 V to 10 V. Depends on the output impedance. Can be set.

#### **property amplitude\_unit**

A string property that controls the units of the amplitude. Valid values are VPP (default), VRMS, and DBM. Can be set.

#### **property arb\_advance**

A string property that selects how the device advances from data point to data point. Can be set to 'TRIG<GER>' or 'SRAT<E>' (default).

#### **property arb\_file**

A string property that selects the arbitrary signal from the volatile memory of the device. String has to match an existing arb signal in volatile memore (set by `data_arb()`). Can be set.

#### **property arb\_filter**

A string property that selects the filter setting for arbitrary signals. Can be set to 'NORM<AL>', 'STEP' and 'OFF'.

#### **property arb\_srate**

An floating point property that sets the sample rate of the currently selected arbitrary signal. Valid values are 1  $\mu$ Sa/s to 250 MSa/s (maximum range, can be lower depending on your device). This can be set.

**beep()**

Causes a system beep.

**property burst\_mode**

A string property that controls the burst mode. Valid values are: TRIG<GERED>, GAT<ED>. This setting can be set.

**property burst\_ncycles**

An integer property that sets the number of cycles to be output when a burst is triggered. Valid values are 1 to 100000. This can be set.

**property burst\_period**

A floating point property that controls the period of subsequent bursts. Has to follow the equation  $\text{burst\_period} > (\text{burst\_ncycles} / \text{frequency}) + 1 \mu\text{s}$ . Valid values are 1  $\mu\text{s}$  to 8000 s. Can be set.

**property burst\_state**

A boolean property that controls whether the burst mode is on (True) or off (False). Can be set.

**clear\_display()**

Removes a text message from the display.

**data\_arb(arb\_name, data\_points, data\_format='DAC')**

Uploads an arbitrary trace into the volatile memory of the device. The data\_points can be given as comma separated 16 bit DAC values (ranging from -32767 to +32767), as comma separated floating point values (ranging from -1.0 to +1.0) or as a binary data stream. Check the manual for more information. The storage depends on the device type and ranges from 8 Sa to 16 MSa (maximum). *TODO: Binary is not yet implemented*

**Parameters**

- **arb\_name** – The name of the trace in the volatile memory. This is used to access the trace.
- **data\_points** – Individual points of the trace. The format depends on the format parameter.  
  
format = 'DAC' (default): Accepts list of integer values ranging from -32767 to +32767. Minimum of 8 a maximum of 65536 points.  
  
format = 'float': Accepts list of floating point values ranging from -1.0 to +1.0. Minimum of 8 a maximum of 65536 points.  
  
format = 'binary': Accepts a binary stream of 8 bit data.
- **data\_format** – Defines the format of data\_points. Can be 'DAC' (default), 'float' or 'binary'. See documentation on parameter data\_points above.

**data\_volatile\_clear()**

Clear all arbitrary signals from the volatile memory. This should be done if the same name is used continuously to load different arbitrary signals into the memory, since an error will occur if a trace is loaded which already exists in the memory.

**property display**

A string property which is displayed on the front panel of the device. Can be set.

**property ext\_trig\_out**

A boolean property that controls whether the trigger out signal is active (True) or not (False). This signal is output from the Ext Trig connector on the rear panel in Burst and Wobbel mode. Can be set.

**property frequency**

A floating point property that controls the frequency of the output waveform in Hz, from 1 uHz to 120 MHz (maximum range, can be lower depending on your device), depending on the specified function. Can be set.

**property id**

Reads the instrument identification

**property offset**

A floating point property that controls the voltage offset of the output waveform in V, from 0 V to 4.995 V, depending on the set voltage amplitude (maximum offset =  $(V_{\text{max}} - \text{voltage}) / 2$ ). Can be set.

**property output**

A boolean property that turns on (True, 'on') or off (False, 'off') the output of the function generator. Can be set.

**property output\_load**

Sets the expected load resistance (should be the load impedance connected to the output). The output impedance is always 50 Ohm, this setting can be used to correct the displayed voltage for loads unmatched to 50 Ohm. Valid values are between 1 and 10 kOhm or INF for high impedance. No validator is used since both numeric and string inputs are accepted, thus a value outside the range will not return an error. Can be set.

**property phase**

A floating point property that controls the phase of the output waveform in degrees, from -360 degrees to 360 degrees. Not available for arbitrary waveforms or noise. Can be set.

**property pulse\_dutycycle**

A floating point property that controls the duty cycle of a pulse waveform function in percent, from 0% to 100%. Can be set.

**property pulse\_hold**

A string property that controls if either the pulse width or the duty cycle is retained when changing the period or frequency of the waveform. Can be set to: WIDT<H> or DCYC<LE>.

**property pulse\_period**

A floating point property that controls the period of a pulse waveform function in seconds, ranging from 33 ns to 1e6 s. Can be set and overwrites the frequency for *all* waveforms. If the period is shorter than the pulse width + the edge time, the edge time and pulse width will be adjusted accordingly.

**property pulse\_transition**

A floating point property that controls the edge time in seconds for both the rising and falling edges. It is defined as the time between the 10% and 90% thresholds of the edge. Valid values are between 8.4 ns to 1 µs. Can be set.

**property pulse\_width**

A floating point property that controls the width of a pulse waveform function in seconds, ranging from 16 ns to 1e6 s, within a set of restrictions depending on the period. Can be set.

**property ramp\_symmetry**

A floating point property that controls the symmetry percentage for the ramp waveform, from 0.0% to 100.0% Can be set.

**property shape**

A string property that controls the output waveform. Can be set to: SIN<USOID>, SQU<ARE>, TRI<ANGLE>, RAMP, PULS<E>, PRBS, NOIS<E>, ARB, DC.

**property square\_dutycycle**

A floating point property that controls the duty cycle of a square waveform function in percent, from 0.01% to 99.98%. The duty cycle is limited by the frequency and the minimal pulse width of 16 ns. See manual for more details. Can be set.

**trigger()**

Send a trigger signal to the function generator.

**property trigger\_source**

A string property that controls the trigger source. Valid values are: IMM<EDIATE> (internal), EXT<ERNAL> (rear input), BUS (via trigger command). This setting can be set.

**property voltage\_high**

A floating point property that controls the upper voltage of the output waveform in V, from -4.990 V to 5 V (must be higher than low voltage by at least 1 mV). Can be set.

**property voltage\_low**

A floating point property that controls the lower voltage of the output waveform in V, from -5 V to 4.990 V (must be lower than high voltage by at least 1 mV). Can be set.

**wait\_for\_trigger**(*timeout=3600*, *should\_stop=<function Agilent33500.<lambda>>>*)

Wait until the triggering has finished or timeout is reached.

**Parameters**

- **timeout** – The maximum time the waiting is allowed to take. If timeout is exceeded, a `TimeoutError` is raised. If timeout is set to zero, no timeout will be used.
- **should\_stop** – Optional function (returning a bool) to allow the waiting to be stopped before its end.

## 7.6.10 Agilent 33521A Function/Arbitrary Waveform Generator

**class** `pymeasure.instruments.agilent.Agilent33521A`(*adapter*, *\*\*kwargs*)

Bases: `pymeasure.instruments.agilent.agilent33500.Agilent33500`

Represents the Agilent 33521A Function/Arbitrary Waveform Generator. This documentation page shows only methods different from the parent class `Agilent33500`.

**property arb\_srate**

An floating point property that sets the sample rate of the currently selected arbitrary signal. Valid values are 1  $\mu$ Sa/s to 250 MSa/s. This can be set.

**property frequency**

A floating point property that controls the frequency of the output waveform in Hz, from 1 uHz to 30 MHz, depending on the specified function. Can be set.

## 7.6.11 Agilent B1500 Semiconductor Parameter Analyzer

**Contents**

- *Agilent B1500 Semiconductor Parameter Analyzer*
  - *General Information*
    - \* *Command Translation*
  - *Examples*
    - \* *Initialization of the Instrument*
    - \* *IV measurement with 4 SMUs*
    - \* *Sampling measurement with 4 SMUs*
  - *Main Classes*

- *Supporting Classes*
- \* *Enumerations*

## General Information

This instrument driver does not support all configuration options of the B1500 mainframe yet. So far, it is possible to interface multiple SMU modules and source/measure currents and voltages, perform sampling and staircase sweep measurements. The implementation of further measurement functionalities is highly encouraged. Meanwhile the model is managed by Keysight, see the corresponding “Programming Guide” for details on the control methods and their parameters

## Command Translation

Alphabetical list of implemented B1500 commands and their corresponding method/attribute names in this instrument driver.

Command	Property/Method
AAD	<i>SMU.adc_type()</i>
AB	<i>abort()</i>
AIT	<i>adc_setup()</i>
AV	<i>adc_averaging()</i>
AZ	<i>adc_auto_zero</i>
BC	<i>clear_buffer()</i>
CL	<i>SMU.disable()</i>
CM	<i>auto_calibration</i>
CMM	<i>SMU.meas_op_mode()</i>
CN	<i>SMU.enable()</i>
DI	<i>SMU.force()</i> mode: 'CURRENT'
DV	<i>SMU.force()</i> mode: 'VOLTAGE'
DZ	<i>force_gnd()</i> , <i>SMU.force_gnd()</i>
ERRX?	<i>check_errors()</i>
FL	<i>SMU.filter</i>
FMT	<i>data_format()</i>
*IDN?	<i>id()</i>
*LRN?	<i>query_learn()</i> , multiple methods to read/format settings directly
MI	<i>SMU.sampling_source()</i> mode: 'CURRENT'
ML	<i>sampling_mode</i>
MM	<i>meas_mode()</i>
MSC	<i>sampling_auto_abort()</i>
MT	<i>sampling_timing()</i>
MV	<i>SMU.sampling_source()</i> mode: 'VOLTAGE'
*OPC?	<i>check_idle()</i>
PA	<i>pause()</i>
PAD	<i>parallel_meas</i>
RI	<i>meas_range_current</i>
RM	<i>SMU.meas_range_current_auto()</i>
*RST	<i>reset()</i>
RV	<i>meas_range_voltage</i>

continues on next page

Table 1 – continued from previous page

Command	Property/Method
SSR	<code>series_resistor</code>
TSC	<code>time_stamp</code>
TSR	<code>clear_timer()</code>
UNT?	<code>query_modules()</code>
WAT	<code>wait_time()</code>
WI	<code>SMU.staircase_sweep_source()</code> mode: 'CURRENT'
WM	<code>sweep_auto_abort()</code>
WSI	<code>SMU.synchronous_sweep_source()</code> mode: 'CURRENT'
WSV	<code>SMU.synchronous_sweep_source()</code> mode: 'VOLTAGE'
WT	<code>sweep_timing()</code>
WV	<code>SMU.staircase_sweep_source()</code> mode: 'VOLTAGE'
XE	<code>send_trigger()</code>

## Examples

### Initialization of the Instrument

```
from pymeasure.instruments.agilent import AgilentB1500

# explicitly define r/w terminations; set sufficiently large timeout in milliseconds or
↳ None.
b1500=AgilentB1500("GPIB0::17::INSTR", read_termination='\r\n', write_termination='\r\n',
↳ timeout=6000000)
# query SMU config from instrument and initialize all SMU instances
b1500.initialize_all_smus()
# set data output format (required!)
b1500.data_format(21, mode=1) #call after SMUs are initialized to get names for the
↳ channels
```

### IV measurement with 4 SMUs

```
# choose measurement mode
b1500.meas_mode('STAIRCASE_SWEEP', *b1500.smu_references) #order in smu_references
↳ determines order of measurement

# settings for individual SMUs
for smu in b1500.smu_references:
    smu.enable() #enable SMU
    smu.adc_type = 'HRADC' #set ADC to high-resolution ADC
    smu.meas_range_current = '1 nA'
    smu.meas_op_mode = 'COMPLIANCE_SIDE' # other choices: Current, Voltage, FORCE_SIDE,
↳ COMPLIANCE_AND_FORCE_SIDE

# General Instrument Settings
# b1500.adc_averaging = 1
# b1500.adc_auto_zero = True
b1500.adc_setup('HRADC', 'AUTO', 6)
```

(continues on next page)



(continued from previous page)

```

#b1500.adc_setup('HRADC','PLC',1)

#Sweep Settings
b1500.sweep_timing(0,5,step_delay=0.1) #hold,delay
b1500.sweep_auto_abort(False,post='STOP') #disable auto abort, set post measurement
↳output condition to stop value of sweep
# Sweep Source
nop = 11
b1500.smu1.staircase_sweep_source('VOLTAGE','LINEAR_DOUBLE','Auto Ranging',0,1,nop,0.
↳001) #type, mode, range, start, stop, steps, compliance
# Synchronous Sweep Source
b1500.smu2.synchronous_sweep_source('VOLTAGE','Auto Ranging',0,1,0.001) #type, range,
↳start, stop, comp
# Constant Output (could also be done using synchronous sweep source with start=stop,
↳but then the output is not ramped up)
b1500.smu3.ramp_source('VOLTAGE','Auto Ranging',-1,stepsize=0.1,pause=20e-3) #output
↳starts immediately! (compared to sweeps)
b1500.smu4.ramp_source('VOLTAGE','Auto Ranging',0,stepsize=0.1,pause=20e-3)

#Start Measurement
b1500.check_errors()
b1500.clear_buffer()
b1500.clear_timer()
b1500.send_trigger()

# read measurement data all at once
b1500.check_idle() #wait until measurement is finished
data = b1500.read_data(2*nop) #Factor 2 beacuse of double sweep

#alternatively: read measurement data live
meas = []
for i in range(nop*2):
    read_data = b1500.read_channels(4+1) # 4 measurement channels, 1 sweep source
↳(returned due to mode=1 of data_format)
    # process live data for plotting etc.
    # data format for every channel (status code, channel name e.g. 'SMU1', data name e.g
↳'Current Measurement (A)', value)
    meas.append(read_data)

#sweep constant sources back to 0V
b1500.smu3.ramp_source('VOLTAGE','Auto Ranging',0,stepsize=0.1,pause=20e-3)
b1500.smu4.ramp_source('VOLTAGE','Auto Ranging',0,stepsize=0.1,pause=20e-3)

```

## Sampling measurement with 4 SMUs

```

# Choose measurement mode
b1500.meas_mode('SAMPLING', *b1500.smu_references) #order in smu_references determines
↳ order of measurement
number_of_channels = len(b1500.smu_references)

# settings for individual SMUs
for smu in b1500.smu_references:
    smu.enable() #enable SMU
    smu.adc_type = 'HSADC' #set ADC to high-speed ADC
    smu.meas_range_current = '1 nA'
    smu.meas_op_mode = 'COMPLIANCE_SIDE' # other choices: Current, Voltage, FORCE_SIDE,
↳ COMPLIANCE_AND_FORCE_SIDE

b1500.sampling_mode = 'LINEAR'
# b1500.adc_averaging = 1
# b1500.adc_auto_zero = True
b1500.adc_setup('HSADC', 'AUTO', 1)
#b1500.adc_setup('HSADC', 'PLC', 1)
nop=11
b1500.sampling_timing(2,0.005,nop) #MT: bias hold time, sampling interval, number of
↳ points
b1500.sampling_auto_abort(False,post='BIAS') #MSC: BASE/BIAS
b1500.time_stamp = True

# Sources
b1500.smu1.sampling_source('VOLTAGE', 'Auto Ranging', 0, 1, 0.001) #MV/MI: type, range, base,
↳ bias, compliance
b1500.smu2.sampling_source('VOLTAGE', 'Auto Ranging', 0, 1, 0.001)
b1500.smu3.ramp_source('VOLTAGE', 'Auto Ranging', -1, stepsize=0.1, pause=20e-3) #output
↳ starts immediately! (compared to sweeps)
b1500.smu4.ramp_source('VOLTAGE', 'Auto Ranging', -1, stepsize=0.1, pause=20e-3)

#Start Measurement
b1500.check_errors()
b1500.clear_buffer()
b1500.clear_timer()
b1500.send_trigger()

meas=[]
for i in range(nop):
    read_data = b1500.read_channels(1+2*number_of_channels) #Sampling Index + (time
↳ stamp + measurement value) * number of channels
    # process live data for plotting etc.
    # data format for every channel (status code, channel name e.g. 'SMU1', data name e.g
↳ 'Current Measurement (A)', value)
    meas.append(read_data)

#sweep constant sources back to 0V
b1500.smu3.ramp_source('VOLTAGE', 'Auto Ranging', 0, stepsize=0.1, pause=20e-3)
b1500.smu4.ramp_source('VOLTAGE', 'Auto Ranging', 0, stepsize=0.1, pause=20e-3)

```

## Main Classes

Classes to communicate with the instrument:

- [\*AgilentB1500\*](#): Main instrument class
- [\*SMU\*](#): Instantiated by main instrument class for every SMU

All *query* commands return a human readable dict of settings. These are intended for debugging/logging/file headers, not for passing to the accompanying setting commands.

**class** `pymeasure.instruments.agilent.agilentB1500.AgilentB1500`(*resourceName*, *\*\*kwargs*)

Bases: [`pymeasure.instruments.instrument.Instrument`](#)

Represents the Agilent B1500 Semiconductor Parameter Analyzer and provides a high-level interface for taking different kinds of measurements.

**property** `smu_references`

Returns all SMU instances.

**property** `smu_names`

Returns all SMU names.

**query\_learn**(*query\_type*)

Queries settings from the instrument (\*LRN?). Returns dict of settings.

**Parameters** `query_type` (*int* or *str*) – Query type (number according to manual)

**query\_learn\_header**(*query\_type*, *\*\*kwargs*)

Queries settings from the instrument (\*LRN?). Returns dict of settings in human readable format for debugging or file headers. For optional arguments check the underlying definition of [`QueryLearn.query\_learn\_header\(\)`](#).

**Parameters** `query_type` (*int* or *str*) – Query type (number according to manual)

**reset**()

Resets the instrument to default settings (\*RST)

**query\_modules**()

Queries module models from the instrument. Returns dictionary of channel and module type.

**Returns** Channel:Module Type

**Return type** dict

**initialize\_smu**(*channel*, *smu\_type*, *name*)

Initializes SMU instance by calling [\*SMU\*](#).

**Parameters**

- **channel** (*int*) – SMU channel
- **smu\_type** (*str*) – SMU type, e.g. 'HRSMU'
- **name** (*str*) – SMU name for pymeasure (data output etc.)

**Returns** SMU instance

**Return type** [\*SMU\*](#)

**initialize\_all\_smus**()

Initialize all SMUs by querying available modules and creating a SMU class instance for each. SMUs are accessible via attributes `.smu1` etc.

**pause**(*pause\_seconds*)

Pauses Command Execution for given time in seconds (PA)

**Parameters** `pause_seconds` (*int*) – Seconds to pause

**abort()**

Aborts the present operation but channels may still output current/voltage (AB)

**force\_gnd()**

Force 0V on all channels immediately. Current Settings can be restored with RZ. (DZ)

**check\_errors()**

Check for errors (ERRX?)

**check\_idle()**

Check if instrument is idle (\*OPC?)

**clear\_buffer()**

Clear output data buffer (BC)

**clear\_timer()**

Clear timer count (TSR)

**send\_trigger()**

Send trigger to start measurement (except High Speed Spot) (XE)

**property auto\_calibration**

Enable/Disable SMU auto-calibration every 30 minutes. (CM)

**Type** bool

**data\_format**(*output\_format*, *mode=0*)

Specifies data output format. Check Documentation for parameters. Should be called once per session to set the data format for interpreting the measurement values read from the instrument. (FMT)

Currently implemented are format 1, 11, and 21.

**Parameters**

- **output\_format** (*str*) – Output format string, e.g. FMT21
- **mode** (*int*, *optional*) – Data output mode, defaults to 0 (only measurement data is returned)

**property parallel\_meas**

**Enable/Disable parallel measurements.** Effective for SMUs using HSADC and measurement modes 1,2,10,18. (PAD)

**Type** bool

**query\_meas\_settings()**

Read settings for TM, AV, CM, FMT and MM commands (31) from the instrument.

**query\_meas\_mode()**

Read settings for MM command (part of 31) from the instrument.

**meas\_mode**(*mode*, *\*args*)

Set Measurement mode of channels. Measurements will be taken in the same order as the SMU references are passed. (MM)

**Parameters**

- **mode** (*MeasMode*) – Measurement mode
  - Spot
  - Staircase Sweep

- Sampling

- **args** (*SMU*) – SMU references

#### **query\_adc\_setup()**

Read ADC settings (55, 56) from the instrument.

#### **adc\_setup**(*adc\_type, mode, N=""*)

Set up operation mode and parameters of ADC for each ADC type. (AIT) Defaults:

- HSADC: Auto N=1, Manual N=1, PLC N=1, Time N=0.000002(s)
- HRADC: Auto N=6, Manual N=3, PLC N=1

#### **Parameters**

- **adc\_type** (*ADCType*) – ADC type
- **mode** (*ADCMode*) – ADC mode
- **N** (*str, optional*) – additional parameter, check documentation, defaults to ''

#### **adc\_averaging**(*number, mode='Auto'*)

Set number of averaging samples of the HSADC. (AV)

Defaults: N=1, Auto

#### **Parameters**

- **number** (*int*) – Number of averages
- **mode** (*AutoManual*, optional) – Mode ('Auto', 'Manual'), defaults to 'Auto'

#### **property adc\_auto\_zero**

Enable/Disable ADC zero function. Halfs the integration time, if off. (AZ)

**Type** bool

#### **property time\_stamp**

Enable/Disable Time Stamp function. (TSC)

**Type** bool

#### **query\_time\_stamp\_setting()**

Read time stamp settings (60) from the instrument.

#### **wait\_time**(*wait\_type, N, offset=0*)

Configure wait time. (WAT)

#### **Parameters**

- **wait\_type** (*WaitTimeType*) – Wait time type
- **N** (*float*) – Coefficient for initial wait time, default: 1
- **offset** (*int, optional*) – Offset for wait time, defaults to 0

#### **query\_staircase\_sweep\_settings()**

Reads Staircase Sweep Measurement settings (33) from the instrument.

#### **sweep\_timing**(*hold, delay, step\_delay=0, step\_trigger\_delay=0, measurement\_trigger\_delay=0*)

Sets Hold Time, Delay Time and Step Delay Time for staircase or multi channel sweep measurement. (WT)

If not set, all parameters are 0.

#### **Parameters**

- **hold** (*float*) – Hold time

- **delay** (*float*) – Delay time
- **step\_delay** (*float*, *optional*) – Step delay time, defaults to 0
- **step\_trigger\_delay** (*float*, *optional*) – Trigger delay time, defaults to 0
- **measurement\_trigger\_delay** (*float*, *optional*) – Measurement trigger delay time, defaults to 0

**sweep\_auto\_abort** (*abort*, *post*='START')

Enables/Disables the automatic abort function. Also sets the post measurement condition. (WM)

**Parameters**

- **abort** (*bool*) – Enable/Disable automatic abort
- **post** (*StaircaseSweepPostOutput*, *optional*) – Output after measurement, defaults to 'Start'

**query\_sampling\_settings**()

Reads Sampling Measurement settings (47) from the instrument.

**property sampling\_mode**

Set linear or logarithmic sampling mode. (ML)

**Type** *SamplingMode*

**sampling\_timing** (*hold\_bias*, *interval*, *number*, *hold\_base*=0)

Sets Timing Parameters for the Sampling Measurement (MT)

**Parameters**

- **hold\_bias** (*float*) – Bias hold time
- **interval** (*float*) – Sampling interval
- **number** (*int*) – Number of Samples
- **hold\_base** (*float*, *optional*) – Base hold time, defaults to 0

**sampling\_auto\_abort** (*abort*, *post*='Bias')

Enables/Disables the automatic abort function. Also sets the post measurement condition. (MSC)

**Parameters**

- **abort** (*bool*) – Enable/Disable automatic abort
- **post** (*SamplingPostOutput*, *optional*) – Output after measurement, defaults to 'Bias'

**read\_data** (*number\_of\_points*)

Reads all data from buffer and returns Pandas DataFrame. Specify number of measurement points for correct splitting of the data list.

**Parameters** **number\_of\_points** (*int*) – Number of measurement points

**Returns** Measurement Data

**Return type** *pd.DataFrame*

**read\_channels** (*nchannels*)

Reads data for 1 measurement point from the buffer. Specify number of measurement channels + sweep sources (depending on data output setting).

**Parameters** **nchannels** (*int*) – Number of channels which return data

**Returns** Measurement data

**Return type** *tuple*

**query\_series\_resistor()**

Read series resistor status (53) for all SMUs.

**query\_meas\_range\_current\_auto()**

Read auto ranging mode status (54) for all SMUs.

**query\_meas\_op\_mode()**

Read SMU measurement operation mode (46) for all SMUs.

**query\_meas\_ranges()**

Read measruement ranging status (32) for all SMUs.

**class** pymeasure.instruments.agilent.agilentB1500.SMU(*parent, channel, smu\_type, name, \*\*kwargs*)

Bases: object

Provides specific methods for the SMUs of the Agilent B1500 mainframe

**Parameters**

- **parent** (*AgilentB1500*) – Instance of the B1500 mainframe class
- **channel** (*int*) – Channel number of the SMU
- **smu\_type** (*str*) – Type of the SMU
- **name** (*str*) – Name of the SMU

**write**(*string*)

Wraps *Instrument.write()* method of B1500.

**ask**(*string*)

Wraps *ask()* method of B1500.

**query\_learn**(*query\_type, command*)

Wraps *query\_learn()* method of B1500.

**check\_errors()**

Wraps *check\_errors()* method of B1500.

**property status**

Query status of the SMU.

**enable()**

Enable Source/Measurement Channel (CN)

**disable()**

Disable Source/Measurement Channel (CL)

**force\_gnd()**

Force 0V immediately. Current Settings can be restored with RZ (not implemented). (DZ)

**property filter**

Enables/Disables SMU Filter. (FL)

**Type** bool

**property series\_resistor**

Enables/Disables 1MOhm series resistor. (SSR)

**Type** bool

**property meas\_op\_mode**

Set SMU measurement operation mode. (CMM)

**Type** *MeasOpMode*

**property adc\_type**

ADC type of individual measurement channel. (AAD)

Type [ADCType](#)

**force**(*source\_type*, *source\_range*, *output*, *comp*="", *comp\_polarity*="", *comp\_range*="")

Applies DC Current or Voltage from SMU immediately. (DI, DV)

**Parameters**

- **source\_type** (*str*) – Source type ('Voltage', 'Current')
- **source\_range** (*int* or *str*) – Output range index or name
- **output** – Source output value in A or V
- **comp** (*float*, *optional*) – Compliance value, defaults to previous setting
- **comp\_polarity** ([CompliancePolarity](#)) – Compliance polarity, defaults to auto
- **comp\_range** (*int* or *str*, *optional*) – Compliance ranging type, defaults to auto

**ramp\_source**(*source\_type*, *source\_range*, *target\_output*, *comp*="", *comp\_polarity*="", *comp\_range*="", *stepsize*=0.001, *pause*=0.02)

Ramps to a target output from the set value with a given step size, each separated by a pause.

**Parameters**

- **source\_type** (*str*) – Source type ('Voltage' or 'Current')
- **target\_output** – Target output voltage or current
- **irange** (*int*) – Output range index
- **comp** (*float*, *optional*) – Compliance, defaults to previous setting
- **comp\_polarity** ([CompliancePolarity](#)) – Compliance polarity, defaults to auto
- **comp\_range** (*int* or *str*, *optional*) – Compliance ranging type, defaults to auto
- **stepsize** – Maximum size of steps
- **pause** – Duration in seconds to wait between steps

Type *target\_output*: float

**property meas\_range\_current**

Current measurement range index. (RI)

Possible settings depend on SMU type, e.g. 0 for Auto Ranging: [SMUCurrentRanging](#)

**property meas\_range\_voltage**

Voltage measurement range index. (RV)

Possible settings depend on SMU type, e.g. 0 for Auto Ranging: [SMUVoltageRanging](#)

**meas\_range\_current\_auto**(*mode*, *rate*=50)

Specifies the auto range operation. Check Documentation. (RM)

**Parameters**

- **mode** (*int*) – Range changing operation mode
- **rate** (*int*, *optional*) – Parameter used to calculate the *current* value, defaults to 50

**staircase\_sweep\_source**(*source\_type*, *mode*, *source\_range*, *start*, *stop*, *steps*, *comp*, *Pcomp*="")

Specifies Staircase Sweep Source (Current or Voltage) and its parameters. (WV or WI)

**Parameters**



- **source\_type** (*str*) – Source type ('Voltage', 'Current')
- **mode** (*SweepMode*) – Sweep mode
- **source\_range** (*int*) – Source range index
- **start** (*float*) – Sweep start value
- **stop** (*float*) – Sweep stop value
- **steps** (*int*) – Number of sweep steps
- **comp** (*float*) – Compliance value
- **Pcomp** (*float*, *optional*) – Power compliance, defaults to not set

**synchronous\_sweep\_source**(*source\_type*, *source\_range*, *start*, *stop*, *comp*, *Pcomp*=")

Specifies Synchronous Staircase Sweep Source (Current or Voltage) and its parameters. (WSV or WSI)

#### Parameters

- **source\_type** (*str*) – Source type ('Voltage', 'Current')
- **source\_range** (*int*) – Source range index
- **start** (*float*) – Sweep start value
- **stop** (*float*) – Sweep stop value
- **comp** (*float*) – Compliance value
- **Pcomp** (*float*, *optional*) – Power compliance, defaults to not set

**sampling\_source**(*source\_type*, *source\_range*, *base*, *bias*, *comp*)

Sets DC Source (Current or Voltage) for sampling measurement. DV/DI commands on the same channel overwrite this setting. (MV or MI)

#### Parameters

- **source\_type** (*str*) – Source type ('Voltage', 'Current')
- **source\_range** (*int*) – Source range index
- **base** (*float*) – Base voltage/current
- **bias** (*float*) – Bias voltage/current
- **comp** (*float*) – Compliance value

## Supporting Classes

Classes that provide additional functionalities:

- *QueryLearn*: Process read out of instrument settings
- *SMUCurrentRanging*, *SMUVoltageRanging*: Allowed ranges for different SMU types and transformation of range names to indices (base: *Ranging*)

**class** pymeasure.instruments.agilent.agilentB1500.QueryLearn

Bases: object

Methods to issue and process \*LRN? (learn) command and response.

**static query\_learn**(*ask*, *query\_type*)

Issues \*LRN? (learn) command to the instrument to read configuration. Returns dictionary of commands and set values.

**Parameters** `query_type` (*int*) – Query type according to the programming guide

**Returns** Dictionary of command and set values

**Return type** dict

**classmethod** `query_learn_header`(*ask, query\_type, smu\_references, single\_command=False*)

Issues \*LRN? (learn) command to the instrument to read configuration. Processes information to human readable values for debugging purposes or file headers.

**Parameters**

- `ask` (*Instrument.ask*) – ask method of the instrument
- `query_type` (*int or str*) – Number according to Programming Guide
- `smu_references` (*dict*) – SMU references by channel
- `single_command` (*str*) – if only a single command should be returned, defaults to False

**Returns** Read configuration

**Return type** dict

**static** `to_dict`(*parameters, names, \*args*)

Takes parameters returned by `query_learn()` and ordered list of corresponding parameter names (optional function) and returns dict of parameters including names.

**Parameters**

- `parameters` (*dict*) – Parameters for one command returned by `query_learn()`
- `names` (*list*) – list of names or (name, function) tuples, ordered

**Returns** Parameter name and (processed) parameter

**Return type** dict

**class** `pymeasure.instruments.agilent.agilentB1500.Ranging`(*supported\_ranges, ranges, fixed\_ranges=False*)

Bases: object

Possible Settings for SMU Current/Voltage Output/Measurement ranges. Transformation of available Voltage/Current Range Names to Index and back.

**Parameters**

- `supported_ranges` (*list*) – Ranges which are supported (list of range indices)
- `ranges` (*dict*) – All range names {Name: Indices}
- `fixed_ranges` – add fixed ranges (negative indices); defaults to False

**\_\_call\_\_**(*input\_value*)

Gives named tuple (name/index) of given Range. Throws error if range is not supported by this SMU.

**Parameters** `input` (*str or int*) – Range name or index

**Returns** named tuple (name/index) of range

**Return type** namedtuple

**class** `pymeasure.instruments.agilent.agilentB1500.SMUCurrentRanging`(*smu\_type*)

Bases: object

Provides Range Name/Index transformation for current measurement/sourcing. Validity of ranges is checked against the type of the SMU.

Omitting the ‘limited auto ranging’/‘range fixed’ specification in the range string for current measurement defaults to ‘limited auto ranging’.

Full specification: ‘1 nA range fixed’ or ‘1 nA limited auto ranging’

‘1 nA’ defaults to ‘1 nA limited auto ranging’

**class** `pymeasure.instruments.agilent.agilentB1500.SMUVoltageRanging(smu_type)`

Bases: `object`

Provides Range Name/Index transformation for voltage measurement/sourcing. Validity of ranges is checked against the type of the SMU.

Omitting the ‘limited auto ranging’/‘range fixed’ specification in the range string for voltage measurement defaults to ‘limited auto ranging’.

Full specification: ‘2 V range fixed’ or ‘2 V limited auto ranging’

‘2 V’ defaults to ‘2 V limited auto ranging’

## Enumerations

Enumerations are used for easy selection of the available parameters (where it is applicable). Methods accept member name or number as input, but name is recommended for readability reasons. The member number is passed to the instrument. Converting an enumeration member into a string gives a title case, whitespace separated string (`__str__()`) which cannot be used to select an enumeration member again. It’s purpose is only logging or documentation.

**class** `pymeasure.instruments.agilent.agilentB1500.CustomIntEnum(value)`

Bases: `enum.IntEnum`

Provides additional methods to `IntEnum`:

- Conversion to string automatically replaces ‘\_’ with ‘ ’ in names and converts to title case
- get classmethod to get enum reference with name or integer

`__str__()`

Gives title case string of enum value

**classmethod** `get(input_value)`

Gives Enum member by specifying name or value.

**Parameters** `input_value` (`str` or `int`) – Enum name or value

**Returns** Enum member

**class** `pymeasure.instruments.agilent.agilentB1500.ADCType(value)`

Bases: `pymeasure.instruments.agilent.agilentB1500.CustomIntEnum`

ADC Type

**HSADC** = 0

High-speed ADC

**HRADC** = 1

High-resolution ADC

**HSADC\_PULSED** = 2

High-resolution ADC for pulsed measurements

**class** `pymeasure.instruments.agilent.agilentB1500.ADCMode(value)`

Bases: `pymeasure.instruments.agilent.agilentB1500.CustomIntEnum`

ADC Mode

**AUTO = 0**

**MANUAL = 1**

**PLC = 2**

**TIME = 3**

**class** pymeasure.instruments.agilent.agilentB1500.**AutoManual**(*value*)  
Bases: [pymeasure.instruments.agilent.agilentB1500.CustomIntEnum](#)  
Auto/Manual selection

**AUTO = 0**

**MANUAL = 1**

**class** pymeasure.instruments.agilent.agilentB1500.**MeasMode**(*value*)  
Bases: [pymeasure.instruments.agilent.agilentB1500.CustomIntEnum](#)  
Measurement Mode

**SPOT = 1**

**STAIRCASE\_SWEEP = 2**

**SAMPLING = 10**

**class** pymeasure.instruments.agilent.agilentB1500.**MeasOpMode**(*value*)  
Bases: [pymeasure.instruments.agilent.agilentB1500.CustomIntEnum](#)  
Measurement Operation Mode

**COMPLIANCE\_SIDE = 0**

**CURRENT = 1**

**VOLTAGE = 2**

**FORCE\_SIDE = 3**

**COMPLIANCE\_AND\_FORCE\_SIDE = 4**

**class** pymeasure.instruments.agilent.agilentB1500.**SweepMode**(*value*)  
Bases: [pymeasure.instruments.agilent.agilentB1500.CustomIntEnum](#)  
Sweep Mode

**LINEAR\_SINGLE = 1**

**LOG\_SINGLE = 2**

**LINEAR\_DOUBLE = 3**

**LOG\_DOUBLE = 4**

**class** pymeasure.instruments.agilent.agilentB1500.**SamplingMode**(*value*)  
Bases: [pymeasure.instruments.agilent.agilentB1500.CustomIntEnum](#)  
Sampling Mode

**LINEAR = 1**

**LOG\_10 = 2**  
Logarithmic 10 data points/decade

**LOG\_25 = 3**  
Logarithmic 25 data points/decade

```

LOG_50 = 4
    Logarithmic 50 data points/decade

LOG_100 = 5
    Logarithmic 100 data points/decade

LOG_250 = 6
    Logarithmic 250 data points/decade

LOG_5000 = 7
    Logarithmic 5000 data points/decade

```

```

class pymeasure.instruments.agilent.agilentB1500.SamplingPostOutput(value)
    Bases: pymeasure.instruments.agilent.agilentB1500.CustomIntEnum

    Output after sampling

    BASE = 1
    BIAS = 2

class pymeasure.instruments.agilent.agilentB1500.StaircaseSweepPostOutput(value)
    Bases: pymeasure.instruments.agilent.agilentB1500.CustomIntEnum

    Output after staircase sweep

    START = 1
    STOP = 2

class pymeasure.instruments.agilent.agilentB1500.CompliancePolarity(value)
    Bases: pymeasure.instruments.agilent.agilentB1500.CustomIntEnum

    Compliance polarity

    AUTO = 0
    MANUAL = 1

class pymeasure.instruments.agilent.agilentB1500.WaitTimeType(value)
    Bases: pymeasure.instruments.agilent.agilentB1500.CustomIntEnum

    Wait time type

    SMU_SOURCE = 1
    SMU_MEASUREMENT = 2
    CMU_MEASUREMENT = 3

```

## 7.7 Ametek

This section contains specific documentation on the Ametek instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.7.1 Ametek 7270 DSP Lockin Amplifier

**class** `pymeasure.instruments.ametek.Ametek7270(resource_name, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

This is the class for the Ametek DSP 7270 lockin amplifier

**property** `adc1`

Reads the input value of ADC1 in Volts

**property** `adc2`

Reads the input value of ADC2 in Volts

**property** `adc3`

Reads the input value of ADC3 in Volts

**property** `adc4`

Reads the input value of ADC4 in Volts

**property** `dac1`

A floating point property that represents the output value on DAC1 in Volts. This property can be set.

**property** `dac2`

A floating point property that represents the output value on DAC2 in Volts. This property can be set.

**property** `dac3`

A floating point property that represents the output value on DAC3 in Volts. This property can be set.

**property** `dac4`

A floating point property that represents the output value on DAC4 in Volts. This property can be set.

**property** `frequency`

A floating point property that represents the lock-in frequency in Hz. This property can be set.

**property** `harmonic`

An integer property that represents the reference harmonic mode control, taking values from 1 to 127. This property can be set.

**property** `id`

Reads the instrument identification

**property** `mag`

Reads the magnitude in Volts

**property** `phase`

A floating point property that represents the reference harmonic phase in degrees. This property can be set.

**property** `sensitivity`

A floating point property that controls the sensitivity range in Volts, which can take discrete values from 2 nV to 1 V. This property can be set.

**set\_channel\_A\_mode()**

Sets instrument to channel A mode – assuming it is in voltage mode

**set\_differential\_mode(line\_filtering=True)**

Sets instrument to differential mode – assuming it is in voltage mode

**set\_voltage\_mode()**

Sets instrument to voltage control mode

**shutdown()**

Ensures the instrument in a safe state

**property slope**

A integer property that controls the filter slope in dB/octave, which can take the values 6, 12, 18, or 24 dB/octave. This property can be set.

**property time\_constant**

A floating point property that controls the time constant in seconds, which takes values from 10 microseconds to 100,000 seconds. This property can be set.

**property voltage**

A floating point property that represents the voltage in Volts. This property can be set.

**property x**

Reads the X value in Volts

**property x1**

Reads the first harmonic X value in Volts

**property x2**

Reads the second harmonic X value in Volts

**property xy**

Reads both the X and Y values in Volts

**property y**

Reads the Y value in Volts

**property y1**

Reads the first harmonic Y value in Volts

**property y2**

Reads the second harmonic Y value in Volts

## 7.8 AMI

This section contains specific documentation on the AMI instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.8.1 AMI 430 Power Supply

**class** `pymeasure.instruments.ami.AMI430(resource_name, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the AMI 430 Power supply and provides a high-level for interacting with the instrument.

```
magnet = AMI430("TCPIP::web.address.com::7180::SOCKET")

magnet.coilconst = 1.182          # kGauss/A
magnet.voltage_limit = 2.2        # Sets the voltage limit in V

magnet.target_current = 10        # Sets the target current to 10 A
magnet.target_field = 1           # Sets target field to 1 kGauss

magnet.ramp_rate_current = 0.0357 # Sets the ramp rate in A/s
magnet.ramp_rate_field = 0.0422   # Sets the ramp rate in kGauss/s
magnet.ramp                       # Initiates the ramping
```

(continues on next page)

(continued from previous page)

```

magnet.pause           # Pauses the ramping
magnet.status          # Returns the status of the magnet

magnet.ramp_to_current(5) # Ramps the current to 5 A

magnet.shutdown()      # Ramps the current to zero and disables_
↪ output

```

**property coilconst**

A floating point property that sets the coil constant in kGauss/A.

**disable\_persistent\_switch()**

Disables the persistent switch.

**enable\_persistent\_switch()**

Enables the persistent switch.

**property field**

Reads the field in kGauss of the magnet.

**has\_persistent\_switch\_enabled()**

Returns a boolean if the persistent switch is enabled.

**property magnet\_current**

Reads the current in Amps of the magnet.

**pause()**

Pauses the ramping of the magnetic field.

**ramp()**

Initiates the ramping of the magnetic field to set current/field with ramping rate previously set.

**property ramp\_rate\_current**

A floating point property that sets the current ramping rate in A/s.

**property ramp\_rate\_field**

A floating point property that sets the field ramping rate in kGauss/s.

**ramp\_to\_current(current, rate)**

Heats up the persistent switch and ramps the current with set ramp rate.

**ramp\_to\_field(field, rate)**

Heats up the persistent switch and ramps the current with set ramp rate.

**shutdown(ramp\_rate=0.0357)**

Turns on the persistent switch, ramps down the current to zero, and turns off the persistent switch.

**property state**

Reads the field in kGauss of the magnet.

**property supply\_current**

Reads the current in Amps of the power supply.

**property target\_current**

A floating point property that sets the target current in A for the magnet.

**property target\_field**

A floating point property that sets the target field in kGauss for the magnet.

**property voltage\_limit**

A floating point property that sets the voltage limit for charging/discharging the magnet.



**wait\_for\_holding**(*should\_stop*=<function AMI430.<lambda>>, *timeout*=800, *interval*=0.1)

**zero**()

Initiates the ramping of the magnetic field to zero current/field with ramping rate previously set.

## 7.9 Anaheim Automation

This section contains specific documentation on the Anaheim Automation instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.9.1 DP-Series Step Motor Controller

The DPSeriesMotorController class implements a base driver class for Anaheim-Automation DP Series stepper motor controllers. There are many controllers sold in this series, all of which implement the same core command set. Some controllers, like the DPY50601, implement additional functionality that is not included in this driver. If these additional features are desired, they should be implemented in a subclass.

```
class pymeasure.instruments.anaheimautomation.DPSeriesMotorController(resourceName,
                                                                    address=0,
                                                                    encoder_enabled=False,
                                                                    **kwargs)
```

Bases: [pymeasure.instruments.instrument.Instrument](#)

Base class to interface with Anaheim Automation DP series stepper motor controllers. This driver has been tested with the DPY50601 and DPE25601 motor controllers.

**property absolute\_position**

Float property representing the value of the motor position measured in absolute units. Note that in DP series motor controller instrument manuals, *absolute position* refers to the 'step\_position' property rather than this property. Also note that use of this property relies on `steps_to_absolute()` and `absolute_to_steps()` being implemented in a subclass. In this way, the user can define the conversion from a motor step position into any desired absolute unit. Absolute units could be the position in meters of a linear stage or the angular position of a gimbal mount, etc. This property can be set.

**absolute\_to\_steps**(*pos*)

Convert an absolute position to a number of steps to move. This must be implemented in subclasses.

**Parameters** *pos* – Absolute position in the units determined by the subclassed `absolute_to_steps()` method.

**property address**

Integer property representing the address that the motor controller uses for serial communications.

**ask**(*command*)

Override the instrument base ask method to add the motor controller's address to the command string.

**Parameters** *command* – command string to be sent to the instrument

**property basespeed**

Integer property that represents the motor controller's starting/homing speed. This property can be set.

**property busy**

Query to see if the controller is currently moving a motor.

**check\_errors**()

Method to read the error codes register and log when an error is detected.

**Return error\_code** one byte with the error codes register contents

**property direction**

A string property that represents the direction in which the stepper motor will rotate upon subsequent step commands. This property can be set. 'CW' corresponds to clockwise rotation and 'CCW' corresponds to counter-clockwise rotation.

**property encoder\_autocorrect**

A boolean property to enable or disable the encoder auto correct function. This property can be set.

**property encoder\_delay**

An integer property that represents the wait time in ms. after a move is finished before the encoder is read for a potential encoder auto-correct action to take place. This property can be set.

**property encoder\_enabled**

A boolean property to represent whether an external encoder is connected and should be used to set the step\_position property.

**property encoder\_motor\_ratio**

An integer property that represents the ratio of the number of encoder pulses per motor step. This property can be set.

**property encoder\_retries**

An integer property that represents the number of times the motor controller will try the encoder auto correct function before setting an error flag. This property can be set.

**property encoder\_window**

An integer property that represents the allowable error in encoder pulses from the desired position before the encoder auto-correct function runs. This property can be set.

**property error\_reg**

Reads the current value of the error codes register.

**home(*home\_mode*)**

Send command to the motor controller to 'home' the motor.

**Parameters** *home\_mode* – 0 or 1 specifying which homing mode to run.

0 will perform a homing operation where the controller moves the motor until a soft limit is reached, then will ramp down to base speed and continue motion until a home limit is reached.

In mode 1, the controller will move the motor until a limit is reached, then will ramp down to base speed, change direction, and run until the limit is released.

**property maxspeed**

Integer property that represents the motor controller's maximum (running) speed. This property can be set.

**move(*direction*)**

**Move the stepper motor continuously in the given direction until a stop command is sent or a limit switch is reached.** This method corresponds to the 'slew' command in the DP series instrument manuals.

**Parameters** *direction* – value to set on the direction property before moving the motor.

**reset\_position()**

Reset the position as counted by the motor controller and an externally connected encoder to 0.

**property step\_position**

Integer property representing the value of the motor position measured in steps counted by the motor controller or, if encoder\_enabled is set, the steps counted by an externally connected encoder. Note that in the

DP series motor controller instrument manuals, this property would be referred to as the ‘absolute position’ while this driver implements a conversion between steps and absolute units for the *absolute position* property. This property can be set.

**steps\_to\_absolute**(*steps*)

Convert a position measured in steps to an absolute position.

**Parameters** **steps** – Position in steps to be converted to an absolute position.

**stop**()

Method that stops all motion on the motor controller.

**values**(*command*, *\*\*kwargs*)

Override the instrument base values method to add the motor controller’s address to the command string.

**Parameters** **command** – command string to be sent to the motor controller.

**wait\_for\_completion**(*interval=0.5*)

Block until the controller is not “busy” (i.e. block until the motor is no longer moving.)

**Parameters** **interval** – (float) seconds between queries to the “busy” flag.

**Returns** None

**write**(*command*)

Override the instrument base write method to add the motor controller’s address to the command string.

**Parameters** **command** – command string to be sent to the motor controller.

## 7.10 Anapico

This section contains specific documentation on the Anapico instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.10.1 Anapico APSIN12G Signal Generator

**class** `pymeasure.instruments.anapico.APSIN12G`(*resourceName*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Anapico APSIN12G Signal Generator with option 9K, HP and GPIB.

**property blanking**

A string property that represents the blanking of output power when frequency is changed. ON makes the output to be blanked (off) while changing frequency. This property can be set.

**disable\_rf**()

Disables the RF output.

**enable\_rf**()

Enables the RF output.

**property frequency**

A floating point property that represents the output frequency in Hz. This property can be set.

**property power**

A floating point property that represents the output power in dBm. This property can be set.

**property reference\_output**

A string property that represents the 10MHz reference output from the synth. This property can be set.

## 7.11 Anritsu

This section contains specific documentation on the Anritsu instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.11.1 Anritsu MG3692C Signal Generator

**class** `pymeasure.instruments.anritsu.AnritsuMG3692C(resource_name, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Anritsu MG3692C Signal Generator

**disable()**

Disables the signal output.

**enable()**

Enables the signal output.

**property frequency**

A floating point property that represents the output frequency in Hz. This property can be set.

**property output**

A boolean property that represents the signal output state. This property can be set to control the output.

**property power**

A floating point property that represents the output power in dBm. This property can be set.

**shutdown()**

Shuts down the instrument, putting it in a safe state.

### 7.11.2 Anritsu MS9710C Optical Spectrum Analyzer

**class** `pymeasure.instruments.anritsu.AnritsuMS9710C(adapter, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Anritsu MS9710C Optical Spectrum Analyzer.

**property analysis**

Analysis Control

**property analysis\_result**

Read back analysis result from current scan.

**property average\_point**

Number of averages to take on each point (2-1000), or OFF

**property average\_sweep**

Number of averages to make on a sweep (2-1000) or OFF

**center\_at\_peak(\*\*kwargs)**

Center the spectrum at the measured peak.

**property data\_memory\_a\_condition**

Returns the data condition of data memory register A. Starting wavelength, and a sampling point (l1, l2, n).

**property data\_memory\_a\_size**

Returns the number of points sampled in data memory register A.

**property data\_memory\_a\_values**

Reads the binary data from memory register A.

**property data\_memory\_b\_condition**

Returns the data condition of data memory register B. Starting wavelength, and a sampling point (l1, l2, n).

**property data\_memory\_b\_size**

Returns the number of points sampled in data memory register B.

**property data\_memory\_b\_values**

Reads the binary data from memory register B.

**property data\_memory\_select**

Memory Data Select.

**property dip\_search**

Dip Search Mode

**property ese2**

Extended Event Status Enable Register 2

**property esr2**

Extended Event Status Register 2

**property level\_lin**

Level Linear Scale (/div)

**property level\_log**

Level Log Scale (/div)

**property level\_opt\_attn**

Optical Attenuation Status (ON/OFF)

**property level\_scale**

Current Level Scale

**property measure\_mode**

Returns the current Measure Mode the OSA is in.

**measure\_peak()**

Measure the peak and return the trace marker.

**property peak\_search**

Peak Search Mode

**read\_memory(slot='A')**

Read the scan saved in a memory slot.

**property resolution**

Resolution (nm)

**property resolution\_actual**

Resolution Actual (ON/OFF)

**property resolution\_vbw**

Video Bandwidth Resolution

**property sampling\_points**

Number of sampling points

**single\_sweep(\*\*kwargs)**

Perform a single sweep and wait for completion.

**property trace\_marker**

Sets the trace marker with a wavelength. Returns the trace wavelength and power.

**property trace\_marker\_center**

Trace Marker at Center. Set to 1 or True to initiate command

**wait**(*n=3, delay=1*)

Query OPC Command and waits for appropriate response.

**wait\_for\_sweep**(*n=20, delay=0.5*)

Wait for a sweep to stop.

This is performed by checking bit 1 of the ESR2.

**property wavelength\_center**

Center Wavelength of Spectrum Scan in nm.

**property wavelength\_marker\_value**

Wavelength Marker Value (wavelength or freq.?)

**property wavelength\_span**

Wavelength Span of Spectrum Scan in nm.

**property wavelength\_start**

Wavelength Start of Spectrum Scan in nm.

**property wavelength\_stop**

Wavelength Stop of Spectrum Scan in nm.

**property wavelength\_value\_in**

Wavelength value in Vacuum or Air

**property wavelengths**

Return a numpy array of the current wavelengths of scans.

## 7.12 Attocube

This section contains specific documentation on the Attocube instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.12.1 Attocube Adapters

**class** `pymeasure.instruments.attocube.adapters.AttocubeConsoleAdapter`(*host, port, passwd, \*\*kwargs*)

Bases: `pymeasure.adapters.telnet.TelnetAdapter`

Adapter class for connecting to the Attocube Standard Console. This console is a Telnet prompt with password authentication.

**Parameters**

- **host** – host address of the instrument
- **port** – TCPIP port
- **passwd** – password required to open the connection
- **kwargs** – Any valid key-word argument for TelnetAdapter

**ask**(*command*)

Writes a command to the instrument and returns the resulting ASCII response

**Parameters** **command** – command string to be sent to the instrument

**Returns** String ASCII response of the instrument

**check\_acknowledgement**(*reply*, *msg=""*)

checks the last reply of the instrument to be 'OK', otherwise a ValueError is raised.

**Parameters**

- **reply** – last reply string of the instrument
- **msg** – optional message for the eventual error

**extract\_value**(*reply*)

preprocess\_reply function for the Attocube console. This function tries to extract <value> from 'name = <value> [unit]'. If <value> can not be identified the original string is returned.

**Parameters** **reply** – reply string

**Returns** string with only the numerical value, or the original string

**read**()

Reads a reply of the instrument which consists of two or more lines. The first ones are the reply to the command while the last one is 'OK' or 'ERROR' to indicate any problem. In case the reply is not OK a ValueError is raised.

**Returns** String ASCII response of the instrument.

**write**(*command*, *check\_ack=True*)

Writes a command to the instrument

**Parameters**

- **command** – command string to be sent to the instrument
- **check\_ack** – boolean flag to decide if the acknowledgement is read back from the instrument. This should be True for set pure commands and False otherwise.

## 7.12.2 Attocube ANC300 Motion Controller

**class** `pymeasure.instruments.attocube.anc300.ANC300Controller`(*host*, *axisnames*, *passwd*, *query\_delay=0.05*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Attocube ANC300 Piezo stage controller with several axes

**Parameters**

- **host** – host address of the instrument
- **axisnames** – a list of axis names which will be used to create properties with these names
- **passwd** – password for the attocube standard console
- **query\_delay** – delay between sending and reading (default 0.05 sec)
- **kwargs** – Any valid key-word argument for TelnetAdapter

**property** `controllerBoardVersion`

Serial number of the controller board

**ground\_all()**

Grounds all axis of the controller.

**stop\_all()**

Stop all movements of the axis.

**property version**

Version number and instrument identification

**class** pymeasure.instruments.attocube.anc300.**Axis**(*controller*, *axis*)

Bases: object

Represents a single open loop axis of the Attocube ANC350

**Parameters**

- **axis** – axis identifier, integer from 1 to 7
- **controller** – ANC300Controller instance used for the communication

**property capacity**

Saved capacity value in nF of the axis.

**property frequency**

Frequency of the stepping motion in Hertz from 1 to 10000 Hz. This property can be set.

**measure\_capacity()**

Obtains a new measurement of the capacity. The mode of the axis returns to 'gnd' after the measurement.

**Returns capacity** the freshly measured capacity in nF.

**property mode**

Axis mode. This can be 'gnd', 'inp', 'cap', 'stp', 'off', 'stp+', 'stp-'. Available modes depend on the actual axis model

**move**(*steps*, *gnd=True*)

Move 'steps' steps in the direction given by the sign of the argument. This method will change the mode of the axis automatically and ground the axis on the end if 'gnd' is True. The method returns only when the movement is finished.

**Parameters**

- **steps** – finite integer value of steps to be performed. A positive sign corresponds to upwards steps, a negative sign to downwards steps.
- **gnd** – bool, flag to decide if the axis should be grounded after completion of the movement

**property offset\_voltage**

Offset voltage in Volts from 0 to 150 V. This property can be set.

**property output\_voltage**

Output voltage in volts.

**property pattern\_down**

step down pattern of the piezo drive. 256 values ranging from 0 to 255 representing the the sequence of output voltages within one step of the piezo drive. This property can be set, the set value needs to be an array with 256 integer values.

**property pattern\_up**

step up pattern of the piezo drive. 256 values ranging from 0 to 255 representing the the sequence of output voltages within one step of the piezo drive. This property can be set, the set value needs to be an array with 256 integer values.



**property serial\_nr**

Serial number of the axis

**property stepd**

Step downwards for N steps. Mode must be 'stp' and N must be positive.

**property stepu**

Step upwards for N steps. Mode must be 'stp' and N must be positive.

**stop()**

Stop any motion of the axis

**property voltage**

Amplitude of the stepping voltage in volts from 0 to 150 V. This property can be set.

## 7.13 BK Precision

This section contains specific documentation on the BK Precision instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.13.1 BK Precision 9130B DC Power Supply

**class** `pymeasure.instruments.bkprecision.BKPrecision9130B(adapter, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the BK Precision 9130B DC Power Supply interface for interacting with the instrument.

**property channel**

An integer property used to control which channel is selected. Can only take values [1, 2, 3].

**check\_errors()**

Read all errors from the instrument.

**Returns** list of error entries

**property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device's Output Queue when all pending selected device operations have been finished.

**property current**

Floating point property used to control current of the selected channel.

**property id**

Requests and returns the identification of the instrument.

**property options**

Requests and returns the device options installed.

**reset()**

Resets the instrument.

**shutdown()**

Brings the instrument to a safe and stable state

**property source\_enabled**

A boolean property that controls whether the source is enabled, takes values True or False.

**property status**

Requests and returns the status byte and Master Summary Status bit.

**property voltage**

Floating point property used to control voltage of the selected channel.

## 7.14 Danfysik

This section contains specific documentation on the Danfysik instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.14.1 Danfysik Serial Adapter

**class** `pymeasure.instruments.danfysik.DanfysikAdapter`(*port*)

Bases: `pymeasure.adapters.serial.SerialAdapter`

Provides a `SerialAdapter` with the specific baudrate and timeout for Danfysik serial communication.

Initiates the adapter to open serial communication over the supplied port.

**Parameters** **port** – A string representing the serial port

**read()**

Overwrites the `SerialAdapter.read` method to automatically raise exceptions if errors are reported by the instrument.

**Returns** String ASCII response of the instrument

**Raises** An Exception if the Danfysik raises an error

**write(command)**

Overwrites the `SerialAdapter.write` method to automatically append a Unix-style linebreak at the end of the command.

**Parameters** **command** – SCPI command string to be sent to the instrument

### 7.14.2 Danfysik 8500 Power Supply

**class** `pymeasure.instruments.danfysik.Danfysik8500`(*port*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Danfysik 8500 Electromagnet Current Supply and provides a high-level interface for interacting with the instrument

To allow user access to the Prolific Technology PL2303 Serial port adapter in Linux, create the file: `/etc/udev/rules.d/50-danfysik.rules`, with contents:

```
SUBSYSTEMS=="usb",ATTRS{idVendor}=="067b",ATTRS{idProduct}=="2303",MODE="0666",
↪SYMLINK+="danfysik"
```

Then reload the udev rules with:

```
sudo udevadm control --reload-rules
sudo udevadm trigger
```

The device will be accessible through the port `/dev/danfysik`.

**add\_ramp\_step**(*current*)

Adds a current step to the ramp set.

**Parameters** *current* – A current in Amps

**clear\_ramp\_set**()

Clears the ramp set.

**clear\_sequence**(*stack*)

Clears the sequence by the stack number.

**Parameters** *stack* – A stack number between 0-15

**property current**

The actual current in Amps. This property can be set through [current\\_ppm](#).

**property current\_ppm**

The current in parts per million. This property can be set.

**property current\_setpoint**

The setpoint for the current, which can deviate from the actual current ([current](#)) while the supply is in the process of setting the value.

**disable**()

Disables the flow of current.

**enable**()

Enables the flow of current.

**property id**

Reads the identification information.

**is\_current\_stable**()

Returns True if the current is within 0.02 A of the setpoint value.

**is\_enabled**()

Returns True if the current supply is enabled.

**is\_ready**()

Returns True if the instrument is in the ready state.

**is\_sequence\_running**(*stack*)

Returns True if a sequence is running with a given stack number

**Parameters** *stack* – A stack number between 0-15

**local**()

Sets the instrument in local mode, where the front panel can be used.

**property polarity**

The polarity of the current supply, being either -1 or 1. This property can be set by supplying one of these values.

**ramp\_to\_current**(*current*, *points*, *delay\_time=1*)

Executes [set\\_ramp\\_to\\_current\(\)](#) and starts the ramp.

**remote**()

Sets the instrument in remote mode, where the front panel is disabled.

**reset\_interlocks**()

Resets the instrument interlocks.

**set\_ramp\_delay**(*time*)

Sets the ramp delay time in seconds.

**Parameters** **time** – The time delay time in seconds

**set\_ramp\_to\_current**(*current, points, delay\_time=1*)

Sets up a linear ramp from the initial current to a different current, with a number of points, and delay time.

**Parameters**

- **current** – The final current in Amps
- **points** – The number of linear points to traverse
- **delay\_time** – A delay time in seconds

**set\_sequence**(*stack, currents, times, multiplier=999999*)

Sets up an arbitrary ramp profile with a list of currents (Amps) and a list of interval times (seconds) on the specified stack number (0-15)

**property** **slew\_rate**

The slew rate of the current sweep.

**start\_ramp**()

Starts the current ramp.

**start\_sequence**(*stack*)

Starts a sequence by the stack number.

**Parameters** **stack** – A stack number between 0-15

**property** **status**

A list of human-readable strings that contain the instrument status information, based on [status\\_hex](#).

**property** **status\_hex**

The status in hexadecimal. This value is parsed in [status](#) into a human-readable list.

**stop\_ramp**()

Stops the current ramp.

**stop\_sequence**()

Stops the currently running sequence.

**sync\_sequence**(*stack, delay=0*)

Arms the ramp sequence to be triggered by a hardware input to pin P33 1&2 (10 to 24 V) or a TS command. If a delay is provided, the sequence will start after the delay.

**Parameters**

- **stack** – A stack number between 0-15
- **delay** – A delay time in seconds

**wait\_for\_current**(*has\_aborted=<function Danfysik8500.<lambda>>, delay=0.01*)

Blocks the process until the current has stabilized. A provided function `has_aborted` can be supplied, which is checked after each delay time (in seconds) in addition to the stability check. This allows an abort feature to be integrated.

**Parameters**

- **has\_aborted** – A function that returns True if the process should stop waiting
- **delay** – The delay time in seconds between each check for stability

**wait\_for\_ready**(*has\_aborted=<function Danfysik8500.<lambda>>, delay=0.01*)

Blocks the process until the instrument is ready. A provided function `has_aborted` can be supplied, which is checked after each delay time (in seconds) in addition to the readiness check. This allows an abort feature to be integrated.

**Parameters**

- **has\_aborted** – A function that returns True if the process should stop waiting
- **delay** – The delay time in seconds between each check for readiness

## 7.15 Delta Elektronika

This section contains specific documentation on the Delta Elektronika instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.15.1 Delta Elektronika SM7045D Power source

**class** pymeasure.instruments.deltalelektronika.**SM7045D**(resourceName, \*\*kwargs)

Bases: *pymeasure.instruments.instrument.Instrument*

This is the class for the SM 70-45 D power supply.

```
source = SM7045D("GPIB::8")

source.ramp_to_zero(1)           # Set output to 0 before enabling
source.enable()                 # Enables the output
source.current = 1              # Sets a current of 1 Amps
```

**property current**

A floating point property that represents the output current of the power supply in Amps. This property can be set.

**disable()**

Enables remote shutdown, hence input will be disabled.

**enable()**

Disable remote shutdown, hence output will be enabled.

**property max\_current**

A floating point property that represents the maximum output current of the power supply in Amps. This property can be set.

**property max\_voltage**

A floating point property that represents the maximum output voltage of the power supply in Volts. This property can be set.

**property measure\_current**

Measures the actual output current of the power supply in Amps.

**property measure\_voltage**

Measures the actual output voltage of the power supply in Volts.

**ramp\_to\_current**(target\_current, current\_step=0.1)

Gradually increase/decrease current to target current.

**Parameters**

- **target\_current** – Float that sets the target current (in A)
- **current\_step** – Optional float that sets the current steps / ramp rate (in A/s)

**ramp\_to\_zero**(*current\_step=0.1*)

Gradually decrease the current to zero.

**Parameters** **current\_step** – Optional float that sets the current steps / ramp rate (in A/s)

**property** **rsd**

Check whether remote shutdown is enabled/disabled and thus if the output of the power supply is disabled/enabled.

**shutdown**()

Set the current to 0 A and disable the output of the power source.

**property** **voltage**

A floating point property that represents the output voltage setting of the power supply in Volts. This property can be set.

## 7.16 Edwards

This section contains specific documentation on the Edwards instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.16.1 Edwards nxds vacuum pump

`pymeasure.instruments.edwards.nxds`

## 7.17 Fluke

This section contains specific documentation on the Fluke instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.17.1 Fluke 7341 Temperature bath

**class** `pymeasure.instruments.fluke.Fluke7341`(*resource\_name*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the compact constant temperature bath from Fluke

**property** **id**

Read the instrument model

**property** **set\_point**

A *float* property to set the bath temperature set-point. Valid values are in the range  $-40$  to  $150$  °C. The unit is as defined in property *unit*. This property can be read

**property** **temperature**

Read the current bath temperature. The unit is as defined in property *unit*.

**property** **unit**

A string property that controls the temperature unit. Possible values are *c* for Celsius and *f* for Fahrenheit.

## 7.18 F.W. Bell

This section contains specific documentation on the F.W. Bell instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.18.1 F.W. Bell 5080 Handheld Gaussmeter

**class** `pymeasure.instruments.fwbell.FWBell5080(port)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the F.W. Bell 5080 Handheld Gaussmeter and provides a high-level interface for interacting with the instrument

**Parameters** `port` – The serial port of the instrument

```
meter = FWBell5080('/dev/ttyUSB0')      # Connects over serial port /dev/ttyUSB0
↳ (Linux)

meter.units = 'gauss'                   # Sets the measurement units to Gauss
meter.range = 3e3                        # Sets the range to 3 kG
print(meter.field)                      # Reads and prints a field measurement in G

fields = meter.fields(100)              # Samples 100 field measurements
print(fields.mean(), fields.std())       # Prints the mean and standard deviation of
↳ the samples
```

**ask(command)**

Overwrites the `Instrument.ask` method to remove the last 2 characters from the output.

**auto\_range()**

Enables the auto range functionality.

**property field**

Reads a floating point value of the field in the appropriate units.

**fields(samples=1)**

Returns a numpy array of field samples for a given sample number.

**Parameters** `samples` – The number of samples to preform

**property id**

Reads the identification information.

**property range**

A floating point property that controls the maximum field range in the active units. This can take the values of 300 G, 3 kG, and 30 kG for Gauss, 30 mT, 300 mT, and 3 T for Tesla, and 23.88 kAm, 238.8 kAm, and 2388 kAm for Amp-meter.

**read()**

Overwrites the `Instrument.read` method to remove the last 2 characters from the output.

**reset()**

Resets the instrument.

**property units**

A string property that controls the field units, which can take the values: 'gauss', 'gauss ac', 'tesla', 'tesla ac', 'amp-meter', and 'amp-meter ac'. The AC versions configure the instrument to measure AC.

**values**(*command*)

Overwrites the `Instrument.values` method to remove the lastv2 characters from the output.

## 7.19 Hewlett Packard

This section contains specific documentation on the Hewlett Packard instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.19.1 HP 33120A Arbitrary Waveform Generator

**class** `pymeasure.instruments.hp.HP33120A(resource_name, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Hewlett Packard 33120A Arbitrary Waveform Generator and provides a high-level interface for interacting with the instrument.

**property amplitude**

A floating point property that controls the voltage amplitude of the output signal. The default units are in peak-to-peak Volts, but can be controlled by `amplitude_units`. The allowed range depends on the waveform shape and can be queried with `max_amplitude` and `min_amplitude`.

**property amplitude\_units**

A string property that controls the units of the amplitude, which can take the values Vpp, Vrms, dBm, and default.

**beep()**

Causes a system beep.

**property frequency**

A floating point property that controls the frequency of the output in Hz. The allowed range depends on the waveform shape and can be queried with `max_frequency` and `min_frequency`.

**property max\_amplitude**

Reads the maximum `amplitude` in Volts for the given shape

**property max\_frequency**

Reads the maximum `frequency` in Hz for the given shape

**property max\_offset**

Reads the maximum `offset` in Volts for the given shape

**property min\_amplitude**

Reads the minimum `amplitude` in Volts for the given shape

**property min\_frequency**

Reads the minimum `frequency` in Hz for the given shape

**property min\_offset**

Reads the minimum `offset` in Volts for the given shape

**property offset**

A floating point property that controls the amplitude voltage offset in Volts. The allowed range depends on the waveform shape and can be queried with `max_offset` and `min_offset`.

**property shape**

A string property that controls the shape of the wave, which can take the values: sinusoid, square, triangle, ramp, noise, dc, and user.



### 7.19.2 HP 34401A Multimeter

**class** `pymeasure.instruments.hp.HP34401A(resourceName, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the HP 34401A instrument.

**property** `current_ac`

AC current, in Amps

**property** `current_dc`

DC current, in Amps

**property** `resistance`

Resistance, in Ohms

**property** `resistance_4w`

Four-wires (remote sensing) resistance, in Ohms

**property** `voltage_ac`

AC voltage, in Volts

**property** `voltage_dc`

DC voltage, in Volts

### 7.19.3 HP 3478A Multimeter

**class** `pymeasure.instruments.hp.HP3478A(resourceName, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Hewlett Packard 3748A 5 1/2 digit multimeter and provides a high-level interface for interacting with the instrument.

**class** `ERRORS(value)`

Bases: `enum.IntFlag`

Enum element for error bit decoding

**method** `GPIB_trigger()`

Initiate trigger via low-level GPIB-command (aka GET - group execute trigger)

**class** `SRQ(value)`

Bases: `enum.IntFlag`

Enum element for SRQ mask bit decoding

**property** `SRQ_mask`

Return current SRQ mask, this property can be set,

bit assignment for SRQ:

Bit (dec)	Description
1	SRQ when Data ready
4	SRQ when Syntax error
8	SRQ when internal error
16	front panel SQR button
32	SRQ by invalid calibration

**property** `active_connectors`

Return selected connectors ("front"/"back"), based on front-panel selector switch

**property auto\_range\_enabled**

Property describing the auto-ranging status

Value	Status
True	auto-range function activated
False	manual range selection / auto-range disabled

The range can be set with the [range](#) property

**property auto\_zero\_enabled**

Return auto-zero status, this property can be set

Value	Status
True	auto-zero active
False	auto-zero disabled

**property calibration\_enabled**

Return calibration enable switch setting, based on front-panel selector switch

Value	Status
True	calbration possible
False	calibration locked

**check\_errors()**

Method to read the error status register

**Return error\_status** one byte with the error status register content

**Rtype error\_status** int

**classmethod decode\_mode(function)**

Method to decode current mode

**Parameters function** – int indicating the measurement function selected

**Return cur\_mode** string with the current measurement mode

**Rtype cur\_mode** str

**classmethod decode\_range(range\_undecoded, function)**

Method to decode current range

**Parameters**

- **range\_undecoded** – int to be decoded
- **function** – int indicating the measurement function selected

**Return cur\_range** float value representing the active measurment range

**Rtype cur\_range** float

**classmethod decode\_status(status\_bytes, field=None)**

Method to handle the decoding of the status bytes into something meaningfull

**Parameters**

- **status\_bytes** – list of bytes to be decoded
- **field** – name of field to be returned

Return `ret_val` int status value

**static** `decode_trigger(status_bytes)`

Method to decode trigger mode

**Parameters** `status_bytes` – list of bytes to be decoded

**Return** `trigger_mode` string with the current trigger mode

**Rtype** `trigger_mode` str

**display\_reset()**

Reset the display of the instrument.

**property** `display_text`

Displays up to 12 upper-case ASCII characters on the display.

**property** `display_text_no_symbol`

Displays up to 12 upper-case ASCII characters on the display and disables all symbols on the display.

**property** `error_status`

Checks the error status register

**get\_status()**

Method to read the status bytes from the instrument :return `current_status`: a byte array representing the instrument status :rtype `current_status`: bytes

**property** `measure_ACI`

Returns the measured value for AC current as a float in A.

**property** `measure_ACV`

Returns the measured value for AC Voltage as a float in V.

**property** `measure_DCI`

Returns the measured value for DC current as a float in A.

**property** `measure_DCV`

Returns the measured value for DC Voltage as a float in V.

**property** `measure_R2W`

Returns the measured value for 2-wire resistance as a float in Ohm.

**property** `measure_R4W`

Returns the measured value for 4-wire resistance as a float in Ohm.

**property** `measure_Rext`

Returns the measured value for extended resistance mode (>30M, 2-wire) resistance as a float in Ohm.

**property** `mode`

Return current selected measurement mode, this property can be set. Allowed values are

Mode	Function
ACI	AC current
ACV	AC voltage
DCI	DC current
DCV	DC voltage
R2W	2-wire resistance
R4W	4-wire resistance
Rext	extended resistance method (requires additional 10 M resistor)

**property** `range`

Returns the current measurement range, this property can be set.

Valid values are :

Mode	Range
ACI	0.3, 3, auto
ACV	0.3, 3, 30, 300, auto
DCI	0.3, 3, auto
DCV	0.03, 0.3, 3, 30, 300, auto
R2W	30, 300, 3000, 3E4, 3E5, 3E6, 3E7, auto
R4W	30, 300, 3000, 3E4, 3E5, 3E6, 3E7, auto
Rext	3E7, auto

#### **reset()**

Initiates a reset (like a power-on reset) of the HP3478A

#### **property resolution**

Returns current selected resolution, this property can be set.

Possible values are 3,4 or 5 (for 3 1/2, 4 1/2 or 5 1/2 digits of resolution)

#### **shutdown()**

provides a way to gracefully close the connection to the HP3478A

#### **property status**

Returns an object representing the current status of the unit.

#### **property trigger**

Return current selected trigger mode, this property can be set

Possible values are:

Value	Meaning
auto	automatic trigger (internal)
internal	automatic trigger (internal)
external	external trigger (connector on back or GET)
hold	holds the measurement
fast	fast trigger for AC measurements

## 7.20 Keithley

This section contains specific documentation on the Keithley instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.20.1 Keithley 2000 Multimeter

**class** `pymeasure.instruments.keithley.Keithley2000`(*adapter*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`, `pymeasure.instruments.keithley.buffer.KeithleyBuffer`

Represents the Keithley 2000 Multimeter and provides a high-level interface for interacting with the instrument.

```
meter = Keithley2000("GPIB::1")
meter.measure_voltage()
print(meter.voltage)
```

**acquire\_reference**(*mode=None*)

Sets the active value as the reference for the active mode, or can set another mode by its name.

**Parameters** *mode* – A valid *mode* name, or None for the active mode

**auto\_range**(*mode=None*)

Sets the active mode to use auto-range, or can set another mode by its name.

**Parameters** *mode* – A valid *mode* name, or None for the active mode

**beep**(*frequency, duration*)

Sounds a system beep.

**Parameters**

- **frequency** – A frequency in Hz between 65 Hz and 2 MHz
- **duration** – A time in seconds between 0 and 7.9 seconds

**property beep\_state**

A string property that enables or disables the system status beeper, which can take the values: :code:'enabled' and :code:'disabled'.

**property buffer\_data**

Returns a numpy array of values from the buffer.

**property buffer\_points**

An integer property that controls the number of buffer points. This does not represent actual points in the buffer, but the configuration value instead.

**check\_errors**()

Read all errors from the instrument.

**Returns** list of error entries

**property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device's Output Queue when all pending selected device operations have been finished.

**config\_buffer**(*points=64, delay=0*)

Configures the measurement buffer for a number of points, to be taken with a specified delay.

**Parameters**

- **points** – The number of points in the buffer.
- **delay** – The delay time in seconds.

**property current**

Reads a DC or AC current measurement in Amps, based on the active *mode*.

**property current\_ac\_bandwidth**

A floating point property that sets the AC current detector bandwidth in Hz, which can take the values 3, 30, and 300 Hz.

**property current\_ac\_digits**

An integer property that controls the number of digits in the AC current readings, which can take values from 4 to 7.

**property current\_ac\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the AC current measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property current\_ac\_range**

A floating point property that controls the AC current range in Amps, which can take values from 0 to 3.1 A. Auto-range is disabled when this property is set.

**property current\_ac\_reference**

A floating point property that controls the AC current reference value in Amps, which can take values from -3.1 to 3.1 A.

**property current\_digits**

An integer property that controls the number of digits in the DC current readings, which can take values from 4 to 7.

**property current\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the DC current measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property current\_range**

A floating point property that controls the DC current range in Amps, which can take values from 0 to 3.1 A. Auto-range is disabled when this property is set.

**property current\_reference**

A floating point property that controls the DC current reference value in Amps, which can take values from -3.1 to 3.1 A.

**disable\_buffer()**

Disables the connection between measurements and the buffer, but does not abort the measurement process.

**disable\_filter(mode=None)**

Disables the averaging filter for the active mode, or can set another mode by its name.

**Parameters** *mode* – A valid *mode* name, or None for the active mode

**disable\_reference(mode=None)**

Disables the reference for the active mode, or can set another mode by its name.

**Parameters** *mode* – A valid *mode* name, or None for the active mode

**enable\_filter(mode=None, type='repeat', count=1)**

Enables the averaging filter for the active mode, or can set another mode by its name.

**Parameters**

- **mode** – A valid *mode* name, or None for the active mode
- **type** – The type of averaging filter, either 'repeat' or 'moving'.
- **count** – A number of averages, which can take values from 1 to 100

**enable\_reference(mode=None)**

Enables the reference for the active mode, or can set another mode by its name.

**Parameters** *mode* – A valid *mode* name, or None for the active mode

**property frequency**

Reads a frequency measurement in Hz, based on the active *mode*.

**property frequency\_aperature**

A floating point property that controls the frequency aperture in seconds, which sets the integration period and measurement speed. Takes values from 0.01 to 1.0 s.

**property frequency\_digits**

An integer property that controls the number of digits in the frequency readings, which can take values from 4 to 7.

**property frequency\_reference**

A floating point property that controls the frequency reference value in Hz, which can take values from 0 to 15 MHz.

**property frequency\_threshold**

A floating point property that controls the voltage signal threshold level in Volts for the frequency measurement, which can take values from 0 to 1010 V.

**property id**

Requests and returns the identification of the instrument.

**is\_buffer\_full()**

Returns True if the buffer is full of measurements.

**local()**

Returns control to the instrument panel, and enables the panel if disabled.

**measure\_continuity()**

Configures the instrument to perform continuity testing.

**measure\_current**(*max\_current=0.01, ac=False*)

Configures the instrument to measure current, based on a maximum current to set the range, and a boolean flag to determine if DC or AC is required.

**Parameters**

- **max\_current** – A current in Volts to set the current range
- **ac** – False for DC current, and True for AC current

**measure\_diode()**

Configures the instrument to perform diode testing.

**measure\_frequency()**

Configures the instrument to measure the frequency.

**measure\_period()**

Configures the instrument to measure the period.

**measure\_resistance**(*max\_resistance=10000000.0, wires=2*)

Configures the instrument to measure voltage, based on a maximum voltage to set the range, and a boolean flag to determine if DC or AC is required.

**Parameters**

- **max\_voltage** – A voltage in Volts to set the voltage range
- **ac** – False for DC voltage, and True for AC voltage

**measure\_temperature()**

Configures the instrument to measure the temperature.

**measure\_voltage**(*max\_voltage=1, ac=False*)

Configures the instrument to measure voltage, based on a maximum voltage to set the range, and a boolean flag to determine if DC or AC is required.

**Parameters**

- **max\_voltage** – A voltage in Volts to set the voltage range
- **ac** – False for DC voltage, and True for AC voltage

**property mode**

A string property that controls the configuration mode for measurements, which can take the values:

:code:'current' (DC), :code:'current ac', :code:'voltage' (DC), :code:'voltage ac', :code:'resistance' (2-wire), :code:'resistance 4W' (4-wire), :code:'period', :code:'frequency', :code:'temperature', :code:'diode', and :code:'frequency'.

**property options**

Requests and returns the device options installed.

**property period**

Reads a period measurement in seconds, based on the active *mode*.

**property period\_aperature**

A floating point property that controls the period aperature in seconds, which sets the integration period and measurement speed. Takes values from 0.01 to 1.0 s.

**property period\_digits**

An integer property that controls the number of digits in the period readings, which can take values from 4 to 7.

**property period\_reference**

A floating point property that controls the period reference value in seconds, which can take values from 0 to 1 s.

**property period\_threshold**

A floating point property that controls the voltage signal threshold level in Volts for the period measurement, which can take values from 0 to 1010 V.

**remote()**

Places the instrument in the remote state, which is does not need to be explicitly called in general.

**remote\_lock()**

Disables and locks the front panel controls to prevent changes during remote operations. This is disabled by calling *local()*.

**reset()**

Resets the instrument state.

**reset\_buffer()**

Resets the buffer.

**property resistance**

Reads a resistance measurement in Ohms for both 2-wire and 4-wire configurations, based on the active *mode*.

**property resistance\_4W\_digits**

An integer property that controls the number of digits in the 4-wire resistance readings, which can take values from 4 to 7.

**property resistance\_4W\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the 4-wire resistance measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property resistance\_4W\_range**

A floating point property that controls the 4-wire resistance range in Ohms, which can take values from 0 to 120 MOhms. Auto-range is disabled when this property is set.

**property resistance\_4W\_reference**

A floating point property that controls the 4-wire resistance reference value in Ohms, which can take values from 0 to 120 MOhms.



**property resistance\_digits**

An integer property that controls the number of digits in the 2-wire resistance readings, which can take values from 4 to 7.

**property resistance\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the 2-wire resistance measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property resistance\_range**

A floating point property that controls the 2-wire resistance range in Ohms, which can take values from 0 to 120 MOhms. Auto-range is disabled when this property is set.

**property resistance\_reference**

A floating point property that controls the 2-wire resistance reference value in Ohms, which can take values from 0 to 120 MOhms.

**shutdown()**

Brings the instrument to a safe and stable state

**start\_buffer()**

Starts the buffer.

**property status**

Requests and returns the status byte and Master Summary Status bit.

**stop\_buffer()**

Aborts the buffering measurement, by stopping the measurement arming and triggering sequence. If possible, a Selected Device Clear (SDC) is used.

**property temperature**

Reads a temperature measurement in Celsius, based on the active *mode*.

**property temperature\_digits**

An integer property that controls the number of digits in the temperature readings, which can take values from 4 to 7.

**property temperature\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the temperature measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property temperature\_reference**

A floating point property that controls the temperature reference value in Celsius, which can take values from -200 to 1372 C.

**property trigger\_count**

An integer property that controls the trigger count, which can take values from 1 to 9,999.

**property trigger\_delay**

A floating point property that controls the trigger delay in seconds, which can take values from 1 to 9,999,999.999 s.

**property voltage**

Reads a DC or AC voltage measurement in Volts, based on the active *mode*.

**property voltage\_ac\_bandwidth**

A floating point property that sets the AC voltage detector bandwidth in Hz, which can take the values 3, 30, and 300 Hz.

**property voltage\_ac\_digits**

An integer property that controls the number of digits in the AC voltage readings, which can take values from 4 to 7.

**property voltage\_ac\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the AC voltage measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property voltage\_ac\_range**

A floating point property that controls the AC voltage range in Volts, which can take values from 0 to 757.5 V. Auto-range is disabled when this property is set.

**property voltage\_ac\_reference**

A floating point property that controls the AC voltage reference value in Volts, which can take values from -757.5 to 757.5 Volts.

**property voltage\_digits**

An integer property that controls the number of digits in the DC voltage readings, which can take values from 4 to 7.

**property voltage\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the DC voltage measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property voltage\_range**

A floating point property that controls the DC voltage range in Volts, which can take values from 0 to 1010 V. Auto-range is disabled when this property is set.

**property voltage\_reference**

A floating point property that controls the DC voltage reference value in Volts, which can take values from -1010 to 1010 V.

**wait\_for\_buffer**(*should\_stop=<function KeithleyBuffer.<lambda>>, timeout=60, interval=0.1*)

Blocks the program, waiting for a full buffer. This function returns early if the `should_stop` function returns True or the timeout is reached before the buffer is full.

**Parameters**

- **should\_stop** – A function that returns True when this function should return early
- **timeout** – A time in seconds after which this function should return early
- **interval** – A time in seconds for how often to check if the buffer is full

## 7.20.2 Keithley 2260B DC Power Supply

```
class pymeasure.instruments.keithley.Keithley2260B(adapter, read_termination='\n', **kwargs)
```

Bases: [pymeasure.instruments.instrument.Instrument](#)

Represents the Keithley 2260B Power Supply (minimal implementation) and provides a high-level interface for interacting with the instrument.

For a connection through tcpip, the device only accepts connections at port 2268, which cannot be configured otherwise. example connection string: 'TCPIP::xxx.xxx.xxx.xxx::2268::SOCKET' the read termination for this interface is

```
source = Keithley2260B("GPIB::1")
source.voltage = 1
print(source.voltage)
print(source.current)
print(source.power)
print(source.applied)
```

**property applied**

Simultaneous control of voltage (volts) and current (amps). Values need to be supplied as tuple of (voltage, current). Depending on whether the instrument is in constant current or constant voltage mode, the values achieved by the instrument will differ from the ones set.

**check\_errors()**

Logs any system errors reported by the instrument.

**property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device's Output Queue when all pending selected device operations have been finished.

**property current**

Reads the current (in Ampere) the dc power supply is putting out.

**property current\_limit**

A floating point property that controls the source current in amps. This is not checked against the allowed range. Depending on whether the instrument is in constant current or constant voltage mode, this might differ from the actual current achieved.

**property enabled**

A boolean property that controls whether the source is enabled, takes values True or False.

**property error**

Returns a tuple of an error code and message from a single error.

**property id**

Requests and returns the identification of the instrument.

**property options**

Requests and returns the device options installed.

**property power**

Reads the power (in Watt) the dc power supply is putting out.

**reset()**

Resets the instrument.

**shutdown()**

Disable output, call parent function

**property status**

Requests and returns the status byte and Master Summary Status bit.

**property voltage**

Reads the voltage (in Volt) the dc power supply is putting out.

**property voltage\_setpoint**

A floating point property that controls the source voltage in volts. This is not checked against the allowed range. Depending on whether the instrument is in constant current or constant voltage mode, this might differ from the actual voltage achieved.

### 7.20.3 Keithley 2400 SourceMeter

**class** `pymeasure.instruments.keithley.Keithley2400`(*adapter*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`, `pymeasure.instruments.keithley.buffer.KeithleyBuffer`

Represents the Keithley 2400 SourceMeter and provides a high-level interface for interacting with the instrument.

```
keithley = Keithley2400("GPIB::1")

keithley.apply_current()           # Sets up to source current
keithley.source_current_range = 10e-3 # Sets the source current range to 10 mA
keithley.compliance_voltage = 10    # Sets the compliance voltage to 10 V
keithley.source_current = 0         # Sets the source current to 0 mA
keithley.enable_source()           # Enables the source output

keithley.measure_voltage()         # Sets up to measure voltage

keithley.ramp_to_current(5e-3)     # Ramps the current to 5 mA
print(keithley.voltage)             # Prints the voltage in Volts

keithley.shutdown()               # Ramps the current to 0 mA and disables_
↪ output
```

**apply\_current**(*current\_range=None*, *compliance\_voltage=0.1*)

Configures the instrument to apply a source current, and uses an auto range unless a current range is specified. The compliance voltage is also set.

#### Parameters

- **compliance\_voltage** – A float in the correct range for a `compliance_voltage`
- **current\_range** – A `current_range` value or None

**apply\_voltage**(*voltage\_range=None*, *compliance\_current=0.1*)

Configures the instrument to apply a source voltage, and uses an auto range unless a voltage range is specified. The compliance current is also set.

#### Parameters

- **compliance\_current** – A float in the correct range for a `compliance_current`
- **voltage\_range** – A `voltage_range` value or None

**property auto\_output\_off**

A boolean property that enables or disables the auto output-off. Valid values are True (output off after measurement) and False (output stays on after measurement).

**auto\_range\_source**()

Configures the source to use an automatic range.

**property auto\_zero**

A property that controls the auto zero option. Valid values are True (enabled) and False (disabled) and 'ONCE' (force immediate).

**beep**(*frequency*, *duration*)

Sounds a system beep.

#### Parameters

- **frequency** – A frequency in Hz between 65 Hz and 2 MHz

- **duration** – A time in seconds between 0 and 7.9 seconds

**property buffer\_data**

Returns a numpy array of values from the buffer.

**property buffer\_points**

An integer property that controls the number of buffer points. This does not represent actual points in the buffer, but the configuration value instead.

**check\_errors()**

Logs any system errors reported by the instrument.

**property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device's Output Queue when all pending selected device operations have been finished.

**property compliance\_current**

A floating point property that controls the compliance current in Amps.

**property compliance\_voltage**

A floating point property that controls the compliance voltage in Volts.

**config\_buffer(*points=64, delay=0*)**

Configures the measurement buffer for a number of points, to be taken with a specified delay.

**Parameters**

- **points** – The number of points in the buffer.
- **delay** – The delay time in seconds.

**property current**

Reads the current in Amps, if configured for this reading.

**property current\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the DC current measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property current\_range**

A floating point property that controls the measurement current range in Amps, which can take values between -1.05 and +1.05 A. Auto-range is disabled when this property is set.

**disable\_buffer()**

Disables the connection between measurements and the buffer, but does not abort the measurement process.

**disable\_output\_trigger()**

Disables the output trigger for the Trigger layer

**disable\_source()**

Disables the source of current or voltage depending on the configuration of the instrument.

**property display\_enabled**

A boolean property that controls whether or not the display of the sourcemeter is enabled. Valid values are True and False.

**enable\_source()**

Enables the source of current or voltage depending on the configuration of the instrument.

**property error**

Returns a tuple of an error code and message from a single error.

**property filter\_count**

A integer property that controls the number of readings that are acquired and stored in the filter buffer for the averaging

**property filter\_state**

A string property that controls if the filter is active.

**property filter\_type**

A String property that controls the filter's type. REP : Repeating filter MOV : Moving filter

**property id**

Requests and returns the identification of the instrument.

**is\_buffer\_full()**

Returns True if the buffer is full of measurements.

**property line\_frequency**

An integer property that controls the line frequency in Hertz. Valid values are 50 and 60.

**property line\_frequency\_auto**

A boolean property that enables or disables auto line frequency. Valid values are True and False.

**property max\_current**

Returns the maximum current from the buffer

**property max\_resistance**

Returns the maximum resistance from the buffer

**property max\_voltage**

Returns the maximum voltage from the buffer

**property maximums**

Returns the calculated maximums for voltage, current, and resistance from the buffer data as a list.

**property mean\_current**

Returns the mean current from the buffer

**property mean\_resistance**

Returns the mean resistance from the buffer

**property mean\_voltage**

Returns the mean voltage from the buffer

**property means**

Reads the calculated means (averages) for voltage, current, and resistance from the buffer data as a list.

**property measure\_concurrent\_functions**

A boolean property that enables or disables the ability to measure more than one function simultaneously. When disabled, volts function is enabled. Valid values are True and False.

**measure\_current**(*nplc=1, current=0.000105, auto\_range=True*)

Configures the measurement of current.

**Parameters**

- **nplc** – Number of power line cycles (NPLC) from 0.01 to 10
- **current** – Upper limit of current in Amps, from -1.05 A to 1.05 A
- **auto\_range** – Enables auto\_range if True, else uses the set current

**measure\_resistance**(*nplc=1, resistance=210000.0, auto\_range=True*)

Configures the measurement of resistance.

**Parameters**

- **nplc** – Number of power line cycles (NPLC) from 0.01 to 10
- **resistance** – Upper limit of resistance in Ohms, from -210 MOhms to 210 MOhms
- **auto\_range** – Enables auto\_range if True, else uses the set resistance

**measure\_voltage**(*nplc=1, voltage=21.0, auto\_range=True*)

Configures the measurement of voltage.

#### Parameters

- **nplc** – Number of power line cycles (NPLC) from 0.01 to 10
- **voltage** – Upper limit of voltage in Volts, from -210 V to 210 V
- **auto\_range** – Enables auto\_range if True, else uses the set voltage

**property min\_current**

Returns the minimum current from the buffer

**property min\_resistance**

Returns the minimum resistance from the buffer

**property min\_voltage**

Returns the minimum voltage from the buffer

**property minimums**

Returns the calculated minimums for voltage, current, and resistance from the buffer data as a list.

**property options**

Requests and returns the device options installed.

**property output\_off\_state**

Select the output-off state of the SourceMeter. HIMP : output relay is open, disconnects external circuitry. NORM : V-Source is selected and set to 0V, Compliance is set to 0.5% full scale of the present current range. ZERO : V-Source is selected and set to 0V, compliance is set to the programmed Source I value or to 0.5% full scale of the present current range, whichever is greater. GUAR : I-Source is selected and set to 0A

**output\_trigger\_on\_external**(*line=1, after='DEL'*)

Configures the output trigger on the specified trigger link line number, with the option of supplying the part of the measurement after which the trigger should be generated (default to delay, which is right before the measurement)

#### Parameters

- **line** – A trigger line from 1 to 4
- **after** – An event string that determines when to trigger

**ramp\_to\_current**(*target\_current, steps=30, pause=0.02*)

Ramps to a target current from the set current value over a certain number of linear steps, each separated by a pause duration.

#### Parameters

- **target\_current** – A current in Amps
- **steps** – An integer number of steps
- **pause** – A pause duration in seconds to wait between steps

**ramp\_to\_voltage**(*target\_voltage, steps=30, pause=0.02*)

Ramps to a target voltage from the set voltage value over a certain number of linear steps, each separated by a pause duration.

**Parameters**

- **target\_voltage** – A voltage in Amps
- **steps** – An integer number of steps
- **pause** – A pause duration in seconds to wait between steps

**reset()**

Resets the instrument and clears the queue.

**reset\_buffer()**

Resets the buffer.

**property resistance**

Reads the resistance in Ohms, if configured for this reading.

**property resistance\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the 2-wire resistance measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property resistance\_range**

A floating point property that controls the resistance range in Ohms, which can take values from 0 to 210 MOhms. Auto-range is disabled when this property is set.

**sample\_continuously()**

Causes the instrument to continuously read samples and turns off any buffer or output triggering

**set\_timed\_arm(interval)**

Sets up the measurement to be taken with the internal trigger at a variable sampling rate defined by the interval in seconds between sampling points

**set\_trigger\_counts(arm, trigger)**

Sets the number of counts for both the sweeps (arm) and the points in those sweeps (trigger), where the total number of points can not exceed 2500

**shutdown()**

Ensures that the current or voltage is turned to zero and disables the output.

**property source\_current**

A floating point property that controls the source current in Amps.

**property source\_current\_range**

A floating point property that controls the source current range in Amps, which can take values between -1.05 and +1.05 A. Auto-range is disabled when this property is set.

**property source\_delay**

A floating point property that sets a manual delay for the source after the output is turned on before a measurement is taken. When this property is set, the auto delay is turned off. Valid values are between 0 [seconds] and 999.9999 [seconds].

**property source\_delay\_auto**

A boolean property that enables or disables auto delay. Valid values are True and False.

**property source\_enabled**

A boolean property that controls whether the source is enabled, takes values True or False. The convenience methods [enable\\_source\(\)](#) and [disable\\_source\(\)](#) can also be used.

**property source\_mode**

A string property that controls the source mode, which can take the values 'current' or 'voltage'. The convenience methods [apply\\_current\(\)](#) and [apply\\_voltage\(\)](#) can also be used.



**property source\_voltage**

A floating point property that controls the source voltage in Volts.

**property source\_voltage\_range**

A floating point property that controls the source voltage range in Volts, which can take values from -210 to 210 V. Auto-range is disabled when this property is set.

**property standard\_devs**

Returns the calculated standard deviations for voltage, current, and resistance from the buffer data as a list.

**start\_buffer()**

Starts the buffer.

**status()**

Requests and returns the status byte and Master Summary Status bit.

**property std\_current**

Returns the current standard deviation from the buffer

**property std\_resistance**

Returns the resistance standard deviation from the buffer

**property std\_voltage**

Returns the voltage standard deviation from the buffer

**stop\_buffer()**

Aborts the buffering measurement, by stopping the measurement arming and triggering sequence. If possible, a Selected Device Clear (SDC) is used.

**triad(*base\_frequency*, *duration*)**

Sounds a musical triad using the system beep.

**Parameters**

- **base\_frequency** – A frequency in Hz between 65 Hz and 1.3 MHz
- **duration** – A time in seconds between 0 and 7.9 seconds

**trigger()**

Executes a bus trigger, which can be used when [trigger\\_on\\_bus\(\)](#) is configured.

**property trigger\_count**

An integer property that controls the trigger count, which can take values from 1 to 9,999.

**property trigger\_delay**

A floating point property that controls the trigger delay in seconds, which can take values from 0 to 999.9999 s.

**trigger\_immediately()**

Configures measurements to be taken with the internal trigger at the maximum sampling rate.

**trigger\_on\_bus()**

Configures the trigger to detect events based on the bus trigger, which can be activated by [trigger\(\)](#).

**trigger\_on\_external(*line=1*)**

Configures the measurement trigger to be taken from a specific line of an external trigger

**Parameters** **line** – A trigger line from 1 to 4

**use\_front\_terminals()**

Enables the front terminals for measurement, and disables the rear terminals.

**use\_rear\_terminals()**

Enables the rear terminals for measurement, and disables the front terminals.

**property voltage**

Reads the voltage in Volts, if configured for this reading.

**property voltage\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the DC voltage measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property voltage\_range**

A floating point property that controls the measurement voltage range in Volts, which can take values from -210 to 210 V. Auto-range is disabled when this property is set.

**wait\_for\_buffer**(*should\_stop*=<function KeithleyBuffer.<lambda>>, *timeout*=60, *interval*=0.1)

Blocks the program, waiting for a full buffer. This function returns early if the *should\_stop* function returns True or the timeout is reached before the buffer is full.

**Parameters**

- **should\_stop** – A function that returns True when this function should return early
- **timeout** – A time in seconds after which this function should return early
- **interval** – A time in seconds for how often to check if the buffer is full

**property wires**

An integer property that controls the number of wires in use for resistance measurements, which can take the value of 2 or 4.

## 7.20.4 Keithley 2450 SourceMeter

**class** pymeasure.instruments.keithley.**Keithley2450**(*adapter*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`, `pymeasure.instruments.keithley.buffer.KeithleyBuffer`

Represents the Keithely 2450 SourceMeter and provides a high-level interface for interacting with the instrument.

```
keithley = Keithley2450("GPIB::1")

keithley.apply_current()           # Sets up to source current
keithley.source_current_range = 10e-3 # Sets the source current range to 10 mA
keithley.compliance_voltage = 10    # Sets the compliance voltage to 10 V
keithley.source_current = 0         # Sets the source current to 0 mA
keithley.enable_source()           # Enables the source output

keithley.measure_voltage()         # Sets up to measure voltage

keithley.ramp_to_current(5e-3)     # Ramps the current to 5 mA
print(keithley.voltage)             # Prints the voltage in Volts

keithley.shutdown()               # Ramps the current to 0 mA and disables
↪ output
```

**apply\_current**(*current\_range*=None, *compliance\_voltage*=0.1)

Configures the instrument to apply a source current, and uses an auto range unless a current range is specified. The compliance voltage is also set.

**Parameters**

- **compliance\_voltage** – A float in the correct range for a `compliance_voltage`

- **current\_range** – A [current\\_range](#) value or None

**apply\_voltage**(*voltage\_range=None, compliance\_current=0.1*)

Configures the instrument to apply a source voltage, and uses an auto range unless a voltage range is specified. The compliance current is also set.

#### Parameters

- **compliance\_current** – A float in the correct range for a [compliance\\_current](#)
- **voltage\_range** – A [voltage\\_range](#) value or None

**auto\_range\_source**()

Configures the source to use an automatic range.

**beep**(*frequency, duration*)

Sounds a system beep.

#### Parameters

- **frequency** – A frequency in Hz between 65 Hz and 2 MHz
- **duration** – A time in seconds between 0 and 7.9 seconds

**property buffer\_data**

Returns a numpy array of values from the buffer.

**property buffer\_points**

An integer property that controls the number of buffer points. This does not represent actual points in the buffer, but the configuration value instead.

**check\_errors**()

Logs any system errors reported by the instrument.

**property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device's Output Queue when all pending selected device operations have been finished.

**property compliance\_current**

A floating point property that controls the compliance current in Amps.

**property compliance\_voltage**

A floating point property that controls the compliance voltage in Volts.

**config\_buffer**(*points=64, delay=0*)

Configures the measurement buffer for a number of points, to be taken with a specified delay.

#### Parameters

- **points** – The number of points in the buffer.
- **delay** – The delay time in seconds.

**property current**

Reads the current in Amps, if configured for this reading.

**property current\_filter\_count**

A integer property that controls the number of readings that are acquired and stored in the filter buffer for the averaging

**property current\_filter\_state**

A string property that controls if the filter is active.

**property current\_filter\_type**

A String property that controls the filter's type for the current. REP : Repeating filter MOV : Moving filter

**property current\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the DC current measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property current\_output\_off\_state**

Select the output-off state of the SourceMeter. HIMP : output relay is open, disconnects external circuitry. NORM : V-Source is selected and set to 0V, Compliance is set to 0.5% full scale of the present current range. ZERO : V-Source is selected and set to 0V, compliance is set to the programmed Source I value or to 0.5% full scale of the present current range, whichever is greater. GUAR : I-Source is selected and set to 0A

**property current\_range**

A floating point property that controls the measurement current range in Amps, which can take values between -1.05 and +1.05 A. Auto-range is disabled when this property is set.

**disable\_buffer()**

Disables the connection between measurements and the buffer, but does not abort the measurement process.

**disable\_source()**

Disables the source of current or voltage depending on the configuration of the instrument.

**enable\_source()**

Enables the source of current or voltage depending on the configuration of the instrument.

**property error**

Returns a tuple of an error code and message from a single error.

**property id**

Requests and returns the identification of the instrument.

**is\_buffer\_full()**

Returns True if the buffer is full of measurements.

**property max\_current**

Returns the maximum current from the buffer

**property max\_resistance**

Returns the maximum resistance from the buffer

**property max\_voltage**

Returns the maximum voltage from the buffer

**property maximums**

Returns the calculated maximums for voltage, current, and resistance from the buffer data as a list.

**property mean\_current**

Returns the mean current from the buffer

**property mean\_resistance**

Returns the mean resistance from the buffer

**property mean\_voltage**

Returns the mean voltage from the buffer

**property means**

Reads the calculated means (averages) for voltage, current, and resistance from the buffer data as a list.

**measure\_current(*nplc=1, current=0.000105, auto\_range=True*)**

Configures the measurement of current.

**Parameters**

- **nplc** – Number of power line cycles (NPLC) from 0.01 to 10
- **current** – Upper limit of current in Amps, from -1.05 A to 1.05 A
- **auto\_range** – Enables auto\_range if True, else uses the set current

**measure\_resistance**(*nplc=1, resistance=210000.0, auto\_range=True*)

Configures the measurement of resistance.

**Parameters**

- **nplc** – Number of power line cycles (NPLC) from 0.01 to 10
- **resistance** – Upper limit of resistance in Ohms, from -210 MOhms to 210 MOhms
- **auto\_range** – Enables auto\_range if True, else uses the set resistance

**measure\_voltage**(*nplc=1, voltage=21.0, auto\_range=True*)

Configures the measurement of voltage.

**Parameters**

- **nplc** – Number of power line cycles (NPLC) from 0.01 to 10
- **voltage** – Upper limit of voltage in Volts, from -210 V to 210 V
- **auto\_range** – Enables auto\_range if True, else uses the set voltage

**property min\_current**

Returns the minimum current from the buffer

**property min\_resistance**

Returns the minimum resistance from the buffer

**property min\_voltage**

Returns the minimum voltage from the buffer

**property minimums**

Returns the calculated minimums for voltage, current, and resistance from the buffer data as a list.

**property options**

Requests and returns the device options installed.

**ramp\_to\_current**(*target\_current, steps=30, pause=0.02*)

Ramps to a target current from the set current value over a certain number of linear steps, each separated by a pause duration.

**Parameters**

- **target\_current** – A current in Amps
- **steps** – An integer number of steps
- **pause** – A pause duration in seconds to wait between steps

**ramp\_to\_voltage**(*target\_voltage, steps=30, pause=0.02*)

Ramps to a target voltage from the set voltage value over a certain number of linear steps, each separated by a pause duration.

**Parameters**

- **target\_voltage** – A voltage in Amps
- **steps** – An integer number of steps
- **pause** – A pause duration in seconds to wait between steps

**reset()**

Resets the instrument and clears the queue.

**reset\_buffer()**

Resets the buffer.

**property resistance**

Reads the resistance in Ohms, if configured for this reading.

**property resistance\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the 2-wire resistance measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property resistance\_range**

A floating point property that controls the resistance range in Ohms, which can take values from 0 to 210 MOhms. Auto-range is disabled when this property is set.

**shutdown()**

Ensures that the current or voltage is turned to zero and disables the output.

**property source\_current**

A floating point property that controls the source current in Amps.

**property source\_current\_delay**

A floating point property that sets a manual delay for the source after the output is turned on before a measurement is taken. When this property is set, the auto delay is turned off. Valid values are between 0 [seconds] and 999.9999 [seconds].

**property source\_current\_delay\_auto**

A boolean property that enables or disables auto delay. Valid values are True and False.

**property source\_current\_range**

A floating point property that controls the source current range in Amps, which can take values between -1.05 and +1.05 A. Auto-range is disabled when this property is set.

**property source\_enabled**

Reads a boolean value that is True if the source is enabled.

**property source\_mode**

A string property that controls the source mode, which can take the values 'current' or 'voltage'. The convenience methods [apply\\_current\(\)](#) and [apply\\_voltage\(\)](#) can also be used.

**property source\_voltage**

A floating point property that controls the source voltage in Volts.

**property source\_voltage\_delay**

A floating point property that sets a manual delay for the source after the output is turned on before a measurement is taken. When this property is set, the auto delay is turned off. Valid values are between 0 [seconds] and 999.9999 [seconds].

**property source\_voltage\_delay\_auto**

A boolean property that enables or disables auto delay. Valid values are True and False.

**property source\_voltage\_range**

A floating point property that controls the source voltage range in Volts, which can take values from -210 to 210 V. Auto-range is disabled when this property is set.

**property standard\_devs**

Returns the calculated standard deviations for voltage, current, and resistance from the buffer data as a list.

**start\_buffer()**

Starts the buffer.

**property status**

Requests and returns the status byte and Master Summary Status bit.

**property std\_current**

Returns the current standard deviation from the buffer

**property std\_resistance**

Returns the resistance standard deviation from the buffer

**property std\_voltage**

Returns the voltage standard deviation from the buffer

**stop\_buffer()**

Aborts the buffering measurement, by stopping the measurement arming and triggering sequence. If possible, a Selected Device Clear (SDC) is used.

**triad(*base\_frequency*, *duration*)**

Sounds a musical triad using the system beep.

**Parameters**

- **base\_frequency** – A frequency in Hz between 65 Hz and 1.3 MHz
- **duration** – A time in seconds between 0 and 7.9 seconds

**trigger()**

Executes a bus trigger.

**use\_front\_terminals()**

Enables the front terminals for measurement, and disables the rear terminals.

**use\_rear\_terminals()**

Enables the rear terminals for measurement, and disables the front terminals.

**property voltage**

Reads the voltage in Volts, if configured for this reading.

**property voltage\_filter\_count**

A integer property that controls the number of readings that are acquired and stored in the filter buffer for the averaging

**property voltage\_filter\_type**

A String property that controls the filter's type for the current. REP : Repeating filter MOV : Moving filter

**property voltage\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the DC voltage measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property voltage\_output\_off\_state**

Select the output-off state of the SourceMeter. HIMP : output relay is open, disconnects external circuitry. NORM : V-Source is selected and set to 0V, Compliance is set to 0.5% full scale of the present current range. ZERO : V-Source is selected and set to 0V, compliance is set to the programmed Source I value or to 0.5% full scale of the present current range, whichever is greater. GUAR : I-Source is selected and set to 0A

**property voltage\_range**

A floating point property that controls the measurement voltage range in Volts, which can take values from -210 to 210 V. Auto-range is disabled when this property is set.

**wait\_for\_buffer**(*should\_stop=<function KeithleyBuffer.<lambda>>, timeout=60, interval=0.1*)

Blocks the program, waiting for a full buffer. This function returns early if the `should_stop` function returns True or the timeout is reached before the buffer is full.

**Parameters**

- **should\_stop** – A function that returns True when this function should return early
- **timeout** – A time in seconds after which this function should return early
- **interval** – A time in seconds for how often to check if the buffer is full

**property wires**

An integer property that controls the number of wires in use for resistance measurements, which can take the value of 2 or 4.

## 7.20.5 Keithley 2700 MultiMeter/Switch System

**class** `pymeasure.instruments.keithley.Keithley2700`(*adapter, \*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`, `pymeasure.instruments.keithley.buffer.KeithleyBuffer`

Represents the Keithley 2700 Multimeter/Switch System and provides a high-level interface for interacting with the instrument.

```
keithley = Keithley2700("GPIB::1")
```

**beep**(*frequency, duration*)

Sounds a system beep.

**Parameters**

- **frequency** – A frequency in Hz between 65 Hz and 2 MHz
- **duration** – A time in seconds between 0 and 7.9 seconds

**property buffer\_data**

Returns a numpy array of values from the buffer.

**property buffer\_points**

An integer property that controls the number of buffer points. This does not represent actual points in the buffer, but the configuration value instead.

**channels\_from\_rows\_columns**(*rows, columns, slot=None*)

Determine the channel numbers between column(s) and row(s) of the 7709 connection matrix. Returns a list of channel numbers. Only one of the parameters ‘rows’ or ‘columns’ can be “all”

**Parameters**

- **rows** – row number or list of numbers; can also be “all”
- **columns** – column number or list of numbers; can also be “all”
- **slot** – slot number (1 or 2) of the 7709 card to be used

**check\_errors**()

Logs any system errors reported by the instrument.

**close\_rows\_to\_columns**(*rows, columns, slot=None*)

Closes (connects) the channels between column(s) and row(s) of the 7709 connection matrix. Only one of the parameters ‘rows’ or ‘columns’ can be “all”

**Parameters**



- **rows** – row number or list of numbers; can also be “all”
- **columns** – column number or list of numbers; can also be “all”
- **slot** – slot number (1 or 2) of the 7709 card to be used

**property closed\_channels**

Parameter that controls the opened and closed channels. All mentioned channels are closed, other channels will be opened.

**property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device’s Output Queue when all pending selected device operations have been finished.

**config\_buffer**(*points=64, delay=0*)

Configures the measurement buffer for a number of points, to be taken with a specified delay.

**Parameters**

- **points** – The number of points in the buffer.
- **delay** – The delay time in seconds.

**determine\_valid\_channels()**

Determine what cards are installed into the Keithley 2700 and from that determine what channels are valid.

**disable\_buffer()**

Disables the connection between measurements and the buffer, but does not abort the measurement process.

**display\_closed\_channels()**

Show the presently closed channels on the display of the Keithley 2700.

**property display\_text**

A string property that controls the text shown on the display of the Keithley 2700. Text can be up to 12 ASCII characters and must be enabled to show.

**property error**

Returns a tuple of an error code and message from a single error.

**get\_state\_of\_channels**(*channels*)

Get the open or closed state of the specified channels

**Parameters** **channels** – a list of channel numbers, or single channel number

**property id**

Requests and returns the identification of the instrument.

**is\_buffer\_full()**

Returns True if the buffer is full of measurements.

**open\_all\_channels()**

Open all channels of the Keithley 2700.

**property open\_channels**

A parameter that opens the specified list of channels. Can only be set.

**open\_rows\_to\_columns**(*rows, columns, slot=None*)

Opens (disconnects) the channels between column(s) and row(s) of the 7709 connection matrix. Only one of the parameters ‘rows’ or ‘columns’ can be “all”

**Parameters**

- **rows** – row number or list of numbers; can also be “all”

- **columns** – column number or list of numbers; can also be “all”
- **slot** – slot number (1 or 2) of the 7709 card to be used

**property options**

Property that lists the installed cards in the Keithley 2700. Returns a dict with the integer card numbers on the position.

**reset()**

Resets the instrument and clears the queue.

**reset\_buffer()**

Resets the buffer.

**shutdown()**

Brings the instrument to a safe and stable state

**start\_buffer()**

Starts the buffer.

**property status**

Requests and returns the status byte and Master Summary Status bit.

**stop\_buffer()**

Aborts the buffering measurement, by stopping the measurement arming and triggering sequence. If possible, a Selected Device Clear (SDC) is used.

**property text\_enabled**

A boolean property that controls whether a text message can be shown on the display of the Keithley 2700.

**triad(*base\_frequency*, *duration*)**

Sounds a musical triad using the system beep.

**Parameters**

- **base\_frequency** – A frequency in Hz between 65 Hz and 1.3 MHz
- **duration** – A time in seconds between 0 and 7.9 seconds

**wait\_for\_buffer(*should\_stop*=<function KeithleyBuffer.<lambda>>, *timeout*=60, *interval*=0.1)**

Blocks the program, waiting for a full buffer. This function returns early if the `should_stop` function returns True or the timeout is reached before the buffer is full.

**Parameters**

- **should\_stop** – A function that returns True when this function should return early
- **timeout** – A time in seconds after which this function should return early
- **interval** – A time in seconds for how often to check if the buffer is full

## 7.20.6 Keithley 6221 AC and DC Current Source

```
class pymeasure.instruments.keithley.Keithley6221(adapter, **kwargs)
```

Bases: `pymeasure.instruments.instrument.Instrument`, `pymeasure.instruments.keithley.buffer.KeithleyBuffer`

Represents the Keithely 6221 AC and DC current source and provides a high-level interface for interacting with the instrument.

```

keithley = Keithley6221("GPIB::1")
keithley.clear()

# Use the keithley as an AC source
keithley.waveform_function = "square" # Set a square waveform
keithley.waveform_amplitude = 0.05    # Set the amplitude in Amps
keithley.waveform_offset = 0          # Set zero offset
keithley.source_compliance = 10       # Set compliance (limit) in V
keithley.waveform_dutycycle = 50      # Set duty cycle of wave in %
keithley.waveform_frequency = 347     # Set the frequency in Hz
keithley.waveform_ranging = "best"    # Set optimal output ranging
keithley.waveform_duration_cycles = 100 # Set duration of the waveform

# Link end of waveform to Service Request status bit
keithley.operation_event_enabled = 128 # OSB listens to end of wave
keithley.srq_event_enabled = 128      # SRQ listens to OSB

keithley.waveform_arm()               # Arm (load) the waveform

keithley.waveform_start()             # Start the waveform

keithley.adapter.wait_for_srq()       # Wait for the pulse to finish

keithley.waveform_abort()            # Disarm (unload) the waveform

keithley.shutdown()                  # Disables output

```

**beep**(frequency, duration)  
Sounds a system beep.

#### Parameters

- **frequency** – A frequency in Hz between 65 Hz and 2 MHz
- **duration** – A time in seconds between 0 and 7.9 seconds

**property buffer\_data**  
Returns a numpy array of values from the buffer.

**property buffer\_points**  
An integer property that controls the number of buffer points. This does not represent actual points in the buffer, but the configuration value instead.

**check\_errors**()  
Logs any system errors reported by the instrument.

**property complete**  
This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device's Output Queue when all pending selected device operations have been finished.

**config\_buffer**(points=64, delay=0)  
Configures the measurement buffer for a number of points, to be taken with a specified delay.

#### Parameters

- **points** – The number of points in the buffer.
- **delay** – The delay time in seconds.

**define\_arbitrary\_waveform**(datapoints, location=1)

Define the data points for the arbitrary waveform and copy the defined waveform into the given storage location.

**Parameters**

- **datapoints** – a list (or numpy array) of the data points; all values have to be between -1 and 1; 100 points maximum.
- **location** – integer storage location to store the waveform in. Value must be in range 1 to 4.

**disable\_buffer**()

Disables the connection between measurements and the buffer, but does not abort the measurement process.

**disable\_output\_trigger**()

Disables the output trigger for the Trigger layer

**disable\_source**()

Disables the source of current or voltage depending on the configuration of the instrument.

**property display\_enabled**

A boolean property that controls whether or not the display of the sourcemeter is enabled. Valid values are True and False.

**enable\_source**()

Enables the source of current or voltage depending on the configuration of the instrument.

**property error**

Returns a tuple of an error code and message from a single error.

**property id**

Requests and returns the identification of the instrument.

**is\_buffer\_full**()

Returns True if the buffer is full of measurements.

**property measurement\_event\_enabled**

An integer value that controls which measurement events are registered in the Measurement Summary Bit (MSB) status bit. Refer to the Model 6220/6221 Reference Manual for more information about programming the status bits.

**property measurement\_events**

An integer value that reads which measurement events have been registered in the Measurement event registers. Refer to the Model 6220/6221 Reference Manual for more information about programming the status bits. Reading this value clears the register.

**property operation\_event\_enabled**

An integer value that controls which operation events are registered in the Operation Summary Bit (OSB) status bit. Refer to the Model 6220/6221 Reference Manual for more information about programming the status bits.

**property operation\_events**

An integer value that reads which operation events have been registered in the Operation event registers. Refer to the Model 6220/6221 Reference Manual for more information about programming the status bits. Reading this value clears the register.

**property options**

Requests and returns the device options installed.

**property output\_low\_grounded**

A boolean property that controls whether the low output of the triax connection is connected to earth ground (True) or is floating (False).

**output\_trigger\_on\_external**(*line=1, after='DEL'*)

Configures the output trigger on the specified trigger link line number, with the option of supplying the part of the measurement after which the trigger should be generated (default to delay, which is right before the measurement)

**Parameters**

- **line** – A trigger line from 1 to 4
- **after** – An event string that determines when to trigger

**property questionable\_event\_enabled**

An integer value that controls which questionable events are registered in the Questionable Summary Bit (QSB) status bit. Refer to the Model 6220/6221 Reference Manual for more information about programming the status bits.

**property questionable\_events**

An integer value that reads which questionable events have been registered in the Questionable event registers. Refer to the Model 6220/6221 Reference Manual for more information about programming the status bits. Reading this value clears the register.

**reset()**

Resets the instrument and clears the queue.

**reset\_buffer()**

Resets the buffer.

**set\_timed\_arm**(*interval*)

Sets up the measurement to be taken with the internal trigger at a variable sampling rate defined by the interval in seconds between sampling points

**shutdown()**

Disables the output.

**property source\_auto\_range**

A boolean property that controls the auto range of the current source. Valid values are True or False.

**property source\_compliance**

A floating point property that controls the compliance of the current source in Volts. valid values are in range 0.1 [V] to 105 [V].

**property source\_current**

A floating point property that controls the source current in Amps.

**property source\_delay**

A floating point property that sets a manual delay for the source after the output is turned on before a measurement is taken. When this property is set, the auto delay is turned off. Valid values are between 1e-3 [seconds] and 999999.999 [seconds].

**property source\_enabled**

A boolean property that controls whether the source is enabled, takes values True or False. The convenience methods [enable\\_source\(\)](#) and [disable\\_source\(\)](#) can also be used.

**property source\_range**

A floating point property that controls the source current range in Amps, which can take values between -0.105 A and +0.105 A. Auto-range is disabled when this property is set.

**property `srq_event_enabled`**

An integer value that controls which event registers trigger the Service Request (SRQ) status bit. Refer to the Model 6220/6221 Reference Manual for more information about programming the status bits.

**property `standard_event_enabled`**

An integer value that controls which standard events are registered in the Event Summary Bit (ESB) status bit. Refer to the Model 6220/6221 Reference Manual for more information about programming the status bits.

**property `standard_events`**

An integer value that reads which standard events have been registered in the Standard event registers. Refer to the Model 6220/6221 Reference Manual for more information about programming the status bits. Reading this value clears the register.

**start\_buffer()**

Starts the buffer.

**property `status`**

Requests and returns the status byte and Master Summary Status bit.

**stop\_buffer()**

Aborts the buffering measurement, by stopping the measurement arming and triggering sequence. If possible, a Selected Device Clear (SDC) is used.

**triad(*base\_frequency*, *duration*)**

Sounds a musical triad using the system beep.

**Parameters**

- **base\_frequency** – A frequency in Hz between 65 Hz and 1.3 MHz
- **duration** – A time in seconds between 0 and 7.9 seconds

**trigger()**

Executes a bus trigger, which can be used when `trigger_on_bus()` is configured.

**trigger\_immediately()**

Configures measurements to be taken with the internal trigger at the maximum sampling rate.

**trigger\_on\_bus()**

Configures the trigger to detect events based on the bus trigger, which can be activated by `trigger()`.

**trigger\_on\_external(*line=1*)**

Configures the measurement trigger to be taken from a specific line of an external trigger

**Parameters** **line** – A trigger line from 1 to 4

**wait\_for\_buffer(*should\_stop=<function KeithleyBuffer.<lambda>>*, *timeout=60*, *interval=0.1*)**

Blocks the program, waiting for a full buffer. This function returns early if the `should_stop` function returns True or the timeout is reached before the buffer is full.

**Parameters**

- **should\_stop** – A function that returns True when this function should return early
- **timeout** – A time in seconds after which this function should return early
- **interval** – A time in seconds for how often to check if the buffer is full

**waveform\_abort()**

Abort the waveform output and disarm the waveform function.

**property waveform\_amplitude**

A floating point property that controls the (peak) amplitude of the waveform in Amps. Valid values are in range 2e-12 to 0.105.

**waveform\_arm()**

Arm the current waveform function.

**property waveform\_duration\_cycles**

A floating point property that controls the duration of the waveform in cycles. Valid values are in range 1e-3 to 99999999900.

**waveform\_duration\_set\_infinity()**

Set the waveform duration to infinity.

**property waveform\_duration\_time**

A floating point property that controls the duration of the waveform in seconds. Valid values are in range 100e-9 to 999999.999.

**property waveform\_dutycycle**

A floating point property that controls the duty-cycle of the waveform in percent for the square and ramp waves. Valid values are in range 0 to 100.

**property waveform\_frequency**

A floating point property that controls the frequency of the waveform in Hertz. Valid values are in range 1e-3 to 1e5.

**property waveform\_function**

A string property that controls the selected wave function. Valid values are “sine”, “ramp”, “square”, “arbitrary1”, “arbitrary2”, “arbitrary3” and “arbitrary4”.

**property waveform\_offset**

A floating point property that controls the offset of the waveform in Amps. Valid values are in range -0.105 to 0.105.

**property waveform\_ranging**

A string property that controls the source ranging of the waveform. Valid values are “best” and “fixed”.

**waveform\_start()**

Start the waveform output. Must already be armed

**property waveform\_use\_phasemarker**

A boolean property that controls whether the phase marker option is turned on or of. Valid values True (on) or False (off). Other settings for the phase marker have not yet been implemented.

## 7.20.7 Keithley 6517B Electrometer

**class** pymeasure.instruments.keithley.**Keithley6517B**(*adapter*, *\*\*kwargs*)

Bases: [pymeasure.instruments.instrument.Instrument](#), [pymeasure.instruments.keithley.buffer.KeithleyBuffer](#)

Represents the Keithley 6517B ElectroMeter and provides a high-level interface for interacting with the instrument.

```
keithley = Keithley6517B("GPIB::1")

keithley.apply_voltage()           # Sets up to source current
keithley.source_voltage_range = 200 # Sets the source voltage
                                   # range to 200 V
```

(continues on next page)

(continued from previous page)

```

keithley.source_voltage = 20      # Sets the source voltage to 20 V
keithley.enable_source()          # Enables the source output

keithley.measure_resistance()     # Sets up to measure resistance

keithley.ramp_to_voltage(50)      # Ramps the voltage to 50 V
print(keithley.resistance)        # Prints the resistance in Ohms

keithley.shutdown()              # Ramps the voltage to 0 V
                                # and disables output

```

**apply\_voltage**(*voltage\_range=None*)

Configures the instrument to apply a source voltage, and uses an auto range unless a voltage range is specified.

**Parameters** **voltage\_range** – A *voltage\_range* value or None (activates auto range)

**auto\_range\_source**()

Configures the source to use an automatic range.

**property buffer\_data**

Returns a numpy array of values from the buffer.

**property buffer\_points**

An integer property that controls the number of buffer points. This does not represent actual points in the buffer, but the configuration value instead.

**check\_errors**()

Logs any system errors reported by the instrument.

**property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device's Output Queue when all pending selected device operations have been finished.

**config\_buffer**(*points=64, delay=0*)

Configures the measurement buffer for a number of points, to be taken with a specified delay.

**Parameters**

- **points** – The number of points in the buffer.
- **delay** – The delay time in seconds.

**property current**

Reads the current in Amps, if configured for this reading.

**property current\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the DC current measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property current\_range**

A floating point property that controls the measurement current range in Amps, which can take values between -20 and +20 mA. Auto-range is disabled when this property is set.

**disable\_buffer**()

Disables the connection between measurements and the buffer, but does not abort the measurement process.



**disable\_source()**

Disables the source of current or voltage depending on the configuration of the instrument.

**enable\_source()**

Enables the source of current or voltage depending on the configuration of the instrument.

**property error**

Returns a tuple of an error code and message from a single error.

**property id**

Requests and returns the identification of the instrument.

**is\_buffer\_full()**

Returns True if the buffer is full of measurements.

**measure\_current(*nplc=1, current=0.000105, auto\_range=True*)**

Configures the measurement of current.

**Parameters**

- **nplc** – Number of power line cycles (NPLC) from 0.01 to 10
- **current** – Upper limit of current in Amps, from -21 mA to 21 mA
- **auto\_range** – Enables auto\_range if True, else uses the current\_range attribut

**measure\_resistance(*nplc=1, resistance=210000.0, auto\_range=True*)**

Configures the measurement of resistance.

**Parameters**

- **nplc** – Number of power line cycles (NPLC) from 0.01 to 10
- **resistance** – Upper limit of resistance in Ohms, from -210 POhms to 210 POhms
- **auto\_range** – Enables auto\_range if True, else uses the resistance\_range attribut

**measure\_voltage(*nplc=1, voltage=21.0, auto\_range=True*)**

Configures the measurement of voltage.

**Parameters**

- **nplc** – Number of power line cycles (NPLC) from 0.01 to 10
- **voltage** – Upper limit of voltage in Volts, from -1000 V to 1000 V
- **auto\_range** – Enables auto\_range if True, else uses the voltage\_range attribut

**property options**

Requests and returns the device options installed.

**ramp\_to\_voltage(*target\_voltage, steps=30, pause=0.02*)**

Ramps to a target voltage from the set voltage value over a certain number of linear steps, each separated by a pause duration.

**Parameters**

- **target\_voltage** – A voltage in Volts
- **steps** – An integer number of steps
- **pause** – A pause duration in seconds to wait between steps

**reset()**

Resets the instrument and clears the queue.

**reset\_buffer()**

Resets the buffer.

**property resistance**

Reads the resistance in Ohms, if configured for this reading.

**property resistance\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the 2-wire resistance measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property resistance\_range**

A floating point property that controls the resistance range in Ohms, which can take values from 0 to 100e18 Ohms. Auto-range is disabled when this property is set.

**shutdown()**

Ensures that the current or voltage is turned to zero and disables the output.

**property source\_current\_resistance\_limit**

Boolean property which enables or disables resistance current limit

**property source\_enabled**

Reads a boolean value that is True if the source is enabled.

**property source\_voltage**

A floating point property that controls the source voltage in Volts.

**property source\_voltage\_range**

A floating point property that controls the source voltage range in Volts, which can take values from -1000 to 1000 V. Auto-range is disabled when this property is set.

**start\_buffer()**

Starts the buffer.

**property status**

Requests and returns the status byte and Master Summary Status bit.

**stop\_buffer()**

Aborts the buffering measurement, by stopping the measurement arming and triggering sequence. If possible, a Selected Device Clear (SDC) is used.

**trigger()**

Executes a bus trigger, which can be used when [trigger\\_on\\_bus\(\)](#) is configured.

**trigger\_immediately()**

Configures measurements to be taken with the internal trigger at the maximum sampling rate.

**trigger\_on\_bus()**

Configures the trigger to detect events based on the bus trigger, which can be activated by [trigger\(\)](#).

**property voltage**

Reads the voltage in Volts, if configured for this reading.

**property voltage\_nplc**

A floating point property that controls the number of power line cycles (NPLC) for the DC voltage measurements, which sets the integration period and measurement speed. Takes values from 0.01 to 10, where 0.1, 1, and 10 are Fast, Medium, and Slow respectively.

**property voltage\_range**

A floating point property that controls the measurement voltage range in Volts, which can take values from -1000 to 1000 V. Auto-range is disabled when this property is set.

**wait\_for\_buffer**(*should\_stop=<function KeithleyBuffer.<lambda>>, timeout=60, interval=0.1*)  
 Blocks the program, waiting for a full buffer. This function returns early if the `should_stop` function returns True or the timeout is reached before the buffer is full.

**Parameters**

- **should\_stop** – A function that returns True when this function should return early
- **timeout** – A time in seconds after which this function should return early
- **interval** – A time in seconds for how often to check if the buffer is full

## 7.20.8 Keithley 2750 Multimeter/Switch System

**class** `pymeasure.instruments.keithley.Keithley2750`(*adapter, \*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Keithley2750 multimeter/switch system and provides a high-level interface for interacting with the instrument.

**check\_errors()**

Read all errors from the instrument.

**Returns** list of error entries

**close**(*channel*)

Closes (connects) the specified channel.

**Parameters** **channel** (*int*) – 3-digit number for the channel

**Returns** None

**property closed\_channels**

Reads the list of closed channels

**property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device's Output Queue when all pending selected device operations have been finished.

**property id**

Requests and returns the identification of the instrument.

**open**(*channel*)

Opens (disconnects) the specified channel.

**Parameters** **channel** (*int*) – 3-digit number for the channel

**Returns** None

**open\_all()**

Opens (disconnects) all the channels on the switch matrix.

**Returns** None

**property options**

Requests and returns the device options installed.

**reset()**

Resets the instrument.

**shutdown()**

Brings the instrument to a safe and stable state

**property status**

Requests and returns the status byte and Master Summary Status bit.

## 7.20.9 Keithley 2600 SourceMeter

**class** `pymeasure.instruments.keithley.Keithley2600(adapter, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Keithley 2600 series (channel A and B) SourceMeter

**check\_errors()**

Logs any system errors reported by the instrument.

**property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device's Output Queue when all pending selected device operations have been finished.

**property error**

Returns a tuple of an error code and message from a single error.

**property id**

Requests and returns the identification of the instrument.

**property options**

Requests and returns the device options installed.

**reset()**

Resets the instrument.

**shutdown()**

Brings the instrument to a safe and stable state

**property status**

Requests and returns the status byte and Master Summary Status bit.

## 7.21 Keysight

This section contains specific documentation on the keysight instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.21.1 Keysight DSOX1102G Oscilloscope

**class** `pymeasure.instruments.keysight.KeysightDSOX1102G(adapter, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Keysight DSOX1102G Oscilloscope interface for interacting with the instrument.

Refer to the Keysight DSOX1102G Oscilloscope Programmer's Guide for further details about using the lower-level methods to interact directly with the scope.

```
scope = KeysightDSOX1102G(resource)
scope.autoscale()
ch1_data_array, ch1_preamble = scope.download_data(source="channel1", points=2000)
# ...
scope.shutdown()
```

Known issues:

- The digitize command will be completed before the operation is. May lead to VI\_ERROR\_TMO (timeout) occurring when sending commands immediately after digitize. Current fix: if deemed necessary, add delay between digitize and follow-up command to scope.

#### **property acquisition\_mode**

A string parameter that sets the acquisition mode. Can be “realtime” or “segmented”.

#### **property acquisition\_type**

A string parameter that sets the type of data acquisition. Can be “normal”, “average”, “hresolution”, or “peak”.

#### **ask(command)**

Writes the command to the instrument through the adapter and returns the read response.

**Parameters** **command** – command string to be sent to the instrument

#### **autoscale()**

Autoscale displayed channels.

#### **check\_errors()**

Read all errors from the instrument.

**Returns** list of error entries

#### **clear()**

Clears the instrument status byte

#### **clear\_status()**

Clear device status.

#### **property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device’s Output Queue when all pending selected device operations have been finished.

**static control**(*get\_command, set\_command, docs, validator=<function Instrument.<lambda>>, values=(), map\_values=False, get\_process=<function Instrument.<lambda>>, set\_process=<function Instrument.<lambda>>, check\_set\_errors=False, check\_get\_errors=False, \*\*kwargs*)

Returns a property for the class based on the supplied commands. This property may be set and read from the instrument.

#### **Parameters**

- **get\_command** – A string command that asks for the value
- **set\_command** – A string command that writes the value
- **docs** – A docstring that will be included in the documentation
- **validator** – A function that takes both a value and a group of valid values and returns a valid value, while it otherwise raises an exception
- **values** – A list, tuple, range, or dictionary of valid values, that can be used as to map values if **map\_values** is True.
- **map\_values** – A boolean flag that determines if the values should be interpreted as a map
- **get\_process** – A function that take a value and allows processing before value mapping, returning the processed value

- **set\_process** – A function that takes a value and allows processing before value mapping, returning the processed value
- **check\_set\_errors** – Toggles checking errors after setting
- **check\_get\_errors** – Toggles checking errors after getting

**default\_setup()**

Default setup, some user settings (like preferences) remain unchanged.

**digitize(source: str)**

Acquire waveforms according to the settings of the :ACQUIRE commands. Ensure a delay between the digitize operation and further commands, as timeout may be reached before digitize has completed. :param source: “channel1”, “channel2”, “function”, “math”, “fft”, “abus”, or “ext”.

**download\_data(source, points=62500)**

Get data from specified source of oscilloscope. Returned objects are a np.ndarray of data values (no temporal axis) and a dict of the waveform preamble, which can be used to build the corresponding time values for all data points.

Multimeter will be stopped for proper acquisition.

**Parameters**

- **source** – measurement source, can be “channel1”, “channel2”, “function”, “fft”, “wmemory1”, “wmemory2”, or “ext”.
- **points** – integer number of points to acquire. Note that oscilloscope may return less points than specified, this is not an issue of this library. Can be 100, 250, 500, 1000, 2000, 5000, 10000, 20000, 50000, or 62500.

**Return data\_ndarray, waveform\_preamble\_dict** see waveform\_preamble property for dict format.

**download\_image(format\_='png', color\_palette='color')**

Get image of oscilloscope screen in bytearray of specified file format.

**Parameters**

- **format** – “bmp”, “bmp8bit”, or “png”
- **color\_palette** – “color” or “grayscale”

**factory\_reset()**

Factory default setup, no user settings remain unchanged.

**property id**

Requests and returns the identification of the instrument.

**static measurement**(get\_command, docs, values=(), map\_values=None, get\_process=<function Instrument.<lambda>>, command\_process=<function Instrument.<lambda>>, check\_get\_errors=False, \*\*kwargs)

Returns a property for the class based on the supplied commands. This is a measurement quantity that may only be read from the instrument, not set.

**Parameters**

- **get\_command** – A string command that asks for the value
- **docs** – A docstring that will be included in the documentation
- **values** – A list, tuple, range, or dictionary of valid values, that can be used as to map values if map\_values is True.
- **map\_values** – A boolean flag that determines if the values should be interpreted as a map

- **get\_process** – A function that take a value and allows processing before value mapping, returning the processed value
- **command\_process** – A function that take a command and allows processing before executing the command, for both getting and setting
- **check\_get\_errors** – Toggles checking errors after getting

**property options**

Requests and returns the device options installed.

**read()**

Reads from the instrument through the adapter and returns the response.

**reset()**

Resets the instrument.

**run()**

Starts repetitive acquisitions. This is the same as pressing the Run key on the front panel.

**static setting**(*set\_command*, *docs*, *validator*=<function Instrument.<lambda>>, *values*=(), *map\_values*=False, *set\_process*=<function Instrument.<lambda>>, *check\_set\_errors*=False, *\*\*kwargs*)

Returns a property for the class based on the supplied commands. This property may be set, but raises an exception when being read from the instrument.

**Parameters**

- **set\_command** – A string command that writes the value
- **docs** – A docstring that will be included in the documentation
- **validator** – A function that takes both a value and a group of valid values and returns a valid value, while it otherwise raises an exception
- **values** – A list, tuple, range, or dictionary of valid values, that can be used as to map values if *map\_values* is True.
- **map\_values** – A boolean flag that determines if the values should be interpreted as a map
- **set\_process** – A function that takes a value and allows processing before value mapping, returning the processed value
- **check\_set\_errors** – Toggles checking errors after setting

**shutdown()**

Brings the instrument to a safe and stable state

**single()**

Causes the instrument to acquire a single trigger of data. This is the same as pressing the Single key on the front panel.

**property status**

Requests and returns the status byte and Master Summary Status bit.

**stop()**

Stops the acquisition. This is the same as pressing the Stop key on the front panel.

**property system\_setup**

A string parameter that sets up the oscilloscope. Must be in IEEE 488.2 format. It is recommended to only set a string previously obtained from this command.

**property timebase**

Read timebase setup as a dict containing the following keys: - “REF”: position on screen of timebase

reference (str) - “MAIN:RANG”: full-scale timebase range (float) - “POS”: interval between trigger and reference point (float) - “MODE”: mode (str)

**property timebase\_mode**

A string parameter that sets the current time base. Can be “main”, “window”, “xy”, or “roll”.

**property timebase\_offset**

A float parameter that sets the time interval in seconds between the trigger event and the reference position (at center of screen by default).

**property timebase\_range**

A float parameter that sets the full-scale horizontal time in seconds for the main window.

**property timebase\_scale**

A float parameter that sets the horizontal scale (units per division) in seconds for the main window.

**timebase\_setup**(*mode=None, offset=None, horizontal\_range=None, scale=None*)

Set up timebase. Unspecified parameters are not modified. Modifying a single parameter might impact other parameters. Refer to oscilloscope documentation and make multiple consecutive calls to `channel_setup` if needed.

**Parameters**

- **mode** – Timebase mode, can be “main”, “window”, “xy”, or “roll”.
- **offset** – Offset in seconds between trigger and center of screen.
- **horizontal\_range** – Full-scale range in seconds.
- **scale** – Units-per-division in seconds.

**values**(*command, \*\*kwargs*)

Reads a set of values from the instrument through the adapter, passing on any key-word arguments.

**property waveform\_data**

Get the binary block of sampled data points transmitted using the IEEE 488.2 arbitrary block data format.

**property waveform\_format**

A string parameter that controls how the data is formatted when sent from the oscilloscope. Can be “ascii”, “word” or “byte”. Words are transmitted in big endian by default.

**property waveform\_points**

An integer parameter that sets the number of waveform points to be transferred with the `waveform_data` method. Can be any of the following values: 100, 250, 500, 1000, 2 000, 5 000, 10 000, 20 000, 50 000, 62 500.

Note that the oscilloscope may provide less than the specified nb of points.

**property waveform\_points\_mode**

A string parameter that sets the data record to be transferred with the `waveform_data` method. Can be “normal”, “maximum”, or “raw”.

**property waveform\_preamble**

Get preamble information for the selected waveform source as a dict with the following keys: - “format”: byte, word, or ascii (str) - “type”: normal, peak detect, or average (str) - “points”: nb of data points transferred (int) - “count”: always 1 (int) - “xincrement”: time difference between data points (float) - “xorigin”: first data point in memory (float) - “xreference”: data point associated with xorigin (int) - “yincrement”: voltage difference between data points (float) - “yorigin”: voltage at center of screen (float) - “yreference”: data point associated with yorigin (int)

**property waveform\_source**

A string parameter that selects the analog channel, function, or reference waveform to be used as the source



for the waveform methods. Can be “channel1”, “channel2”, “function”, “fft”, “wmemory1”, “wmemory2”, or “ext”.

**write**(*command*)

Writes the command to the instrument through the adapter.

**Parameters** **command** – command string to be sent to the instrument

## 7.21.2 Keysight N5767A Power Supply

**class** `pymeasure.instruments.keysight.KeysightN5767A(adapter, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Keysight N5767A Power supply interface for interacting with the instrument.

**check\_errors**()

Read all errors from the instrument.

**Returns** list of error entries

**property complete**

This property allows synchronization between a controller and a device. The Operation Complete query places an ASCII character 1 into the device’s Output Queue when all pending selected device operations have been finished.

**property current**

Reads a setting current in Amps.

**property current\_range**

A floating point property that controls the DC current range in Amps, which can take values from 0 to 25 A. Auto-range is disabled when this property is set.

**disable**()

Disables the flow of current.

**enable**()

Enables the flow of current.

**property id**

Requests and returns the identification of the instrument.

**is\_enabled**()

Returns True if the current supply is enabled.

**property options**

Requests and returns the device options installed.

**reset**()

Resets the instrument.

**shutdown**()

Brings the instrument to a safe and stable state

**property status**

Requests and returns the status byte and Master Summary Status bit.

**property voltage**

Reads a DC voltage measurement in Volts.

**property voltage\_range**

A floating point property that controls the DC voltage range in Volts, which can take values from 0 to 60 V. Auto-range is disabled when this property is set.

## 7.22 Lake Shore Cryogenics

This section contains specific documentation on the Lake Shore Cryogenics instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.22.1 Lake Shore Adapters

**class** `pymeasure.instruments.lakeshore.LakeShoreUSBAdapter(port)`

Bases: `pymeasure.adapters.serial.SerialAdapter`

Provides a `SerialAdapter` with the specific baudrate, timeout, parity, and byte size for LakeShore USB communication.

Initiates the adapter to open serial communication over the supplied port.

**Parameters** `port` – A string representing the serial port

**`_format_binary_values(values, datatype='f', is_big_endian=False, header_fmt='ieee')`**

Format values in binary format, used internally in `write_binary_values()`.

**Parameters**

- **values** – data to be written to the device.
- **datatype** – the format string for a single element. See struct module.
- **is\_big\_endian** – boolean indicating endianness.
- **header\_fmt** – Format of the header prefixing the data (“ieee”, “hp”, “empty”).

**Returns** binary string.

**Return type** bytes

**`ask(command)`**

Writes the command to the instrument and returns the resulting ASCII response

**Parameters** `command` – SCPI command string to be sent to the instrument

**Returns** String ASCII response of the instrument

**`binary_values(command, header_bytes=0, dtype=<class 'numpy.float32'>)`**

Returns a numpy array from a query for binary data

**Parameters**

- **command** – SCPI command to be sent to the instrument
- **header\_bytes** – Integer number of bytes to ignore in header
- **dtype** – The NumPy data type to format the values with

**Returns** NumPy array of values

**`read()`**

Reads until the buffer is empty and returns the resulting ASCII response

**Returns** String ASCII response of the instrument.

**values**(*command*, *separator*=' ', *cast*=<class 'float'>, *preprocess\_reply*=None)

Writes a command to the instrument and returns a list of formatted values from the result

#### Parameters

- **command** – SCPI command to be sent to the instrument
- **separator** – A separator character to split the string into a list
- **cast** – A type to cast the result
- **preprocess\_reply** – optional callable used to preprocess values received from the instrument. The callable returns the processed string. If not specified, the Adapter default is used if available, otherwise no preprocessing is done.

**Returns** A list of the desired type, or strings where the casting fails

**write**(*command*)

Overwrites the [SerialAdapter.write](#) method to automatically append a Unix-style linebreak at the end of the command.

**Parameters** **command** – SCPI command string to be sent to the instrument

**write\_binary\_values**(*command*, *values*, *\*\*kwargs*)

Write binary data to the instrument, e.g. waveform for signal generators

#### Parameters

- **command** – SCPI command to be sent to the instrument
- **values** – iterable representing the binary values
- **kwargs** – Key-word arguments to pass onto [\\_format\\_binary\\_values\(\)](#)

**Returns** number of bytes written

## 7.22.2 Lake Shore 331 Temperature Controller

**class** `pymeasure.instruments.lakeshore.LakeShore331`(*adapter*, *\*\*kwargs*)

Bases: [pymeasure.instruments.instrument.Instrument](#)

Represents the Lake Shore 331 Temperature Controller and provides a high-level interface for interacting with the instrument.

```
controller = LakeShore331("GPIB::1")

print(controller.setpoint_1)      # Print the current setpoint for loop 1
controller.setpoint_1 = 50        # Change the setpoint to 50 K
controller.heater_range = 'low'   # Change the heater range to Low
controller.wait_for_temperature() # Wait for the temperature to stabilize
print(controller.temperature_A)    # Print the temperature at sensor A
```

**disable\_heater**()

Turns the [heater\\_range](#) to off to disable the heater.

**property heater\_range**

A string property that controls the heater range, which can take the values: off, low, medium, and high. These values correlate to 0, 0.5, 5 and 50 W respectively.

**property setpoint\_1**

A floating point property that controls the setpoint temperature in Kelvin for Loop 1.

**property setpoint\_2**

A floating point property that controls the setpoint temperature in Kelvin for Loop 2.

**property temperature\_A**

Reads the temperature of the sensor A in Kelvin.

**property temperature\_B**

Reads the temperature of the sensor B in Kelvin.

**wait\_for\_temperature**(*accuracy=0.1, interval=0.1, sensor='A', setpoint=1, timeout=360,*  
*should\_stop=<function LakeShore331.<lambda>>>)*

Blocks the program, waiting for the temperature to reach the setpoint within the accuracy (%), checking this each interval time in seconds.

**Parameters**

- **accuracy** – An acceptable percentage deviation between the setpoint and temperature
- **interval** – A time in seconds that controls the refresh rate
- **sensor** – The desired sensor to read, either A or B
- **setpoint** – The desired setpoint loop to read, either 1 or 2
- **timeout** – A timeout in seconds after which an exception is raised
- **should\_stop** – A function that returns True if waiting should stop, by default this always returns False

### 7.22.3 Lake Shore 425 Gaussmeter

**class** `pymeasure.instruments.lakeshore.LakeShore425`(*port*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the LakeShore 425 Gaussmeter and provides a high-level interface for interacting with the instrument

To allow user access to the LakeShore 425 Gaussmeter in Linux, create the file: `/etc/udev/rules.d/52-lakeshore425.rules`, with contents:

```
SUBSYSTEMS=="usb",ATTRS{idVendor}=="1fb9",ATTRS{idProduct}=="0401",MODE="0666",
↳SYMLINK+="lakeshore425"
```

Then reload the udev rules with:

```
sudo udevadm control --reload-rules
sudo udevadm trigger
```

The device will be accessible through `/dev/lakeshore425`.

**ac\_mode**(*wideband=True*)

Sets up a measurement of an oscillating (AC) field

**auto\_range**()

Sets the field range to automatically adjust

**dc\_mode**(*wideband=True*)

Sets up a steady-state (DC) measurement of the field

**property field**

Returns the field in the current units

**measure**(*points*, *has\_aborted*=<function LakeShore425.<lambda>>, *delay*=0.001)

Returns the mean and standard deviation of a given number of points while blocking

**property range**

A floating point property that controls the field range in units of Gauss, which can take the values 35, 350, 3500, and 35,000 G.

**property unit**

A string property that controls the units of the instrument, which can take the values of G, T, Oe, or A/m.

**zero\_probe()**

Initiates the zero field sequence to calibrate the probe

## 7.23 Newport

This section contains specific documentation on the Newport instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.23.1 ESP 300 Motion Controller

**class** `pymeasure.instruments.newport.ESP300`(*resourceName*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Newport ESP 300 Motion Controller and provides a high-level for interacting with the instrument.

By default this instrument is constructed with x, y, and phi attributes that represent axes 1, 2, and 3. Custom implementations can overwrite this depending on the available axes. Axes are controlled through an [Axis](#) class.

**property axes**

A list of the [Axis](#) objects that are present.

**clear\_errors()**

Clears the error messages by checking until a 0 code is received.

**disable()**

Disables all of the axes associated with this controller.

**enable()**

Enables all of the axes associated with this controller.

**property error**

Reads an error code from the motion controller.

**property errors**

Returns a list of error Exceptions that can be later raised, or used to diagnose the situation.

**shutdown()**

Shuts down the controller by disabling all of the axes.

**class** `pymeasure.instruments.newport.esp300.Axis`(*axis*, *controller*)

Bases: `object`

Represents an axis of the Newport ESP300 Motor Controller, which can have independent parameters from the other axes.

**define\_position**(*position*)

Overwrites the value of the current position with the given value.

**disable()**

Disables motion for the axis.

**enable()**

Enables motion for the axis.

**property enabled**

Returns a boolean value that is True if the motion for this axis is enabled.

**home(*type=1*)**

Drives the axis to the home position, which may be the negative hardware limit for some actuators (e.g. LTA-HS). *type* can take integer values from 0 to 6.

**property left\_limit**

A floating point property that controls the left software limit of the axis.

**property motion\_done**

Returns a boolean that is True if the motion is finished.

**property position**

A floating point property that controls the position of the axis. The units are defined based on the actuator. Use the [`wait\_for\_stop\(\)`](#) method to ensure the position is stable.

**property right\_limit**

A floating point property that controls the right software limit of the axis.

**property units**

A string property that controls the displacement units of the axis, which can take values of: encoder count, motor step, millimeter, micrometer, inches, milli-inches, micro-inches, degree, gradient, radian, milliradian, and microradian.

**wait\_for\_stop(*delay=0, interval=0.05*)**

Blocks the program until the motion is completed. A further delay can be specified in seconds.

**zero()**

Resets the axis position to be zero at the current position.

**class** `pymeasure.instruments.newport.esp300.AxisError(code)`

Bases: `Exception`

Raised when a particular axis causes an error for the Newport ESP300.

**class** `pymeasure.instruments.newport.esp300.GeneralError(code)`

Bases: `Exception`

Raised when the Newport ESP300 has a general error.

## 7.24 National Instruments

This section contains specific documentation on the National Instruments instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

## 7.24.1 NI Virtual Bench

### General Information

The `armstrap/pyvirtualbench` Python wrapper for the VirtualBench C-API is required. This Instrument driver only interfaces the `pyvirtualbench` Python wrapper.

### Examples

To be documented. Check the examples in the `pyvirtualbench` repository to get an idea.

Simple Example to switch digital lines of the DIO module.

```
from pymeasure.instruments.ni import VirtualBench

vb = VirtualBench(device_name='VB8012-3057E1C')
line = 'dig/2' # may be list of lines
# initialize DIO module -> available via vb.dio
vb.acquire_digital_input_output(line, reset=False)

vb.dio.write(self.line, {True})
sleep(1000)
vb.dio.write(self.line, {False})

vb.shutdown()
```

### Instrument Class

`class pymeasure.instruments.ni.virtualbench.VirtualBench(device_name="", name='VirtualBench')`  
 Bases: object

Represents National Instruments Virtual Bench main frame.

Subclasses implement the functionalities of the different modules:

- Mixed-Signal-Oscilloscope (MSO)
- Digital Input Output (DIO)
- Function Generator (FGEN)
- Power Supply (PS)
- Serial Peripheral Interface (SPI) -> not implemented for pymeasure yet
- Inter Integrated Circuit (I2C) -> not implemented for pymeasure yet

For every module exist methods to save/load the configuration to file. These methods are not wrapped so far, checkout the `pyvirtualbench` file.

All calibration methods and classes are not wrapped so far, since these are not required on a very regular basis. Also the connections via network are not yet implemented. Check the `pyvirtualbench` file, if you need the functionality.

#### Parameters

- **device\_name** (*str*) – Full unique device name
- **name** (*str*) – Name for display in pymeasure

**class DigitalInputOutput**(*virtualbench, lines, reset, vb\_name=""*)

Bases: `pymeasure.instruments.ni.virtualbench.VirtualBench.VirtualBenchInstrument`

Represents Digital Input Output (DIO) Module of Virtual Bench device. Allows to read/write digital channels and/or set channels to export the start signal of FGGEN module or trigger of MSO module.

**export\_signal**(*line, digitalSignalSource*)

Exports a signal to the specified line.

**Parameters**

- **line** (*str*) – Line string
- **digitalSignalSource** (*int*) – 0 for FGGEN start or 1 for MSO trigger

**query\_export\_signal**(*line*)

Indicates the signal being exported on the specified line.

**Parameters** **line** (*str*) – Line string

**Returns** Exported signal (FGGEN start or MSO trigger)

**Return type** enum

**query\_line\_configuration**()

Indicates the current line configurations. Tristate Lines, Static Lines, and Export Lines contain comma-separated range\_data and/or colon-delimited lists of all acquired lines

**read**(*lines*)

Reads the current state of the specified lines.

**Parameters** **lines** (*str*) – Line string, requires full name specification e.g. 'VB8012-xxxxxxx/dig/0:7' since instrument\_handle is not required (only library\_handle)

**Returns** List of line states (HIGH/LOW)

**Return type** list

**reset\_instrument**()

Resets the session configuration to default values, and resets the device and driver software to a known state.

**shutdown**()

Removes the session and deallocates any resources acquired during the session. If output is enabled on any channels, they remain in their current state.

**tristate\_lines**(*lines*)

Sets all specified lines to a high-impedance state. (Default)

**validate\_lines**(*lines, return\_single\_lines=False, validate\_init=False*)

Validate lines string Allowed patterns (case sensitive):

- 'VBxxxx-xxxxxxx/dig/0:7'
- 'VBxxxx-xxxxxxx/dig/0'
- 'dig/0'
- 'VBxxxx-xxxxxxx/trig'
- 'trig'

Allowed Line Numbers: 0-7 or trig

**Parameters**

- **lines** (*str*) – Line string to test
- **return\_single\_lines** (*bool, optional*) – Return list of line numbers as well, defaults to False
- **validate\_init** (*bool, optional*) – Check if lines are initialized (in self.\_line\_numbers), defaults to False

**Returns** Line string, optional list of single line numbers

**Return type** str, optional (str, list)



**write**(*lines, data*)

Writes data to the specified lines.

**Parameters**

- **lines** (*str*) – Line string
- **data** (*list or tuple*) – List of data, (True = High, False = Low)

**class DigitalMultimeter**(*virtualbench, reset, vb\_name=""*)

Bases: `pymeasure.instruments.ni.virtualbench.VirtualBench.VirtualBenchInstrument`

Represents Digital Multimeter (DMM) Module of Virtual Bench device. Allows to measure either DC/AC voltage or current, Resistance or Diodes.

**configure\_ac\_current**(*auto\_range\_terminal*)

Configure auto range terminal for AC current measurement

**Parameters** **auto\_range\_terminal** – Terminal to perform auto ranging ('LOW' or 'HIGH')

**configure\_dc\_current**(*auto\_range\_terminal*)

Configure auto range terminal for DC current measurement

**Parameters** **auto\_range\_terminal** – Terminal to perform auto ranging ('LOW' or 'HIGH')

**configure\_dc\_voltage**(*dmm\_input\_resistance*)

Configure DC voltage input resistance

**Parameters** **dmm\_input\_resistance** (*int or str*) – Input resistance ('TEN\_MEGA\_OHM' or 'TEN\_GIGA\_OHM')

**configure\_measurement**(*dmm\_function, auto\_range=True, manual\_range=1.0*)

Configure Instrument to take a DMM measurement

**Parameters**

- **function index or name** (*dmm\_function:DMM*) –
  - 'DC\_VOLTS', 'AC\_VOLTS'
  - 'DC\_CURRENT', 'AC\_CURRENT'
  - 'RESISTANCE'
  - 'DIODE'
- **auto\_range** (*bool*) – Enable/Disable auto ranging
- **manual\_range** (*float*) – Manually set measurement range

**query\_ac\_current**()

Indicates auto range terminal for AC current measurement

**query\_dc\_current**()

Indicates auto range terminal for DC current measurement

**query\_dc\_voltage**()

Indicates input resistance setting for DC voltage measurement

**query\_measurement**()

Query DMM measurement settings from the instrument

**Returns** Auto range, range data

**Return type** (bool, float)

**read**()

Read measurement value from the instrument

**Returns** Measurement value

**Return type** float

**reset\_instrument()**

Reset the DMM module to defaults

**shutdown()**

Removes the session and deallocates any resources acquired during the session. If output is enabled on any channels, they remain in their current state.

**validate\_auto\_range\_terminal**(*auto\_range\_terminal*)

Check value for choosing the auto range terminal for DC current measurement

**Parameters** **auto\_range\_terminal** (*int* or *str*) – Terminal to perform auto ranging ('LOW' or 'HIGH')

**Returns** Auto range terminal to pass to the instrument

**Return type** *int*

**validate\_dmm\_function**(*dmm\_function*)

Check if DMM function *dmm\_function* exists

**Parameters** **dmm\_function** (*int* or *str*) – DMM function index or name:

- 'DC\_VOLTS', 'AC\_VOLTS'
- 'DC\_CURRENT', 'AC\_CURRENT'
- 'RESISTANCE'
- 'DIODE'

**Returns** DMM function index to pass to the instrument

**Return type** *int*

**static validate\_range**(*dmm\_function*, *range*)

Checks if *range* is valid for the chosen *dmm\_function*

**Parameters**

- **dmm\_function** (*int*) – DMM Function
- **range** (*int* or *float*) – Range value, e.g. maximum value to measure

**Returns** Range value to pass to instrument

**Return type** *int*

**class FunctionGenerator**(*virtualbench*, *reset*, *vb\_name=""*)

Bases: `pymeasure.instruments.ni.virtualbench.VirtualBench.VirtualBenchInstrument`

Represents Function Generator (FGEN) Module of Virtual Bench device.

**configure\_arbitrary\_waveform**(*waveform*, *sample\_period*)

Configures the instrument to output a waveform. The waveform is output either after the end of the current waveform if output is enabled, or immediately after output is enabled.

**Parameters**

- **waveform** (*list*) – Waveform as list of values
- **sample\_period** (*float*) – Time between two waveform points (maximum of 125MS/s, which equals 80ns)

**configure\_arbitrary\_waveform\_gain\_and\_offset**(*gain*, *dc\_offset*)

Configures the instrument to output an arbitrary waveform with a specified gain and offset value. The waveform is output either after the end of the current waveform if output is enabled, or immediately after output is enabled.

**Parameters**

- **gain** (*float*) – Gain, multiplier of waveform values
- **dc\_offset** (*float*) – DC offset in volts

**configure\_standard\_waveform**(*waveform\_function*, *amplitude*, *dc\_offset*, *frequency*, *duty\_cycle*)

Configures the instrument to output a standard waveform. Check instrument manual for maximum ratings which depend on load.

**Parameters**

- **waveform\_function** (*int* or *str*) – Waveform function ("SINE", "SQUARE", "TRIANGLE/RAMP", "DC")
- **amplitude** (*float*) – Amplitude in volts
- **dc\_offset** (*float*) – DC offset in volts
- **frequency** (*float*) – Frequency in Hz
- **duty\_cycle** (*int*) – Duty cycle in %

**property filter**

Enables or disables the filter on the instrument.

**Parameters** **enable\_filter** (*bool*) – Enable/Disable filter

**query\_arbitrary\_waveform()**

Returns the samples per second for arbitrary waveform generation.

**Returns** Samples per second

**Return type** int

**query\_arbitrary\_waveform\_gain\_and\_offset()**

Returns the settings for arbitrary waveform generation that includes gain and offset settings.

**Returns** Gain, DC offset

**Return type** (float, float)

**query\_generation\_status()**

Returns the status of waveform generation on the instrument.

**Returns** Status

**Return type** enum

**query\_standard\_waveform()**

Returns the settings for a standard waveform generation.

**Returns** Waveform function, amplitude, dc\_offset, frequency, duty\_cycle

**Return type** (enum, float, float, float, int)

**query\_waveform\_mode()**

Indicates whether the waveform output by the instrument is a standard or arbitrary waveform.

**Returns** Waveform mode

**Return type** enum

**reset\_instrument()**

Resets the session configuration to default values, and resets the device and driver software to a known state.

**run()**

Transitions the session from the Stopped state to the Running state.

**self\_calibrate()**

Performs offset nulling calibration on the device. You must run FGEN Initialize prior to running this method.

**shutdown()**

Removes the session and deallocates any resources acquired during the session. If output is enabled on any channels, they remain in their current state.

**stop()**

Transitions the acquisition from either the Triggered or Running state to the Stopped state.

**class MixedSignalOscilloscope**(*virtualbench, reset, vb\_name=""*)

Bases: `pymeasure.instruments.ni.virtualbench.VirtualBench.VirtualBenchInstrument`

Represents Mixed Signal Oscilloscope (MSO) Module of Virtual Bench device. Allows to measure oscilloscope data from analog and digital channels.

Methods from pyvirtualbench not implemented in pymeasure yet:

- `enable_digital_channels`
- `configure_digital_threshold`
- `configure_advanced_digital_timing`
- `configure_state_mode`
- `configure_digital_edge_trigger`
- `configure_digital_pattern_trigger`
- `configure_digital_glitch_trigger`
- `configure_digital_pulse_width_trigger`
- `query_digital_channel`
- `query_enabled_digital_channels`
- `query_digital_threshold`
- `query_advanced_digital_timing`
- `query_state_mode`
- `query_digital_edge_trigger`
- `query_digital_pattern_trigger`
- `query_digital_glitch_trigger`
- `query_digital_pulse_width_trigger`
- `read_digital_u64`

**auto\_setup()**

Automatically configure the instrument

**configure\_analog\_channel**(*channel, enable\_channel, vertical\_range, vertical\_offset, probe\_attenuation, vertical\_coupling*)

Configure analog measurement channel

**Parameters**

- **channel** (*str*) – Channel string

- **enable\_channel** (*bool*) – Enable/Disable channel
- **vertical\_range** (*float*) – Vertical measurement range (0V - 20V), the instrument discretizes to these ranges: [20, 10, 5, 2, 1, 0.5, 0.2, 0.1, 0.05] which are 5x the values shown in the native UI.
- **vertical\_offset** (*float*) – Vertical offset to correct for (inverted compared to VB native UI, -20V - +20V, resolution 0.1mV)
- **probe\_attenuation** (*int or str*) – Probe attenuation ('ATTENUATION\_10X' or 'ATTENUATION\_1X')
- **vertical\_coupling** (*int or str*) – Vertical coupling ('AC' or 'DC')

**configure\_analog\_channel\_characteristics**(*channel, input\_impedance, bandwidth\_limit*)  
Configure electrical characteristics of the specified channel

#### Parameters

- **channel** (*str*) – Channel string
- **input\_impedance** (*int or str*) – Input Impedance ('ONE\_MEGA\_OHM' or 'FIFTY\_OHMS')
- **bandwidth\_limit** (*int*) – Bandwidth limit (100MHz or 20MHz)

**configure\_analog\_edge\_trigger**(*trigger\_source, trigger\_slope, trigger\_level, trigger\_hysteresis, trigger\_instance*)

Configures a trigger to activate on the specified source when the analog edge reaches the specified levels.

#### Parameters

- **trigger\_source** (*str*) – Channel string
- **trigger\_slope** (*int or str*) – Trigger slope ('RISING', 'FALLING' or 'EITHER')
- **trigger\_level** (*float*) – Trigger level
- **trigger\_hysteresis** (*float*) – Trigger hysteresis
- **trigger\_instance** (*int or str*) – Trigger instance

**configure\_analog\_pulse\_width\_trigger**(*trigger\_source, trigger\_polarity, trigger\_level, comparison\_mode, lower\_limit, upper\_limit, trigger\_instance*)

Configures a trigger to activate on the specified source when the analog edge reaches the specified levels within a specified window of time.

#### Parameters

- **trigger\_source** (*str*) – Channel string
- **trigger\_polarity** (*int or str*) – Trigger slope ('POSITIVE' or 'NEGATIVE')
- **trigger\_level** (*float*) – Trigger level
- **comparison\_mode** (*int or str*) – Mode of comparison ('GREATER\_THAN\_UPPER\_LIMIT', 'LESS\_THAN\_LOWER\_LIMIT', 'INSIDE\_LIMITS' or 'OUTSIDE\_LIMITS')
- **lower\_limit** (*float*) – Lower limit
- **upper\_limit** (*float*) – Upper limit

- **trigger\_instance** (*int or str*) – Trigger instance

**configure\_immediate\_trigger()**

Configures a trigger to immediately activate on the specified channels after the pretrigger time has expired.

**configure\_timing**(*sample\_rate, acquisition\_time, pretrigger\_time, sampling\_mode*)

Configure timing settings of the MSO

**Parameters**

- **sample\_rate** (*int*) – Sample rate (15.26kS - 1GS)
- **acquisition\_time** (*float*) – Acquisition time (1ns - 68.711s)
- **pretrigger\_time** (*float*) – Pretrigger time (0s - 10s)
- **sampling\_mode** – Sampling mode ('SAMPLE' or 'PEAK\_DETECT')

**configure\_trigger\_delay**(*trigger\_delay*)

Configures the amount of time to wait after a trigger condition is met before triggering.

**param float trigger\_delay** Trigger delay (0s - 17.1799s)

**force\_trigger()**

Causes a software-timed trigger to occur after the pretrigger time has expired.

**query\_acquisition\_status()**

Returns the status of a completed or ongoing acquisition.

**query\_analog\_channel**(*channel*)

Indicates the vertical configuration of the specified channel.

**Returns** Channel enabled, vertical range, vertical offset, probe attenuation, vertical coupling

**Return type** (bool, float, float, enum, enum)

**query\_analog\_channel\_characteristics**(*channel*)

Indicates the properties that control the electrical characteristics of the specified channel. This method returns an error if too much power is applied to the channel.

**return** Input impedance, bandwidth limit

**rtype** (enum, float)

**query\_analog\_edge\_trigger**(*trigger\_instance*)

Indicates the analog edge trigger configuration of the specified instance.

**Returns** Trigger source, trigger slope, trigger level, trigger hysteresis

**Return type** (str, enum, float, float)

**query\_analog\_pulse\_width\_trigger**(*trigger\_instance*)

Indicates the analog pulse width trigger configuration of the specified instance.

**Returns** Trigger source, trigger polarity, trigger level, comparison mode, lower limit, upper limit

**Return type** (str, enum, float, enum, float, float)

**query\_enabled\_analog\_channels()**

Returns String of enabled analog channels.

**Returns** Enabled analog channels

**Return type** str

**query\_timing()**

Indicates the timing configuration of the MSO. Call directly before measurement to read the actual timing configuration and write it to the corresponding class variables. Necessary to interpret the measurement data, since it contains no time information.

**Returns** Sample rate, acquisition time, pretrigger time, sampling mode

**Return type** (float, float, float, enum)

**query\_trigger\_delay()**

Indicates the trigger delay setting of the MSO.

**Returns** Trigger delay

**Return type** float

**query\_trigger\_type(trigger\_instance)**

Indicates the trigger type of the specified instance.

**Parameters** **trigger\_instance** – Trigger instance ('A' or 'B')

**Returns** Trigger type

**Return type** str

**read\_analog\_digital\_dataframe()**

Transfers data from the instrument and returns a pandas dataframe of the analog measurement data, including time coordinates

**Returns** Dataframe with time and measurement data

**Return type** pd.DataFrame

**read\_analog\_digital\_u64()**

Transfers data from the instrument as long as the acquisition state is Acquisition Complete. If the state is either Running or Triggered, this method will wait until the state transitions to Acquisition Complete. If the state is Stopped, this method returns an error.

**Returns** Analog data out, analog data stride, analog t0, digital data out, digital timestamps out, digital t0, trigger timestamp, trigger reason

**Return type** (list, int, pyvb.Timestamp, list, list, pyvb.Timestamp, pyvb.Timestamp, enum)

**reset\_instrument()**

Resets the session configuration to default values, and resets the device and driver software to a known state.

**run(autoTrigger=True)**

Transitions the acquisition from the Stopped state to the Running state. If the current state is Triggered, the acquisition is first transitioned to the Stopped state before transitioning to the Running state. This method returns an error if too much power is applied to any enabled channel.

**Parameters** **autoTrigger** (*bool*) – Enable/Disable auto triggering

**shutdown()**

Removes the session and deallocates any resources acquired during the session. If output is enabled on any channels, they remain in their current state.

**stop()**

Transitions the acquisition from either the Triggered or Running state to the Stopped state.

**validate\_channel(channel)**

Check if `channel` is a correct specification

**Parameters** `channel` (*str*) – Channel string

**Returns** Channel string

**Return type** `str`

**static** `validate_trigger_instance(trigger_instance)`

Check if `trigger_instance` is a valid choice

**Parameters** `trigger_instance` (*int or str*) – Trigger instance ('A' or 'B')

**Returns** Trigger instance

**Return type** `int`

**class** `PowerSupply`(*virtualbench, reset, vb\_name=""*)

Bases: `pymeasure.instruments.ni.virtualbench.VirtualBench.VirtualBenchInstrument`

Represents Power Supply (PS) Module of Virtual Bench device

**configure\_current\_output**(*channel, current\_level, voltage\_limit*)

Configures a current output on the specified channel. This method should be called once for every channel you want to configure to output current.

**configure\_voltage\_output**(*channel, voltage\_level, current\_limit*)

Configures a voltage output on the specified channel. This method should be called once for every channel you want to configure to output voltage.

**property** `outputs_enabled`

Enables or disables all outputs on all channels of the instrument.

**Parameters** `enable_outputs` (*bool*) – Enable/Disable outputs

**query\_current\_output**(*channel*)

Indicates the current output settings on the specified channel.

**query\_voltage\_output**(*channel*)

Indicates the voltage output settings on the specified channel.

**read\_output**(*channel*)

Reads the voltage and current levels and outout mode of the specified channel.

**reset\_instrument**()

Resets the session configuration to default values, and resets the device and driver software to a known state.

**shutdown**()

Removes the session and deallocates any resources acquired during the session. If output is enabled on any channels, they remain in their current state.

**property** `tracking`

Enables or disables tracking between the positive and negative 25V channels. If enabled, any configuration change on the positive 25V channel is mirrored to the negative 25V channel, and any writes to the negative 25V channel are ignored.

**Parameters** `enable_tracking` (*bool*) – Enable/Disable tracking

**validate\_channel**(*channel, current=False, voltage=False*)

Check if channel string is valid and if output current/voltage are within the output ranges of the channel

**Parameters**

- **channel** (*str*) – Channel string ("ps/+6V", "ps/+25V", "ps/-25V")



- **current** (*bool*, *optional*) – Current output, defaults to False
- **voltage** (*bool*, *optional*) – Voltage output, defaults to False

**Returns** channel or channel, current & voltage

**Return type** str or (str, float, float)

**acquire\_digital\_input\_output**(*lines*, *reset=False*)

Establishes communication with the DIO module. This method should be called once per session.

**Parameters**

- **lines** (*str*) – Lines to acquire, reading is possible on all lines
- **reset** (*bool*, *optional*) – Reset DIO module, defaults to False

**acquire\_digital\_multimeter**(*reset=False*)

Establishes communication with the DMM module. This method should be called once per session.

**Parameters** **reset** (*bool*, *optional*) – Reset the DMM module, defaults to False

**acquire\_function\_generator**(*reset=False*)

Establishes communication with the FGEN module. This method should be called once per session.

**Parameters** **reset** (*bool*, *optional*) – Reset the FGEN module, defaults to False

**acquire\_mixed\_signal\_oscilloscope**(*reset=False*)

Establishes communication with the MSO module. This method should be called once per session.

**Parameters** **reset** (*bool*, *optional*) – Reset the MSO module, defaults to False

**acquire\_power\_supply**(*reset=False*)

Establishes communication with the PS module. This method should be called once per session.

**Parameters** **reset** (*bool*, *optional*) – Reset the PS module, defaults to False

**collapse\_channel\_string**(*names\_in*)

Collapses a channel string into a comma and colon-delimited equivalent. Last element is the number of channels.

**Parameters** **names\_in** (*str*) – Channel string

**Returns** Channel string with colon notation where possible, number of channels

**Return type** (str, int)

**convert\_timestamp\_to\_values**(*timestamp*)

Converts a timestamp to seconds and fractional seconds

**Parameters** **timestamp** (*pyvb.Timestamp*) – VirtualBench timestamp

**Returns** (seconds\_since\_1970, fractional seconds)

**Return type** (int, float)

**convert\_values\_to\_datetime**(*timestamp*)

Converts timestamp to datetime object

**Parameters** **timestamp** (*pyvb.Timestamp*) – VirtualBench timestamp

**Returns** Timestamp as DateTime object

**Return type** DateTime

**convert\_values\_to\_timestamp**(*seconds\_since\_1970, fractional\_seconds*)

Converts seconds and fractional seconds to a timestamp

**Parameters**

- **seconds\_since\_1970** (*int*) – Date/Time in seconds since 1970
- **fractional\_seconds** (*float*) – Fractional seconds

**Returns** VirtualBench timestamp

**Return type** pyvb.Timestamp

**expand\_channel\_string**(*names\_in*)

Expands a channel string into a comma-delimited (no colon) equivalent. Last element is the number of channels. 'dig/0:2' -> ('dig/0, dig/1, dig/2',3)

**Parameters** **names\_in** (*str*) – Channel string

**Returns** Channel string with all channels separated by comma, number of channels

**Return type** (*str*, *int*)

**get\_calibration\_information**()

Returns calibration information for the specified device, including the last calibration date and calibration interval.

**Returns** Calibration date, recommended calibration interval in months, calibration interval in months

**Return type** (pyvb.Timestamp, *int*, *int*)

**get\_library\_version**()

Return the version of the VirtualBench runtime library

**shutdown**()

Finalize the VirtualBench library.

**class** pymeasure.instruments.ni.virtualbench.**VirtualBench\_Direct**(\*args: Any, \*\*kwargs: Any)  
Bases: pyvirtualbench.

Represents National Instruments Virtual Bench main frame. This class provides direct access to the arm-strap/pyvirtualbench Python wrapper.

## 7.25 Oxford Instruments

This section contains specific documentation on the Oxford Instruments instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.25.1 Oxford Instruments Intelligent Temperature Controller 503

**class** pymeasure.instruments.oxfordinstruments.**ITC503**(resourceName, clear\_buffer=True,  
max\_temperature=301, min\_temperature=0,  
\*\*kwargs)

Bases: [pymeasure.instruments.instrument.Instrument](#)

Represents the Oxford Intelligent Temperature Controller 503.

```
itc = ITC503("GPIB::24")           # Default channel for the ITC503

itc.control_mode = "RU"             # Set the control mode to remote
itc.heater_gas_mode = "AUTO"        # Turn on auto heater and flow
```

(continues on next page)

(continued from previous page)

```

itc.auto_pid = True          # Turn on auto-pid

print(itc.temperature_setpoint) # Print the current set-point
itc.temperature_setpoint = 300 # Change the set-point to 300 K
itc.wait_for_temperature()    # Wait for the temperature to stabilize
print(itc.temperature_1)      # Print the temperature at sensor 1

```

**property auto\_pid**

A boolean property that sets the Auto-PID mode on (True) or off (False).

**property control\_mode**

A string property that sets the ITC in LOCAL or REMOTE and LOCKES, or UNLOCKES, the LOC/REM button. Allowed values are: LL: LOCAL & LOCKED RL: REMOTE & LOCKED LU: LOCAL & UNLOCKED RU: REMOTE & UNLOCKED.

**property gasflow**

A floating point property that controls gas flow when in manual mode. The value is expressed as a percentage of the maximum gas flow. Valid values are in range 0 [off] to 99.9 [%].

**property heater**

A floating point property that sets the required heater output when in manual mode. The parameter is expressed as a percentage of the maximum voltage. Valid values are in range 0 [off] to 99.9 [%].

**property heater\_gas\_mode**

A string property that sets the heater and gas flow control to AUTO or MANUAL. Allowed values are: MANUAL: HEATER MANUAL, GAS MANUAL AM: HEATER AUTO, GAS MANUAL MA: HEATER MANUAL, GAS AUTO AUTO: HEATER AUTO, GAS AUTO.

**program\_sweep**(temperatures, sweep\_time, hold\_time, steps=None)

Program a temperature sweep in the controller. Stops any running sweep. After programming the sweep, it can be started using OxfordITC503.sweep\_status = 1.

**Parameters**

- **temperatures** – An array containing the temperatures for the sweep
- **sweep\_time** – The time (or an array of times) to sweep to a set-point in minutes (between 0 and 1339.9).
- **hold\_time** – The time (or an array of times) to hold at a set-point in minutes (between 0 and 1339.9).
- **steps** – The number of steps in the sweep, if given, the temperatures, sweep\_time and hold\_time will be interpolated into (approximately) equal segments

**property sweep\_status**

An integer property that sets the sweep status. Values are: 0: Sweep not running 1: Start sweep / sweeping to first set-point 2P - 1: Sweeping to set-point P 2P: Holding at set-point P.

**property sweep\_table**

A property that sets values in the sweep table. Relies on the xpointer and ypointer to point at the location in the table that is to be set.

**property temperature\_1**

Reads the temperature of the sensor 1 in Kelvin.

**property temperature\_2**

Reads the temperature of the sensor 2 in Kelvin.

**property temperature\_3**

Reads the temperature of the sensor 3 in Kelvin.

**property temperature\_error**

Reads the difference between the set-point and the measured temperature in Kelvin. Positive when set-point is larger than measured.

**property temperature\_setpoint**

A floating point property that controls the temperature set-point of the ITC in kelvin.

**wait\_for\_temperature**(*error=0.01, timeout=3600, check\_interval=0.5, stability\_interval=10, thermalize\_interval=300, should\_stop=<function ITC503.<lambda>>, max\_comm\_errors=None*)

Wait for the ITC to reach the set-point temperature.

**Parameters**

- **error** – The maximum error in Kelvin under which the temperature is considered at set-point
- **timeout** – The maximum time the waiting is allowed to take. If timeout is exceeded, a `TimeoutError` is raised. If timeout is set to zero, no timeout will be used.
- **check\_interval** – The time between temperature queries to the ITC.
- **stability\_interval** – The time over which the temperature\_error is to be below error to be considered stable.
- **thermalize\_interval** – The time to wait after stabilizing for the system to thermalize.
- **should\_stop** – Optional function (returning a bool) to allow the waiting to be stopped before its end.
- **max\_comm\_errors** – The maximum number of communication errors that are allowed before the wait is stopped. if set to `None` (default), no maximum will be used.

**property xpointer**

An integer property to set pointers into tables for loading and examining values in the table. For programming the sweep table values from 1 to 16 are allowed, corresponding to the maximum number of steps.

**property ypointer**

An integer property to set pointers into tables for loading and examining values in the table. For programming the sweep table the allowed values are: 1: Setpoint temperature, 2: Sweep-time to set-point, 3: Hold-time at set-point.

## 7.26 Parker

This section contains specific documentation on the Parker instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.26.1 Parker GV6 Servo Motor Controller

**class** `pymeasure.instruments.parker.ParkerGV6(port)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Parker Gemini GV6 Servo Motor Controller and provides a high-level interface for interacting with the instrument

**property** `angle`

Returns the angle in degrees based on the position and whether relative or absolute positioning is enabled, returning None on error

**property** `angle_error`

Returns the angle error in degrees based on the position error, or returns None on error

**disable()**

Disables the motor from moving

**echo(*enable=False*)**

Enables (True) or disables (False) the echoing of all commands that are sent to the instrument

**enable()**

Enables the motor to move

**is\_moving()**

Returns True if the motor is currently moving

**kill()**

Stops the motor

**move()**

Initiates the motor to move to the setpoint

**property** `position`

Returns an integer number of counts that correspond to the angular position where 1 revolution equals 4000 counts

**property** `position_error`

Returns the error in the number of counts that corresponds to the error in the angular position where 1 revolution equals 4000 counts

**read()**

Overwrites the Instrument.read command to provide the correct functionality

**reset()**

Resets the motor controller while blocking and (CAUTION) resets the absolute position value of the motor

**set\_defaults()**

Sets up the default values for the motor, which is run upon construction

**set\_hardware\_limits(*positive=True, negative=True*)**

Enables (True) or disables (False) the hardware limits for the motor

**set\_software\_limits(*positive, negative*)**

Sets the software limits for motion based on the count unit where 4000 counts is 1 revolution

**property** `status`

Returns a list of the motor status in readable format

**stop()**

Stops the motor during movement

**use\_absolute\_position()**

Sets the motor to accept setpoints from an absolute zero position

**use\_relative\_position()**

Sets the motor to accept setpoints that are relative to the last position

**write(*command*)**

Overwrites the Instrumnet.write command to provide the correct line break syntax

## 7.27 Pendulum

This section contains specific documentation on the Pendulum instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.27.1 Pendulum CNT91 frequency counter

**class** pymeasure.instruments.pendulum.cnt91.CNT91(*resourceName*, **\*\*kwargs**)

Bases: [pymeasure.instruments.instrument.Instrument](#)

Represents a Pendulum CNT-91 frequency counter.

**property batch\_size**

Maximum number of buffer entries that can be transmitted at once.

**buffer\_frequency\_time\_series(*channel*, *n\_samples*, *sample\_rate*, *trigger\_source=None*)**

Record a time series to the buffer and read it out after completion.

#### Parameters

- **channel** – Channel that should be used
- **n\_samples** – The number of samples
- **sample\_rate** – Sample rate in Hz
- **trigger\_source** – Optionally specify a trigger source to start the measurement

**configure\_frequency\_array\_measurement(*n\_samples*, *channel*)**

Configure the counter for an array of measurements.

#### Parameters

- **n\_samples** – The number of samples
- **channel** – Measurment channel (A, B, C, E, INTREF)

**property continuous**

Controls whether to perform continuous measurements.

**property external\_arming\_start\_slope**

Set slope for the start arming condition.

**property external\_start\_arming\_source**

Select arming input or switch off the start arming function. Options are 'A', 'B' and 'E' (rear). 'IMM' turns trigger off.

**property format**

Reponse format (ASCII or REAL).

**property interpolator\_autocalibrated**

Controls if interpolators should be calibrated automatically.

**property measurement\_time**

Gate time for one measurement in s.

**read\_buffer**(*expected\_length=0*)

Read out the entire buffer.

**Parameters** **expected\_length** – The expected length of the buffer. If more data is read, values at the end are removed. Defaults to 0, which means that the entire buffer is returned independent of its length.

**Returns** Frequency values from the buffer.

## 7.28 Razorbill

This section contains specific documentation on the Razorbill instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.28.1 Razorbill RP100 custrom power supply for Razorbill Instrums stress & strain cells

**class** `pymeasure.instruments.razorbill.razorbillRP100`(*adapter*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents Razorbill RP100 strain cell controller

```
scontrol = razorbillRP100("ASRL/dev/ttyACM0::INSTR")

scontrol.output_1 = True      # turns output on
scontrol.slew_rate_1 = 1     # sets slew rate to 1V/s
scontrol.voltage_1 = 10      # sets voltage on output 1 to 10V
```

**property contact\_current\_1**

Returns the current in amps present at the front panel output of channel 1

**property contact\_current\_2**

Returns the current in amps present at the front panel output of channel 2

**property contact\_voltage\_1**

Returns the Voltage in volts present at the front panel output of channel 1

**property contact\_voltage\_2**

Returns the Voltage in volts present at the front panel output of channel 2

**property instant\_voltage\_1**

Returns the instantaneous output of source one in volts

**property instant\_voltage\_2**

Returns the instanteneous output of source two in volts

**property output\_1**

Turns output of channel 1 on or off

**property output\_2**

Turns output of channel 2 on or off

**property slew\_rate\_1**

Sets or queries the source slew rate in volts/sec of channel 1

**property** `slew_rate_2`

Sets or queries the source slew rate in volts/sec of channel 2

**property** `voltage_1`

Sets or queries the output voltage of channel 1

**property** `voltage_2`

Sets or queries the output voltage of channel 2

## 7.29 Rohde & Schwarz

This section contains specific documentation on the Rohde & Schwarz instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.29.1 R&S SFM TV test transmitter

**class** `pymeasure.instruments.rohdeschwarz.sfm.SFM(resource_name, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Rohde&Schwarz SFM TV test transmitter interface for interacting with the instrument.

---

**Note:** The current implementation only works with the first system in this unit.

Further source extension for system 2-6 would be required.

The intermodulation subsystem is also not yet implemented.

---

**property** `R75_out`

A bool property that controls the use of the 75R output (if installed)

Value	Meaning
False	50R output active (N)
True	75R output active (BNC)

refer also to chapter 3.6.5 of the manual

**property** `TV_country`

A string property that controls the country specifics of the video/sound system to be used

Possible values are:



Value	Meaning
BG_G	BG General
DK_G	DK General
I_G	I General
L_G	L General
GERM	Germany
BELG	Belgium
NETH	Netherlands
FIN	Finland
AUST	Australia
BG_T	BG Th
DENM	Denmark
NORW	Norway
SWED	Sweden
GUS	Russia
POL1	Poland
POL2	Poland
HUNG	Hungary
CHEC	Czech Republic
CHINA1	China
CHINA2	China
GRE	Great Britain
SAFR	South Africa
FRAN	France
USA	United States
KOR	Korea
JAP	Japan
CAN	Canada
SAM	South America

Please confirm with the manual about the details for these settings.

#### **property TV\_standard**

A string property that controls the type of video standard

Possible values are:

Value	Lines	System
BG	625	PAL
DK	625	SECAM
I	625	PAL
K1	625	SECAM
L	625	SECAM
M	525	NTSC
N	625	NTSC

Please confirm with the manual about the details for these settings.

#### **property basic\_info**

A String property containing information about the hardware modules installed in the unit

#### **property beeper\_enabled**

A bool property that controls the beeper status,

refer also to chapter 3.6.8 of the manual

**calibration**(*number=1, subsystem=None*)

Function to either calibrate the whole modulator, when subsystem parameter is omitted, or calibrate a subsystem of the modulator.

Valid subsystem selections: “NICam, VISION, SOUND1, SOUND2, CODer”

**channel\_down\_relative()**

Decreases the output frequency to the next low channel/special channel based on the current country settings

**property channel\_sweep\_start**

A float property controlling the start frequency for channel sweep in Hz

- Minimum 5 MHz
- Maximum 1 GHz

**property channel\_sweep\_step**

A float property controlling the start frequency for channel sweep in Hz

- Minimum 5 MHz
- Maximum 1 GHz

**property channel\_sweep\_stop**

A float property controlling the start frequency for channel sweep in Hz

- Minimum 5 MHz
- Maximum 1 GHz

**property channel\_table**

A string property controlling which channel table is used

Possible selections are:

Value	Meaning
DEF	Default channel table
USR1	User table No. 1
USR2	User table No. 2
USR3	User table No. 3
USR4	User table No. 4
USR5	User table No. 5

refer also to chapter 3.6.6.1 of the manual

**channel\_up\_relative()**

Increases the output frequency to the next higher channel/special channel based on the current country settings

**coder\_adjust()**

Starts the automatic setting of the differential deviation

refer also to chapter 3.6.6.4 of the manual

**property coder\_id\_frequency**

A int property that controls the frequency of the identification of the coder

valid range 0 .. 200 Hz

**property coder\_modulation\_degree**

A float property that controls the modulation degree of the identification of the coder

valid range: 0 .. 0.9

**property coder\_pilot\_deviation**

A int property that controls deviation of the pilot frequency of the coder

valid range: 1 .. 4 kHz

**property coder\_pilot\_frequency**

A int property that controls the pilot frequency of the coder

valid range: 40 .. 60 kHz

**property cw\_frequency**

A float property controlling the CW-frequency in Hz

- Minimum 5 MHz
- Maximum 1 GHz

**property date**

A list property for the date of the RTC in the unit

**property event\_reg**

Content of the event register of the Status Operation Register refer also to chapter 3.6.7 of the manual

**property ext\_ref\_base\_unit**

A bool property for the external reference for the basic unit

Value	Meaning
False	Internal 10 MHz is used
True	External 10 MHz is used

**property ext\_ref\_extension**

A bool property for the external reference for the extension frame

Value	Meaning
False	Internal 10 MHz is used
True	External 10 MHz is used

**property ext\_vid\_connector**

A string property controlling which connector is used as the input of the video source

Possible selections are:

Value	Meaning
HIGH	Front connector - Hi-Z
LOW	Front connector - 75R
REAR1	Rear connector 1
REAR2	Rear connector 2
AUTO	Automatic assignment

**property external\_modulation\_frequency**

A int property that controls the setting for the external modulator frequency

valid range: 32 .. 46 MHz

**property external\_modulation\_power**

A int property that controls the setting for the external modulator output power

valid range: -7..0 dBm

refer also to chapter 3.6.6.5 of the manual

**property external\_modulation\_source**

A bool property for the modulation source selection

refer also to chapter 3.6.6.8 of the manual

**property frequency**

A float property controlling the frequency in Hz

- Minimum 5 MHz
- Maximum 1 GHz

**property frequency\_mode**

A string property controlling which the unit is used in

Possible selections are:

Value	Meaning
CW	Continuous wave mode
FIXED	fixed frequency mode
CHSW	Channel sweep
RFSW	Frequency sweep

---

**Note:** selecting the sweep mode, will start the sweep immediately!

---

**property gpib\_address**

A int property that controls the GPIB address of the unit

valid range: 0..30

**property high\_frequency\_resolution**

A property that controls the frequency resolution,

Possible selections are:

Value	Meaning
False	Low resolution (1000Hz)
True	High resolution (1Hz)

**property level**

A float property controlling the output level in dBm,

- Minimum -99dBm
- Maximum 10dBm (depending on output mode)

refer also to chapter 3.6.6.2 of the manual

**property level\_mode**

A string property controlling the output attenuator and linearity mode

Possible selections are:

Value	Meaning	max. output level
NORM	Normal mode	+6 dBm
LOWN	low noise mode	+10 dBm
CONT	continous mode	+10 dBm
LOWD	low distortion mode	+0 dBm

Continous mode allows up to 14 dB of level setting without use of the mechanical attenuator.

**property lower\_sideband\_enabled**

A bool property that controls the use of the lower sideband

refer also to chapter 3.6.6.10 of the manual

**property modulation\_enabled**

A bool property that controls the modulation status

**property nicam\_IQ\_inverted**

A bool property that controls if the NICAM IQ signals are inverted or not

Value	Meaning
False	normal (IQ)
True	inverted (QI)

**property nicam\_additional\_bits**

A int property that controls the additional data in the NICAM modulator

valid range: 0 .. 2047

**property nicam\_audio\_frequency**

A int property that controls the frequency of the internal sound generator

valid range: 0 Hz .. 15 kHz

**property nicam\_audio\_volume**

A float property that controls the audio volume in the NICAM modulator in dB

valid range: 0..60 dB

**property nicam\_bit\_error\_enabled**

A bool property that controls the status of an artifical bit error rate to be applied

**property nicam\_bit\_error\_rate**

A float property that controls the artifical bit error rate.

valid range: 1.2E-7 .. 2E-3

**property nicam\_carrier\_enabled**

A bool property that controls if the NICAM carrier is switched on or off

**property nicam\_carrier\_frequency**

A float property that controls the frequency of the NICAM carrier

valid range: 33.05 MHz +/- 0.2 Mhz

**property nicam\_carrier\_level**

A float property that controls the value of the NICAM carrier

valid range: -40 .. -13 dB

**property nicam\_control\_bits**

A int property that controls the additional data in the NICAM modulator

valid range: 0 .. 3

**property nicam\_data**

A int property that controls the data in the NICAM modulator

valid range: 0 .. 2047

**property nicam\_intercarrier\_frequency**

A float property that controls the inter-carrier frequency of the NICAM carrier

valid range: 5 .. 9 MHz

**property nicam\_mode**

A string property that controls the signal type to be sent via NICAM

Possible values are:

Value	Meaning
MON	Mono sound + NICAM data
STER	Stereo sound
DUAL	Dual channel sound
DATA	NICAM data only

refer also to chapter 3.6.6.6 of the manual

**property nicam\_preemphasis\_enabled**

A bool property that controls the status of the J17 preemphasis

**property nicam\_source**

A string property that controls the signal source for NICAM

Possible values are:

Value	Meaning
INT	Internal audio generator(s)
EXT	External audio source
CW	Continuous wave signal
RAND	Random data stream
TEST	Test signal

**property nicam\_test\_signal**

A int property that controls the selection of the test signal applied

Value	Meaning
1	Test signal 1 (91 kHz square wave, I&Q 90deg apart)
2	Test signal 2 (45.5 kHz square wave, I&Q 90deg apart)
3	Test signal 3 (182 kHz sine wave, I&Q in phase)

**property normal\_channel**

A int property controlling the current selected regular/normal channel number valid selections are based on the country settings.

**property operation\_enable\_reg**

Content of the enable register of the Status Operation Register

Valid range: 0...32767

**property output\_voltage**

A float property controlling the output level in Volt,

Minimum 2.50891e-6, Maximum 0.707068 (depending on output mode) refer also to chapter 3.6.6.12 of the manual

**property questionable\_event\_reg**

Content of the event register of the Status Questionable Operation Register

**property questionable\_operation\_enable\_reg**

Content of the enable register of the Status Questionable Operation Register

Valid range 0...32767

**property questionable\_status\_reg**

Content of the condition register of the Status Questionable Operation Register

**property remote\_interfaces**

A string property controlling the selection of interfaces for remote control

Possible selections are:

Value	Meaning
OFF	no remote control
GPIB	GPIB only enabled
SER	RS232 only enabled
BOTH	GPIB & RS232 enabled

**property rf\_out\_enabled**

A bool property that controls the status of the RF-output

**property rf\_sweep\_center**

A float property controlling the center frequency for sweep in Hz

- Minimum 5 MHz
- Maximum 1 GHz

**property rf\_sweep\_span**

A float property controlling the sweep span in Hz,

- Minimum 1 kHz
- Maximum 1 GHz

**property rf\_sweep\_start**

A float property controlling the start frequency for sweep in Hz

- Minimum 5 MHz
- Maximum 1 GHz

**property rf\_sweep\_step**

A float property controlling the stepwidth for sweep in Hz,

- Minimum 1 kHz
- Maximum 1 GHz

**property rf\_sweep\_stop**

A float property controlling the stop frequency for sweep in Hz

- Minimum 5 MHz
- Maximum 1 GHz

**property scale\_volt**

A string property that controls the unit to be used for voltage entries on the unit

Possible values are: AV,FV, PV, NV, UV, MV, V, KV, MAV, GV, TV, PEV, EV, DBAV, DBFV, DBPV, DBNV, DBUV, DBMV, DBV, DBKV, DBMAV, DBGV, DBTV, DBPEv, DBEV

refer also to chapter 3.6.9 of the manual

**property serial\_baud**

A int property that controls the serial communication speed ,

Possible values are: 110,300,600,1200,4800,9600,19200

**property serial\_bits**

A int property that controls the number of bits used in serial communication

Possible values are: 7 or 8

**property serial\_flowcontrol**

A string property that controls the serial handshake type used in serial communication

Possible values are:

Value	Meaning
NONE	no flow-control/handshake
XON	XON/XOFF flow-control
ACK	hardware handshake with RTS&CTS

**property serial\_parity**

A string property that controls the parity type used for serial communication

Possible values are:

Value	Meaning
NONE	no parity
EVEN	even parity
ODD	odd parity
ONE	parity bit fixed to 1
ZERO	parity bit fixed to 0

**property serial\_stopbits**

A int property that controls the number of stop-bits used in serial communication,

Possible values are: 1 or 2

**property sound\_mode**

A string property that controls the type of audio signal

Possible values are:



Value	Meaning
MONO	MOnoaural sound
PIL	pilot-carrier + mono
BTSC	BTSC + mono
STER	Stereo sound
DUAL	Dual channel sound
NIC	NICAM + Mono

**property special\_channel**

A int property controlling the current selected special channel number valid selections are based on the country settings.

**property status\_info\_shown**

A bool property that controls if the display shows information during remote control

**status\_preset()**

partly resets the SCPI status reporting structures

**property status\_reg**

Content of the condition register of the Status Operation Register

**property subsystem\_info**

A String property containing information about the system configuration

**property system\_number**

A int property for the selected systems (if more than 1 available)

- Minimum 1
- Maximum 6

**property time**

A list property for the time of the RTC in the unit

**property vision\_average\_enabled**

A bool property that controls the average mode for the vision system

**property vision\_balance**

A float property that controls the balance of the vision modulator

valid range: -0.5 .. 0.5

**property vision\_carrier\_enabled**

A bool property that controls the vision carrier status

refer also to chapter 3.6.6.9 of the manual

**property vision\_carrier\_frequency**

A float property that controls the frequency of the vision carrier

valid range: 32 .. 46 MHz

**property vision\_clamping\_average**

A float property that controls the operation point of the vision modulator

valid range: -0.5 .. 0.5

**property vision\_clamping\_enabled**

A bool property that controls the clamping behavior of the vision modulator

**property vision\_clamping\_mode**

A string property that controls the clamping mode of the vision modulator

Possible selections are HARD or SOFT

**property vision\_precorrection\_enabled**

A bool property that controls the precorrection behavior of the vision modulator

**property vision\_residual\_carrier\_level**

A float property that controls the value of the residual carrier

valid range: 0 .. 0.3 (30%)

**property vision\_sideband\_filter\_enabled**

A bool property that controls the use of the VSBF (vestigial sideband filter) in the vision modulator

**property vision\_videosignal\_enabled**

A bool property that controls if the video signal is switched on or off

**class** `pymessage.instruments.rohdeschwarz.sfm.Sound_Channel`(*instrument, number*)

Bases: object

Class object for the two sound channels

refere also to chapter 3.6.6.7 of the user manual

**property carrier\_enabled**

A bool property that controls if the audio carrier is switched on or off

**property carrier\_frequency**

A float property that controls the frequency of the sound carrier

valid range: 32 .. 46 MHz

**property carrier\_level**

A float property that controls the level of the audio carrier in dB relative to the vision carrier (0dB)

valid range: -34 .. -6 dB

**property deviation**

A int property that controls deviation of the selected audio signal

valid range: 0 .. 110 kHz

**property frequency**

A int property that controls the frequency of the internal sound generator

valid range: 300 Hz .. 15 kHz

**property modulation\_degree**

A float property that controls the modulation depth for the audio signal (Note: only for the use of AM in Standard L)

valid range: 0 .. 1 (100%)

**property modulation\_enabled**

A bool property that controls the audio modulation status

Value	Meaning
False	modulation disabled
True	modulation enabled

**property preemphasis\_enabled**

A bool property that controls if the preemphasis for the audio is switched on or off

**property preemphasis\_time**

A int property that controls if the mode of the preemphasis for the audio signal

Value	Meaning
50	50 us preemphasis
75	75 us preemphasis

**property use\_external\_source**

A bool property for the audio source selection

Value	Meaning
False	Internal audio generator(s)
True	External signal source

**values(command, \*\*kwargs)**

Reads a set of values from the instrument through the adapter, passing on any keyword arguments.

## 7.29.2 R&S FSL spectrum analyzer

### Connecting to the instrument via network

Once connected to the network, the instrument's IP address can be found by clicking the "Setup" button and navigating to "General Settings" -> "Network Address".

It can then be connected like this:

```
from pymeasure.instruments.rohdeschwarz import FSL
fsl = FSL("TCPIP::192.168.1.123::INSTR")
```

### Getting and setting parameters

Most parameters are implemented as properties, which means they can be read and written (getting and setting) in a consistent and simple way. If numerical values are provided, base units are used (s, Hz, dB, ...). Alternatively, the values can also be provided with a unit, e.g. "1.5 GHz" or "1.5GHz". Return values are always numerical.

```
# Getting the current center frequency
fsl.freq_center

9000000000.0
```

```
# Changing it to 10 MHz by providing the numerical value
fsl.freq_center = 10e6
```

```
# Verifying:
fsl.freq_center

10000000.0
```

```
# Changing it to 9 GHz by providing a string and verifying the result
fsl.freq_center = '9GHz'
fsl.freq_center

9000000000.0
```

```
# Setting the span to maximum
fsl.freq_span = '7 GHz'
```

## Reading a trace

We will read the current trace

```
x, y = fsl.read_trace()
```

## Markers

Markers are implemented as their own class. You can create them like this:

```
m1 = fsl.create_marker()
```

Set peak excursion:

```
m1.peak_excursion = 3
```

Set marker to a specific position:

```
m1.x = 10e9
```

Find the next peak to the left and get the level:

```
m1.to_next_peak('left')
m1.y

-34.9349060059
```

## Delta markers

Delta markers can be created by setting the appropriate keyword.

```
d2 = fsl.create_marker(is_delta_marker=True)
d2.name

'DELT2'
```

## Example program

Here is an example of a simple script for recording the peak of a signal.

```
m1 = fsl.create_marker() # create marker 1

# Set standard settings, set to full span
fsl.continuous_sweep = False
fsl.freq_span = '18 GHz'
fsl.res_bandwidth = "AUTO"
fsl.video_bandwidth = "AUTO"
fsl.sweep_time = "AUTO"

# Perform a sweep on full span, set the marker to the peak and some to that marker
fsl.single_sweep()
m1.to_peak()
m1.zoom('20 MHz')

# take data from the zoomed-in region
fsl.single_sweep()
x, y = fsl.read_trace()
```

**class** `pymeasure.instruments.rohdeschwarz.fsl.FSL(resourceName, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents a Rohde&Schwarz FSL spectrum analyzer.

All physical values that can be set can either be as a string of a value and a unit (e.g. “1.2 GHz”) or as a float value in the base units (Hz, dBm, etc.).

**property** `attenuation`

Attenuation in dB.

**continue\_single\_sweep()**

Continue with single sweep with synchronization.

**property** `continuous_sweep`

Continuous (True) or single sweep (False)

**create\_marker(num=1, is\_delta\_marker=False)**

Create a marker.

**Parameters**

- **num** – The marker number (1-4)
- **is\_delta\_marker** – True if the marker is a delta marker, default is False.

**Returns** The marker object.

**property** `freq_center`

Center frequency in Hz.

**property** `freq_span`

Frequency span in Hz.

**property** `freq_start`

Start frequency in Hz.

**property** `freq_stop`

Stop frequency in Hz.

**read\_trace**(*n\_trace=1*)

Read trace data.

**Parameters** *n\_trace* – The trace number (1-6). Default is 1.

**Returns** 2d numpy array of the trace data, [[frequency], [amplitude]].

**property res\_bandwidth**

Resolution bandwidth in Hz. Can be set to 'AUTO'

**single\_sweep**()

Perform a single sweep with synchronization.

**property sweep\_time**

Sweep time in s. Can be set to 'AUTO'.

**property trace\_mode**

Trace mode ('WRIT', 'MAXH', 'MINH', 'AVER' or 'VIEW')

**property video\_bandwidth**

Video bandwidth in Hz. Can be set to 'AUTO'

## 7.30 Signal Recovery

This section contains specific documentation on the Signal Recovery instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.30.1 DSP 7265 Lock-in Amplifier

**class** pymeasure.instruments.signalrecovery.DSP7265(*resourceName, \*\*kwargs*)

Bases: *pymeasure.instruments.instrument.Instrument*

This is the class for the DSP 7265 lockin amplifier

**property adc1**

Reads the input value of ADC1 in Volts

**property adc2**

Reads the input value of ADC2 in Volts

**buffer\_to\_float**(*buffer\_data, sensitivity=None, sensitivity2=None, raise\_error=True*)

Method that converts the fixed-point buffer data to floating point data. The provided data is converted as much as possible, but there are some requirements to the data if all provided columns are to be converted; if a key in the provided data cannot be converted it will be omitted in the returned data or an exception will be raised, depending on the value of *raise\_error*.

The requirements for converting the data are as follows:

- Converting X, Y, magnitude and noise requires sensitivity data, which can either be part of the provided data or can be provided via the *sensitivity* argument
- The same holds for X2, Y2, magnitude2 and noise2 with *sensitivity2*.
- Converting the frequency requires both 'frequency part 1' and 'frequency part 2'.

**Parameters**

- **buffer\_data** (*dict*) – The data to be converted. Must be in the format as returned by the *get\_buffer* method: a dict of numpy arrays.

- **sensitivity** – If provided, the sensitivity used to convert X, Y, magnitude and noise. Can be provided as a float or as an array that matches the length elements in the `buffer_data`. If both a sensitivity is provided and present in the `buffer_data`, the provided value is used for the conversion, but the sensitivity in the `buffer_data` is stored in the returned dict.
- **sensitivity2** – Same as the first sensitivity argument, but for X2, Y2, magnitude2 and noise2.
- **raise\_error** (*bool*) – Boolean that determines whether an exception is raised in case not all keys provided in `buffer_data` can be converted. If False, the columns that cannot be converted are omitted in the returned dict.

#### **property curve\_buffer\_bits**

An integer property that controls which data outputs are stored in the curve buffer. Valid values are values between 1 and 65,535 (or 2,097,151 in dual reference mode).

#### **property curve\_buffer\_interval**

An integer property that controls Sets the time interval between successive points being acquired in the curve buffer. The time interval is specified in ms with a resolution of 5 ms; input values are rounded up to a multiple of 5. Valid values are values between 0 and 1,000,000,000 (corresponding to 12 days). The interval may be set to 0, which sets the rate of data storage to the curve buffer to 1.25 ms/point (800 Hz). However this only allows storage of the X and Y channel outputs. There is no need to issue a CBD 3 command to set this up since it happens automatically when acquisition starts.

#### **property curve\_buffer\_length**

An integer property that controls the length of the curve buffer. Valid values are values between 1 and 32,768, but the actual maximum amount of points is determined by the amount of curves that are stored, as set via the `curve_buffer_bits` property ( $32,768 / n$ )

#### **property curve\_buffer\_status**

A property that represents the status of the curve buffer acquisition with four values: the first value represents the status with 5 possibilities (0: no activity, 1: acquisition via TD command running, 2: acquisition by a TDC command running, 5: acquisition via TD command halted, 6: acquisition by a TDC command halted); the second value is the number of sweeps that is acquired; the third value is the decimal representation of the status byte (the same response as the ST command); the fourth value is the number of points acquired in the curve buffer.

#### **property dac1**

A floating point property that represents the output value on DAC1 in Volts. This property can be set.

#### **property dac2**

A floating point property that represents the output value on DAC2 in Volts. This property can be set.

#### **property dac3**

A floating point property that represents the output value on DAC3 in Volts. This property can be set.

#### **property dac4**

A floating point property that represents the output value on DAC4 in Volts. This property can be set.

#### **property frequency**

A floating point property that represents the lock-in frequency in Hz. This property can be set.

#### **get\_buffer** (*quantity=None, convert\_to\_float=True, wait\_for\_buffer=True*)

Method that retrieves the buffer after it has been filled. The data retrieved from the lock-in is in a fixed-point format, which requires translation before it can be interpreted as meaningful data. When *convert\_to\_float* is True the conversion is performed (if possible) before returning the data.

#### **Parameters**

- **quantity** (*str*) – If provided, names the quantity that is to be retrieved from the curve buffer; can be any of: ‘x’, ‘y’, ‘magnitude’, ‘phase’, ‘sensitivity’, ‘adc1’, ‘adc2’, ‘adc3’, ‘dac1’, ‘dac2’, ‘noise’, ‘ratio’, ‘log ratio’, ‘event’, ‘frequency part 1’ and ‘frequency part 2’; for both dual modes, additional options are: ‘x2’, ‘y2’, ‘magnitude2’, ‘phase2’, ‘sensitivity2’. If no quantity is provided, all available data is retrieved.
- **convert\_to\_float** (*bool*) – Bool that determines whether to convert the fixed-point buffer-data to meaningful floating point values via the *buffer\_to\_float* method. If True, this method tries to convert all the available data to meaningful values; if this is not possible, an exception will be raised. If False, this conversion is not performed and the raw buffer-data is returned.
- **wait\_for\_buffer** (*bool*) – Bool that determines whether to wait for the data acquisition to finished if this method is called before the acquisition is finished. If True, the method waits until the buffer is filled before continuing; if False, the method raises an exception if the acquisition is not finished when the method is called.

**property harmonic**

An integer property that represents the reference harmonic mode control, taking values from 1 to 65535. This property can be set.

**property id**

Reads the instrument identification

**property imode**

Property that controls the voltage/current mode. can be ‘voltage mode’, ‘current mode’, or ‘low noise current mode’

**init\_curve\_buffer()**

Initializes the curve storage memory and status variables. All record of previously taken curves is removed.

**property log\_ratio**

Reads the log ratio output, defined as  $\log(X/ADC1)$

**property mag**

Reads the magnitude in Volts

**property phase**

Reads the phase in degrees

**property ratio**

Reads the ratio output, defined as  $X/ADC1$

**property reference**

Controls the oscillator reference. Can be “internal”, “external rear” or “external front”

**property reference\_phase**

A floating point property that represents the reference harmonic phase in degrees. This property can be set.

**property sensitivity**

A floating point property that controls the sensitivity range in Volts (for voltage mode) or Amps (for current modes). When in Volts it takes discrete values from 2 nV to 1 V. When in Amps it takes discrete values from 2 fA to 1  $\mu$ A (for normal current mode) or up to 10 nA (for low noise current mode). This property can be set.

**setDifferentialMode(lineFiltering=True)**

Sets lockin to differential mode, measuring A-B

**set\_buffer(points, quantities=None, interval=0.01)**

Method that prepares the curve buffer for a measurement.



**Parameters**

- **points** (*int*) – Number of points to be recorded in the curve buffer
- **quantities** (*list*) – List containing the quantities (strings) that are to be recorded in the curve buffer, can be any of: 'x', 'y', 'magnitude', 'phase', 'sensitivity', 'adc1', 'adc2', 'adc3', 'dac1', 'dac2', 'noise', 'ratio', 'log ratio', 'event', 'frequency' (or 'frequency part 1' and 'frequency part 2'); for both dual modes, additional options are: 'x2', 'y2', 'magnitude2', 'phase2', 'sensitivity2'. Default is 'x' and 'y'.
- **interval** (*float*) – The interval between two subsequent points stored in the curve buffer in s. Default is 10 ms.

**shutdown()**

Brings the instrument to a safe and stable state

**property slope**

A integer property that controls the filter slope in dB/octave, which can take the values 6, 12, 18, or 24 dB/octave. This property can be set.

**start\_buffer()**

Initiates data acquisition. Acquisition starts at the current position in the curve buffer and continues at the rate set by the STR command until the buffer is full.

**property time\_constant**

A floating point property that controls the time constant in seconds, which takes values from 10 microseconds to 50,000 seconds. This property can be set.

**property voltage**

A floating point property that represents the voltage in Volts. This property can be set.

**wait\_for\_buffer**(*timeout=None, delay=0.1*)

Method that waits until the curve buffer is filled

**property x**

Reads the X value in Volts

**property xy**

Reads both the X and Y values in Volts

**property y**

Reads the Y value in Volts

## 7.31 Stanford Research Systems

This section contains specific documentation on the Stanford Research Systems (SRS) instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.31.1 SR510 Lock-in Amplifier

**class** `pymeasure.instruments.srs.SR510(resourceName, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

**property frequency**

A float property representing the SR510 input reference frequency

**property output**

A float property that represents the SR510 output voltage in Volts.

**property phase**

A float property that represents the SR510 reference to input phase offset in degrees. Queries return values between -180 and 180 degrees. This property can be set with a range of values between -999 to 999 degrees. Set values are mapped internal in the lockin to -180 and 180 degrees.

**property sensitivity**

A float property that represents the SR510 sensitivity value. This property can be set.

**property status**

A string property representing the bits set within the SR510 status byte

**property time\_constant**

A float property that represents the SR510 PRE filter time constant. This property can be set.

### 7.31.2 SR570 Lock-in Amplifier

**class** `pymeasure.instruments.srs.SR570(resourceName, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

**property bias\_enabled**

Boolean that turns the bias on or off. Allowed values are: True (bias on) and False (bias off)

**property bias\_level**

A floating point value in V that sets the bias voltage level of the amplifier, in the [-5V,+5V] limits. The values are up to 1 mV precision level.

**blank\_front()**

“Blanks the frontend output of the device

**clear\_overload()**

“Reset the filter capacitors to clear an overload condition

**disable\_bias()**

Turns the bias voltage off

**disable\_offset\_current()**

“Disables the offset current

**enable\_bias()**

Turns the bias voltage on

**enable\_offset\_current()**

“Enables the offset current

**property filter\_type**

A string that sets the filter type. Allowed values are: ['6dB Highpass', '12dB Highpass', '6dB Bandpass', '6dB Lowpass', '12dB Lowpass', 'none']

**property front\_blanked**

Boolean that blanks(True) or un-blanks (False) the front panel

**property gain\_mode**

A string that sets the gain mode. Allowed values are: ['Low Noise', 'High Bandwidth', 'Low Drift']

**property high\_freq**

A floating point value that sets the highpass frequency of the amplifier, which takes a discrete value in a 1-3 sequence. Values are truncated to the closest allowed value if not exact. Allowed values range from 0.03 Hz to 1 MHz.

**property invert\_signal\_sign**

An boolean sets the signal invert sense. Allowed values are: True (inverted) and False (not inverted).

**property low\_freq**

A floating point value that sets the lowpass frequency of the amplifier, which takes a discrete value in a 1-3 sequence. Values are truncated to the closest allowed value if not exact. Allowed values range from 0.03 Hz to 1 MHz.

**property offset\_current**

A floating point value in A that sets the absolute value of the offset current of the amplifier, in the [1pA,5mA] limits. The offset current takes discrete values in a 1-2-5 sequence. Values are truncated to the closest allowed value if not exact.

**property offset\_current\_enabled**

Boolean that turns the offset current on or off. Allowed values are: True (current on) and False (current off).

**property offset\_current\_sign**

An string that sets the offset current sign. Allowed values are: 'positive' and 'negative'.

**property sensitivity**

A floating point value that sets the sensitivity of the amplifier, which takes discrete values in a 1-2-5 sequence. Values are truncated to the closest allowed value if not exact. Allowed values range from 1 pA/V to 1 mA/V.

**property signal\_inverted**

Boolean that inverts the signal if True

**unblank\_front()**

Un-blanks the frontend output of the device

### 7.31.3 SR830 Lock-in Amplifier

**class** `pymeasure.instruments.srs.SR830(resource_name, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

**property adc1**

Reads the Aux input 1 value in Volts with 1/3 mV resolution.

**property adc2**

Reads the Aux input 2 value in Volts with 1/3 mV resolution.

**property adc3**

Reads the Aux input 3 value in Volts with 1/3 mV resolution.

**property adc4**

Reads the Aux input 4 value in Volts with 1/3 mV resolution.

**auto\_offset(channel)**

Offsets the channel (X, Y, or R) to zero

**property aux\_in\_1**

Reads the Aux input 1 value in Volts with 1/3 mV resolution.

**property aux\_in\_2**

Reads the Aux input 2 value in Volts with 1/3 mV resolution.

**property aux\_in\_3**

Reads the Aux input 3 value in Volts with 1/3 mV resolution.

**property aux\_in\_4**

Reads the Aux input 4 value in Volts with 1/3 mV resolution.

**property aux\_out\_1**

A floating point property that controls the output of Aux output 1 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property aux\_out\_2**

A floating point property that controls the output of Aux output 2 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property aux\_out\_3**

A floating point property that controls the output of Aux output 3 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property aux\_out\_4**

A floating point property that controls the output of Aux output 4 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property channel1**

A string property that represents the type of Channel 1, taking the values X, R, X Noise, Aux In 1, or Aux In 2. This property can be set.

**property channel2**

A string property that represents the type of Channel 2, taking the values Y, Theta, Y Noise, Aux In 3, or Aux In 4. This property can be set.

**property dac1**

A floating point property that controls the output of Aux output 1 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property dac2**

A floating point property that controls the output of Aux output 2 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property dac3**

A floating point property that controls the output of Aux output 3 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property dac4**

A floating point property that controls the output of Aux output 4 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property err\_status**

Reads the value of the lockin error (ERR) status byte. Returns an IntFlag type with positions within the string corresponding to different error flags: bit 0: unused bit 1: backup error bit 2: RAM error bit 3: unused bit 4: ROM error bit 5: GPIB error bit 6: DSP error bit 7: Math error

**property filter\_slope**

An integer property that controls the filter slope, which can take on the values 6, 12, 18, and 24 dB/octave. Values are truncated to the next highest level if they are not exact.

**property frequency**

A floating point property that represents the lock-in frequency in Hz. This property can be set.

**get\_buffer(channel=1, start=0, end=None)**

Acquires the 32 bit floating point data through binary transfer

**get\_scaling(channel)**

Returns the offset present and the expansion term that are used to scale the channel in question

**property harmonic**

An integer property that controls the harmonic that is measured. Allowed values are 1 to 19999. Can be set.

**property input\_config**

An string property that controls the input configuration. Allowed values are: ['A', 'A - B', 'I (1 MOhm)', 'I (100 MOhm)']

**property input\_coupling**

An string property that controls the input coupling. Allowed values are: ['AC', 'DC']

**property input\_grounding**

An string property that controls the input shield grounding. Allowed values are: ['Float', 'Ground']

**property input\_notch\_config**

An string property that controls the input line notch filter status. Allowed values are: ['None', 'Line', '2 x Line', 'Both']

**is\_out\_of\_range()**

Returns True if the magnitude is out of range

**property lia\_status**

Reads the value of the lockin amplifier (LIA) status byte. Returns a binary string with positions within the string corresponding to different status flags: bit 0: Input/Amplifier overload bit 1: Time constant filter overload bit 2: Output overload bit 3: Reference unlock bit 4: Detection frequency range switched bit 5: Time constant changed indirectly bit 6: Data storage triggered bit 7: unused

**property magnitude**

Reads the magnitude in Volts.

**output\_conversion(channel)**

Returns a function that can be used to determine the signal from the channel output (X, Y, or R)

**property phase**

A floating point property that represents the lock-in phase in degrees. This property can be set.

**quick\_range()**

While the magnitude is out of range, increase the sensitivity by one setting

**property reference\_source**

An string property that controls the reference source. Allowed values are: ['External', 'Internal']

**property sample\_frequency**

Gets the sample frequency in Hz

**property sensitivity**

A floating point property that controls the sensitivity in Volts, which can take discrete values from 2 nV to 1 V. Values are truncated to the next highest level if they are not exact.

**set\_scaling(channel, percent, expand=0)**

Sets the offset of a channel (X=1, Y=2, R=3) to a certain percent (-105% to 105%) of the signal, with an optional expansion term (0, 10=1, 100=2)

**property sine\_voltage**

A floating point property that represents the reference sine-wave voltage in Volts. This property can be set.

**snap**(*val1*='X', *val2*='Y', \**vals*)

Method that records and retrieves 2 to 6 parameters at a single instant. The parameters can be one of: X, Y, R, Theta, Aux In 1, Aux In 2, Aux In 3, Aux In 4, Frequency, CH1, CH2. Default is “X” and “Y”.

**Parameters**

- **val1** – first parameter to retrieve
- **val2** – second parameter to retrieve
- **vals** – other parameters to retrieve (optional)

**property theta**

Reads the theta value in degrees.

**property time\_constant**

A floating point property that controls the time constant in seconds, which can take discrete values from 10 microseconds to 30,000 seconds. Values are truncated to the next highest level if they are not exact.

**wait\_for\_buffer**(*count*, *has\_aborted*=<function SR830.<lambda>>, *timeout*=60, *timestep*=0.01)

Wait for the buffer to fill a certain count

**property x**

Reads the X value in Volts.

**property xy**

Reads the X and Y values in Volts.

**property y**

Reads the Y value in Volts.

## 7.31.4 SR860 Lock-in Amplifier

**class** `pymeasure.instruments.srs.SR860`(*resourceName*, \*\**kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

**property adc1**

Reads the Aux input 1 value in Volts with 1/3 mV resolution.

**property adc2**

Reads the Aux input 2 value in Volts with 1/3 mV resolution.

**property adc3**

Reads the Aux input 3 value in Volts with 1/3 mV resolution.

**property adc4**

Reads the Aux input 4 value in Volts with 1/3 mV resolution.

**property aux\_in\_1**

Reads the Aux input 1 value in Volts with 1/3 mV resolution.

**property aux\_in\_2**

Reads the Aux input 2 value in Volts with 1/3 mV resolution.

**property aux\_in\_3**

Reads the Aux input 3 value in Volts with 1/3 mV resolution.

**property aux\_in\_4**

Reads the Aux input 4 value in Volts with 1/3 mV resolution.

**property aux\_out\_1**

A floating point property that controls the output of Aux output 1 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property aux\_out\_2**

A floating point property that controls the output of Aux output 2 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property aux\_out\_3**

A floating point property that controls the output of Aux output 3 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property aux\_out\_4**

A floating point property that controls the output of Aux output 4 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property dac1**

A floating point property that controls the output of Aux output 1 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property dac2**

A floating point property that controls the output of Aux output 2 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property dac3**

A floating point property that controls the output of Aux output 3 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property dac4**

A floating point property that controls the output of Aux output 4 in Volts, taking values between -10.5 V and +10.5 V. This property can be set.

**property dcmode**

A string property that represents the sine out dc mode. This property can be set. Allowed values are: ['COM', 'DIF', 'common', 'difference']

**property detectedfrequency**

Returns the actual detected frequency in HZ.

**property extfrequency**

Returns the external frequency in Hz.

**property filer\_synchronous**

A string property that represents the synchronous filter. This property can be set. Allowed values are: ['Off', 'On']

**property filter\_advanced**

A string property that represents the advanced filter. This property can be set. Allowed values are: ['Off', 'On']

**property filter\_slope**

A integer property that sets the filter slope to 6 dB/oct(i=0), 12 DB/oct(i=1), 18 dB/oct(i=2), 24 dB/oct(i=3).

**property frequency**

A floating point property that represents the lock-in frequency in Hz. This property can be set.

**property frequencypreset1**

A floating point property that represents the preset frequency for the F1 preset button. This property can be set.

**property frequencypreset2**

A floating point property that represents the preset frequency for the F2 preset button. This property can be set.

**property frequencypreset3**

A floating point property that represents the preset frequency for the F3 preset button. This property can be set.

**property frequencypreset4**

A floating point property that represents the preset frequency for the F4 preset button. This property can be set.

**property front\_panel**

Turns the front panel blanking on(i=0) or off(i=1).

**property get\_noise\_bandwidth**

Returns the equivalent noise bandwidth, in hertz.

**property get\_signal\_strength\_indicator**

Returns the signal strength indicator.

**property gettimebase**

Returns the current 10 MHz timebase source.

**property harmonic**

An integer property that controls the harmonic that is measured. Allowed values are 1 to 99. Can be set.

**property harmonicdual**

An integer property that controls the harmonic in dual reference mode that is measured. Allowed values are 1 to 99. Can be set.

**property horizontal\_time\_div**

A integer property that sets the horizontal time/div according to the following table: ['0=0.5s', '1=1s', '2=2s', '3=5s', '4=10s', '5=30s', '6=1min', '7=2min', '8=5min', '9=10min', '10=30min', '11=1hour', '12=2hour', '13=6hour', '14=12hour', '15=1day', '16=2days']

**property input\_coupling**

A string property that represents the input coupling. This property can be set. Allowed values are: ['AC', 'DC']

**property input\_current\_gain**

A string property that represents the current input gain. This property can be set. Allowed values are: ['1MEG', '100MEG']

**property input\_range**

A string property that represents the input range. This property can be set. Allowed values are: ['1V', '300M', '100M', '30M', '10M']

**property input\_shields**

A string property that represents the input shield grounding. This property can be set. Allowed values are: ['Float', 'Ground']

**property input\_signal**

A string property that represents the signal input. This property can be set. Allowed values are: ['VOLT', 'CURR', 'voltage', 'current']



**property input\_voltage\_mode**

A string property that represents the voltage input mode. This property can be set. Allowed values are: ['A', 'A-B']

**property internalfrequency**

A floating property that represents the internal lock-in frequency in Hz This property can be set.

**property magnitude**

Reads the magnitude in Volts.

**property parameter\_DAT1**

A integer property that assigns a parameter to data channel 1(green). This parameters can be set. Allowed values are: ['i=', '0=Xoutput', '1=Youtput', '2=Routput', 'Thetaoutput', '4=Aux IN1', '5=Aux IN2', '6=Aux IN3', '7=Aux IN4', '8=Xnoise', '9=Ynoise', '10=AUXOut1', '11=AuxOut2', '12=Phase', '13=Sine Out amplitude', '14=DCLlevel', '15I=nt.referenceFreq', '16=Ext.referenceFreq']

**property parameter\_DAT2**

A integer property that assigns a parameter to data channel 2(blue). This parameters can be set. Allowed values are: ['i=', '0=Xoutput', '1=Youtput', '2=Routput', 'Thetaoutput', '4=Aux IN1', '5=Aux IN2', '6=Aux IN3', '7=Aux IN4', '8=Xnoise', '9=Ynoise', '10=AUXOut1', '11=AuxOut2', '12=Phase', '13=Sine Out amplitude', '14=DCLlevel', '15I=nt.referenceFreq', '16=Ext.referenceFreq']

**property parameter\_DAT3**

A integer property that assigns a parameter to data channel 3(yellow). This parameters can be set. Allowed values are: ['i=', '0=Xoutput', '1=Youtput', '2=Routput', 'Thetaoutput', '4=Aux IN1', '5=Aux IN2', '6=Aux IN3', '7=Aux IN4', '8=Xnoise', '9=Ynoise', '10=AUXOut1', '11=AuxOut2', '12=Phase', '13=Sine Out amplitude', '14=DCLlevel', '15I=nt.referenceFreq', '16=Ext.referenceFreq']

**property parameter\_DAT4**

A integer property that assigns a parameter to data channel 3(orange). This parameters can be set. Allowed values are: ['i=', '0=Xoutput', '1=Youtput', '2=Routput', 'Thetaoutput', '4=Aux IN1', '5=Aux IN2', '6=Aux IN3', '7=Aux IN4', '8=Xnoise', '9=Ynoise', '10=AUXOut1', '11=AuxOut2', '12=Phase', '13=Sine Out amplitude', '14=DCLlevel', '15I=nt.referenceFreq', '16=Ext.referenceFreq']

**property phase**

A floating point property that represents the lock-in phase in degrees. This property can be set.

**property reference\_externalinput**

A string property that represents the external reference input. This property can be set. Allowed values are: ['50OHMS', '1MEG']

**property reference\_source**

A string property that represents the reference source. This property can be set. Allowed values are: ['INT', 'EXT', 'DUAL', 'CHOP']

**property reference\_triggermode**

A string property that represents the external reference trigger mode. This property can be set. Allowed values are: ['SIN', 'POS', 'NEG', 'POSTTL', 'NEGTL']

**property screen\_layout**

A integer property that Sets the screen layout to trend(i=0), full strip chart history(i=1), half strip chart history(i=2), full FFT(i=3), half FFT(i=4) or big numerical(i=5).

**screenshot()**

Take screenshot on device The DCAP command saves a screenshot to a USB memory stick. This command is the same as pressing the [Screen Shot] key. A USB memory stick must be present in the front panel USB port.

**property sensitivity**

A floating point property that controls the sensitivity in Volts, which can take discrete values from 2 nV

to 1 V. Values are truncated to the next highest level if they are not exact.

**property sine\_amplitudepreset1**

A floating point property that represents the preset sine out amplitude, for the A1 preset button. This property can be set.

**property sine\_amplitudepreset2**

A floating point property that represents the preset sine out amplitude, for the A2 preset button. This property can be set.

**property sine\_amplitudepreset3**

A floating point property that represents the preset sine out amplitude, for the A3 preset button. This property can be set.

**property sine\_amplitudepreset4**

A floating point property that represents the preset sine out amplitude, for the A3 preset button. This property can be set.

**property sine\_dclevelpreset1**

A floating point property that represents the preset sine out dc level for the L1 button. This property can be set.

**property sine\_dclevelpreset2**

A floating point property that represents the preset sine out dc level for the L2 button. This property can be set.

**property sine\_dclevelpreset3**

A floating point property that represents the preset sine out dc level for the L3 button. This property can be set.

**property sine\_dclevelpreset4**

A floating point property that represents the preset sine out dc level for the L4 button. This property can be set.

**property sine\_voltage**

A floating point property that represents the reference sine-wave voltage in Volts. This property can be set.

**snap**(val1='X', val2='Y', val3=None)

retrieve 2 or 3 parameters at once parameters can be chosen by index, or enumeration as follows:

j enumeration parameter j enumeration parameter

0 X X output 9 YNOise Ynoise 1 Y Youtput 10 OUT1 Aux Out1 2 R R output 11 OUT2 Aux Out2 3  
THeta output 12 PHase Reference Phase 4 IN1 Aux In1 13 SAMp Sine Out Amplitude 5 IN2 Aux In2 14  
LEVel DC Level 6 IN3 Aux In3 15 FInt Int. Ref. Frequency 7 IN4 Aux In4 16 FExt Ext. Ref. Frequency  
8 XNOise Xnoise

**Parameters**

- **val1** – parameter enumeration/index
- **val2** – parameter enumeration/index
- **val3** – parameter enumeration/index (optional)

**Defaults:** val1 = “X” val2 = “Y” val3 = None

**property strip\_chart\_dat1**

A integer property that turns the strip chart graph of data channel 1 off(i=0) or on(i=1).

**property strip\_chart\_dat2**

A integer property that turns the strip chart graph of data channel 2 off(*i*=0) or on(*i*=1).

**property strip\_chart\_dat3**

A integer property that turns the strip chart graph of data channel 1 off(*i*=0) or on(*i*=1).

**property strip\_chart\_dat4**

A integer property that turns the strip chart graph of data channel 4 off(*i*=0) or on(*i*=1).

**property theta**

Reads the theta value in degrees.

**property time\_constant**

A floating point property that controls the time constant in seconds, which can take discrete values from 10 microseconds to 30,000 seconds. Values are truncated to the next highest level if they are not exact.

**property timebase**

Sets the external 10 MHz timebase to auto(*i*=0) or internal(*i*=1).

**property x**

Reads the X value in Volts

**property y**

Reads the Y value in Volts

## 7.32 Tektronix

This section contains specific documentation on the Tektronix instruments that are implemented. If you are interested in an instrument not included, please consider [adding the instrument](#).

### 7.32.1 TDS2000 Oscilloscope

**class** `pymeasure.instruments.tektronix.TDS2000`(*resourceName*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Tektronix TDS 2000 Oscilloscope and provides a high-level for interacting with the instrument

### 7.32.2 AFG3152C Arbitrary function generator

**class** `pymeasure.instruments.tektronix.AFG3152C`(*adapter*, *\*\*kwargs*)

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Tektronix AFG 3000 series (one or two channels) arbitrary function generator and provides a high-level for interacting with the instrument.

```
afg=AFG3152C("GPIB::1") # AFG on GPIB 1
afg.reset() # Reset to default
afg.ch1.shape='sinusoidal' # Sinusoidal shape
afg.ch1.unit='VPP' # Sets CH1 unit to VPP
afg.ch1.amp_vpp=1 # Sets the CH1 level to 1 VPP
afg.ch1.frequency=1e3 # Sets the CH1 frequency to 1KHz
afg.ch1.enable() # Enables the output from CH1
```

## 7.33 Thorlabs

This section contains specific documentation on the Thorlabs instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.33.1 Thorlabs PM100USB Powermeter

**class** `pymeasure.instruments.thorlabs.ThorlabsPM100USB(adapter, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents Thorlabs PM100USB powermeter.

**property** `energy`

Energy in J.

**property** `power`

Power in W.

**property** `wavelength`

Wavelength in nm.

**property** `wavelength_max`

Maximum wavelength, in nm

**property** `wavelength_min`

Minimum wavelength, in nm

### 7.33.2 Thorlabs Pro 8000 modular laser driver

**class** `pymeasure.instruments.thorlabs.ThorlabsPro8000(resourceName, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents Thorlabs Pro 8000 modular laser driver

**property** `LDCCurrent`

Laser current.

**property** `LDCCurrentLimit`

Set Software current Limit (value must be lower than hardware current limit).

**property** `LDCPolarity`

Set laser diode polarity. Allowed values are: ['AG', 'CG']

**property** `LDCStatus`

Set laser diode status. Allowed values are: ['ON', 'OFF']

**property** `TEDSetTemperature`

Set TEC temperature

**property** `TEDStatus`

Set TEC status. Allowed values are: ['ON', 'OFF']

**property** `slot`

Slot selection. Allowed values are: range(1, 9)

## 7.34 Yokogawa

This section contains specific documentation on the Yokogawa instruments that are implemented. If you are interested in an instrument not included, please consider *adding the instrument*.

### 7.34.1 Yokogawa 7651 Programmable Supply

**class** `pymeasure.instruments.yokogawa.Yokogawa7651(adapter, **kwargs)`

Bases: `pymeasure.instruments.instrument.Instrument`

Represents the Yokogawa 7651 Programmable DC Source and provides a high-level for interacting with the instrument.

```
yoko = Yokogawa7651("GPIB::1")

yoko.apply_current()           # Sets up to source current
yoko.source_current_range = 10e-3 # Sets the current range to 10 mA
yoko.compliance_voltage = 10    # Sets the compliance voltage to 10 V
yoko.source_current = 0         # Sets the source current to 0 mA

yoko.enable_source()          # Enables the current output
yoko.ramp_to_current(5e-3)     # Ramps the current to 5 mA

yoko.shutdown()               # Ramps the current to 0 mA and disables output
```

**apply\_current**(*max\_current=0.001, compliance\_voltage=1*)

Configures the instrument to apply a source current, which can take optional parameters that defer to the `source_current_range` and `compliance_voltage` properties.

**apply\_voltage**(*max\_voltage=1, compliance\_current=0.01*)

Configures the instrument to apply a source voltage, which can take optional parameters that defer to the `source_voltage_range` and `compliance_current` properties.

**property compliance\_current**

A floating point property that sets the compliance current in Amps, which can take values from 5 to 120 mA.

**property compliance\_voltage**

A floating point property that sets the compliance voltage in Volts, which can take values between 1 and 30 V.

**disable\_source()**

Disables the source of current or voltage depending on the configuration of the instrument.

**enable\_source()**

Enables the source of current or voltage depending on the configuration of the instrument.

**property id**

Returns the identification of the instrument

**ramp\_to\_current**(*current, steps=25, duration=0.5*)

Ramps the current to a value in Amps by traversing a linear spacing of current steps over a duration, defined in seconds.

#### Parameters

- **steps** – A number of linear steps to traverse

- **duration** – A time in seconds over which to ramp

**ramp\_to\_voltage**(*voltage*, *steps*=25, *duration*=0.5)

Ramps the voltage to a value in Volts by traversing a linear spacing of voltage steps over a duration, defined in seconds.

**Parameters**

- **steps** – A number of linear steps to traverse
- **duration** – A time in seconds over which to ramp

**shutdown**()

Shuts down the instrument, and ramps the current or voltage to zero before disabling the source.

**property source\_current**

A floating point property that controls the source current in Amps, if that mode is active.

**property source\_current\_range**

A floating point property that sets the current voltage range in Amps, which can take values: 1 mA, 10 mA, and 100 mA. Currents are truncated to an appropriate value if needed.

**property source\_enabled**

Reads a boolean value that is True if the source is enabled, determined by checking if the 5th bit of the OC flag is a binary 1.

**property source\_mode**

A string property that controls the source mode, which can take the values ‘current’ or ‘voltage’. The convenience methods [apply\\_current\(\)](#) and [apply\\_voltage\(\)](#) can also be used.

**property source\_voltage**

A floating point property that controls the source voltage in Volts, if that mode is active.

**property source\_voltage\_range**

A floating point property that sets the source voltage range in Volts, which can take values: 10 mV, 100 mV, 1 V, 10 V, and 30 V. Voltages are truncated to an appropriate value if needed.

## 7.34.2 Yokogawa GS200 Source

**class** `pymeasure.instruments.yokogawa.YokogawaGS200`(*adapter*, *\*\*kwargs*)

Bases: [pymeasure.instruments.instrument.Instrument](#)

Represents the Yokogawa GS200 source and provides a high-level interface for interacting with the instrument.

**property current\_limit**

Floating point number that controls the current limit. “Limit” refers to maximum value of the electrical value that is conjugate to the mode (current is conjugate to voltage, and vice versa). Thus, current limit is only applicable when in ‘voltage’ mode

**property source\_enabled**

A boolean property that controls whether the source is enabled, takes values True or False.

**property source\_level**

Floating point number that controls the output level, either a voltage or a current, depending on the source mode.

**property source\_mode**

String property that controls the source mode. Can be either ‘current’ or ‘voltage’.

**property source\_range**

Floating point number that controls the range (either in voltage or current) of the output. “Range” refers to the maximum source level.

**trigger\_ramp\_to\_level**(*level*, *ramp\_time*)

Ramp the output level from its current value to “level” in time “ramp\_time”. This method will NOT wait until the ramp is finished (thus, it will not block further code evaluation).

**Parameters**

- **level** (*float*) – final output level
- **ramp\_time** (*float*) – time in seconds to ramp

**Returns** None**property voltage\_limit**

Floating point number that controls the voltage limit. “Limit” refers to maximum value of the electrical value that is conjugate to the mode (current is conjugate to voltage, and vice versa). Thus, voltage limit is only applicable when in ‘current’ mode





## CONTRIBUTING

Contributions to the instrument repository and the main code base are highly encouraged. This section outlines the basic work-flow for new contributors.

### 8.1 Using the development version

New features are added to the development version of PyMeasure, hosted on [GitHub](#). We use [Git version control](#) to track and manage changes to the source code. On Windows, we recommend using [GitHub Desktop](#). Make sure you have an appropriate version of Git (or GitHub Desktop) installed and that you have a GitHub account.

In order to add your feature, you need to first [fork](#) PyMeasure. This will create a copy of the repository under your GitHub account.

The instructions below assume that you have set up Anaconda, as described in the [Quick Start guide](#) and describe the terminal commands necessary. If you are using GitHub Desktop, take a look through [their documentation](#) to understand the corresponding steps.

Clone your fork of PyMeasure `your-github-username/pymasure`. In the following terminal commands replace your desired path and GitHub username.

```
cd /path/for/code
git clone https://github.com/your-github-username/pymasure.git
```

If you had already installed PyMeasure using `pip`, make sure to uninstall it before continuing.

```
pip uninstall pymasure
```

Install PyMeasure in the editable mode.

```
cd /path/for/code/pymasure
pip install -e .
```

This will allow you to edit the files of PyMeasure and see the changes reflected. Make sure to reset your notebook kernel or Python console when doing so. Now you have your own copy of the development version of PyMeasure installed!

## 8.2 Working on a new feature

We use branches in Git to allow multiple features to be worked on simultaneously, without causing conflicts. The master branch contains the stable development version. Instead of working on the master branch, you will create your own branch off the master and merge it back into the master when you are finished.

Create a new branch for your feature before editing the code. For example, if you want to add the new instrument “Extreme 5000” you will make a new branch “dev/extreme-5000”.

```
git branch dev/extreme-5000
```

You can also [make a new branch](#) on GitHub. If you do so, you will have to fetch these changes before the branch will show up on your local computer.

```
git fetch
```

Once you have created the branch, change your current branch to match the new one.

```
git checkout dev/extreme-5000
```

Now you are ready to write your new feature and make changes to the code. To ensure consistency, please follow the [coding standards for PyMeasure](#). Use `git status` to check on the files that have been changed. As you go, commit your changes and push them to your fork.

```
git add file-that-changed.py
git commit -m "A short description about what changed"
git push
```

## 8.3 Making a pull-request

While you are working, its helpful to start a pull-request (PR) on the master branch of `pymasure/pymasure`. This will allow you to discuss your feature with other contributors. We encourage you to start this pull-request after your first commit.

[Start a pull-request on the PyMeasure GitHub page](#).

Your pull-request will be merged by the PyMeasure maintainers once it meets the coding standards and passes unit tests. You will notice that your pull-request is automatically checked with the unit tests.

## 8.4 Unit testing

Unit tests are run each time a new commit is made to a branch. The purpose is to catch changes that break the current functionality, by testing each feature unit. PyMeasure relies on `pytest` to preform these tests, which are run on TravisCI and Appveyor for Linux/macOS and Windows respectively.

Running the unit tests while you develop is highly encouraged. This will ensure that you have a working contribution when you create a pull request.

```
python setup.py test
```

If your feature can be tested, unit tests are required. This will ensure that your features keep working as new features are added.

Now you are familiar with all the pieces of the PyMeasure development work-flow. We look forward to seeing your pull-request!



## **REPORTING AN ERROR**

Please report all errors to the [Issues section](#) of the PyMeasure GitHub repository. Use the search function to determine if there is an existing or resolved issued before posting.



## ADDING INSTRUMENTS

You can make a significant contribution to PyMeasure by adding a new instrument to the `pymeasure.instruments` package. Even adding an instrument with a few features can help get the ball rolling, since its likely that others are interested in the same instrument.

Before getting started, become familiar with the [contributing work-flow](#) for PyMeasure, which steps through the process of adding a new feature (like an instrument) to the development version of the source code. This section will describe how to lay out your instrument code.

### 10.1 File structure

Your new instrument should be placed in the directory corresponding to the manufacturer of the instrument. For example, if you are going to add an “Extreme 5000” instrument you should add the following files assuming “Extreme” is the manufacturer. Use lowercase for all filenames to distinguish packages from CamelCase Python classes.

```
pymeasure/pymeasure/instruments/extreme/  
|--> __init__.py  
|--> extreme5000.py
```

#### 10.1.1 Updating the init file

The `__init__.py` file in the manufacturer directory should import all of the instruments that correspond to the manufacturer, to allow the files to be easily imported. For a new manufacturer, the manufacturer should also be added to `pymeasure/pymeasure/instruments/__init__.py`.

#### 10.1.2 Adding documentation

Documentation for each instrument is required, and helps others understand the features you have implemented. Add a new reStructuredText file to the documentation.

```
pymeasure/docs/api/instruments/extreme/  
|--> index.rst  
|--> extreme5000.rst
```

Copy an existing instrument documentation file, which will automatically generate the documentation for the instrument. The `index.rst` file should link to the `extreme5000` file. For a new manufacturer, the manufacturer should be also linked in `pymeasure/docs/api/instruments/index.rst`.

## 10.2 Instrument file

All standard instruments should be child class of *Instrument*. This provides the basic functionality for working with *Adapters*, which perform the actual communication.

The most basic instrument, for our “Extreme 5000” example starts like this:

```
#
# This file is part of the PyMeasure package.
#
# Copyright (c) 2013-2021 PyMeasure Developers
#
# Permission is hereby granted, free of charge, to any person obtaining a copy
# of this software and associated documentation files (the "Software"), to deal
# in the Software without restriction, including without limitation the rights
# to use, copy, modify, merge, publish, distribute, sublicense, and/or sell
# copies of the Software, and to permit persons to whom the Software is
# furnished to do so, subject to the following conditions:
#
# The above copyright notice and this permission notice shall be included in
# all copies or substantial portions of the Software.
#
# THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR
# IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY,
# FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE
# AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER
# LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM,
# OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN
# THE SOFTWARE.
#
# from pymeasure.instruments import Instrument
```

This is a minimal instrument definition:

```
class Extreme5000(Instrument):
    """ Represents the imaginary Extreme 5000 instrument.
    """

    def __init__(self, resourceName, **kwargs):
        super().__init__(
            resourceName,
            "Extreme 5000",
            **kwargs
        )
```

Make sure to include the PyMeasure license to each file, and add yourself as an author to the `AUTHORS.txt` file.

In principle you are free to write any functions that are necessary for interacting with the instrument. When doing so, make sure to use the `self.ask(command)`, `self.write(command)`, and `self.read()` methods to issue command instead of calling the adapter directly.

In practice, we have developed a number of convenience functions for making instruments easy to write and maintain. The following sections detail these conveniences and are highly encouraged.



## 10.3 Writing properties

In PyMeasure, [Python properties](#) are the preferred method for dealing with variables that are read or set. PyMeasure comes with two convenience functions for making properties for classes. The `Instrument.measurement` function returns a property that issues a GPIB/SCPI requests when the value is used. For example, if our “Extreme 5000” has the `*IDN?` command we can write the following property to be added above the `def __init__` line in our above example class, or added to the class after the fact as in the code here:

```
Extreme5000.id = Instrument.measurement(
    "*IDN?", """ Reads the instrument identification """
)
```

You will notice that a documentation string is required, and should be descriptive and specific.

When we use this property we will get the identification information.

```
>>> extreme = Extreme5000("GPIB::1")
>>> extreme.id           # Reads "*IDN?"
'Extreme 5000 identification from instrument'
```

The `Instrument.control` function extends this behavior by creating a property that you can read and set. For example, if our “Extreme 5000” has the `:VOLT?` and `:VOLT <float>` commands that are in Volts, we can write the following property.

```
Extreme5000.voltage = Instrument.control(
    ":VOLT?", ":VOLT %g",
    """ A floating point property that controls the voltage
    in Volts. This property can be set.
    """
)
```

You will notice that we use the [Python string format](#) `%g` to pass through the floating point.

We can use this property to set the voltage to 100 mV, which will execute the command and then request the current voltage.

```
>>> extreme = Extreme5000("GPIB::1")
>>> extreme.voltage = 0.1      # Executes ":VOLT 0.1"
>>> extreme.voltage          # Reads ":VOLT?"
0.1
```

Using `Instrument.control` and `Instrument.measurement` functions, you can create a number of properties for basic measurements and controls.

The `Instrument.control` function can be used with multiple values at once, passed as a tuple. Say, we may set voltages and frequencies in our “Extreme 5000”, and the the commands for this are `:VOLT:FREQ?` and `:VOLT:FREQ <float>,<float>`, we could use the following property:

```
Extreme5000.combination = Instrument.control(
    ":VOLT:FREQ?", ":VOLT:FREQ %g,%g",
    """ A floating point property that simultaneously controls the voltage
    in Volts and the frequency in Hertz. This property can be set by a tuple.
    """
)
```

In use, we could set the voltage to 200 mV, and the Frequency to 931 Hz, and read both values immediately afterwards.

```
>>> extreme = Extreme5000("GPIB::1")
>>> extreme.combination = (0.2, 931)      # Executes ":VOLT:FREQ 0.2,931"
>>> extreme.combination                  # Reads ":VOLT:FREQ?"
[0.2, 931.0]
```

The next section details additional features of `Instrument.control` that allow you to write properties that cover specific ranges, or have to map between a real value to one used in the command. Furthermore it is shown how to perform more complex processing of return values from your device.

## 10.4 Advanced properties

Many GPIB/SCPI commands are more restrictive than our basic examples above. The `Instrument.control` function has the ability to encode these restrictions using *validators*. A validator is a function that takes a value and a set of values, and returns a valid value or raises an exception. There are a number of pre-defined validators in `pymeasure.instruments.validators` that should cover most situations. We will cover the four basic types here.

In the examples below we assume you have imported the validators.

In many situations you will also need to process the return string in order to extract the wanted quantity or process a value before sending it to the device. The `Instrument.control`, `Instrument.measurement` and `Instrument.setting` function also provide means to achieve this.

### 10.4.1 In a restricted range

If you have a property with a restricted range, you can use the `strict_range` and `truncated_range` functions.

For example, if our “Extreme 5000” can only support voltages from -1 V to 1 V, we can modify our previous example to use a strict validator over this range.

```
Extreme5000.voltage = Instrument.control(
    ":VOLT?", ":VOLT %g",
    """ A floating point property that controls the voltage
    in Volts, from -1 to 1 V. This property can be set. """ ,
    validator=strict_range,
    values=[-1, 1]
)
```

Now our voltage will raise a `ValueError` if the value is out of the range.

```
>>> extreme = Extreme5000("GPIB::1")
>>> extreme.voltage = 100
Traceback (most recent call last):
...
ValueError: Value of 100 is not in range [-1,1]
```

This is useful if you want to alert the programmer that they are using an invalid value. However, sometimes it can be nicer to truncate the value to be within the range.

```
Extreme5000.voltage = Instrument.control(
    ":VOLT?", ":VOLT %g",
    """ A floating point property that controls the voltage
    in Volts, from -1 to 1 V. Invalid voltages are truncated.
```

(continues on next page)

(continued from previous page)

```

    This property can be set. """,
    validator=truncated_range,
    values=[-1, 1]
)

```

Now our voltage will not raise an error, and will truncate the value to the range bounds.

```

>>> extreme = Extreme5000("GPIB::1")
>>> extreme.voltage = 100          # Executes ":VOLT 1"
>>> extreme.voltage
1.0

```

### 10.4.2 In a discrete set

Often a control property should only take a few discrete values. You can use the `strict_discrete_set` and `truncated_discrete_set` functions to handle these situations. The strict version raises an error if the value is not in the set, as in the range examples above.

For example, if our “Extreme 5000” has a `:RANG <float>` command that sets the voltage range that can take values of 10 mV, 100 mV, and 1 V in Volts, then we can write a control as follows.

```

Extreme5000.voltage = Instrument.control(
    ":RANG?", ":RANG %g",
    """ A floating point property that controls the voltage
    range in Volts. This property can be set.
    """,
    validator=truncated_discrete_set,
    values=[10e-3, 100e-3, 1]
)

```

Now we can set the voltage range, which will automatically truncate to an appropriate value.

```

>>> extreme = Extreme5000("GPIB::1")
>>> extreme.voltage = 0.08
>>> extreme.voltage
0.1

```

### 10.4.3 Using maps

Now that you are familiar with the validators, you can additionally use maps to satisfy instruments which require non-physical values. The `map_values` argument of `Instrument.control` enables this feature.

If your set of values is a list, then the command will use the index of the list. For example, if our “Extreme 5000” instead has a `:RANG <integer>`, where 0, 1, and 2 correspond to 10 mV, 100 mV, and 1 V, then we can use the following control.

```

Extreme5000.voltage = Instrument.control(
    ":RANG?", ":RANG %d",
    """ A floating point property that controls the voltage
    range in Volts, which takes values of 10 mV, 100 mV and 1 V.
    This property can be set. """,

```

(continues on next page)

(continued from previous page)

```

        validator=truncated_discrete_set,
        values=[10e-3, 100e-3, 1],
        map_values=True
    )

```

Now the actual GPIB/SCIP command is “:RANG 1” for a value of 100 mV, since the index of 100 mV in the values list is 1.

```

>>> extreme = Extreme5000("GPIB::1")
>>> extreme.voltage = 100e-3
>>> extreme.read()
'1'
>>> extreme.voltage = 1
>>> extreme.voltage
1

```

Dictionaries provide a more flexible method for mapping between real-values and those required by the instrument. If instead the :RANG <integer> took 1, 2, and 3 to correspond to 10 mV, 100 mV, and 1 V, then we can replace our previous control with the following.

```

Extreme5000.voltage = Instrument.control(
    ":RANG?", ":RANG %d",
    """ A floating point property that controls the voltage
    range in Volts, which takes values of 10 mV, 100 mV and 1 V.
    This property can be set. """ ,
    validator=truncated_discrete_set,
    values={10e-3:1, 100e-3:2, 1:3},
    map_values=True
)

```

```

>>> extreme = Extreme5000("GPIB::1")
>>> extreme.voltage = 10e-3
>>> extreme.read()
'1'
>>> extreme.voltage = 100e-3
>>> extreme.voltage
0.1

```

The dictionary now maps the keys to specific values. The values and keys can be any type, so this can support properties that use strings:

```

Extreme5000.channel = Instrument.control(
    ":CHAN?", ":CHAN %d",
    """ A string property that controls the measurement channel,
    which can take the values X, Y, or Z.
    """ ,
    validator=strict_discrete_set,
    values={'X':1, 'Y':2, 'Z':3},
    map_values=True
)

```

```
>>> extreme = Extreme5000("GPIB::1")
>>> extreme.channel = 'X'
>>> extreme.read()
'1'
>>> extreme.channel = 'Y'
>>> extreme.channel
'Y'
```

As you have seen, the `Instrument.control` function can be significantly extended by using validators and maps.

#### 10.4.4 Processing of set values

The `Instrument.control`, and `Instrument.setting` allow a keyword argument `set_process` which must be a function that takes a value after validation and performs processing before value mapping. This function must return the processed value. This can be typically used for unit conversions as in the following example:

```
Extreme5000.current = Instrument.setting(
    ":CURR %g",
    """ A floating point property that takes the measurement current in A
    """,
    validator=strict_range,
    values=[0, 10],
    set_process=lambda v: 1e3*v, # convert current to mA
)
```

```
>>> extreme = Extreme5000("GPIB::1")
>>> extreme.current = 1 # set current to 1000 mA
```

#### 10.4.5 Processing of return values

Similar to `set_process` the `Instrument.control`, and `Instrument.measurement` functions allow a `get_process` argument which if specified must be a function that takes a value and performs processing before value mapping. The function must return the processed value. In analogy to the example above this can be used for example for unit conversion:

```
Extreme5000.current = Instrument.control(
    ":CURR?", ":CURR %g",
    """ A floating point property representing the measurement current in A
    """,
    validator=strict_range,
    values=[0, 10],
    set_process=lambda v: 1e3*v, # convert to mA
    get_process=lambda v: 1e-3*v, # convert to A
)
```

```
>>> extreme = Extreme5000("GPIB::1")
>>> extreme.current = 3.1
>>> extreme.current
3.1
```

`get_process` can also be used to perform string processing. Let's say your instrument returns a value with its unit which has to be removed. This could be achieved by the following code:

```
Extreme5000.capacity = Instrument.measurement(  
    ":CAP?",  
    """ A measurement returning a capacity in nF in the format '<cap> nF'  
    """,  
    get_process=lambda v: float(v.replace('nF', ''))  
)
```

The same can be also achieved by the `preprocess_reply` keyword argument to `Instrument.control` or `Instrument.measurement`. This function is forwarded to `Adapter.values` and runs directly after receiving the reply from the device. One can therefore take advantage of the built in casting abilities and simplify the code accordingly:

```
Extreme5000.capacity = Instrument.measurement(  
    ":CAP?",  
    """ A measurement returning a capacity in nF in the format '<cap> nF'  
    """,  
    preprocess_reply=lambda v: v.replace('nF', '')  
    # notice how we don't need to cast to float anymore  
)
```

The real purpose of `preprocess_reply` is, however, for instruments where many/all properties need similar reply processing. `preprocess_reply` can be applied to all `Instrument.control` or `Instrument.measurement` properties, for example if all quantities are returned with a unit as in the example above. To avoid running into troubles for other properties this `preprocess_reply` should be clever enough to skip the processing in case it is not appropriate, for example if some identification string is returned. Typically this can be achieved by regular expression matching. In case of no match the reply is returned unchanged:

```
import re  
_reg_value = re.compile(r"([+-]?[0-9]*\.[0-9]*)\s+[a-zA-Z]+" )  
  
def extract_value(reply):  
    """ extract numerical value from reply. If none can be found the reply  
    is returned unchanged.  
  
    :param reply: reply string  
    :returns: string with only the numerical value  
    """  
    r = _reg_value.search(reply)  
    if r:  
        return r.groups()[0]  
    else:  
        return reply  
  
class Extreme5001(Instrument):  
    """ Represents the imaginary Extreme 5001 instrument. This instrument  
    sends numerical values including their units in an format "<value>  
    <unit>".  
    """  
    capacity = Instrument.measurement(  
        ":CAP?",  
        """ A measurement returning a capacity in nF in the format '<cap> nF'  
        """)
```

(continues on next page)

(continued from previous page)

```

)

voltage = Instrument.measurement(
    ":VOLT?",
    """ A measurement returning a voltage in V in the format '<volt> V'
    """
)

id = Instrument.measurement(
    "*idn?",
    """ The identification of the instrument.
    """
)

def __init__(self, resourceName, **kwargs):
    super().__init__(
        resourceName,
        "Extreme 5000",
        preprocess_reply=extract_value,
        **kwargs,
    )

```

In cases where the general `preprocess_reply` function should not run it can be also overwritten in the property definition:

```

Extreme5001.channel = Instrument.control(
    ":CHAN?", ":CHAN %d",
    """ A string property that controls the measurement channel,
    which can take the values X, Y, or Z.
    """,
    validator=strict_discrete_set,
    values=[1,2,3],
    preprocess_reply=lambda v: v,
)

```

Using a combination of the decribed abilities also complex communication schemes can be achieved.





## CODING STANDARDS

In order to maintain consistency across the different instruments in the PyMeasure repository, we enforce the following standards.

### 11.1 Python style guides

The [PEP8 style guide](#) and [PEP257 docstring conventions](#) should be followed.

Function and variable names should be lower case with underscores as needed to separate words. CamelCase should only be used for class names, unless working with Qt, where its use is common.

There are no plans to support type hinting in PyMeasure code. This adds a lot of additional code to manage, without a clear advantage for this project. Type documentation should be placed in the docstring where not clear from the variable name.

### 11.2 Documentation

PyMeasure documents code using reStructuredText and the [Sphinx documentation generator](#). All functions, classes, and methods should be documented in the code using a [docstring](#).

### 11.3 Usage of getter and setter functions

Getter and setter functions are discouraged, since properties provide a more fluid experience. Given the extensive tools available for defining properties, detailed in the [Advanced properties](#) section, these types of properties are preferred.



## AUTHORS

PyMeasure was started in 2013 by Colin Jermain and Graham Rowlands at Cornell University, when it became apparent that both were working on similar Python packages for scientific measurements. PyMeasure combined these efforts and continues to gain valuable contributions from other scientists who are interested in advancing measurement software.

The following developers have contributed to the PyMeasure package:

Colin Jermain  
Graham Rowlands  
Minh-Hai Nguyen  
Guen Prawiro-Atmodjo  
Tim van Bortel  
Davide Spirito  
Marcos Guimaraes  
Ghislain Antony Vaillant  
Ben Feinstein  
Neal Reynolds  
Christoph Buchner  
Julian Dlugosch  
Sylvain Karlen  
Joseph Mittelstaedt  
Troy Fox  
Vikram Sekar  
Casper Schippers  
Sumatran Tiger  
Michael Schneider  
Dennis Feng  
Stefano Pirotta  
Moritz Jung  
Richard Schlitz  
Manuel Zahn  
Mikhaël Myara  
Domenic Prete  
Mathieu Jeannin  
Paul Goulain  
John McMaster  
Dominik Kriegner

Jonathan Larochelle  
Dominic Caron  
Mathieu Plante  
Michele Sardo  
Steven Siegl  
Benjamin Klebel-Knobloch  
Demetra Adrahtas  
Dan McDonald  
Hud Wahab  
Nicola Corna  
Robert Eckelmann  
Sam Condon  
Andreas Maeder  
Bastian Leykauf

**LICENSE**

Copyright (c) 2013-2021 PyMeasure Developers

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the “Software”), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED “AS IS”, WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.



## PYTHON MODULE INDEX

### p

- `pymeasure.display.browser`, 49
- `pymeasure.display.curves`, 49
- `pymeasure.display.inputs`, 50
- `pymeasure.display.listeners`, 52
- `pymeasure.display.log`, 52
- `pymeasure.display.manager`, 52
- `pymeasure.display.plotter`, 53
- `pymeasure.display.thread`, 54
- `pymeasure.display.widgets`, 54
- `pymeasure.display.windows`, 57
- `pymeasure.experiment.experiment`, 41
- `pymeasure.experiment.listeners`, 42
- `pymeasure.experiment.parameters`, 44
- `pymeasure.experiment.procedure`, 43
- `pymeasure.experiment.results`, 47
- `pymeasure.experiment.workers`, 46
- `pymeasure.instruments`, 59
- `pymeasure.instruments.advantest`, 65
- `pymeasure.instruments.advantest.advantestR3767CG`, 66
- `pymeasure.instruments.agilent`, 66
- `pymeasure.instruments.agilent.agilent4156`, 76
- `pymeasure.instruments.agilent.agilentB1500`, 101
- `pymeasure.instruments.ametek`, 103
- `pymeasure.instruments.ami`, 105
- `pymeasure.instruments.anaheimautomation`, 107
- `pymeasure.instruments.anapico`, 109
- `pymeasure.instruments.anritsu`, 109
- `pymeasure.instruments.attocube`, 112
- `pymeasure.instruments.bkprecision`, 115
- `pymeasure.instruments.comedi`, 65
- `pymeasure.instruments.danfysik`, 116
- `pymeasure.instruments.deltalelektronika`, 119
- `pymeasure.instruments.edwards`, 120
- `pymeasure.instruments.fluke`, 120
- `pymeasure.instruments.fwbell`, 120
- `pymeasure.instruments.hp`, 122
- `pymeasure.instruments.keithley`, 126
- `pymeasure.instruments.keysight`, 158
- `pymeasure.instruments.lakeshore`, 164
- `pymeasure.instruments.newport`, 167
- `pymeasure.instruments.ni`, 168
- `pymeasure.instruments.oxfordinstruments`, 180
- `pymeasure.instruments.parker`, 182
- `pymeasure.instruments.pendulum`, 184
- `pymeasure.instruments.razorbill`, 185
- `pymeasure.instruments.rohdeschwarz`, 186
- `pymeasure.instruments.signalrecovery`, 200
- `pymeasure.instruments.srs`, 203
- `pymeasure.instruments.tektronix`, 213
- `pymeasure.instruments.thorlabs`, 213
- `pymeasure.instruments.validators`, 63
- `pymeasure.instruments.yokogawa`, 214





## Symbols

- `__call__()` (*pymeasure.instruments.agilent.agilentB1500.Range* method), 100
- `__str__()` (*pymeasure.instruments.agilent.agilentB1500.CustomIntEnum* method), 101
- `_format_binary_values()` (*pymeasure.adapters.PrologixAdapter* method), 35
- `_format_binary_values()` (*pymeasure.adapters.SerialAdapter* method), 33
- `_format_binary_values()` (*pymeasure.instruments.lakeshore.LakeShoreUSBAdapter* method), 164
- A**
- `abort()` (*pymeasure.display.manager.Manager* method), 53
- `abort()` (*pymeasure.instruments.agilent.agilentB1500.AgilentB1500* method), 94
- `absolute_position` (*pymeasure.instruments.anaheimautomation.DPSeriesMotorController* property), 107
- `absolute_to_steps()` (*pymeasure.instruments.anaheimautomation.DPSeriesMotorController* method), 107
- `ac_current` (*pymeasure.instruments.agilent.AgilentE4980* property), 71
- `ac_mode()` (*pymeasure.instruments.lakeshore.LakeShore425* method), 166
- `ac_voltage` (*pymeasure.instruments.agilent.AgilentE4980* property), 71
- `acquire_digital_input_output()` (*pymeasure.instruments.ni.virtualbench.VirtualBench* method), 179
- `acquire_digital_multimeter()` (*pymeasure.instruments.ni.virtualbench.VirtualBench* method), 179
- `acquire_function_generator()` (*pymeasure.instruments.ni.virtualbench.VirtualBench* method), 179
- `acquire_mixed_signal_oscilloscope()` (*pymeasure.instruments.ni.virtualbench.VirtualBench* method), 179
- `acquire_power_supply()` (*pymeasure.instruments.ni.virtualbench.VirtualBench* method), 179
- `acquire_reference()` (*pymeasure.instruments.keithley.Keithley2000* method), 126
- `acquisition_mode` (*pymeasure.instruments.keysight.KeysightDSOX1102G* property), 159
- `acquisition_type` (*pymeasure.instruments.keysight.KeysightDSOX1102G* property), 159
- `active_connectors` (*pymeasure.instruments.hp.HP3478A* property), 123
- `Adapter` (class in *pymeasure.adapters*), 31
- `adc1` (*pymeasure.instruments.ametek.Ametek7270* property), 104
- `adc1` (*pymeasure.instruments.signalrecovery.DSP7265* property), 200
- `adc1` (*pymeasure.instruments.srs.SR830* property), 205
- `adc1` (*pymeasure.instruments.srs.SR860* property), 208
- `adc2` (*pymeasure.instruments.ametek.Ametek7270* property), 104
- `adc2` (*pymeasure.instruments.signalrecovery.DSP7265* property), 200
- `adc2` (*pymeasure.instruments.srs.SR830* property), 205
- `adc2` (*pymeasure.instruments.srs.SR860* property), 208
- `adc3` (*pymeasure.instruments.ametek.Ametek7270* property), 104
- `adc3` (*pymeasure.instruments.srs.SR830* property), 205
- `adc3` (*pymeasure.instruments.srs.SR860* property), 208
- `adc4` (*pymeasure.instruments.ametek.Ametek7270* property), 104
- `adc4` (*pymeasure.instruments.srs.SR830* property), 205
- `adc4` (*pymeasure.instruments.srs.SR860* property), 208
- `adc_auto_zero` (*pymeasure.instruments.agilent.agilentB1500.AgilentB1500* property), 95
- `adc_averaging()` (*pymeasure.instruments.agilent.agilentB1500.AgilentB1500* property), 95

- method*), 95
- `adc_setup()` (`pymeasure.instruments.agilent.agilentB1500.AgilentB1500` *property*), 66
- `adc_type` (`pymeasure.instruments.agilent.agilentB1500.SMA100` *property*), 97
- `ADCMode` (class in `pymeasure.instruments.agilent.agilentB1500`), 101
- `ADCType` (class in `pymeasure.instruments.agilent.agilentB1500`), 101
- `add()` (`pymeasure.display.browser.Browser` *method*), 49
- `add_ramp_step()` (`pymeasure.instruments.danfysik.Danfysik8500` *method*), 116
- `address` (`pymeasure.instruments.anaheimautomation.DPSeriesMotorControl` *property*), 107
- `AdvantestR3767CG` (class in `pymeasure.instruments.advantest.advantestR3767CG`), 66
- `AFG3152C` (class in `pymeasure.instruments.tektronix`), 213
- `Agilent33220A` (class in `pymeasure.instruments.agilent`), 83
- `Agilent33500` (class in `pymeasure.instruments.agilent`), 85
- `Agilent33521A` (class in `pymeasure.instruments.agilent`), 88
- `Agilent34410A` (class in `pymeasure.instruments.agilent`), 72
- `Agilent34450A` (class in `pymeasure.instruments.agilent`), 73
- `Agilent4156` (class in `pymeasure.instruments.agilent.agilent4156`), 76
- `Agilent8257D` (class in `pymeasure.instruments.agilent`), 66
- `Agilent8722ES` (class in `pymeasure.instruments.agilent`), 69
- `AgilentB1500` (class in `pymeasure.instruments.agilent.agilentB1500`), 93
- `AgilentE4408B` (class in `pymeasure.instruments.agilent`), 70
- `AgilentE4980` (class in `pymeasure.instruments.agilent`), 71
- `Ametek7270` (class in `pymeasure.instruments.ametek`), 104
- `AMI430` (class in `pymeasure.instruments.ami`), 105
- `amplitude` (`pymeasure.instruments.agilent.Agilent33220A` *property*), 83
- `amplitude` (`pymeasure.instruments.agilent.Agilent33500` *property*), 85
- `amplitude` (`pymeasure.instruments.hp.HP33120A` *property*), 122
- `amplitude_depth` (`pymeasure.instruments.agilent.Agilent8257D` *property*), 66
- `amplitude_source` (`pymeasure.instruments.agilent.Agilent8257D` *property*), 66
- `amplitude_unit` (`pymeasure.instruments.agilent.Agilent33220A` *property*), 83
- `amplitude_unit` (`pymeasure.instruments.agilent.Agilent33500` *property*), 85
- `amplitude_units` (`pymeasure.instruments.hp.HP33120A` *property*), 122
- `analysis` (`pymeasure.instruments.anritsu.AnritsuMS9710C` *property*), 110
- `analysis_result` (`pymeasure.instruments.anritsu.AnritsuMS9710C` *property*), 110
- `analyzer_mode` (`pymeasure.instruments.agilent.agilent4156.Agilent4156` *property*), 78
- `ANC300Controller` (class in `pymeasure.instruments.attocube.anc300`), 113
- `angle` (`pymeasure.instruments.parker.ParkerGV6` *property*), 183
- `angle_error` (`pymeasure.instruments.parker.ParkerGV6` *property*), 183
- `AnritsuMG3692C` (class in `pymeasure.instruments.anritsu`), 110
- `AnritsuMS9710C` (class in `pymeasure.instruments.anritsu`), 110
- `aperture()` (`pymeasure.instruments.agilent.AgilentE4980` *method*), 71
- `append()` (`pymeasure.display.curves.BufferCurve` *method*), 49
- `applied` (`pymeasure.instruments.keithley.Keithley2260B` *property*), 133
- `apply_current()` (`pymeasure.instruments.keithley.Keithley2400` *method*), 134
- `apply_current()` (`pymeasure.instruments.keithley.Keithley2450` *method*), 140
- `apply_current()` (`pymeasure.instruments.yokogawa.Yokogawa7651` *method*), 215
- `apply_voltage()` (`pymeasure.instruments.keithley.Keithley2400` *method*), 134
- `apply_voltage()` (`pymeasure.instruments.keithley.Keithley2450` *method*), 141
- `apply_voltage()` (`pymeasure.instruments.keithley.Keithley6517B` *method*), 154

apply\_voltage() (pymeasure.instruments.yokogawa.Yokogawa7651 method), 215  
 APSIN12G (class in pymeasure.instruments.anapico), 109  
 arb\_advance (pymeasure.instruments.agilent.Agilent33500 property), 85  
 arb\_file (pymeasure.instruments.agilent.Agilent33500 property), 85  
 arb\_filter (pymeasure.instruments.agilent.Agilent33500 property), 85  
 arb\_srate (pymeasure.instruments.agilent.Agilent33500 property), 85  
 arb\_srate (pymeasure.instruments.agilent.Agilent33521A property), 88  
 ask() (pymeasure.adapters.Adapter method), 31  
 ask() (pymeasure.adapters.FakeAdapter method), 32  
 ask() (pymeasure.adapters.PrologixAdapter method), 35  
 ask() (pymeasure.adapters.SerialAdapter method), 33  
 ask() (pymeasure.adapters.TelnetAdapter method), 39  
 ask() (pymeasure.adapters.VISAAdapter method), 37  
 ask() (pymeasure.adapters.VXIIAdapter method), 38  
 ask() (pymeasure.instruments.agilent.agilentB1500.SMU method), 97  
 ask() (pymeasure.instruments.anaheimautomation.DPSeriesMotorController method), 107  
 ask() (pymeasure.instruments.attocube.adapters.AttocubeConsoleAdapter method), 112  
 ask() (pymeasure.instruments.fwbell.FWBell5080 method), 121  
 ask() (pymeasure.instruments.Instrument method), 61  
 ask() (pymeasure.instruments.keysight.KeysightDSOX1102G method), 159  
 ask() (pymeasure.instruments.lakeshore.LakeShoreUSBAdapter method), 164  
 ask\_raw() (pymeasure.adapters.VXIIAdapter method), 38  
 ask\_values() (pymeasure.adapters.VISAAdapter method), 37  
 attenuation (pymeasure.instruments.rohdeschwarz.fsl.FSL property), 199  
 AttocubeConsoleAdapter (class in pymeasure.instruments.attocube.adapters), 112  
 AUTO (pymeasure.instruments.agilent.agilentB1500.ADCMode attribute), 101  
 AUTO (pymeasure.instruments.agilent.agilentB1500.AutoManual attribute), 102  
 AUTO (pymeasure.instruments.agilent.agilentB1500.CompliancePolarity attribute), 103  
 auto\_calibration (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 property), 94  
 auto\_offset() (pymeasure.instruments.srs.SR830 method), 205  
 auto\_output\_off (pymeasure.instruments.keithley.Keithley2400 property), 134  
 auto\_pid (pymeasure.instruments.oxfordinstruments.ITC503 property), 181  
 auto\_range() (pymeasure.instruments.fwbell.FWBell5080 method), 121  
 auto\_range() (pymeasure.instruments.keithley.Keithley2000 method), 127  
 auto\_range() (pymeasure.instruments.lakeshore.LakeShore425 method), 166  
 auto\_range\_enabled (pymeasure.instruments.hp.HP3478A property), 124  
 auto\_range\_source() (pymeasure.instruments.keithley.Keithley2400 method), 134  
 auto\_range\_source() (pymeasure.instruments.keithley.Keithley2450 method), 141  
 auto\_range\_source() (pymeasure.instruments.keithley.Keithley6517B method), 154  
 auto\_setup() (pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope method), 174  
 auto\_zero (pymeasure.instruments.keithley.Keithley2400 property), 134  
 auto\_zero\_enabled (pymeasure.instruments.hp.HP3478A property), 124  
 AutoManual (class in pymeasure.instruments.agilent.agilentB1500), 102  
 autoscale() (pymeasure.instruments.keysight.KeysightDSOX1102G method), 159  
 aux\_in\_1 (pymeasure.instruments.srs.SR830 property), 205  
 aux\_in\_1 (pymeasure.instruments.srs.SR860 property), 208  
 aux\_in\_2 (pymeasure.instruments.srs.SR830 property), 206  
 aux\_in\_2 (pymeasure.instruments.srs.SR860 property), 206  
 aux\_in\_3 (pymeasure.instruments.srs.SR830 property), 206  
 aux\_in\_3 (pymeasure.instruments.srs.SR860 property), 208  
 aux\_in\_4 (pymeasure.instruments.srs.SR830 property), 206  
 aux\_in\_4 (pymeasure.instruments.srs.SR860 property), 208  
 aux\_out\_1 (pymeasure.instruments.srs.SR830 property), 208

- 206  
[aux\\_out\\_1](#) (*pymeasure.instruments.srs.SR860* property), 209  
[aux\\_out\\_2](#) (*pymeasure.instruments.srs.SR830* property), 206  
[aux\\_out\\_2](#) (*pymeasure.instruments.srs.SR860* property), 209  
[aux\\_out\\_3](#) (*pymeasure.instruments.srs.SR830* property), 206  
[aux\\_out\\_3](#) (*pymeasure.instruments.srs.SR860* property), 209  
[aux\\_out\\_4](#) (*pymeasure.instruments.srs.SR830* property), 206  
[aux\\_out\\_4](#) (*pymeasure.instruments.srs.SR860* property), 209  
[average\\_point](#) (*pymeasure.instruments.anritsu.AnritsuMS9710C* property), 110  
[average\\_sweep](#) (*pymeasure.instruments.anritsu.AnritsuMS9710C* property), 110  
[averages](#) (*pymeasure.instruments.agilent.Agilent8722ES* property), 69  
[averaging\\_enabled](#) (*pymeasure.instruments.agilent.Agilent8722ES* property), 69  
[axes](#) (*pymeasure.instruments.newport.ESP300* property), 167  
[Axis](#) (class in *pymeasure.instruments.attocube.anc300*), 114  
[Axis](#) (class in *pymeasure.instruments.newport.esp300*), 167  
[AxisError](#) (class in *pymeasure.instruments.newport.esp300*), 168
- ## B
- [BASE](#) (*pymeasure.instruments.agilent.agilentB1500.SamplingPostOutput* attribute), 103  
[basespeed](#) (*pymeasure.instruments.anaheimautomation.DPSeriesMotorController* property), 107  
[basic\\_info](#) (*pymeasure.instruments.rohdeschwarz.sfm.SFM9130B* property), 187  
[batch\\_size](#) (*pymeasure.instruments.pendulum.cnr91.CNT91* property), 184  
[beep\(\)](#) (*pymeasure.instruments.agilent.Agilent33220A* method), 83  
[beep\(\)](#) (*pymeasure.instruments.agilent.Agilent33500* method), 85  
[beep\(\)](#) (*pymeasure.instruments.agilent.Agilent34450A* method), 73  
[beep\(\)](#) (*pymeasure.instruments.hp.HP33120A* method), 122  
[beep\(\)](#) (*pymeasure.instruments.keithley.Keithley2000* method), 127  
[beep\(\)](#) (*pymeasure.instruments.keithley.Keithley2400* method), 134  
[beep\(\)](#) (*pymeasure.instruments.keithley.Keithley2450* method), 141  
[beep\(\)](#) (*pymeasure.instruments.keithley.Keithley2700* method), 146  
[beep\(\)](#) (*pymeasure.instruments.keithley.Keithley6221* method), 149  
[beep\\_state](#) (*pymeasure.instruments.keithley.Keithley2000* property), 127  
[beeper\\_enabled](#) (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 187  
[beeper\\_state](#) (*pymeasure.instruments.agilent.Agilent33220A* property), 83  
[BIAS](#) (*pymeasure.instruments.agilent.agilentB1500.SamplingPostOutput* attribute), 103  
[bias\\_enabled](#) (*pymeasure.instruments.srs.SR570* property), 204  
[bias\\_level](#) (*pymeasure.instruments.srs.SR570* property), 204  
[binary\\_values\(\)](#) (*pymeasure.adapters.Adapter* method), 31  
[binary\\_values\(\)](#) (*pymeasure.adapters.FakeAdapter* method), 32  
[binary\\_values\(\)](#) (*pymeasure.adapters.PrologixAdapter* method), 35  
[binary\\_values\(\)](#) (*pymeasure.adapters.SerialAdapter* method), 33  
[binary\\_values\(\)](#) (*pymeasure.adapters.TelnetAdapter* method), 40  
[binary\\_values\(\)](#) (*pymeasure.adapters.VISAAdapter* method), 37  
[binary\\_values\(\)](#) (*pymeasure.adapters.VXI11Adapter* method), 38  
[binary\\_values\(\)](#) (*pymeasure.instruments.lakeshore.LakeShoreUSBAdapter* method), 164  
[BKPrecision9130B](#) (class in *pymeasure.instruments.bkprecision*), 115  
[blank\\_front\(\)](#) (*pymeasure.instruments.srs.SR570* method), 204  
[blanking](#) (*pymeasure.instruments.anapico.APSIN12G* property), 109  
[BooleanInput](#) (class in *pymeasure.display.inputs*), 50  
[BooleanParameter](#) (class in *pymeasure.experiment.parameters*), 44  
[Browser](#) (class in *pymeasure.display.browser*), 49  
[BrowserItem](#) (class in *pymeasure.display.browser*), 49  
[BrowserWidget](#) (class in *pymeasure.display.widgets*), 54  
[buffer\\_data](#) (*pymeasure.instruments.keithley.Keithley2000* property), 127



[buffer\\_data \(pymeasure.instruments.keithley.Keithley2400C property\), 135](#)  
[buffer\\_data \(pymeasure.instruments.keithley.Keithley2450 property\), 141](#)  
[buffer\\_data \(pymeasure.instruments.keithley.Keithley2700 property\), 146](#)  
[buffer\\_data \(pymeasure.instruments.keithley.Keithley6221 property\), 149](#)  
[buffer\\_data \(pymeasure.instruments.keithley.Keithley6517B property\), 154](#)  
[buffer\\_frequency\\_time\\_series\(\) \(pymeasure.instruments.pendulum.cnt91.CNT91 method\), 184](#)  
[buffer\\_points \(pymeasure.instruments.keithley.Keithley2000 property\), 127](#)  
[buffer\\_points \(pymeasure.instruments.keithley.Keithley2400 property\), 135](#)  
[buffer\\_points \(pymeasure.instruments.keithley.Keithley2450 property\), 141](#)  
[buffer\\_points \(pymeasure.instruments.keithley.Keithley2700 property\), 146](#)  
[buffer\\_points \(pymeasure.instruments.keithley.Keithley6221 property\), 149](#)  
[buffer\\_points \(pymeasure.instruments.keithley.Keithley6517B property\), 154](#)  
[buffer\\_to\\_float\(\) \(pymeasure.instruments.signalrecovery.DSP7265 method\), 200](#)  
[BufferCurve \(class in pymeasure.display.curves\), 49](#)  
[burst\\_mode \(pymeasure.instruments.agilent.Agilent33220A property\), 83](#)  
[burst\\_mode \(pymeasure.instruments.agilent.Agilent33500 property\), 86](#)  
[burst\\_ncycles \(pymeasure.instruments.agilent.Agilent33220A property\), 83](#)  
[burst\\_ncycles \(pymeasure.instruments.agilent.Agilent33500 property\), 86](#)  
[burst\\_period \(pymeasure.instruments.agilent.Agilent33500 property\), 86](#)  
[burst\\_state \(pymeasure.instruments.agilent.Agilent33220A property\), 83](#)  
[burst\\_state \(pymeasure.instruments.agilent.Agilent33500 property\), 86](#)  
[busy \(pymeasure.instruments.anaheimautomation.DPSeriesMotorController property\), 107](#)  
[calibration\(\) \(pymeasure.instruments.rohdeschwarz.sfm.SFM method\), 188](#)  
[calibration\\_enabled \(pymeasure.instruments.hp.HP3478A property\), 124](#)  
[capacitance \(pymeasure.instruments.agilent.Agilent34450A property\), 73](#)  
[capacitance\\_auto\\_range \(pymeasure.instruments.agilent.Agilent34450A property\), 73](#)  
[capacitance\\_range \(pymeasure.instruments.agilent.Agilent34450A property\), 73](#)  
[capacity \(pymeasure.instruments.attocube.anc300.Axis property\), 114](#)  
[carrier\\_enabled \(pymeasure.instruments.rohdeschwarz.sfm.Sound\\_Channel property\), 196](#)  
[carrier\\_frequency \(pymeasure.instruments.rohdeschwarz.sfm.Sound\\_Channel property\), 196](#)  
[carrier\\_level \(pymeasure.instruments.rohdeschwarz.sfm.Sound\\_Channel property\), 196](#)  
[center\\_at\\_peak\(\) \(pymeasure.instruments.anritsu.AnritsuMS9710C method\), 110](#)  
[center\\_frequency \(pymeasure.instruments.advantest.advantestR3767CG.AdvantestR3767CG property\), 66](#)  
[center\\_frequency \(pymeasure.instruments.agilent.Agilent8257D property\), 67](#)  
[center\\_frequency \(pymeasure.instruments.agilent.AgilentE4408B property\), 70](#)  
[channel \(pymeasure.instruments.bkprecision.BKPrecision9130B property\), 115](#)  
[channel1 \(pymeasure.instruments.srs.SR830 property\), 206](#)  
[channel2 \(pymeasure.instruments.srs.SR830 property\), 206](#)  
[channel\\_down\\_relative\(\) \(pymeasure.instruments.rohdeschwarz.sfm.SFM method\), 188](#)  
[channel\\_function \(pymeasure.instruments.agilent.agilent4156.SMU property\), 79](#)  
[channel\\_function \(pymeasure.instruments.agilent.agilent4156.VSU property\), 82](#)  
[channel\\_mode \(pymeasure.instruments.agilent.agilent4156.VSU property\), 82](#)

- `sure.instruments.agilent.agilent4156.SMU` (pymea-  
property), 79
- `channel_mode` (pymea-  
`sure.instruments.agilent.agilent4156.VMU`  
property), 82
- `channel_mode` (pymea-  
`sure.instruments.agilent.agilent4156.VSU`  
property), 82
- `channel_sweep_start` (pymea-  
`sure.instruments.rohdeschwarz.sfm.SFM`  
property), 188
- `channel_sweep_step` (pymea-  
`sure.instruments.rohdeschwarz.sfm.SFM`  
property), 188
- `channel_sweep_stop` (pymea-  
`sure.instruments.rohdeschwarz.sfm.SFM`  
property), 188
- `channel_table` (pymea-  
`sure.instruments.rohdeschwarz.sfm.SFM`  
property), 188
- `channel_up_relative()` (pymea-  
`sure.instruments.rohdeschwarz.sfm.SFM`  
method), 188
- `channels_from_rows_columns()` (pymea-  
`sure.instruments.keithley.Keithley2700`  
method), 146
- `check_acknowledgement()` (pymea-  
`sure.instruments.attocube.adapters.AttocubeConsolidator`  
method), 113
- `check_errors()` (pymea-  
`sure.instruments.agilent.agilentB1500.AgilentB1500`  
method), 94
- `check_errors()` (pymea-  
`sure.instruments.agilent.agilentB1500.SMU`  
method), 97
- `check_errors()` (pymea-  
`sure.instruments.anaheimautomation.DPSeriesModule`  
method), 107
- `check_errors()` (pymea-  
`sure.instruments.bkprecision.BKPrecision9130B`  
method), 115
- `check_errors()` (pymea-  
`sure.instruments.hp.HP3478A`  
method), 124
- `check_errors()` (pymea-  
`sure.instruments.Instrument`  
method), 61
- `check_errors()` (pymea-  
`sure.instruments.keithley.Keithley2000`  
method), 127
- `check_errors()` (pymea-  
`sure.instruments.keithley.Keithley2260B`  
method), 133
- `check_errors()` (pymea-  
`sure.instruments.keithley.Keithley2400`  
method), 135
- `check_errors()` (pymea-  
`sure.instruments.keithley.Keithley2450`  
method), 141
- `check_errors()` (pymea-  
`sure.instruments.keithley.Keithley2600`  
method), 158
- `check_errors()` (pymea-  
`sure.instruments.keithley.Keithley2700`  
method), 146
- `check_errors()` (pymea-  
`sure.instruments.keithley.Keithley2750`  
method), 157
- `check_errors()` (pymea-  
`sure.instruments.keithley.Keithley6221`  
method), 149
- `check_errors()` (pymea-  
`sure.instruments.keithley.Keithley6517B`  
method), 154
- `check_errors()` (pymea-  
`sure.instruments.keysight.KeysightDSOX1102G`  
method), 159
- `check_errors()` (pymea-  
`sure.instruments.keysight.KeysightN5767A`  
method), 163
- `check_get_estimates_signature()` (pymea-  
`sure.display.widgets.EstimatorWidget` method),  
54
- `check_idle()` (pymea-  
`sure.instruments.agilent.agilentB1500.AgilentB1500`  
method), 94
- `check_parameters()` (pymea-  
`sure.experiment.procedure.Procedure` method),  
43
- `check_stop()` (pymea-  
`sure.display.windows.PlotterWindow` method),  
59
- `check_timeout()` (pymea-  
`sure.experiment.parameters.ListParameter`  
property), 45
- `clear()` (pymea-  
`sure.display.manager.Manager` method),  
53
- `clear()` (pymea-  
`sure.instruments.Instrument` method), 61
- `clear()` (pymea-  
`sure.instruments.keysight.KeysightDSOX1102G`  
method), 159
- `clear_buffer()` (pymea-  
`sure.instruments.agilent.agilentB1500.AgilentB1500`  
method), 94
- `clear_display()` (pymea-  
`sure.instruments.agilent.Agilent33500` method),  
86
- `clear_errors()` (pymea-  
`sure.instruments.newport.ESP300` method),  
167
- `clear_overload()` (pymea-  
`sure.instruments.srs.SR570` method), 204

- `clear_plot()` (*pymeasure.experiment.experiment.Experiment* method), 41
- `clear_ramp_set()` (*pymeasure.instruments.danfysik.Danfysik8500* method), 117
- `clear_sequence()` (*pymeasure.instruments.danfysik.Danfysik8500* method), 117
- `clear_status()` (*pymeasure.instruments.keysight.KeysightDSOX1102G* method), 159
- `clear_timer()` (*pymeasure.instruments.agilent.agilentB1500.AgilentB1500* method), 94
- `close()` (*pymeasure.instruments.keithley.Keithley2750* method), 157
- `close_rows_to_columns()` (*pymeasure.instruments.keithley.Keithley2700* method), 146
- `closed_channels` (*pymeasure.instruments.keithley.Keithley2700* property), 147
- `closed_channels` (*pymeasure.instruments.keithley.Keithley2750* property), 157
- `CMU_MEASUREMENT` (*pymeasure.instruments.agilent.agilentB1500.WaitTimeType* attribute), 103
- `CNT91` (class in *pymeasure.instruments.pendulum.cnt91*), 184
- `coder_adjust()` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* method), 188
- `coder_id_frequency` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 188
- `coder_modulation_degree` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 188
- `coder_pilot_deviation` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 189
- `coder_pilot_frequency` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 189
- `coilconst` (*pymeasure.instruments.ami.AMI430* property), 106
- `collapse_channel_string()` (*pymeasure.instruments.ni.virtualbench.VirtualBench* method), 179
- `complete` (*pymeasure.instruments.bkprecision.BKPrecision9130B* property), 115
- `complete` (*pymeasure.instruments.Instrument* property), 61
- `complete` (*pymeasure.instruments.keithley.Keithley2000* property), 127
- `complete` (*pymeasure.instruments.keithley.Keithley2260B* property), 133
- `complete` (*pymeasure.instruments.keithley.Keithley2400* property), 135
- `complete` (*pymeasure.instruments.keithley.Keithley2450* property), 141
- `complete` (*pymeasure.instruments.keithley.Keithley2600* property), 158
- `complete` (*pymeasure.instruments.keithley.Keithley2700* property), 147
- `complete` (*pymeasure.instruments.keithley.Keithley2750* property), 157
- `complete` (*pymeasure.instruments.keithley.Keithley6221* property), 149
- `complete` (*pymeasure.instruments.keithley.Keithley6517B* property), 154
- `complete` (*pymeasure.instruments.keysight.KeysightDSOX1102G* property), 159
- `complete` (*pymeasure.instruments.keysight.KeysightN5767A* property), 163
- `compliance` (*pymeasure.instruments.agilent.agilent4156.SMU* property), 80
- `compliance` (*pymeasure.instruments.agilent.agilent4156.VARD* property), 81
- `compliance` (*pymeasure.instruments.agilent.agilent4156.VARX* property), 81
- `COMPLIANCE_AND_FORCE_SIDE` (*pymeasure.instruments.agilent.agilentB1500.MeasOpMode* attribute), 102
- `compliance_current` (*pymeasure.instruments.keithley.Keithley2400* property), 135
- `compliance_current` (*pymeasure.instruments.keithley.Keithley2450* property), 141
- `compliance_current` (*pymeasure.instruments.yokogawa.Yokogawa7651* property), 215
- `COMPLIANCE_SIDE` (*pymeasure.instruments.agilent.agilentB1500.MeasOpMode* attribute), 102
- `compliance_voltage` (*pymeasure.instruments.keithley.Keithley2400* property), 135
- `compliance_voltage` (*pymeasure.instruments.keithley.Keithley2450* property), 141
- `compliance_voltage` (*pymeasure.instruments.yokogawa.Yokogawa7651* property), 215
- `CompliancePolarity` (class in *pymeasure*

- `sure.instruments.agilent.agilentB1500`), 103
- `config_amplitude_modulation()` (`pymea-  
sure.instruments.agilent.Agilent8257D`  
method), 67
- `config_buffer()` (`pymea-  
sure.instruments.keithley.Keithley2000`  
method), 127
- `config_buffer()` (`pymea-  
sure.instruments.keithley.Keithley2400`  
method), 135
- `config_buffer()` (`pymea-  
sure.instruments.keithley.Keithley2450`  
method), 141
- `config_buffer()` (`pymea-  
sure.instruments.keithley.Keithley2700`  
method), 147
- `config_buffer()` (`pymea-  
sure.instruments.keithley.Keithley6221`  
method), 149
- `config_buffer()` (`pymea-  
sure.instruments.keithley.Keithley6517B`  
method), 154
- `config_low_freq_out()` (`pymea-  
sure.instruments.agilent.Agilent8257D`  
method), 67
- `config_pulse_modulation()` (`pymea-  
sure.instruments.agilent.Agilent8257D`  
method), 67
- `config_step_sweep()` (`pymea-  
sure.instruments.agilent.Agilent8257D`  
method), 67
- `configure()` (`pymea-  
sure.instruments.agilent.agilent4156`  
method), 78
- `configure_ac_current()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter`  
method), 171
- `configure_analog_channel()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope`  
method), 174
- `configure_analog_channel_characteristics()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope`  
method), 175
- `configure_analog_edge_trigger()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope`  
method), 175
- `configure_analog_pulse_width_trigger()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope`  
method), 175
- `configure_arbitrary_waveform()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.FunctionGenerator`  
method), 172
- `configure_arbitrary_waveform_gain_and_offset()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.FunctionGenerator`  
method), 172
- `configure_capacitance()` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
method), 73
- `configure_continuity()` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
method), 73
- `configure_current()` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
method), 73
- `configure_current_output()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.PowerSupply`  
method), 178
- `configure_dc_current()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter`  
method), 171
- `configure_dc_voltage()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter`  
method), 171
- `configure_diode()` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
method), 73
- `configure_frequency()` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
method), 73
- `configure_frequency_array_measurement()` (`pymea-  
sure.instruments.pendulum.cnt91.CNT91`  
method), 184
- `configure_immediate_trigger()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope`  
method), 176
- `configure_measurement()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter`  
method), 171
- `configure_resistance()` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
method), 74
- `configure_standard_waveform()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.FunctionGenerator`  
method), 173
- `configure_temperature()` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
method), 74
- `configure_timing()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope`  
method), 176
- `configure_trigger_delay()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope`  
method), 176
- `configure_voltage()` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
method), 74
- `configure_voltage_output()` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.PowerSupply`  
method), 178



- `constant_value` (`pymeasure.instruments.agilent.agilent4156.SMU` property), 80
- `constant_value` (`pymeasure.instruments.agilent.agilent4156.VSU` property), 82
- `contact_current_1` (`pymeasure.instruments.razorbill.razorbillRP100` property), 185
- `contact_current_2` (`pymeasure.instruments.razorbill.razorbillRP100` property), 185
- `contact_voltage_1` (`pymeasure.instruments.razorbill.razorbillRP100` property), 185
- `contact_voltage_2` (`pymeasure.instruments.razorbill.razorbillRP100` property), 185
- `continue_single_sweep()` (`pymeasure.instruments.rohdeschwarz.fsl.FSL` method), 199
- `continuity` (`pymeasure.instruments.agilent.Agilent34450A` property), 74
- `continuous` (`pymeasure.instruments.pendulum.cnt91.CNT91` property), 184
- `continuous_sweep` (`pymeasure.instruments.rohdeschwarz.fsl.FSL` property), 199
- `control()` (`pymeasure.instruments.Instrument` static method), 61
- `control()` (`pymeasure.instruments.keysight.KeysightDSOX1102G` static method), 159
- `control_mode` (`pymeasure.instruments.oxfordinstruments.ITC503` property), 181
- `controllerBoardVersion` (`pymeasure.instruments.attocube.anc300.ANC300Controller` property), 113
- `convert_timestamp_to_values()` (`pymeasure.instruments.ni.virtualbench.VirtualBench` method), 179
- `convert_values_to_datetime()` (`pymeasure.instruments.ni.virtualbench.VirtualBench` method), 179
- `convert_values_to_timestamp()` (`pymeasure.instruments.ni.virtualbench.VirtualBench` method), 179
- `create_filename()` (in module `pymeasure.experiment.experiment`), 42
- `create_marker()` (`pymeasure.instruments.rohdeschwarz.fsl.FSL` method), 199
- `Crosshairs` (class in `pymeasure.display.curves`), 49
- `CSVFormatter` (class in `pymeasure.experiment.results`), 47
- `current` (`pymeasure.instruments.agilent.Agilent34450A` property), 74
- `CURRENT` (`pymeasure.instruments.agilent.agilentB1500.MeasOpMode` attribute), 102
- `current` (`pymeasure.instruments.bkprecision.BKPrecision9130B` property), 115
- `current` (`pymeasure.instruments.danfysik.Danfysik8500` property), 117
- `current` (`pymeasure.instruments.deltaelektronika.SM7045D` property), 119
- `current` (`pymeasure.instruments.keithley.Keithley2000` property), 127
- `current` (`pymeasure.instruments.keithley.Keithley2260B` property), 133
- `current` (`pymeasure.instruments.keithley.Keithley2400` property), 135
- `current` (`pymeasure.instruments.keithley.Keithley2450` property), 141
- `current` (`pymeasure.instruments.keithley.Keithley6517B` property), 154
- `current` (`pymeasure.instruments.keysight.KeysightN5767A` property), 163
- `current_ac` (`pymeasure.instruments.agilent.Agilent34410A` property), 72
- `current_ac` (`pymeasure.instruments.agilent.Agilent34450A` property), 74
- `current_ac` (`pymeasure.instruments.hp.HP34401A` property), 123
- `current_ac_auto_range` (`pymeasure.instruments.agilent.Agilent34450A` property), 74
- `current_ac_bandwidth` (`pymeasure.instruments.keithley.Keithley2000` property), 127
- `current_ac_digits` (`pymeasure.instruments.keithley.Keithley2000` property), 127
- `current_ac_nplc` (`pymeasure.instruments.keithley.Keithley2000` property), 127
- `current_ac_range` (`pymeasure.instruments.agilent.Agilent34450A` property), 74
- `current_ac_range` (`pymeasure.instruments.keithley.Keithley2000` property), 127
- `current_ac_reference` (`pymeasure.instruments.keithley.Keithley2000` property), 128
- `current_ac_resolution` (`pymeasure.instruments.agilent.Agilent34450A` property), 74
- `current_auto_range` (`pymeasure` property), 74

- `sure.instruments.agilent.Agilent34450A` (property), 74
- `current_dc` (`pymeasure.instruments.agilent.Agilent34410A` property), 72
- `current_dc` (`pymeasure.instruments.hp.HP34401A` property), 123
- `current_digits` (`pymeasure.instruments.keithley.Keithley2000` property), 128
- `current_filter_count` (`pymeasure.instruments.keithley.Keithley2450` property), 141
- `current_filter_state` (`pymeasure.instruments.keithley.Keithley2450` property), 141
- `current_filter_type` (`pymeasure.instruments.keithley.Keithley2450` property), 141
- `current_limit` (`pymeasure.instruments.keithley.Keithley2260B` property), 133
- `current_limit` (`pymeasure.instruments.yokogawa.YokogawaGS200` property), 216
- `current_name` (`pymeasure.instruments.agilent.agilent4156.SMU` property), 80
- `current_nplc` (`pymeasure.instruments.keithley.Keithley2000` property), 128
- `current_nplc` (`pymeasure.instruments.keithley.Keithley2400` property), 135
- `current_nplc` (`pymeasure.instruments.keithley.Keithley2450` property), 142
- `current_nplc` (`pymeasure.instruments.keithley.Keithley6517B` property), 154
- `current_output_off_state` (`pymeasure.instruments.keithley.Keithley2450` property), 142
- `current_ppm` (`pymeasure.instruments.danfysik.Danfysik8500` property), 117
- `current_range` (`pymeasure.instruments.agilent.Agilent34450A` property), 74
- `current_range` (`pymeasure.instruments.keithley.Keithley2000` property), 128
- `current_range` (`pymeasure.instruments.keithley.Keithley2400` property), 135
- `current_range` (`pymeasure.instruments.keithley.Keithley2450` property), 142
- `current_range` (`pymeasure.instruments.keithley.Keithley6517B` property), 154
- `current_range` (`pymeasure.instruments.keysight.KeysightN5767A` property), 163
- `current_reference` (`pymeasure.instruments.keithley.Keithley2000` property), 128
- `current_resolution` (`pymeasure.instruments.agilent.Agilent34450A` property), 75
- `current_setpoint` (`pymeasure.instruments.danfysik.Danfysik8500` property), 117
- `curve_buffer_bits` (`pymeasure.instruments.signalrecovery.DSP7265` property), 201
- `curve_buffer_interval` (`pymeasure.instruments.signalrecovery.DSP7265` property), 201
- `curve_buffer_length` (`pymeasure.instruments.signalrecovery.DSP7265` property), 201
- `curve_buffer_status` (`pymeasure.instruments.signalrecovery.DSP7265` property), 201
- `CustomIntEnum` (class in `pymeasure.instruments.agilent.agilentB1500`), 101
- `cw_frequency` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 189
- ## D
- `dac1` (`pymeasure.instruments.ametek.Ametek7270` property), 104
- `dac1` (`pymeasure.instruments.signalrecovery.DSP7265` property), 201
- `dac1` (`pymeasure.instruments.srs.SR830` property), 206
- `dac1` (`pymeasure.instruments.srs.SR860` property), 209
- `dac2` (`pymeasure.instruments.ametek.Ametek7270` property), 104
- `dac2` (`pymeasure.instruments.signalrecovery.DSP7265` property), 201
- `dac2` (`pymeasure.instruments.srs.SR830` property), 206
- `dac2` (`pymeasure.instruments.srs.SR860` property), 209
- `dac3` (`pymeasure.instruments.ametek.Ametek7270` property), 104
- `dac3` (`pymeasure.instruments.signalrecovery.DSP7265` property), 201
- `dac3` (`pymeasure.instruments.srs.SR830` property), 206
- `dac3` (`pymeasure.instruments.srs.SR860` property), 209

- `dac4` (`pymeasure.instruments.ametek.Ametek7270` property), 104
- `dac4` (`pymeasure.instruments.signalrecovery.DSP7265` property), 201
- `dac4` (`pymeasure.instruments.srs.SR830` property), 206
- `dac4` (`pymeasure.instruments.srs.SR860` property), 209
- `Danfysik8500` (class in `pymeasure.instruments.danfysik`), 116
- `DanfysikAdapter` (class in `pymeasure.instruments.danfysik`), 116
- `data` (`pymeasure.experiment.experiment.Experiment` property), 41
- `data` (`pymeasure.instruments.agilent.Agilent8722ES` property), 69
- `data_arb()` (`pymeasure.instruments.agilent.Agilent33500` method), 86
- `data_complex` (`pymeasure.instruments.agilent.Agilent8722ES` property), 69
- `data_format()` (`pymeasure.instruments.agilent.agilentB1500.AgilentB1500` method), 94
- `data_log_magnitude` (`pymeasure.instruments.agilent.Agilent8722ES` property), 69
- `data_magnitude` (`pymeasure.instruments.agilent.Agilent8722ES` property), 69
- `data_memory_a_condition` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 110
- `data_memory_a_size` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 110
- `data_memory_a_values` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 110
- `data_memory_b_condition` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 111
- `data_memory_b_size` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 111
- `data_memory_b_values` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 111
- `data_memory_select` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 111
- `data_phase` (`pymeasure.instruments.agilent.Agilent8722ES` property), 69
- `data_variables` (`pymeasure.instruments.agilent.agilent4156.Agilent4156` property), 78
- `data_volatile_clear()` (`pymeasure.instruments.agilent.Agilent33500` method), 86
- `date` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 189
- `dc_mode()` (`pymeasure.instruments.lakeshore.LakeShore425` method), 166
- `dcmode` (`pymeasure.instruments.srs.SR860` property), 209
- `decode_mode()` (`pymeasure.instruments.hp.HP3478A` class method), 124
- `decode_range()` (`pymeasure.instruments.hp.HP3478A` class method), 124
- `decode_status()` (`pymeasure.instruments.hp.HP3478A` class method), 124
- `decode_trigger()` (`pymeasure.instruments.hp.HP3478A` static method), 125
- `default_setup()` (`pymeasure.instruments.keysight.KeysightDSOX1102G` method), 160
- `define_arbitrary_waveform()` (`pymeasure.instruments.keithley.Keithley6221` method), 149
- `define_position()` (`pymeasure.instruments.newport.esp300.Axis` method), 167
- `delay_time` (`pymeasure.instruments.agilent.agilent4156.Agilent4156` property), 78
- `detectedfrequency` (`pymeasure.instruments.srs.SR860` property), 209
- `determine_valid_channels()` (`pymeasure.instruments.keithley.Keithley2700` method), 147
- `deviation` (`pymeasure.instruments.rohdeschwarz.sfm.Sound_Channel` property), 196
- `digitize()` (`pymeasure.instruments.keysight.KeysightDSOX1102G` method), 160
- `diode` (`pymeasure.instruments.agilent.Agilent34450A` property), 75
- `dip_search` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 111
- `direction` (`pymeasure.instruments.anaheimautomation.DPSeriesMotorCo` property), 107
- `DirectoryLineEdit` (class in `pymeasure.display.widgets`), 54
- `disable` (`pymeasure.instruments.agilent.agilent4156.SMU` property), 80
- `disable` (`pymeasure.instruments.agilent.agilent4156.VMU` property), 82
- `disable` (`pymeasure.instruments.agilent.agilent4156.VSU` property), 82
- `disable()` (`pymeasure.instruments.agilent.Agilent8257D` method), 82

- method*), 67
- `disable()` (*pymeasure.instruments.agilent.agilentB1500.SMW* *method*), 97
- `disable()` (*pymeasure.instruments.anritsu.AnritsuMG3692C* *method*), 110
- `disable()` (*pymeasure.instruments.danfysik.Danfysik8500* *method*), 117
- `disable()` (*pymeasure.instruments.deltaelektronika.SM7045D* *method*), 119
- `disable()` (*pymeasure.instruments.keysight.KeysightN5767A* *method*), 163
- `disable()` (*pymeasure.instruments.newport.ESP300* *method*), 167
- `disable()` (*pymeasure.instruments.newport.esp300.Axis* *method*), 167
- `disable()` (*pymeasure.instruments.parker.ParkerGV6* *method*), 183
- `disable_all()` (*pymeasure.instruments.agilent.agilent4156.Agilent4156* *method*), 78
- `disable_amplitude_modulation()` (*pymeasure.instruments.agilent.Agilent8257D* *method*), 67
- `disable_averaging()` (*pymeasure.instruments.agilent.Agilent8722ES* *method*), 69
- `disable_bias()` (*pymeasure.instruments.srs.SR570* *method*), 204
- `disable_buffer()` (*pymeasure.instruments.keithley.Keithley2000* *method*), 128
- `disable_buffer()` (*pymeasure.instruments.keithley.Keithley2400* *method*), 135
- `disable_buffer()` (*pymeasure.instruments.keithley.Keithley2450* *method*), 142
- `disable_buffer()` (*pymeasure.instruments.keithley.Keithley2700* *method*), 147
- `disable_buffer()` (*pymeasure.instruments.keithley.Keithley6221* *method*), 150
- `disable_buffer()` (*pymeasure.instruments.keithley.Keithley6517B* *method*), 154
- `disable_filter()` (*pymeasure.instruments.keithley.Keithley2000* *method*), 128
- `disable_heater()` (*pymeasure.instruments.lakeshore.LakeShore331* *method*), 165
- `disable_low_freq_out()` (*pymeasure.instruments.agilent.Agilent8257D* *method*), 67
- `disable_modulation()` (*pymeasure.instruments.agilent.Agilent8257D* *method*), 67
- `disable_offset_current()` (*pymeasure.instruments.srs.SR570* *method*), 204
- `disable_output_trigger()` (*pymeasure.instruments.keithley.Keithley2400* *method*), 135
- `disable_output_trigger()` (*pymeasure.instruments.keithley.Keithley6221* *method*), 150
- `disable_persistent_switch()` (*pymeasure.instruments.ami.AMI430* *method*), 106
- `disable_pulse_modulation()` (*pymeasure.instruments.agilent.Agilent8257D* *method*), 67
- `disable_reference()` (*pymeasure.instruments.keithley.Keithley2000* *method*), 128
- `disable_rf()` (*pymeasure.instruments.anapico.APSIN12G* *method*), 109
- `disable_source()` (*pymeasure.instruments.keithley.Keithley2400* *method*), 135
- `disable_source()` (*pymeasure.instruments.keithley.Keithley2450* *method*), 142
- `disable_source()` (*pymeasure.instruments.keithley.Keithley6221* *method*), 150
- `disable_source()` (*pymeasure.instruments.keithley.Keithley6517B* *method*), 154
- `disable_source()` (*pymeasure.instruments.yokogawa.Yokogawa7651* *method*), 215
- `display` (*pymeasure.instruments.agilent.Agilent33500* *property*), 86
- `display_closed_channels()` (*pymeasure.instruments.keithley.Keithley2700* *method*), 147
- `display_enabled` (*pymeasure.instruments.keithley.Keithley2400* *property*), 135
- `display_enabled` (*pymeasure.instruments.keithley.Keithley6221* *property*), 150
- `display_estimates()` (*pymeasure.display.widgets.EstimatorWidget* *method*), 54
- `display_reset()` (*pymeasure.instruments.hp.HP3478A* *method*), 125

- [display\\_text](#) ([pymeasure.instruments.hp.HP3478A](#) property), 125  
[display\\_text](#) ([pymeasure.instruments.keithley.Keithley2700](#) property), 147  
[display\\_text\\_no\\_symbol](#) ([pymeasure.instruments.hp.HP3478A](#) property), 125  
[download\\_data\(\)](#) ([pymeasure.instruments.keysight.KeysightDSOX1102G](#) method), 160  
[download\\_image\(\)](#) ([pymeasure.instruments.keysight.KeysightDSOX1102G](#) method), 160  
[DPSeriesMotorController](#) (class in [pymeasure.instruments.anaheimautomation](#)), 107  
[DSP7265](#) (class in [pymeasure.instruments.signalrecovery](#)), 200  
[dwell\\_time](#) ([pymeasure.instruments.agilent.Agilent8257D](#) property), 67
- ## E
- [echo\(\)](#) ([pymeasure.instruments.parker.ParkerGV6](#) method), 183  
[emit\(\)](#) ([pymeasure.display.log.LogHandler](#) method), 52  
[emit\(\)](#) ([pymeasure.experiment.workers.Worker](#) method), 46  
[enable\(\)](#) ([pymeasure.instruments.agilent.Agilent8257D](#) method), 67  
[enable\(\)](#) ([pymeasure.instruments.agilent.agilentB1500.SMU](#) method), 97  
[enable\(\)](#) ([pymeasure.instruments.anritsu.AnritsuMG3692C](#) method), 110  
[enable\(\)](#) ([pymeasure.instruments.danfysik.Danfysik8500](#) method), 117  
[enable\(\)](#) ([pymeasure.instruments.deltaelektronika.SM7045B](#) method), 119  
[enable\(\)](#) ([pymeasure.instruments.keysight.KeysightN5767A](#) method), 163  
[enable\(\)](#) ([pymeasure.instruments.newport.ESP300](#) method), 167  
[enable\(\)](#) ([pymeasure.instruments.newport.esp300.Axis](#) method), 168  
[enable\(\)](#) ([pymeasure.instruments.parker.ParkerGV6](#) method), 183  
[enable\\_amplitude\\_modulation\(\)](#) ([pymeasure.instruments.agilent.Agilent8257D](#) method), 67  
[enable\\_averaging\(\)](#) ([pymeasure.instruments.agilent.Agilent8722ES](#) method), 69  
[enable\\_bias\(\)](#) ([pymeasure.instruments.srs.SR570](#) method), 204  
[enable\\_filter\(\)](#) ([pymeasure.instruments.keithley.Keithley2000](#) method), 128  
[enable\\_low\\_freq\\_out\(\)](#) ([pymeasure.instruments.agilent.Agilent8257D](#) method), 67  
[enable\\_offset\\_current\(\)](#) ([pymeasure.instruments.srs.SR570](#) method), 204  
[enable\\_persistent\\_switch\(\)](#) ([pymeasure.instruments.ami.AMI430](#) method), 106  
[enable\\_pulse\\_modulation\(\)](#) ([pymeasure.instruments.agilent.Agilent8257D](#) method), 67  
[enable\\_reference\(\)](#) ([pymeasure.instruments.keithley.Keithley2000](#) method), 128  
[enable\\_rf\(\)](#) ([pymeasure.instruments.anapico.APSIN12G](#) method), 109  
[enable\\_source\(\)](#) ([pymeasure.instruments.keithley.Keithley2400](#) method), 135  
[enable\\_source\(\)](#) ([pymeasure.instruments.keithley.Keithley2450](#) method), 142  
[enable\\_source\(\)](#) ([pymeasure.instruments.keithley.Keithley6221](#) method), 150  
[enable\\_source\(\)](#) ([pymeasure.instruments.keithley.Keithley6517B](#) method), 155  
[enable\\_source\(\)](#) ([pymeasure.instruments.yokogawa.Yokogawa7651](#) method), 215  
[enabled](#) ([pymeasure.instruments.keithley.Keithley2260B](#) property), 133  
[enabled](#) ([pymeasure.instruments.newport.esp300.Axis](#) property), 168  
[encoder\\_autocorrect](#) ([pymeasure.instruments.anaheimautomation.DPSeriesMotorController](#) property), 108  
[encoder\\_delay](#) ([pymeasure.instruments.anaheimautomation.DPSeriesMotorController](#) property), 108  
[encoder\\_enabled](#) ([pymeasure.instruments.anaheimautomation.DPSeriesMotorController](#) property), 108  
[encoder\\_motor\\_ratio](#) ([pymeasure.instruments.anaheimautomation.DPSeriesMotorController](#) property), 108  
[encoder\\_retries](#) ([pymeasure.instruments.anaheimautomation.DPSeriesMotorController](#) property), 108  
[encoder\\_window](#) ([pymeasure.instruments.anaheimautomation.DPSeriesMotorController](#) property), 108



- property*), 108
  - `energy` (*pymeasure.instruments.thorlabs.ThorlabsPM100USB* *property*), 214
  - `err_status` (*pymeasure.instruments.srs.SR830* *property*), 206
  - `error` (*pymeasure.instruments.keithley.Keithley2260B* *property*), 133
  - `error` (*pymeasure.instruments.keithley.Keithley2400* *property*), 135
  - `error` (*pymeasure.instruments.keithley.Keithley2450* *property*), 142
  - `error` (*pymeasure.instruments.keithley.Keithley2600* *property*), 158
  - `error` (*pymeasure.instruments.keithley.Keithley2700* *property*), 147
  - `error` (*pymeasure.instruments.keithley.Keithley6221* *property*), 150
  - `error` (*pymeasure.instruments.keithley.Keithley6517B* *property*), 155
  - `error` (*pymeasure.instruments.newport.ESP300* *property*), 167
  - `error_reg` (*pymeasure.instruments.anaheimautomation.DPSeries4400* *property*), 108
  - `error_status` (*pymeasure.instruments.hp.HP3478A* *property*), 125
  - `errors` (*pymeasure.instruments.newport.ESP300* *property*), 167
  - `ese2` (*pymeasure.instruments.anritsu.AnritsuMS9710C* *property*), 111
  - `ESP300` (*class in pymeasure.instruments.newport*), 167
  - `esr2` (*pymeasure.instruments.anritsu.AnritsuMS9710C* *property*), 111
  - `EstimatorThread` (*class in pymeasure.display.widgets*), 54
  - `EstimatorWidget` (*class in pymeasure.display.widgets*), 54
  - `eval_string()` (*pymeasure.display.widgets.SequencerWidget* *static method*), 56
  - `event_reg` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* *property*), 189
  - `execute()` (*pymeasure.experiment.procedure.Procedure* *method*), 43
  - `expand_channel_string()` (*pymeasure.instruments.ni.virtualbench.VirtualBench* *method*), 180
  - `Experiment` (*class in pymeasure.display.manager*), 52
  - `Experiment` (*class in pymeasure.experiment.experiment*), 41
  - `ExperimentQueue` (*class in pymeasure.display.manager*), 52
  - `export_signal()` (*pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalInputQuantum* *method*), 170
  - `ext_ref_base_unit` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* *property*), 189
  - `ext_ref_extension` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* *property*), 189
  - `ext_trig_out` (*pymeasure.instruments.agilent.Agilent33500* *property*), 86
  - `ext_vid_connector` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* *property*), 189
  - `external_arming_start_slope` (*pymeasure.instruments.pendulum.cnt91.CNT91* *property*), 184
  - `external_modulation_frequency` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* *property*), 189
  - `external_modulation_power` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* *property*), 189
  - `external_modulation_source` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* *property*), 190
  - `external_start_arming_source` (*pymeasure.instruments.pendulum.cnt91.CNT91* *property*), 184
  - `extfrequency` (*pymeasure.instruments.srs.SR860* *property*), 209
  - `extract_value()` (*pymeasure.instruments.attocube.adapters.AttocubeConsoleAdapter* *method*), 113
- ## F
- `factory_reset()` (*pymeasure.instruments.keysight.KeysightDSOX1102G* *method*), 160
  - `FakeAdapter` (*class in pymeasure.adapters*), 32
  - `field` (*pymeasure.instruments.ami.AMI430* *property*), 106
  - `field` (*pymeasure.instruments.fwbell.FWBell5080* *property*), 121
  - `field` (*pymeasure.instruments.lakeshore.LakeShore425* *property*), 166
  - `fields()` (*pymeasure.instruments.fwbell.FWBell5080* *method*), 121
  - `filer_synchronous` (*pymeasure.instruments.srs.SR860* *property*), 209
  - `filter` (*pymeasure.instruments.agilent.agilentB1500.SMU* *property*), 97
  - `filter` (*pymeasure.instruments.ni.virtualbench.VirtualBench.FunctionGeneratorAdvanced* *property*), 173
  - `filter_advanced` (*pymeasure.instruments.srs.SR860* *property*), 209

[filter\\_count](#) (`pymeasure.instruments.keithley.Keithley2400` property), 135  
[filter\\_slope](#) (`pymeasure.instruments.srs.SR830` property), 206  
[filter\\_slope](#) (`pymeasure.instruments.srs.SR860` property), 209  
[filter\\_state](#) (`pymeasure.instruments.keithley.Keithley2400` property), 136  
[filter\\_type](#) (`pymeasure.instruments.keithley.Keithley2400` property), 136  
[filter\\_type](#) (`pymeasure.instruments.srs.SR570` property), 204  
[find\\_img\\_index\(\)](#) (`pymeasure.display.curves.ResultsImage` method), 50  
[FloatInput](#) (class in `pymeasure.display.inputs`), 50  
[FloatParameter](#) (class in `pymeasure.experiment.parameters`), 44  
[Fluke7341](#) (class in `pymeasure.instruments.fluke`), 120  
[force\(\)](#) (`pymeasure.instruments.agilent.agilentB1500.SMU` method), 98  
[force\\_gnd\(\)](#) (`pymeasure.instruments.agilent.agilentB1500.SMU` method), 94  
[force\\_gnd\(\)](#) (`pymeasure.instruments.agilent.agilentB1500.SMU` method), 97  
[FORCE\\_SIDE](#) (`pymeasure.instruments.agilent.agilentB1500` attribute), 102  
[force\\_trigger\(\)](#) (`pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope` method), 176  
[format](#) (`pymeasure.instruments.pendulum.cnt91.CNT91` property), 184  
[format\(\)](#) (`pymeasure.experiment.results.CSVFormatter` method), 47  
[format\(\)](#) (`pymeasure.experiment.results.Results` method), 47  
[freq\\_center](#) (`pymeasure.instruments.rohdeschwarz.fsl.FSL` property), 199  
[freq\\_span](#) (`pymeasure.instruments.rohdeschwarz.fsl.FSL` property), 199  
[freq\\_start](#) (`pymeasure.instruments.rohdeschwarz.fsl.FSL` property), 199  
[freq\\_stop](#) (`pymeasure.instruments.rohdeschwarz.fsl.FSL` property), 199  
[freq\\_sweep\(\)](#) (`pymeasure.instruments.agilent.AgilentE4980` method), 71  
[frequencies](#) (`pymeasure.instruments.agilent.Agilent8722ES` property), 69  
[frequencies](#) (`pymeasure.instruments.agilent.AgilentE4408B` property), 70  
[frequency](#) (`pymeasure.instruments.agilent.Agilent33220A` property), 83  
[frequency](#) (`pymeasure.instruments.agilent.Agilent33500` property), 86  
[frequency](#) (`pymeasure.instruments.agilent.Agilent33521A` property), 88  
[frequency](#) (`pymeasure.instruments.agilent.Agilent34450A` property), 75  
[frequency](#) (`pymeasure.instruments.agilent.Agilent8257D` property), 68  
[frequency](#) (`pymeasure.instruments.agilent.AgilentE4980` property), 71  
[frequency](#) (`pymeasure.instruments.ametek.Ametek7270` property), 104  
[frequency](#) (`pymeasure.instruments.anapico.APSIN12G` property), 109  
[frequency](#) (`pymeasure.instruments.anritsu.AnritsuMG3692C` property), 110  
[frequency](#) (`pymeasure.instruments.attocube.anc300.Axis` property), 114  
[frequency](#) (`pymeasure.instruments.hp.HP33120A` property), 122  
[frequency](#) (`pymeasure.instruments.keithley.Keithley2000` property), 128  
[frequency](#) (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 190  
[frequency](#) (`pymeasure.instruments.rohdeschwarz.sfm.Sound_Channel` property), 196  
[frequency](#) (`pymeasure.instruments.signalrecovery.DSP7265` property), 201  
[frequency](#) (`pymeasure.instruments.srs.SR510` property), 206  
[frequency](#) (`pymeasure.instruments.srs.SR830` property), 206  
[frequency](#) (`pymeasure.instruments.srs.SR860` property), 209  
[frequency\\_aperature](#) (`pymeasure.instruments.keithley.Keithley2000` property), 128  
[frequency\\_aperture](#) (`pymeasure.instruments.agilent.Agilent34450A` property), 75  
[frequency\\_current\\_auto\\_range](#) (`pymeasure.instruments.agilent.Agilent34450A` property), 75  
[frequency\\_current\\_range](#) (`pymeasure.instruments.agilent.Agilent34450A` property), 75  
[frequency\\_digits](#) (`pymeasure.instruments.keithley.Keithley2000` property), 128  
[frequency\\_mode](#) (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 190  
[frequency\\_points](#) (`pymeasure.instruments.agilent.Agilent33220A` property), 83

- sure.instruments.agilent.AgilentE4408B* property), 70
- frequency\_reference** (*pymea-  
sure.instruments.keithley.Keithley2000* prop-  
erty), 128
- frequency\_step** (*pymea-  
sure.instruments.agilent.AgilentE4408B*  
property), 70
- frequency\_threshold** (*pymea-  
sure.instruments.keithley.Keithley2000* prop-  
erty), 129
- frequency\_voltage\_auto\_range** (*pymea-  
sure.instruments.agilent.Agilent34450A*  
property), 75
- frequency\_voltage\_range** (*pymea-  
sure.instruments.agilent.Agilent34450A*  
property), 75
- frequencypreset1** (*pymea-  
sure.instruments.srs.SR860* property), 209
- frequencypreset2** (*pymea-  
sure.instruments.srs.SR860* property), 210
- frequencypreset3** (*pymea-  
sure.instruments.srs.SR860* property), 210
- frequencypreset4** (*pymea-  
sure.instruments.srs.SR860* property), 210
- front\_blanked** (*pymea-  
sure.instruments.srs.SR570* property), 204
- front\_panel** (*pymea-  
sure.instruments.srs.SR860* prop-  
erty), 210
- FSL** (class in *pymea-  
sure.instruments.rohdeschwarz.fsl*), 199
- FWBell15080** (class in *pymea-  
sure.instruments.fwbell*), 121
- ## G
- gain\_mode** (*pymea-  
sure.instruments.srs.SR570* property), 204
- gasflow** (*pymea-  
sure.instruments.oxfordinstruments.ITC503* property), 181
- gen\_measurement()** (*pymea-  
sure.experiment.procedure.Procedure* method), 43
- GeneralError** (class in *pymea-  
sure.instruments.newport.esp300*), 168
- get()** (*pymea-  
sure.instruments.agilent.agilentB1500.Custom* class method), 101
- get\_array()** (in module *pymea-  
sure.experiment.experiment*), 42
- get\_array\_steps()** (in module *pymea-  
sure.experiment.experiment*), 42
- get\_array\_zero()** (in module *pymea-  
sure.experiment.experiment*), 42
- get\_buffer()** (*pymea-  
sure.instruments.signalrecovery.DSP7265* method), 201
- get\_buffer()** (*pymea-  
sure.instruments.srs.SR830* method), 207
- get\_calibration\_information()** (*pymea-  
sure.instruments.ni.virtualbench.VirtualBench* method), 180
- get\_data()** (*pymea-  
sure.instruments.agilent.agilent4156.Agilent4156* method), 78
- get\_estimates()** (*pymea-  
sure.display.widgets.EstimatorWidget* method), 54
- get\_estimates()** (*pymea-  
sure.experiment.procedure.Procedure* method), 43
- get\_library\_version()** (*pymea-  
sure.instruments.ni.virtualbench.VirtualBench* method), 180
- get\_noise\_bandwidth** (*pymea-  
sure.instruments.srs.SR860* property), 210
- get\_procedure()** (*pymea-  
sure.display.widgets.InputsWidget* method), 55
- get\_scaling()** (*pymea-  
sure.instruments.srs.SR830* method), 207
- get\_sequence\_from\_tree()** (*pymea-  
sure.display.widgets.SequencerWidget* method), 56
- get\_signal\_strength\_indicator** (*pymea-  
sure.instruments.srs.SR860* property), 210
- get\_state\_of\_channels()** (*pymea-  
sure.instruments.keithley.Keithley2700* method), 147
- get\_status()** (*pymea-  
sure.instruments.hp.HP3478A* method), 125
- get\_time()** (*pymea-  
sure.instruments.Mock* method), 63
- get\_voltage()** (*pymea-  
sure.instruments.Mock* method), 63
- get\_wave()** (*pymea-  
sure.instruments.Mock* method), 63
- getAI()** (in module *pymea-  
sure.instruments.comedi*), 65
- getAO()** (in module *pymea-  
sure.instruments.comedi*), 65
- gettimebase** (*pymea-  
sure.instruments.srs.SR860* prop-  
erty), 210
- gpib()** (*pymea-  
sure.adapters.PrologixAdapter* method), 35
- gpib\_address** (*pymea-  
sure.instruments.rohdeschwarz.sfm.SFM* property), 190
- GPiB\_trigger()** (*pymea-  
sure.instruments.hp.HP3478A* method), 123
- ground\_all()** (*pymea-  
sure.instruments.attocube.anc300.ANC300Controller* method), 113



## H

- `handle_abort()` (*pymea-*  
*sure.experiment.workers.Worker*  
method), 46
- `handle_error()` (*pymea-*  
*sure.experiment.workers.Worker*  
method), 46
- `harmonic` (*pymea-*  
*sure.instruments.ametek.Ametek7270*  
property), 104
- `harmonic` (*pymea-*  
*sure.instruments.signalrecovery.DSP7265*  
property), 202
- `harmonic` (*pymea-*  
*sure.instruments.srs.SR830* property), 207
- `harmonic` (*pymea-*  
*sure.instruments.srs.SR860* property), 210
- `harmonicdual` (*pymea-*  
*sure.instruments.srs.SR860* prop-  
erty), 210
- `has_amplitude_modulation` (*pymea-*  
*sure.instruments.agilent.Agilent8257D* prop-  
erty), 68
- `has_modulation` (*pymea-*  
*sure.instruments.agilent.Agilent8257D* prop-  
erty), 68
- `has_next()` (*pymea-*  
*sure.display.manager.ExperimentQueue*  
method), 52
- `has_persistent_switch_enabled()` (*pymea-*  
*sure.instruments.ami.AMI430* method), 106
- `has_pulse_modulation` (*pymea-*  
*sure.instruments.agilent.Agilent8257D* prop-  
erty), 68
- `has_supported_version()` (*pymea-*  
*sure.adapters.VISAAdapter* static method), 37
- `header()` (*pymea-*  
*sure.experiment.results.Results*  
method), 47
- `heater` (*pymea-*  
*sure.instruments.oxfordinstruments.ITC503*  
property), 181
- `heater_gas_mode` (*pymea-*  
*sure.instruments.oxfordinstruments.ITC503*  
property), 181
- `heater_range` (*pymea-*  
*sure.instruments.lakeshore.LakeShore331*  
property), 165
- `high_freq` (*pymea-*  
*sure.instruments.srs.SR570* property), 205
- `high_frequency_resolution` (*pymea-*  
*sure.instruments.rohdeschwarz.sfm.SFM*  
property), 190
- `hold_time` (*pymea-*  
*sure.instruments.agilent.agilent4156.Agilent4156*  
property), 78
- `home()` (*pymea-*  
*sure.instruments.anaheimautomation.DPSeriesMotionController*  
method), 108
- `home()` (*pymea-*  
*sure.instruments.newport.esp300.Axis*  
method), 168
- `horizontal_time_div` (*pymea-*  
*sure.instruments.srs.SR860* property), 210
- `HP33120A` (class in *pymea-*  
*sure.instruments.hp*), 122
- `HP34401A` (class in *pymea-*  
*sure.instruments.hp*), 123
- `HP3478A` (class in *pymea-*  
*sure.instruments.hp*), 123
- `HP3478A.ERRORS` (class in *pymea-*  
*sure.instruments.hp*), 123
- `HP3478A.SRQ` (class in *pymea-*  
*sure.instruments.hp*), 123
- `HRADC` (*pymea-*  
*sure.instruments.agilent.agilentB1500.ADCType*  
attribute), 101
- `HSADC` (*pymea-*  
*sure.instruments.agilent.agilentB1500.ADCType*  
attribute), 101
- `HSADC_PULSED` (*pymea-*  
*sure.instruments.agilent.agilentB1500.ADCType*  
attribute), 101
- `id` (*pymea-*  
*sure.instruments.advantest.advantestR3767CG.AdvantestR3767CG*  
property), 66
- `id` (*pymea-*  
*sure.instruments.agilent.Agilent33500* prop-  
erty), 86
- `id` (*pymea-*  
*sure.instruments.ametek.Ametek7270* prop-  
erty), 104
- `id` (*pymea-*  
*sure.instruments.bkprecision.BKPrecision9130B*  
property), 115
- `id` (*pymea-*  
*sure.instruments.danfysik.Danfysik8500* prop-  
erty), 117
- `id` (*pymea-*  
*sure.instruments.fluke.Fluke7341* property), 120
- `id` (*pymea-*  
*sure.instruments.fwbell.FWBell5080* property), 121
- `id` (*pymea-*  
*sure.instruments.Instrument* property), 62
- `id` (*pymea-*  
*sure.instruments.keithley.Keithley2000* prop-  
erty), 129
- `id` (*pymea-*  
*sure.instruments.keithley.Keithley2260B* prop-  
erty), 133
- `id` (*pymea-*  
*sure.instruments.keithley.Keithley2400* prop-  
erty), 136
- `id` (*pymea-*  
*sure.instruments.keithley.Keithley2450* prop-  
erty), 142
- `id` (*pymea-*  
*sure.instruments.keithley.Keithley2600* prop-  
erty), 158
- `id` (*pymea-*  
*sure.instruments.keithley.Keithley2700* prop-  
erty), 147
- `id` (*pymea-*  
*sure.instruments.keithley.Keithley2750* prop-  
erty), 157
- `id` (*pymea-*  
*sure.instruments.keithley.Keithley6221* prop-  
erty), 150
- `id` (*pymea-*  
*sure.instruments.keithley.Keithley6517B* prop-  
erty), 155
- `id` (*pymea-*  
*sure.instruments.keysight.KeysightDSOX1102G*  
property), 160
- `id` (*pymea-*  
*sure.instruments.keysight.KeysightN5767A*  
property), 163

<code>id</code> ( <code>pymeasure.instruments.signalrecovery.DSP7265</code> property), 202	<code>internal_frequency</code> ( <code>pymeasure.instruments.agilent.Agilent8257D</code> property), 68
<code>id</code> ( <code>pymeasure.instruments.yokogawa.Yokogawa7651</code> property), 215	<code>internal_shape</code> ( <code>pymeasure.instruments.agilent.Agilent8257D</code> property), 68
<code>ImageFrame</code> (class in <code>pymeasure.display.widgets</code> ), 55	<code>internalfrequency</code> ( <code>pymeasure.instruments.srs.SR860</code> property), 211
<code>ImageWidget</code> (class in <code>pymeasure.display.widgets</code> ), 55	<code>interpolator_autocalibrated</code> ( <code>pymeasure.instruments.pendulum.cnt91.CNT91</code> property), 184
<code>imode</code> ( <code>pymeasure.instruments.signalrecovery.DSP7265</code> property), 202	<code>invert_signal_sign</code> ( <code>pymeasure.instruments.srs.SR570</code> property), 205
<code>impedance</code> ( <code>pymeasure.instruments.agilent.AgilentE4980</code> property), 71	<code>is_averaging()</code> ( <code>pymeasure.instruments.agilent.Agilent8722ES</code> method), 69
<code>init_curve_buffer()</code> ( <code>pymeasure.instruments.signalrecovery.DSP7265</code> method), 202	<code>is_buffer_full()</code> ( <code>pymeasure.instruments.keithley.Keithley2000</code> method), 129
<code>initialize_all_smus()</code> ( <code>pymeasure.instruments.agilent.agilentB1500.AgilentB1500</code> method), 93	<code>is_buffer_full()</code> ( <code>pymeasure.instruments.keithley.Keithley2400</code> method), 136
<code>initialize_smu()</code> ( <code>pymeasure.instruments.agilent.agilentB1500.AgilentB1500</code> method), 93	<code>is_buffer_full()</code> ( <code>pymeasure.instruments.keithley.Keithley2450</code> method), 142
<code>Input</code> (class in <code>pymeasure.display.inputs</code> ), 50	<code>is_buffer_full()</code> ( <code>pymeasure.instruments.keithley.Keithley2700</code> method), 147
<code>input_config</code> ( <code>pymeasure.instruments.srs.SR830</code> property), 207	<code>is_buffer_full()</code> ( <code>pymeasure.instruments.keithley.Keithley6221</code> method), 150
<code>input_coupling</code> ( <code>pymeasure.instruments.srs.SR830</code> property), 207	<code>is_buffer_full()</code> ( <code>pymeasure.instruments.keithley.Keithley6517B</code> method), 155
<code>input_coupling</code> ( <code>pymeasure.instruments.srs.SR860</code> property), 210	<code>is_current_stable()</code> ( <code>pymeasure.instruments.danfysik.Danfysik8500</code> method), 117
<code>input_current_gain</code> ( <code>pymeasure.instruments.srs.SR860</code> property), 210	<code>is_enabled</code> ( <code>pymeasure.instruments.agilent.Agilent8257D</code> property), 68
<code>input_grounding</code> ( <code>pymeasure.instruments.srs.SR830</code> property), 207	<code>is_enabled()</code> ( <code>pymeasure.instruments.danfysik.Danfysik8500</code> method), 117
<code>input_notch_config</code> ( <code>pymeasure.instruments.srs.SR830</code> property), 207	<code>is_enabled()</code> ( <code>pymeasure.instruments.keysight.KeysightN5767A</code> method), 163
<code>input_range</code> ( <code>pymeasure.instruments.srs.SR860</code> property), 210	<code>is_moving()</code> ( <code>pymeasure.instruments.parker.ParkerGV6</code> method), 183
<code>input_shields</code> ( <code>pymeasure.instruments.srs.SR860</code> property), 210	<code>is_out_of_range()</code> ( <code>pymeasure.instruments.srs.SR830</code> method), 207
<code>input_signal</code> ( <code>pymeasure.instruments.srs.SR860</code> property), 210	<code>is_ready()</code> ( <code>pymeasure.instruments.danfysik.Danfysik8500</code> method), 117
<code>input_voltage_mode</code> ( <code>pymeasure.instruments.srs.SR860</code> property), 210	<code>is_running()</code> ( <code>pymeasure.display.manager.Manager</code> method), 53
<code>InputsWidget</code> (class in <code>pymeasure.display.widgets</code> ), 55	<code>is_sequence_running()</code> ( <code>pymeasure</code>
<code>instant_voltage_1</code> ( <code>pymeasure.instruments.razorbill.razorbillRP100</code> property), 185	
<code>instant_voltage_2</code> ( <code>pymeasure.instruments.razorbill.razorbillRP100</code> property), 185	
<code>Instrument</code> (class in <code>pymeasure.instruments</code> ), 61	
<code>IntegerInput</code> (class in <code>pymeasure.display.inputs</code> ), 51	
<code>IntegerParameter</code> (class in <code>pymeasure.experiment.parameters</code> ), 44	
<code>integration_time</code> ( <code>pymeasure.instruments.agilent.agilent4156.Agilent4156</code> property), 79	

- sure.instruments.danfysik.Danfysik8500* method), 117
- `is_set()` (*pymasure.experiment.parameters.Parameter* method), 46
- ITC503 (class in *pymasure.instruments.oxfordinstruments*), 180
- ## J
- `join()` (*pymasure.display.thread.StoppableQThread* method), 54
- `join()` (*pymasure.experiment.workers.Worker* method), 46
- ## K
- Keithley2000 (class in *pymasure.instruments.keithley*), 126
- Keithley2260B (class in *pymasure.instruments.keithley*), 132
- Keithley2400 (class in *pymasure.instruments.keithley*), 134
- Keithley2450 (class in *pymasure.instruments.keithley*), 140
- Keithley2600 (class in *pymasure.instruments.keithley*), 158
- Keithley2700 (class in *pymasure.instruments.keithley*), 146
- Keithley2750 (class in *pymasure.instruments.keithley*), 157
- Keithley6221 (class in *pymasure.instruments.keithley*), 148
- Keithley6517B (class in *pymasure.instruments.keithley*), 153
- KeysightDSOX1102G (class in *pymasure.instruments.keysight*), 158
- KeysightN5767A (class in *pymasure.instruments.keysight*), 163
- `kill()` (*pymasure.instruments.parker.ParkerGV6* method), 183
- ## L
- `labels()` (*pymasure.experiment.results.Results* method), 47
- LakeShore331 (class in *pymasure.instruments.lakeshore*), 165
- LakeShore425 (class in *pymasure.instruments.lakeshore*), 166
- LakeShoreUSBAdapter (class in *pymasure.instruments.lakeshore*), 164
- LDCCurrent (*pymasure.instruments.thorlabs.ThorlabsPro8000* property), 214
- LDCCurrentLimit (*pymasure.instruments.thorlabs.ThorlabsPro8000* property), 214
- LDCPolarity (*pymasure.instruments.thorlabs.ThorlabsPro8000* property), 214
- LDCStatus (*pymasure.instruments.thorlabs.ThorlabsPro8000* property), 214
- `left_limit` (*pymasure.instruments.newport.esp300.Axis* property), 168
- `level` (*pymasure.instruments.rohdeschwarz.sfm.SFM* property), 190
- `level_lin` (*pymasure.instruments.anritsu.AnritsuMS9710C* property), 111
- `level_log` (*pymasure.instruments.anritsu.AnritsuMS9710C* property), 111
- `level_mode` (*pymasure.instruments.rohdeschwarz.sfm.SFM* property), 190
- `level_opt_attn` (*pymasure.instruments.anritsu.AnritsuMS9710C* property), 111
- `level_scale` (*pymasure.instruments.anritsu.AnritsuMS9710C* property), 111
- `lia_status` (*pymasure.instruments.srs.SR830* property), 207
- `line_frequency` (*pymasure.instruments.keithley.Keithley2400* property), 136
- `line_frequency_auto` (*pymasure.instruments.keithley.Keithley2400* property), 136
- LINEAR (*pymasure.instruments.agilent.agilentB1500.SamplingMode* attribute), 102
- LINEAR\_DOUBLE (*pymasure.instruments.agilent.agilentB1500.SweepMode* attribute), 102
- LINEAR\_SINGLE (*pymasure.instruments.agilent.agilentB1500.SweepMode* attribute), 102
- `list_resources()` (in module *pymasure.instruments*), 65
- `Listener` (class in *pymasure.experiment.listeners*), 42
- `ListInput` (class in *pymasure.display.inputs*), 51
- `ListParameter` (class in *pymasure.experiment.parameters*), 45
- `load()` (*pymasure.display.manager.Manager* method), 53
- `load()` (*pymasure.display.widgets.ImageWidget* method), 55
- `load()` (*pymasure.display.widgets.PlotWidget* method), 55
- `load()` (*pymasure.display.widgets.TabWidget* method), 56
- `load()` (*pymasure.experiment.results.Results* static method), 47
- `load_sequence()` (*pymasure.display.widgets.SequencerWidget* method), 56

- `local()` (`pymeasure.instruments.danfysik.Danfysik8500` method), 117
- `local()` (`pymeasure.instruments.keithley.Keithley2000` method), 129
- `LOG_10` (`pymeasure.instruments.agilent.agilentB1500.SamplingMode` attribute), 102
- `LOG_100` (`pymeasure.instruments.agilent.agilentB1500.SamplingMode` attribute), 103
- `LOG_25` (`pymeasure.instruments.agilent.agilentB1500.SamplingMode` attribute), 102
- `LOG_250` (`pymeasure.instruments.agilent.agilentB1500.SamplingMode` attribute), 103
- `LOG_50` (`pymeasure.instruments.agilent.agilentB1500.SamplingMode` attribute), 103
- `LOG_5000` (`pymeasure.instruments.agilent.agilentB1500.SamplingMode` attribute), 103
- `LOG_DOUBLE` (`pymeasure.instruments.agilent.agilentB1500.SweepMode` attribute), 102
- `log_magnitude()` (`pymeasure.instruments.agilent.Agilent8722ES` method), 69
- `log_ratio` (`pymeasure.instruments.signalrecovery.DSP7265` property), 202
- `LOG_SINGLE` (`pymeasure.instruments.agilent.agilentB1500.SamplingMode` attribute), 102
- `LogHandler` (class in `pymeasure.display.log`), 52
- `LogHandler.Emitter` (class in `pymeasure.display.log`), 52
- `LogWidget` (class in `pymeasure.display.widgets`), 55
- `low_freq` (`pymeasure.instruments.srs.SR570` property), 205
- `low_freq_out_amplitude` (`pymeasure.instruments.agilent.Agilent8257D` property), 68
- `low_freq_out_source` (`pymeasure.instruments.agilent.Agilent8257D` property), 68
- `lower_sideband_enabled` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 191
- M**
- `mag` (`pymeasure.instruments.ametek.Ametek7270` property), 104
- `mag` (`pymeasure.instruments.signalrecovery.DSP7265` property), 202
- `magnet_current` (`pymeasure.instruments.ami.AMI430` property), 106
- `magnitude` (`pymeasure.instruments.srs.SR830` property), 207
- `magnitude` (`pymeasure.instruments.srs.SR860` property), 211
- `magnitude()` (`pymeasure.instruments.agilent.Agilent8722ES` method), 69
- `ManagedImageWindow` (class in `pymeasure.display.windows`), 57
- `ManagedWindow` (class in `pymeasure.display.windows`), 57
- `ManagedWindowBase` (class in `pymeasure.display.windows`), 57
- `Manager` (class in `pymeasure.display.manager`), 52
- `MANUAL` (`pymeasure.instruments.agilent.agilentB1500.ADCMode` attribute), 102
- `MANUAL` (`pymeasure.instruments.agilent.agilentB1500.AutoManual` attribute), 102
- `MANUAL` (`pymeasure.instruments.agilent.agilentB1500.CompliancePolarity` attribute), 102
- `max_amplitude` (`pymeasure.instruments.hp.HP33120A` property), 122
- `max_current` (`pymeasure.instruments.deltaelektronika.SM7045D` property), 119
- `max_current` (`pymeasure.instruments.keithley.Keithley2400` property), 136
- `max_current` (`pymeasure.instruments.keithley.Keithley2450` property), 142
- `max_frequency` (`pymeasure.instruments.hp.HP33120A` property), 122
- `max_offset` (`pymeasure.instruments.hp.HP33120A` property), 122
- `max_resistance` (`pymeasure.instruments.keithley.Keithley2400` property), 136
- `max_resistance` (`pymeasure.instruments.keithley.Keithley2450` property), 142
- `max_voltage` (`pymeasure.instruments.deltaelektronika.SM7045D` property), 119
- `max_voltage` (`pymeasure.instruments.keithley.Keithley2400` property), 136
- `max_voltage` (`pymeasure.instruments.keithley.Keithley2450` property), 142
- `maximums` (`pymeasure.instruments.keithley.Keithley2400` property), 136
- `maximums` (`pymeasure.instruments.keithley.Keithley2450` property), 142
- `maxspeed` (`pymeasure.instruments.anaheimautomation.DPSeriesMotorControl` property), 108
- `mean_current` (`pymeasure.instruments.keithley.Keithley2400` property), 136
- `mean_current` (`pymeasure.instruments.keithley.Keithley2450` property), 142
- `mean_resistance` (`pymeasure.instruments.keithley.Keithley2400` property), 136
- `mean_resistance` (`pymeasure.instruments.keithley.Keithley2450` property), 136



- erty), 142
- mean\_voltage (pymeasure.instruments.keithley.Keithley2400 property), 136
- mean\_voltage (pymeasure.instruments.keithley.Keithley2450 property), 142
- means (pymeasure.instruments.keithley.Keithley2400 property), 136
- means (pymeasure.instruments.keithley.Keithley2450 property), 142
- meas\_mode() (pymeasure.instruments.agilent.agilentB1500.SMU method), 94
- meas\_op\_mode (pymeasure.instruments.agilent.agilentB1500.SMU property), 97
- meas\_range\_current (pymeasure.instruments.agilent.agilentB1500.SMU property), 98
- meas\_range\_current\_auto() (pymeasure.instruments.agilent.agilentB1500.SMU method), 98
- meas\_range\_voltage (pymeasure.instruments.agilent.agilentB1500.SMU property), 98
- MeasMode (class in pymeasure.instruments.agilent.agilentB1500), 102
- MeasOpMode (class in pymeasure.instruments.agilent.agilentB1500), 102
- Measurable (class in pymeasure.experiment.parameters), 45
- measure() (pymeasure.instruments.agilent.agilent4156.Agilent4156 property), 125
- measure() (pymeasure.instruments.lakeshore.LakeShore425 method), 166
- measure\_ACI (pymeasure.instruments.hp.HP3478A property), 125
- measure\_ACV (pymeasure.instruments.hp.HP3478A property), 125
- measure\_capacity() (pymeasure.instruments.attocube.anc300.Axis method), 114
- measure\_concurrent\_functions (pymeasure.instruments.keithley.Keithley2400 property), 136
- measure\_continuity() (pymeasure.instruments.keithley.Keithley2000 method), 129
- measure\_current (pymeasure.instruments.deltaelektronika.SM7045D property), 119
- measure\_current() (pymeasure.instruments.keithley.Keithley2000 method), 129
- measure\_current() (pymeasure.instruments.keithley.Keithley2400 method), 136
- measure\_current() (pymeasure.instruments.keithley.Keithley2450 method), 142
- measure\_current() (pymeasure.instruments.keithley.Keithley6517B method), 155
- measure\_DCI (pymeasure.instruments.hp.HP3478A property), 125
- measure\_DMM (pymeasure.instruments.hp.HP3478A property), 125
- measure\_diode() (pymeasure.instruments.keithley.Keithley2000 method), 129
- measure\_frequency() (pymeasure.instruments.keithley.Keithley2000 method), 129
- measure\_mode (pymeasure.instruments.anritsu.AnritsuMS9710C property), 111
- measure\_peak() (pymeasure.instruments.anritsu.AnritsuMS9710C method), 111
- measure\_period() (pymeasure.instruments.keithley.Keithley2000 method), 129
- measure\_R2W (pymeasure.instruments.hp.HP3478A property), 125
- measure\_R4W (pymeasure.instruments.hp.HP3478A property), 125
- measure\_resistance() (pymeasure.instruments.keithley.Keithley2000 method), 129
- measure\_resistance() (pymeasure.instruments.keithley.Keithley2400 method), 136
- measure\_resistance() (pymeasure.instruments.keithley.Keithley2450 method), 143
- measure\_resistance() (pymeasure.instruments.keithley.Keithley6517B method), 155
- measure\_Rext (pymeasure.instruments.hp.HP3478A property), 125
- measure\_temperature() (pymeasure.instruments.keithley.Keithley2000 method), 129
- measure\_voltage (pymeasure.instruments.deltaelektronika.SM7045D property), 119
- measure\_voltage() (pymeasure.instruments.keithley.Keithley2000 method), 129

- method*), 129
- `measure_voltage()` (*pymeasure.instruments.keithley.Keithley2400 method*), 137
- `measure_voltage()` (*pymeasure.instruments.keithley.Keithley2450 method*), 143
- `measure_voltage()` (*pymeasure.instruments.keithley.Keithley6517B method*), 155
- `measurement()` (*pymeasure.instruments.Instrument static method*), 62
- `measurement()` (*pymeasure.instruments.keysight.KeysightDSOX1102G static method*), 160
- `measurement_event_enabled` (*pymeasure.instruments.keithley.Keithley6221 property*), 150
- `measurement_events` (*pymeasure.instruments.keithley.Keithley6221 property*), 150
- `measurement_time` (*pymeasure.instruments.pendulum.cnt91.CNT91 property*), 184
- `message_waiting()` (*pymeasure.experiment.listeners.Listener method*), 42
- `min_amplitude` (*pymeasure.instruments.hp.HP33120A property*), 122
- `min_current` (*pymeasure.instruments.keithley.Keithley2400 property*), 137
- `min_current` (*pymeasure.instruments.keithley.Keithley2450 property*), 143
- `min_frequency` (*pymeasure.instruments.hp.HP33120A property*), 122
- `min_offset` (*pymeasure.instruments.hp.HP33120A property*), 122
- `min_resistance` (*pymeasure.instruments.keithley.Keithley2400 property*), 137
- `min_resistance` (*pymeasure.instruments.keithley.Keithley2450 property*), 143
- `min_voltage` (*pymeasure.instruments.keithley.Keithley2400 property*), 137
- `min_voltage` (*pymeasure.instruments.keithley.Keithley2450 property*), 143
- `minimums` (*pymeasure.instruments.keithley.Keithley2400 property*), 137
- `minimums` (*pymeasure.instruments.keithley.Keithley2450 property*), 143
- `Mock` (*class in pymeasure.instruments*), 63
- `mode` (*pymeasure.instruments.agilent.AgilentE4980 property*), 71
- `mode` (*pymeasure.instruments.attocube.anc300.Axis property*), 114
- `mode` (*pymeasure.instruments.hp.HP3478A property*), 125
- `mode` (*pymeasure.instruments.keithley.Keithley2000 property*), 129
- `modulation_degree` (*pymeasure.instruments.rohdeschwarz.sfm.Sound\_Channel property*), 196
- `modulation_enabled` (*pymeasure.instruments.rohdeschwarz.sfm.SFM property*), 191
- `modulation_enabled` (*pymeasure.instruments.rohdeschwarz.sfm.Sound\_Channel property*), 196
- `module`
- `pymeasure.display.browser`, 49
  - `pymeasure.display.curves`, 49
  - `pymeasure.display.inputs`, 50
  - `pymeasure.display.listeners`, 52
  - `pymeasure.display.log`, 52
  - `pymeasure.display.manager`, 52
  - `pymeasure.display.plotter`, 53
  - `pymeasure.display.thread`, 54
  - `pymeasure.display.widgets`, 54
  - `pymeasure.display.windows`, 57
  - `pymeasure.experiment.experiment`, 41
  - `pymeasure.experiment.listeners`, 42
  - `pymeasure.experiment.parameters`, 44
  - `pymeasure.experiment.procedure`, 43
  - `pymeasure.experiment.results`, 47
  - `pymeasure.experiment.workers`, 46
  - `pymeasure.instruments`, 59
  - `pymeasure.instruments.advantest`, 65
  - `pymeasure.instruments.advantest.advantestR3767CG`, 66
  - `pymeasure.instruments.agilent`, 66
  - `pymeasure.instruments.agilent.agilent4156`, 76
  - `pymeasure.instruments.agilent.agilentB1500`, 101
  - `pymeasure.instruments.ametek`, 103
  - `pymeasure.instruments.ami`, 105
  - `pymeasure.instruments.anaheimautomation`, 107
  - `pymeasure.instruments.anapico`, 109
  - `pymeasure.instruments.anritsu`, 109
  - `pymeasure.instruments.attocube`, 112
  - `pymeasure.instruments.bkprecision`, 115
  - `pymeasure.instruments.comedi`, 65
  - `pymeasure.instruments.danfysik`, 116
  - `pymeasure.instruments.deltaelektronika`, 119
  - `pymeasure.instruments.edwards`, 120

- `pymeasure.instruments.fluke`, 120
- `pymeasure.instruments.fwbell`, 120
- `pymeasure.instruments.hp`, 122
- `pymeasure.instruments.keithley`, 126
- `pymeasure.instruments.keysight`, 158
- `pymeasure.instruments.lakeshore`, 164
- `pymeasure.instruments.newport`, 167
- `pymeasure.instruments.ni`, 168
- `pymeasure.instruments.oxfordinstruments`, 180
- `pymeasure.instruments.parker`, 182
- `pymeasure.instruments.pendulum`, 184
- `pymeasure.instruments.razorbill`, 185
- `pymeasure.instruments.rohdeschwarz`, 186
- `pymeasure.instruments.signalrecovery`, 200
- `pymeasure.instruments.srs`, 203
- `pymeasure.instruments.tektronix`, 213
- `pymeasure.instruments.thorlabs`, 213
- `pymeasure.instruments.validators`, 63
- `pymeasure.instruments.yokogawa`, 214
- `Monitor` (class in `pymeasure.display.listeners`), 52
- `Monitor` (class in `pymeasure.experiment.listeners`), 42
- `motion_done` (`pymeasure.instruments.newport.esp300.Axis` property), 168
- `mouseMoved()` (`pymeasure.display.curves.Crosshairs` method), 49
- `move()` (`pymeasure.instruments.anaheimautomation.DPSeriesMotionController` method), 108
- `move()` (`pymeasure.instruments.attocube.anc300.Axis` method), 114
- `move()` (`pymeasure.instruments.parker.ParkerGV6` method), 183
- N**
- `new_curve()` (`pymeasure.display.widgets.ImageWidget` method), 55
- `new_curve()` (`pymeasure.display.widgets.PlotWidget` method), 55
- `new_curve()` (`pymeasure.display.widgets.TabWidget` method), 56
- `next()` (`pymeasure.display.manager.ExperimentQueue` method), 52
- `next()` (`pymeasure.display.manager.Manager` method), 53
- `nicam_additional_bits` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 191
- `nicam_audio_frequency` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 191
- `nicam_audio_volume` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 191
- `nicam_bit_error_enabled` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 191
- `nicam_bit_error_rate` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 191
- `nicam_carrier_enabled` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 191
- `nicam_carrier_frequency` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 191
- `nicam_carrier_level` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 191
- `nicam_control_bits` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 191
- `nicam_data` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 192
- `nicam_intercarrier_frequency` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 192
- `nicam_IQ_inverted` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 191
- `nicam_preemphasis_enabled` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 192
- `nicam_source` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 192
- `nicam_test_signal` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 192
- `normal_channel` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 192
- `nxds` (in module `pymeasure.instruments.edwards`), 120
- O**
- `offset` (`pymeasure.instruments.agilent.Agilent33220A` property), 83
- `offset` (`pymeasure.instruments.agilent.Agilent33500` property), 87
- `offset` (`pymeasure.instruments.agilent.agilent4156.VARD` property), 81
- `offset` (`pymeasure.instruments.hp.HP33120A` property), 122
- `offset_current` (`pymeasure.instruments.srs.SR570` property), 205

- `offset_current_enabled` (`pymeasure.instruments.srs.SR570` property), 205
  - `offset_current_sign` (`pymeasure.instruments.srs.SR570` property), 205
  - `offset_voltage` (`pymeasure.instruments.attocube.anc300.Axis` property), 114
  - `open()` (`pymeasure.instruments.keithley.Keithley2750` method), 157
  - `open_all()` (`pymeasure.instruments.keithley.Keithley2750` method), 157
  - `open_all_channels()` (`pymeasure.instruments.keithley.Keithley2700` method), 147
  - `open_channels` (`pymeasure.instruments.keithley.Keithley2700` property), 147
  - `open_file_externally()` (`pymeasure.display.windows.ManagedWindowBase` method), 58
  - `open_rows_to_columns()` (`pymeasure.instruments.keithley.Keithley2700` method), 147
  - `operation_enable_reg` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 192
  - `operation_event_enabled` (`pymeasure.instruments.keithley.Keithley6221` property), 150
  - `operation_events` (`pymeasure.instruments.keithley.Keithley6221` property), 150
  - `options` (`pymeasure.instruments.bkprecision.BKPrecision9030B` property), 115
  - `options` (`pymeasure.instruments.Instrument` property), 62
  - `options` (`pymeasure.instruments.keithley.Keithley2000` property), 130
  - `options` (`pymeasure.instruments.keithley.Keithley2260B` property), 133
  - `options` (`pymeasure.instruments.keithley.Keithley2400` property), 137
  - `options` (`pymeasure.instruments.keithley.Keithley2450` property), 143
  - `options` (`pymeasure.instruments.keithley.Keithley2600` property), 158
  - `options` (`pymeasure.instruments.keithley.Keithley2700` property), 148
  - `options` (`pymeasure.instruments.keithley.Keithley2750` property), 157
  - `options` (`pymeasure.instruments.keithley.Keithley6221` property), 150
  - `options` (`pymeasure.instruments.keithley.Keithley6517B` property), 155
  - `options` (`pymeasure.instruments.keysight.KeysightDSOX1102G` property), 161
  - `options` (`pymeasure.instruments.keysight.KeysightN5767A` property), 163
  - `output` (`pymeasure.instruments.agilent.Agilent33220A` property), 84
  - `output` (`pymeasure.instruments.agilent.Agilent33500` property), 87
  - `output` (`pymeasure.instruments.anritsu.AnritsuMG3692C` property), 110
  - `output` (`pymeasure.instruments.srs.SR510` property), 204
  - `output_1` (`pymeasure.instruments.razorbill.razorbillRP100` property), 185
  - `output_2` (`pymeasure.instruments.razorbill.razorbillRP100` property), 185
  - `output_conversion()` (`pymeasure.instruments.srs.SR830` method), 207
  - `output_load` (`pymeasure.instruments.agilent.Agilent33500` property), 87
  - `output_low_grounded` (`pymeasure.instruments.keithley.Keithley6221` property), 150
  - `output_off_state` (`pymeasure.instruments.keithley.Keithley2400` property), 137
  - `output_trigger_on_external()` (`pymeasure.instruments.keithley.Keithley2400` method), 137
  - `output_trigger_on_external()` (`pymeasure.instruments.keithley.Keithley6221` method), 151
  - `output_voltage` (`pymeasure.instruments.attocube.anc300.Axis` property), 114
  - `output_voltage` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 193
  - `outputs_enabled` (`pymeasure.instruments.ni.virtualbench.VirtualBench.PowerSupply` property), 178
- ## P
- `parallel_meas` (`pymeasure.instruments.agilent.agilentB1500.AgilentB1500` property), 94
  - `Parameter` (class in `pymeasure.experiment.parameters`), 45
  - `parameter` (`pymeasure.display.inputs.Input` property), 51
  - `parameter_DAT1` (`pymeasure.instruments.srs.SR860` property), 211
  - `parameter_DAT2` (`pymeasure.instruments.srs.SR860` property), 211



- `parameter_DAT3` (`pymeasure.instruments.srs.SR860` property), 211
- `parameter_DAT4` (`pymeasure.instruments.srs.SR860` property), 211
- `parameter_objects()` (`pymeasure.experiment.procedure.Procedure` method), 43
- `parameter_values()` (`pymeasure.experiment.procedure.Procedure` method), 43
- `parameters_are_set()` (`pymeasure.experiment.procedure.Procedure` method), 43
- `ParkerGV6` (class in `pymeasure.instruments.parker`), 183
- `parse()` (`pymeasure.experiment.results.Results` method), 47
- `parse_axis()` (`pymeasure.display.widgets.PlotFrame` method), 55
- `parse_header()` (`pymeasure.experiment.results.Results` static method), 47
- `pattern_down` (`pymeasure.instruments.attocube.anc300.Axis` property), 114
- `pattern_up` (`pymeasure.instruments.attocube.anc300.Axis` property), 114
- `pause()` (`pymeasure.instruments.agilent.agilentB1500.Agilent8722ES` method), 93
- `pause()` (`pymeasure.instruments.ami.AMI430` method), 106
- `pcolor()` (`pymeasure.experiment.experiment.Experiment` method), 42
- `peak_search` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 111
- `period` (`pymeasure.instruments.keithley.Keithley2000` property), 130
- `period_aperature` (`pymeasure.instruments.keithley.Keithley2000` property), 130
- `period_digits` (`pymeasure.instruments.keithley.Keithley2000` property), 130
- `period_reference` (`pymeasure.instruments.keithley.Keithley2000` property), 130
- `period_threshold` (`pymeasure.instruments.keithley.Keithley2000` property), 130
- `phase` (`pymeasure.instruments.agilent.Agilent33500` property), 87
- `phase` (`pymeasure.instruments.ametek.Ametek7270` property), 104
- `phase` (`pymeasure.instruments.signalrecovery.DSP7265` property), 202
- `phase` (`pymeasure.instruments.srs.SR510` property), 204
- `phase` (`pymeasure.instruments.srs.SR830` property), 207
- `phase` (`pymeasure.instruments.srs.SR860` property), 211
- `phase()` (`pymeasure.instruments.agilent.Agilent8722ES` method), 69
- `PhysicalParameter` (class in `pymeasure.experiment.parameters`), 46
- `PLC` (`pymeasure.instruments.agilent.agilentB1500.ADCMode` attribute), 102
- `plot()` (`pymeasure.experiment.experiment.Experiment` method), 42
- `plot_live()` (`pymeasure.experiment.experiment.Experiment` method), 42
- `PlotFrame` (class in `pymeasure.display.widgets`), 55
- `Plotter` (class in `pymeasure.display.plotter`), 53
- `PlotterWindow` (class in `pymeasure.display.windows`), 59
- `PlotWidget` (class in `pymeasure.display.widgets`), 55
- `points` (`pymeasure.instruments.agilent.agilent4156.VAR2` property), 81
- `polarity` (`pymeasure.instruments.danfysik.Danfysik8500` property), 117
- `position` (`pymeasure.instruments.newport.esp300.Axis` property), 168
- `position` (`pymeasure.instruments.parker.ParkerGV6` property), 183
- `position_error` (`pymeasure.instruments.parker.ParkerGV6` property), 183
- `power` (`pymeasure.instruments.agilent.Agilent8257D` property), 68
- `power` (`pymeasure.instruments.anapico.APSIN12G` property), 109
- `power` (`pymeasure.instruments.anritsu.AnritsuMG3692C` property), 110
- `power` (`pymeasure.instruments.keithley.Keithley2260B` property), 133
- `power` (`pymeasure.instruments.thorlabs.ThorlabsPM100USB` property), 214
- `preemphasis_enabled` (`pymeasure.instruments.rohdeschwarz.sfm.Sound_Channel` property), 196
- `preemphasis_time` (`pymeasure.instruments.rohdeschwarz.sfm.Sound_Channel` property), 196
- `prepare()` (`pymeasure.display.curves.BufferCurve` method), 49
- `Procedure` (class in `pymeasure.experiment.procedure`), 43
- `program_sweep()` (`pymeasure.instruments.oxfordinstruments.ITC503` method), 181
- `PrologixAdapter` (class in `pymeasure.adapters`), 34
- `pulse_dutycycle` (`pymeasure.instruments.agilent.Agilent33220A` property), 181

[property](#)), 84

[pulse\\_dutycycle](#) ([pymea-](#)  
[sure.instruments.agilent.Agilent33500](#) [prop-](#)  
[erty](#)), 87

[pulse\\_frequency](#) ([pymea-](#)  
[sure.instruments.agilent.Agilent8257D](#) [prop-](#)  
[erty](#)), 68

[pulse\\_hold](#) ([pymea-](#)[sure.instruments.agilent.Agilent33220A](#)  
[property](#)), 84

[pulse\\_hold](#) ([pymea-](#)[sure.instruments.agilent.Agilent33500](#)  
[property](#)), 87

[pulse\\_input](#) ([pymea-](#)[sure.instruments.agilent.Agilent8257D](#)  
[property](#)), 68

[pulse\\_period](#) ([pymea-](#)  
[sure.instruments.agilent.Agilent33220A](#)  
[property](#)), 84

[pulse\\_period](#) ([pymea-](#)  
[sure.instruments.agilent.Agilent33500](#) [prop-](#)  
[erty](#)), 87

[pulse\\_source](#) ([pymea-](#)  
[sure.instruments.agilent.Agilent8257D](#) [prop-](#)  
[erty](#)), 68

[pulse\\_transition](#) ([pymea-](#)  
[sure.instruments.agilent.Agilent33220A](#)  
[property](#)), 84

[pulse\\_transition](#) ([pymea-](#)  
[sure.instruments.agilent.Agilent33500](#) [prop-](#)  
[erty](#)), 87

[pulse\\_width](#) ([pymea-](#)[sure.instruments.agilent.Agilent33220A](#)  
[property](#)), 84

[pulse\\_width](#) ([pymea-](#)[sure.instruments.agilent.Agilent33500](#)  
[property](#)), 87

[pymea-](#)[sure.display.browser](#)  
[module](#), 49

[pymea-](#)[sure.display.curves](#)  
[module](#), 49

[pymea-](#)[sure.display.inputs](#)  
[module](#), 50

[pymea-](#)[sure.display.listeners](#)  
[module](#), 52

[pymea-](#)[sure.display.log](#)  
[module](#), 52

[pymea-](#)[sure.display.manager](#)  
[module](#), 52

[pymea-](#)[sure.display.plotter](#)  
[module](#), 53

[pymea-](#)[sure.display.thread](#)  
[module](#), 54

[pymea-](#)[sure.display.widgets](#)  
[module](#), 54

[pymea-](#)[sure.display.windows](#)  
[module](#), 57

[pymea-](#)[sure.experiment.experiment](#)  
[module](#), 41

[pymea-](#)[sure.experiment.listeners](#)  
[module](#), 42

[pymea-](#)[sure.experiment.parameters](#)  
[module](#), 44

[pymea-](#)[sure.experiment.procedure](#)  
[module](#), 43

[pymea-](#)[sure.experiment.results](#)  
[module](#), 47

[pymea-](#)[sure.experiment.workers](#)  
[module](#), 46

[pymea-](#)[sure.instruments](#)  
[module](#), 59

[pymea-](#)[sure.instruments.advantest](#)  
[module](#), 65

[pymea-](#)[sure.instruments.advantest.advantestR3767CG](#)  
[module](#), 66

[pymea-](#)[sure.instruments.agilent](#)  
[module](#), 66

[pymea-](#)[sure.instruments.agilent.agilent4156](#)  
[module](#), 76

[pymea-](#)[sure.instruments.agilent.agilentB1500](#)  
[module](#), 101

[pymea-](#)[sure.instruments.ametek](#)  
[module](#), 103

[pymea-](#)[sure.instruments.ami](#)  
[module](#), 105

[pymea-](#)[sure.instruments.anaheimautomation](#)  
[module](#), 107

[pymea-](#)[sure.instruments.anapico](#)  
[module](#), 109

[pymea-](#)[sure.instruments.anritsu](#)  
[module](#), 109

[pymea-](#)[sure.instruments.attocube](#)  
[module](#), 112

[pymea-](#)[sure.instruments.bkprecision](#)  
[module](#), 115

[pymea-](#)[sure.instruments.comedi](#)  
[module](#), 65

[pymea-](#)[sure.instruments.danfysik](#)  
[module](#), 116

[pymea-](#)[sure.instruments.deltaelektronika](#)  
[module](#), 119

[pymea-](#)[sure.instruments.edwards](#)  
[module](#), 120

[pymea-](#)[sure.instruments.fluke](#)  
[module](#), 120

[pymea-](#)[sure.instruments.fwbell](#)  
[module](#), 120

[pymea-](#)[sure.instruments.hp](#)  
[module](#), 122

[pymea-](#)[sure.instruments.keithley](#)  
[module](#), 126

[pymea-](#)[sure.instruments.keysight](#)  
[module](#), 158

`pymeasure.instruments.lakeshore`  
module, 164

`pymeasure.instruments.newport`  
module, 167

`pymeasure.instruments.ni`  
module, 168

`pymeasure.instruments.oxfordinstruments`  
module, 180

`pymeasure.instruments.parker`  
module, 182

`pymeasure.instruments.pendulum`  
module, 184

`pymeasure.instruments.razorbill`  
module, 185

`pymeasure.instruments.rohdeschwarz`  
module, 186

`pymeasure.instruments.signalrecovery`  
module, 200

`pymeasure.instruments.srs`  
module, 203

`pymeasure.instruments.tektronix`  
module, 213

`pymeasure.instruments.thorlabs`  
module, 213

`pymeasure.instruments.validators`  
module, 63

`pymeasure.instruments.yokogawa`  
module, 214

## Q

`QListener` (class in `pymeasure.display.listeners`), 52

`query_ac_current()` (pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter method), 171

`query_acquisition_status()` (pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope method), 176

`query_adc_setup()` (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 95

`query_analog_channel()` (pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope method), 176

`query_analog_channel_characteristics()` (pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope method), 176

`query_analog_edge_trigger()` (pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope method), 176

`query_analog_pulse_width_trigger()` (pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope method), 176

`query_arbitrary_waveform()` (pymeasure.instruments.ni.virtualbench.VirtualBench.FunctionGenerator method), 173

`query_arbitrary_waveform_gain_and_offset()` (pymeasure.instruments.ni.virtualbench.VirtualBench.FunctionGenerator method), 173

`query_current_output()` (pymeasure.instruments.ni.virtualbench.VirtualBench.PowerSupply method), 178

`query_dc_current()` (pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter method), 171

`query_dc_voltage()` (pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter method), 171

`query_enabled_analog_channels()` (pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope method), 176

`query_export_signal()` (pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalInputOutput method), 170

`query_generation_status()` (pymeasure.instruments.ni.virtualbench.VirtualBench.FunctionGenerator method), 173

`query_learn()` (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 93

`query_learn()` (pymeasure.instruments.agilent.agilentB1500.QueryLearn static method), 99

`query_learn()` (pymeasure.instruments.agilent.agilentB1500.SMU method), 97

`query_learn_header()` (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 93

`query_learn_header()` (pymeasure.instruments.agilent.agilentB1500.QueryLearn class method), 100

`query_line_configuration()` (pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalInputOutput method), 170

`query_meas_mode()` (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 94

`query_meas_op_mode()` (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 97

`query_meas_range_current_auto()` (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 97

`query_meas_ranges()` (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 97

`query_meas_settings()` (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 97

method), 94

query\_measurement() (pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter method), 171

query\_modules() (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 93

query\_sampling\_settings() (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 96

query\_series\_resistor() (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 96

query\_staircase\_sweep\_settings() (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 95

query\_standard\_waveform() (pymeasure.instruments.ni.virtualbench.VirtualBench.FunctionGenerator method), 173

query\_time\_stamp\_setting() (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 95

query\_timing() (pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope method), 176

query\_trigger\_delay() (pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope method), 177

query\_trigger\_type() (pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope method), 177

query\_voltage\_output() (pymeasure.instruments.ni.virtualbench.VirtualBench.PowerSupply method), 178

query\_waveform\_mode() (pymeasure.instruments.ni.virtualbench.VirtualBench.FunctionGenerator method), 173

QueryLearn (class in pymeasure.instruments.agilent.agilentB1500), 99

questionable\_event\_enabled (pymeasure.instruments.keithley.Keithley6221 property), 151

questionable\_event\_reg (pymeasure.instruments.rohdeschwarz.sfm.SFM property), 193

questionable\_events (pymeasure.instruments.keithley.Keithley6221 property), 151

questionable\_operation\_enable\_reg (pymeasure.instruments.rohdeschwarz.sfm.SFM property), 193

questionable\_status\_reg (pymeasure.instruments.rohdeschwarz.sfm.SFM property), 193

queue() (pymeasure.display.manager.Manager method), 53

queue() (pymeasure.display.windows.ManagedWindowBase method), 58

queue\_sequence() (pymeasure.display.widgets.SequencerWidget method), 56

quick\_range() (pymeasure.instruments.srs.SR830 method), 207

## R

R75\_out (pymeasure.instruments.rohdeschwarz.sfm.SFM property), 186

ramp() (pymeasure.instruments.ami.AMI430 method), 106

ramp\_rate\_current (pymeasure.instruments.ami.AMI430 property), 106

ramp\_rate\_field (pymeasure.instruments.ami.AMI430 property), 106

ramp\_source() (pymeasure.instruments.agilent.agilentB1500.SMU method), 98

ramps\_symmetry (pymeasure.instruments.agilent.Agilent33220A property), 84

ramps\_symmetry (pymeasure.instruments.agilent.Agilent33500 property), 87

ramp\_to\_current() (pymeasure.instruments.ami.AMI430 method), 106

ramp\_to\_current() (pymeasure.instruments.danfysik.Danfysik8500 method), 117

ramp\_to\_current() (pymeasure.instruments.deltaelektronika.SM7045D method), 119

ramp\_to\_current() (pymeasure.instruments.keithley.Keithley2400 method), 137

ramp\_to\_current() (pymeasure.instruments.keithley.Keithley2450 method), 143

ramp\_to\_current() (pymeasure.instruments.yokogawa.Yokogawa7651 method), 215

ramp\_to\_field() (pymeasure.instruments.ami.AMI430 method), 106

ramp\_to\_voltage() (pymeasure.instruments.keithley.Keithley2400 method), 137

ramp\_to\_voltage() (pymeasure.instruments.keithley.Keithley2450 method), 143

- `ramp_to_voltage()` (*pymeasure.instruments.keithley.Keithley6517B* method), 155
- `ramp_to_voltage()` (*pymeasure.instruments.yokogawa.Yokogawa7651* method), 216
- `ramp_to_zero()` (*pymeasure.instruments.deltalelektronika.SM7045D* method), 119
- `range` (*pymeasure.instruments.fwbell.FWBell5080* property), 121
- `range` (*pymeasure.instruments.hp.HP3478A* property), 125
- `range` (*pymeasure.instruments.lakeshore.LakeShore425* property), 167
- `Ranging` (class in *pymeasure.instruments.agilent.agilentB1500*), 100
- `ratio` (*pymeasure.instruments.agilent.agilent4156.VARD* property), 81
- `ratio` (*pymeasure.instruments.signalrecovery.DSP7265* property), 202
- `razorbillRP100` (class in *pymeasure.instruments.razorbill*), 185
- `read()` (*pymeasure.adapters.Adapter* method), 31
- `read()` (*pymeasure.adapters.FakeAdapter* method), 32
- `read()` (*pymeasure.adapters.PrologixAdapter* method), 35
- `read()` (*pymeasure.adapters.SerialAdapter* method), 34
- `read()` (*pymeasure.adapters.TelnetAdapter* method), 40
- `read()` (*pymeasure.adapters.VISAAAdapter* method), 37
- `read()` (*pymeasure.adapters.VXIIAdapter* method), 39
- `read()` (*pymeasure.instruments.attocube.adapters.AttocubeReferenceOutput* method), 113
- `read()` (*pymeasure.instruments.danfysik.DanfysikAdapter* method), 116
- `read()` (*pymeasure.instruments.fwbell.FWBell5080* method), 121
- `read()` (*pymeasure.instruments.Instrument* method), 62
- `read()` (*pymeasure.instruments.keysight.KeysightDSOX1102G* method), 161
- `read()` (*pymeasure.instruments.lakeshore.LakeShoreUSBAdapter* method), 164
- `read()` (*pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalInputOutput* method), 170
- `read()` (*pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalMultiexperiment.procedure.Procedure* method), 171
- `read()` (*pymeasure.instruments.parker.ParkerGV6* method), 183
- `read_analog_digital_dataframe()` (*pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope* method), 177
- `read_analog_digital_u64()` (*pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalInterfaces* method), 177
- `read_buffer()` (*pymeasure.instruments.pendulum.cnt91.CNT91* method), 185
- `read_bytes()` (*pymeasure.adapters.VISAAAdapter* method), 37
- `read_channels()` (*pymeasure.instruments.agilent.agilentB1500.AgilentB1500* method), 96
- `read_data()` (*pymeasure.instruments.agilent.agilentB1500.AgilentB1500* method), 96
- `read_memory()` (*pymeasure.instruments.anritsu.AnritsuMS9710C* method), 111
- `read_output()` (*pymeasure.instruments.ni.virtualbench.VirtualBench.PowerSupply* method), 178
- `read_raw()` (*pymeasure.adapters.VXIIAdapter* method), 39
- `read_trace()` (*pymeasure.instruments.rohdeschwarz.fsl.FSL* method), 199
- `readAI()` (in module *pymeasure.instruments.comedi*), 65
- `receive()` (*pymeasure.experiment.listeners.Listener* method), 42
- `record` (*pymeasure.display.log.LogHandler.Emitter* attribute), 52
- `Recorder` (class in *pymeasure.experiment.listeners*), 42
- `reference` (*pymeasure.instruments.signalrecovery.DSP7265* property), 202
- `reference_externalinput` (*pymeasure.instruments.srs.SR860* property), 211
- `reference_output` (*pymeasure.instruments.anapico.APSIN12G* property), 109
- `reference_phase` (*pymeasure.instruments.signalrecovery.DSP7265* property), 202
- `reference_source` (*pymeasure.instruments.srs.SR830* property), 207
- `reference_source` (*pymeasure.instruments.srs.SR860* property), 211
- `reference_triggermode` (*pymeasure.instruments.srs.SR860* property), 211
- `refresh_parameters()` (*pymeasure.experiment.procedure.Procedure* method), 43
- `reload()` (*pymeasure.experiment.results.Results* method), 47
- `remote()` (*pymeasure.instruments.danfysik.Danfysik8500* method), 167
- `remote()` (*pymeasure.instruments.keithley.Keithley2000* method), 130
- `remote_interfaces` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* method), 167



- property), 193
- remote\_local\_state (pymeasure.instruments.agilent.Agilent33220A property), 84
- remote\_lock() (pymeasure.instruments.keithley.Keithley2000 method), 130
- remove() (pymeasure.display.manager.Manager method), 53
- remove() (pymeasure.display.widgets.ImageWidget method), 55
- remove() (pymeasure.display.widgets.PlotWidget method), 55
- remove() (pymeasure.display.widgets.TabWidget method), 56
- replace\_placeholders() (in module pymeasure.experiment.results), 47
- res\_bandwidth (pymeasure.instruments.rohdeschwarz.fsl.FSL property), 200
- reset() (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 93
- reset() (pymeasure.instruments.bkprecision.BKPrecision9130B method), 115
- reset() (pymeasure.instruments.fwbell.FWBell5080 method), 121
- reset() (pymeasure.instruments.hp.HP3478A method), 126
- reset() (pymeasure.instruments.Instrument method), 62
- reset() (pymeasure.instruments.keithley.Keithley2000 method), 130
- reset() (pymeasure.instruments.keithley.Keithley2260B method), 133
- reset() (pymeasure.instruments.keithley.Keithley2400 method), 138
- reset() (pymeasure.instruments.keithley.Keithley2450 method), 143
- reset() (pymeasure.instruments.keithley.Keithley2600 method), 158
- reset() (pymeasure.instruments.keithley.Keithley2700 method), 148
- reset() (pymeasure.instruments.keithley.Keithley2750 method), 157
- reset() (pymeasure.instruments.keithley.Keithley6221 method), 151
- reset() (pymeasure.instruments.keithley.Keithley6517B method), 155
- reset() (pymeasure.instruments.keysight.KeysightDSOX1102G method), 161
- reset() (pymeasure.instruments.keysight.KeysightN5767A method), 163
- reset() (pymeasure.instruments.parker.ParkerGV6 method), 183
- reset\_buffer() (pymeasure.instruments.keithley.Keithley2000 method), 130
- reset\_buffer() (pymeasure.instruments.keithley.Keithley2400 method), 138
- reset\_buffer() (pymeasure.instruments.keithley.Keithley2450 method), 144
- reset\_buffer() (pymeasure.instruments.keithley.Keithley2700 method), 148
- reset\_buffer() (pymeasure.instruments.keithley.Keithley6221 method), 151
- reset\_buffer() (pymeasure.instruments.keithley.Keithley6517B method), 155
- reset\_instrument() (pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalInputOutput method), 170
- reset\_instrument() (pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter method), 171
- reset\_instrument() (pymeasure.instruments.ni.virtualbench.VirtualBench.FunctionGenerator method), 173
- reset\_instrument() (pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope method), 177
- reset\_instrument() (pymeasure.instruments.ni.virtualbench.VirtualBench.PowerSupply method), 178
- reset\_interlocks() (pymeasure.instruments.danfysik.Danfysik8500 method), 117
- reset\_position() (pymeasure.instruments.anaheimautomation.DPSeriesMotorController method), 108
- reset\_time() (pymeasure.instruments.Mock method), 63
- resistance (pymeasure.instruments.agilent.Agilent34410A property), 72
- resistance (pymeasure.instruments.agilent.Agilent34450A property), 75
- resistance (pymeasure.instruments.hp.HP34401A property), 123
- resistance (pymeasure.instruments.keithley.Keithley2000 property), 130
- resistance (pymeasure.instruments.keithley.Keithley2400 property), 138
- resistance (pymeasure.instruments.keithley.Keithley2450 property), 144
- resistance (pymeasure.instruments.keithley.Keithley6517B property), 156

- `resistance_4w` (`pymeasure.instruments.agilent.Agilent34410A` property), 72
- `resistance_4w` (`pymeasure.instruments.agilent.Agilent34450A` property), 75
- `resistance_4w` (`pymeasure.instruments.hp.HP34401A` property), 123
- `resistance_4w_auto_range` (`pymeasure.instruments.agilent.Agilent34450A` property), 75
- `resistance_4W_digits` (`pymeasure.instruments.keithley.Keithley2000` property), 130
- `resistance_4W_nplc` (`pymeasure.instruments.keithley.Keithley2000` property), 130
- `resistance_4w_range` (`pymeasure.instruments.agilent.Agilent34450A` property), 75
- `resistance_4W_range` (`pymeasure.instruments.keithley.Keithley2000` property), 130
- `resistance_4W_reference` (`pymeasure.instruments.keithley.Keithley2000` property), 130
- `resistance_4w_resolution` (`pymeasure.instruments.agilent.Agilent34450A` property), 75
- `resistance_auto_range` (`pymeasure.instruments.agilent.Agilent34450A` property), 75
- `resistance_digits` (`pymeasure.instruments.keithley.Keithley2000` property), 130
- `resistance_nplc` (`pymeasure.instruments.keithley.Keithley2000` property), 131
- `resistance_nplc` (`pymeasure.instruments.keithley.Keithley2400` property), 138
- `resistance_nplc` (`pymeasure.instruments.keithley.Keithley2450` property), 144
- `resistance_nplc` (`pymeasure.instruments.keithley.Keithley6517B` property), 156
- `resistance_range` (`pymeasure.instruments.agilent.Agilent34450A` property), 75
- `resistance_range` (`pymeasure.instruments.keithley.Keithley2000` property), 131
- `resistance_range` (`pymeasure.instruments.keithley.Keithley2400` property), 138
- `resistance_range` (`pymeasure.instruments.keithley.Keithley2450` property), 144
- `resistance_range` (`pymeasure.instruments.keithley.Keithley6517B` property), 156
- `resistance_reference` (`pymeasure.instruments.keithley.Keithley2000` property), 131
- `resistance_resolution` (`pymeasure.instruments.agilent.Agilent34450A` property), 75
- `resolution` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 111
- `resolution` (`pymeasure.instruments.hp.HP3478A` property), 126
- `resolution_actual` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 111
- `resolution_vbw` (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 111
- `Results` (class in `pymeasure.experiment.results`), 47
- `ResultsClass` (`pymeasure.display.widgets.ImageFrame` attribute), 55
- `ResultsClass` (`pymeasure.display.widgets.PlotFrame` attribute), 55
- `ResultsCurve` (class in `pymeasure.display.curves`), 50
- `ResultsDialog` (class in `pymeasure.display.widgets`), 56
- `ResultsImage` (class in `pymeasure.display.curves`), 50
- `resume()` (`pymeasure.display.manager.Manager` method), 53
- `rf_out_enabled` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 193
- `rf_sweep_center` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 193
- `rf_sweep_span` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 193
- `rf_sweep_start` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 193
- `rf_sweep_step` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 193
- `rf_sweep_stop` (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 193
- `right_limit` (`pymeasure.instruments.newport.esp300.Axis` property), 168

[round\\_up\(\)](#) (`pymeasure.display.curves.ResultsImage` method), 50  
[rsd](#) (`pymeasure.instruments.deltaelektronika.SM7045D` property), 120  
[run\(\)](#) (`pymeasure.display.plotter.Plotter` method), 53  
[run\(\)](#) (`pymeasure.experiment.workers.Worker` method), 46  
[run\(\)](#) (`pymeasure.instruments.keysight.KeysightDSOX11026` method), 161  
[run\(\)](#) (`pymeasure.instruments.ni.virtualbench.VirtualBenchFunctionGenerator` method), 173  
[run\(\)](#) (`pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope` method), 177  
**S**  
[sample\\_continuously\(\)](#) (`pymeasure.instruments.keithley.Keithley2400` method), 138  
[sample\\_frequency](#) (`pymeasure.instruments.srs.SR830` property), 207  
[SAMPLING](#) (`pymeasure.instruments.agilent.agilentB1500.MeasMode` attribute), 102  
[sampling\\_auto\\_abort\(\)](#) (`pymeasure.instruments.agilent.agilentB1500.AgilentB1500` method), 96  
[sampling\\_mode](#) (`pymeasure.instruments.agilent.agilentB1500.AgilentB1500` property), 96  
[sampling\\_points](#) (`pymeasure.instruments.anritsu.AnritsuMS9710C` property), 111  
[sampling\\_source\(\)](#) (`pymeasure.instruments.agilent.agilentB1500.SMU` method), 99  
[sampling\\_timing\(\)](#) (`pymeasure.instruments.agilent.agilentB1500.AgilentB1500` method), 96  
[SamplingMode](#) (class in `pymeasure.instruments.agilent.agilentB1500`), 102  
[SamplingPostOutput](#) (class in `pymeasure.instruments.agilent.agilentB1500`), 103  
[save\(\)](#) (`pymeasure.instruments.agilent.agilent4156.Agilent4156` method), 79  
[save\\_var\(\)](#) (`pymeasure.instruments.agilent.agilent4156.Agilent4156` method), 79  
[scale\\_volt](#) (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 194  
[scan\(\)](#) (`pymeasure.instruments.agilent.Agilent8722ES` method), 69  
[scan\\_continuous\(\)](#) (`pymeasure.instruments.agilent.Agilent8722ES` method), 69  
[scan\\_points](#) (`pymeasure.instruments.agilent.Agilent8722ES` property), 70  
[scan\\_single\(\)](#) (`pymeasure.instruments.agilent.Agilent8722ES` method), 70  
[ScientificInput](#) (class in `pymeasure.display.inputs`), 51  
[screen\\_layout](#) (`pymeasure.instruments.srs.SR860` property), 211  
[screenshot\(\)](#) (`pymeasure.instruments.srs.SR860` method), 211  
[SelfCalibrate\(\)](#) (`pymeasure.instruments.ni.virtualbench.VirtualBench.FunctionGenerator` method), 173  
[send\\_trigger\(\)](#) (`pymeasure.instruments.agilent.agilentB1500.AgilentB1500` method), 94  
[sensitivity](#) (`pymeasure.instruments.ametek.Ametek7270` property), 104  
[sensitivity](#) (`pymeasure.instruments.signalrecovery.DSP7265` property), 202  
[sensitivity](#) (`pymeasure.instruments.srs.SR510` property), 204  
[sensitivity](#) (`pymeasure.instruments.srs.SR570` property), 205  
[sensitivity](#) (`pymeasure.instruments.srs.SR830` property), 207  
[sensitivity](#) (`pymeasure.instruments.srs.SR860` property), 211  
[SequenceEvaluationException](#), 56  
[SequencerWidget](#) (class in `pymeasure.display.widgets`), 56  
[serial\\_baud](#) (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 194  
[serial\\_bits](#) (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 194  
[serial\\_flowcontrol](#) (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 194  
[serial\\_nr](#) (`pymeasure.instruments.attocube.anc300.Axis` property), 114  
[serial\\_parity](#) (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 194  
[serial\\_stopbits](#) (`pymeasure.instruments.rohdeschwarz.sfm.SFM` property), 194  
[SerialAdapter](#) (class in `pymeasure.adapters`), 33  
[series\\_resistance](#) (`pymeasure.instruments.agilent.agilent4156.SMU` property), 80  
[series\\_resistor](#) (`pymeasure.instruments.agilent.agilentB1500.SMU` property), 97  
[set\\_averaging\(\)](#) (`pymeasure.instruments.agilent.Agilent8722ES` method), 69



- method*), 70
- `set_buffer()` (*pymeasure.instruments.signalrecovery.DSP7265 method*), 202
- `set_channel_A_mode()` (*pymeasure.instruments.ametek.Ametek7270 method*), 104
- `set_color()` (*pymeasure.display.widgets.PlotWidget method*), 56
- `set_color()` (*pymeasure.display.widgets.TabWidget method*), 56
- `set_defaults()` (*pymeasure.adapters.PrologixAdapter method*), 36
- `set_defaults()` (*pymeasure.instruments.parker.ParkerGV6 method*), 183
- `set_differential_mode()` (*pymeasure.instruments.ametek.Ametek7270 method*), 104
- `set_fixed_frequency()` (*pymeasure.instruments.agilent.Agilent8722ES method*), 70
- `set_hardware_limits()` (*pymeasure.instruments.parker.ParkerGV6 method*), 183
- `set_IF_bandwidth()` (*pymeasure.instruments.agilent.Agilent8722ES method*), 70
- `set_output_voltage()` (*pymeasure.instruments.Mock method*), 63
- `set_parameter()` (*pymeasure.display.inputs.BooleanInput method*), 50
- `set_parameter()` (*pymeasure.display.inputs.FloatInput method*), 50
- `set_parameter()` (*pymeasure.display.inputs.Input method*), 51
- `set_parameter()` (*pymeasure.display.inputs.IntegerInput method*), 51
- `set_parameter()` (*pymeasure.display.inputs.ListInput method*), 51
- `set_parameter()` (*pymeasure.display.inputs.ScientificInput method*), 51
- `set_parameters()` (*pymeasure.display.windows.ManagedWindowBase method*), 59
- `set_parameters()` (*pymeasure.experiment.procedure.Procedure method*), 43
- `set_point` (*pymeasure.instruments.fluke.Fluke7341 property*), 120
- `set_ramp_delay()` (*pymeasure.instruments.danfysik.Danfysik8500 method*), 117
- `set_ramp_to_current()` (*pymeasure.instruments.danfysik.Danfysik8500 method*), 118
- `set_scaling()` (*pymeasure.instruments.srs.SR830 method*), 207
- `set_sequence()` (*pymeasure.instruments.danfysik.Danfysik8500 method*), 118
- `set_software_limits()` (*pymeasure.instruments.parker.ParkerGV6 method*), 183
- `set_time()` (*pymeasure.instruments.Mock method*), 63
- `set_timed_arm()` (*pymeasure.instruments.keithley.Keithley2400 method*), 138
- `set_timed_arm()` (*pymeasure.instruments.keithley.Keithley6221 method*), 151
- `set_trigger_counts()` (*pymeasure.instruments.keithley.Keithley2400 method*), 138
- `set_voltage_mode()` (*pymeasure.instruments.ametek.Ametek7270 method*), 104
- `setDifferentialMode()` (*pymeasure.instruments.signalrecovery.DSP7265 method*), 202
- `setpoint_1` (*pymeasure.instruments.lakeshore.LakeShore331 property*), 165
- `setpoint_2` (*pymeasure.instruments.lakeshore.LakeShore331 property*), 165
- `setting()` (*pymeasure.instruments.Instrument static method*), 62
- `setting()` (*pymeasure.instruments.keysight.KeysightDSOX1102G static method*), 161
- `setup_plot()` (*pymeasure.display.plotter.Plotter method*), 53
- `SFM` (*class in pymeasure.instruments.rohdeschwarz.sfm*), 186
- `shape` (*pymeasure.instruments.agilent.Agilent33220A property*), 84
- `shape` (*pymeasure.instruments.agilent.Agilent33500 property*), 87
- `shape` (*pymeasure.instruments.hp.HP33120A property*), 122
- `shutdown()` (*pymeasure.experiment.procedure.Procedure method*), 44
- `shutdown()` (*pymeasure.experiment.workers.Worker method*), 47
- `shutdown()` (*pymeasure.instruments.agilent.Agilent8257D method*), 68
- `shutdown()` (*pymeasure.instruments.ametek.Ametek7270*

method), 104

shutdown() (pymeasure.instruments.ami.AMI430 method), 106

shutdown() (pymeasure.instruments.anritsu.AnritsuMG3692C method), 110

shutdown() (pymeasure.instruments.bkprecision.BKPrecision4300 method), 115

shutdown() (pymeasure.instruments.deltaelektronika.SM7045 method), 120

shutdown() (pymeasure.instruments.hp.HP3478A method), 126

shutdown() (pymeasure.instruments.Instrument method), 63

shutdown() (pymeasure.instruments.keithley.Keithley2000 method), 131

shutdown() (pymeasure.instruments.keithley.Keithley2260B method), 133

shutdown() (pymeasure.instruments.keithley.Keithley2400 method), 138

shutdown() (pymeasure.instruments.keithley.Keithley2450 method), 144

shutdown() (pymeasure.instruments.keithley.Keithley2600 method), 158

shutdown() (pymeasure.instruments.keithley.Keithley2700 method), 148

shutdown() (pymeasure.instruments.keithley.Keithley2750 method), 157

shutdown() (pymeasure.instruments.keithley.Keithley6221 method), 151

shutdown() (pymeasure.instruments.keithley.Keithley6517B method), 156

shutdown() (pymeasure.instruments.keysight.KeysightDSOX1102G method), 161

shutdown() (pymeasure.instruments.keysight.KeysightN5761A method), 163

shutdown() (pymeasure.instruments.newport.ESP300 method), 167

shutdown() (pymeasure.instruments.ni.virtualbench.VirtualBench method), 180

shutdown() (pymeasure.instruments.ni.virtualbench.VirtualBench method), 170

shutdown() (pymeasure.instruments.ni.virtualbench.VirtualBench method), 172

shutdown() (pymeasure.instruments.ni.virtualbench.VirtualBench method), 174

shutdown() (pymeasure.instruments.ni.virtualbench.VirtualBench method), 177

shutdown() (pymeasure.instruments.ni.virtualbench.VirtualBench method), 178

shutdown() (pymeasure.instruments.signalrecovery.DSP7265 method), 203

shutdown() (pymeasure.instruments.yokogawa.Yokogawa7651 method), 216

signal\_inverted (pymeasure.instruments.srs.SR570 property), 205

sine\_amplitudepreset1 (pymeasure.instruments.srs.SR860 property), 212

sine\_amplitudepreset2 (pymeasure.instruments.srs.SR860 property), 212

sine\_amplitudepreset3 (pymeasure.instruments.srs.SR860 property), 212

sine\_amplitudepreset4 (pymeasure.instruments.srs.SR860 property), 212

sine\_dclevelpreset1 (pymeasure.instruments.srs.SR860 property), 212

sine\_dclevelpreset2 (pymeasure.instruments.srs.SR860 property), 212

sine\_dclevelpreset3 (pymeasure.instruments.srs.SR860 property), 212

sine\_dclevelpreset4 (pymeasure.instruments.srs.SR860 property), 212

sine\_voltage (pymeasure.instruments.srs.SR830 property), 207

sine\_voltage (pymeasure.instruments.srs.SR860 property), 212

single() (pymeasure.instruments.keysight.KeysightDSOX1102G method), 161

single\_sweep() (pymeasure.instruments.anritsu.AnritsuMS9710C method), 111

single\_sweep() (pymeasure.instruments.rohdeschwarz.fsl.FSL method), 200

blew\_rate (pymeasure.instruments.danfysik.Danfysik8500 property), 118

blew\_rate\_1 (pymeasure.instruments.razorbill.razorbillRP100 property), 185

blew\_rate\_2 (pymeasure.instruments.razorbill.razorbillRP100 property), 185

slope (pymeasure.instruments.ametek.Ametek7270 property), 104

SRope (pymeasure.instruments.signalrecovery.DSP7265 property), 203

SRock (pymeasure.instruments.thorlabs.ThorlabsPro8000 property), 214

SM7045DigitalMultimeter (class in pymeasure.instruments.deltaelektronika), 119

SMU.FunctionGenerator (class in pymeasure.instruments.agilent.agilent4156), 79

SMU.MixedSignalOscilloscope (class in pymeasure.instruments.agilent.agilentB1500), 97

SMU\_MEASUREMENT (class in pymeasure.instruments.agilent.agilentB1500.WaitTimeType attribute), 103

smu\_names (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 property), 93

smu\_references (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 property), 93

- property), 93
- SMU\_SOURCE (pymeasure.instruments.agilent.agilentB1500.  
attribute), 103
- SMUCurrentRanging (class in pymeasure.instruments.agilent.agilentB1500), 100
- SMUVoltageRanging (class in pymeasure.instruments.agilent.agilentB1500), 101
- snap() (pymeasure.instruments.srs.SR830 method), 208
- snap() (pymeasure.instruments.srs.SR860 method), 212
- Sound\_Channel (class in pymeasure.instruments.rohdeschwarz.sfm), 196
- sound\_mode (pymeasure.instruments.rohdeschwarz.sfm.SFM  
property), 194
- source\_auto\_range (pymeasure.instruments.keithley.Keithley6221  
property), 151
- source\_compliance (pymeasure.instruments.keithley.Keithley6221  
property), 151
- source\_current (pymeasure.instruments.keithley.Keithley2400  
property), 138
- source\_current (pymeasure.instruments.keithley.Keithley2450  
property), 144
- source\_current (pymeasure.instruments.keithley.Keithley6221  
property), 151
- source\_current (pymeasure.instruments.yokogawa.Yokogawa7651  
property), 216
- source\_current\_delay (pymeasure.instruments.keithley.Keithley2450  
property), 144
- source\_current\_delay\_auto (pymeasure.instruments.keithley.Keithley2450  
property), 144
- source\_current\_range (pymeasure.instruments.keithley.Keithley2400  
property), 138
- source\_current\_range (pymeasure.instruments.keithley.Keithley2450  
property), 144
- source\_current\_range (pymeasure.instruments.yokogawa.Yokogawa7651  
property), 216
- source\_current\_resistance\_limit (pymeasure.instruments.keithley.Keithley6517B  
property), 156
- source\_delay (pymeasure.instruments.keithley.Keithley2400  
property), 138
- source\_delay (pymeasure.instruments.keithley.Keithley6221  
property), 151
- source\_delay\_auto (pymeasure.instruments.keithley.Keithley2400  
property), 138
- source\_enabled (pymeasure.instruments.bkprecision.BKPrecision9130B  
property), 115
- source\_enabled (pymeasure.instruments.keithley.Keithley2400  
property), 138
- source\_enabled (pymeasure.instruments.keithley.Keithley2450  
property), 144
- source\_enabled (pymeasure.instruments.keithley.Keithley6221  
property), 151
- source\_enabled (pymeasure.instruments.keithley.Keithley6517B  
property), 156
- source\_enabled (pymeasure.instruments.yokogawa.Yokogawa7651  
property), 216
- source\_enabled (pymeasure.instruments.yokogawa.YokogawaGS200  
property), 216
- source\_level (pymeasure.instruments.yokogawa.YokogawaGS200  
property), 216
- source\_mode (pymeasure.instruments.keithley.Keithley2400  
property), 138
- source\_mode (pymeasure.instruments.keithley.Keithley2450  
property), 144
- source\_mode (pymeasure.instruments.yokogawa.Yokogawa7651  
property), 216
- source\_mode (pymeasure.instruments.yokogawa.YokogawaGS200  
property), 216
- source\_range (pymeasure.instruments.keithley.Keithley6221  
property), 151
- source\_range (pymeasure.instruments.yokogawa.YokogawaGS200  
property), 216
- source\_voltage (pymeasure.instruments.keithley.Keithley2400  
property), 138
- source\_voltage (pymeasure.instruments.keithley.Keithley2450  
property), 144
- source\_voltage (pymeasure.instruments.keithley.Keithley6517B  
property), 156
- source\_voltage (pymeasure.instruments.yokogawa.Yokogawa7651  
property), 216

<code>source_voltage_delay</code>	( <code>pymeasure.instruments.keithley.Keithley2450</code> property), 144	<code>standard_devs</code>	( <code>pymeasure.instruments.keithley.Keithley2450</code> property), 144
<code>source_voltage_delay_auto</code>	( <code>pymeasure.instruments.keithley.Keithley2450</code> property), 144	<code>standard_event_enabled</code>	( <code>pymeasure.instruments.keithley.Keithley6221</code> property), 152
<code>source_voltage_range</code>	( <code>pymeasure.instruments.keithley.Keithley2400</code> property), 139	<code>standard_events</code>	( <code>pymeasure.instruments.keithley.Keithley6221</code> property), 152
<code>source_voltage_range</code>	( <code>pymeasure.instruments.keithley.Keithley2450</code> property), 144	<code>start</code>	( <code>pymeasure.instruments.agilent.agilent4156.VARX</code> property), 81
<code>source_voltage_range</code>	( <code>pymeasure.instruments.keithley.Keithley6517B</code> property), 156	<code>START</code>	( <code>pymeasure.instruments.agilent.agilentB1500.StaircaseSweepPostOutput</code> attribute), 103
<code>source_voltage_range</code>	( <code>pymeasure.instruments.yokogawa.Yokogawa7651</code> property), 216	<code>start()</code>	( <code>pymeasure.experiment.experiment.Experiment</code> method), 42
<code>spacing</code>	( <code>pymeasure.instruments.agilent.agilent4156.VAR1</code> property), 80	<code>start_buffer()</code>	( <code>pymeasure.instruments.keithley.Keithley2000</code> method), 131
<code>span_frequency</code>	( <code>pymeasure.instruments.advantest.advantestR3767CG.AdvantestR3767CG</code> property), 66	<code>start_buffer()</code>	( <code>pymeasure.instruments.keithley.Keithley2400</code> method), 139
<code>special_channel</code>	( <code>pymeasure.instruments.rohdeschwarz.sfm.SFM</code> property), 195	<code>start_buffer()</code>	( <code>pymeasure.instruments.keithley.Keithley2450</code> method), 144
<code>SPOT</code>	( <code>pymeasure.instruments.agilent.agilentB1500.MeasMode</code> attribute), 102	<code>start_buffer()</code>	( <code>pymeasure.instruments.keithley.Keithley2700</code> method), 148
<code>square_dutycycle</code>	( <code>pymeasure.instruments.agilent.Agilent33220A</code> property), 84	<code>start_buffer()</code>	( <code>pymeasure.instruments.keithley.Keithley6221</code> method), 152
<code>square_dutycycle</code>	( <code>pymeasure.instruments.agilent.Agilent33500</code> property), 87	<code>start_buffer()</code>	( <code>pymeasure.instruments.keithley.Keithley6517B</code> method), 156
<code>SR510</code>	(class in <code>pymeasure.instruments.srs</code> ), 204	<code>start_buffer()</code>	( <code>pymeasure.instruments.signalrecovery.DSP7265</code> method), 203
<code>SR570</code>	(class in <code>pymeasure.instruments.srs</code> ), 204	<code>start_frequency</code>	( <code>pymeasure.instruments.advantest.advantestR3767CG.AdvantestR3767CG</code> property), 66
<code>SR830</code>	(class in <code>pymeasure.instruments.srs</code> ), 205	<code>start_frequency</code>	( <code>pymeasure.instruments.agilent.Agilent8257D</code> property), 68
<code>SR860</code>	(class in <code>pymeasure.instruments.srs</code> ), 208	<code>start_frequency</code>	( <code>pymeasure.instruments.agilent.Agilent8722ES</code> property), 70
<code>srq_event_enabled</code>	( <code>pymeasure.instruments.keithley.Keithley6221</code> property), 151	<code>start_frequency</code>	( <code>pymeasure.instruments.agilent.AgilentE4408B</code> property), 70
<code>SRQ_mask</code>	( <code>pymeasure.instruments.hp.HP3478A</code> property), 123	<code>start_power</code>	( <code>pymeasure.instruments.agilent.Agilent8257D</code> property), 68
<code>STAIRCASE_SWEEP</code>	( <code>pymeasure.instruments.agilent.agilentB1500.MeasMode</code> attribute), 102	<code>start_ramp()</code>	( <code>pymeasure.instruments.danfysik.Danfysik8500</code> method), 118
<code>staircase_sweep_source()</code>	( <code>pymeasure.instruments.agilent.agilentB1500.SMU</code> method), 98	<code>start_sequence()</code>	( <code>pymeasure.instruments.keithley.Keithley2400</code> property), 139
<code>StaircaseSweepPostOutput</code>	(class in <code>pymeasure.instruments.agilent.agilentB1500</code> ), 103		



- sure.instruments.danfysik.Danfysik8500* method), 118
- `start_step_sweep()` (*pymea-  
sure.instruments.agilent.Agilent8257D*  
method), 68
- `startup()` (*pymea-  
sure.instruments.procedure.Procedure*  
method), 44
- `startup()` (*pymea-  
sure.instruments.procedure.UnknownProcedure*  
method), 44
- `state` (*pymea-  
sure.instruments.ami.AMI430* property), 106
- `status` (*pymea-  
sure.instruments.agilent.agilentB1500.SMU*  
property), 97
- `status` (*pymea-  
sure.instruments.bkprecision.BKPrecision9130B*  
property), 115
- `status` (*pymea-  
sure.instruments.danfysik.Danfysik8500*  
property), 118
- `status` (*pymea-  
sure.instruments.hp.HP3478A* property), 126
- `status` (*pymea-  
sure.instruments.Instrument* property), 63
- `status` (*pymea-  
sure.instruments.keithley.Keithley2000*  
property), 131
- `status` (*pymea-  
sure.instruments.keithley.Keithley2260B*  
property), 133
- `status` (*pymea-  
sure.instruments.keithley.Keithley2450*  
property), 145
- `status` (*pymea-  
sure.instruments.keithley.Keithley2600*  
property), 158
- `status` (*pymea-  
sure.instruments.keithley.Keithley2700*  
property), 148
- `status` (*pymea-  
sure.instruments.keithley.Keithley2750*  
property), 157
- `status` (*pymea-  
sure.instruments.keithley.Keithley6221*  
property), 152
- `status` (*pymea-  
sure.instruments.keithley.Keithley6517B*  
property), 156
- `status` (*pymea-  
sure.instruments.keysight.KeysightDSOX1102G*  
property), 161
- `status` (*pymea-  
sure.instruments.keysight.KeysightN5767A*  
property), 163
- `status` (*pymea-  
sure.instruments.parker.ParkerGV6* prop-  
erty), 183
- `status` (*pymea-  
sure.instruments.srs.SR510* property), 204
- `status()` (*pymea-  
sure.instruments.keithley.Keithley2400*  
method), 139
- `status_hex` (*pymea-  
sure.instruments.danfysik.Danfysik8500*  
property), 118
- `status_info_shown` (*pymea-  
sure.instruments.rohdeschwarz.sfm.SFM*  
property), 195
- `status_preset()` (*pymea-  
sure.instruments.rohdeschwarz.sfm.SFM*  
method), 195
- `status_reg` (*pymea-  
sure.instruments.rohdeschwarz.sfm.SFM*  
property), 195
- `std_current` (*pymea-  
sure.instruments.keithley.Keithley2400*  
property), 139
- `std_current` (*pymea-  
sure.instruments.keithley.Keithley2450*  
property), 145
- `std_resistance` (*pymea-  
sure.instruments.keithley.Keithley2400* prop-  
erty), 139
- `std_resistance` (*pymea-  
sure.instruments.keithley.Keithley2450* prop-  
erty), 145
- `std_voltage` (*pymea-  
sure.instruments.keithley.Keithley2400*  
property), 139
- `std_voltage` (*pymea-  
sure.instruments.keithley.Keithley2450*  
property), 145
- `step` (*pymea-  
sure.instruments.agilent.agilent4156.VARX*  
property), 81
- `step_points` (*pymea-  
sure.instruments.agilent.Agilent8257D*  
property), 68
- `step_position` (*pymea-  
sure.instruments.anaheimautomation.DPSeriesMotorController*  
property), 108
- `stepd` (*pymea-  
sure.instruments.attocube.anc300.Axis*  
property), 115
- `steps_to_absolute()` (*pymea-  
sure.instruments.anaheimautomation.DPSeriesMotorController*  
method), 109
- `stepu` (*pymea-  
sure.instruments.attocube.anc300.Axis*  
property), 115
- `stop` (*pymea-  
sure.instruments.agilent.agilent4156.VARX*  
property), 82
- `STOP` (*pymea-  
sure.instruments.agilent.agilentB1500.StaircaseSweepPostOut*  
attribute), 103
- `stop()` (*pymea-  
sure.instruments.experiment.listeners.Recorder*  
method), 43
- `stop()` (*pymea-  
sure.instruments.agilent.agilent4156.Agilent4156*  
method), 79
- `stop()` (*pymea-  
sure.instruments.anaheimautomation.DPSeriesMotorContro*  
method), 109
- `stop()` (*pymea-  
sure.instruments.attocube.anc300.Axis*  
method), 115
- `stop()` (*pymea-  
sure.instruments.keysight.KeysightDSOX1102G*  
method), 161
- `stop()` (*pymea-  
sure.instruments.ni.virtualbench.VirtualBench.FunctionGen*  
method), 174
- `stop()` (*pymea-  
sure.instruments.ni.virtualbench.VirtualBench.MixedSignal*  
method), 177
- `stop()` (*pymea-  
sure.instruments.parker.ParkerGV6*  
method), 183
- `stop_all()` (*pymea-  
sure.instruments.attocube.anc300.ANC300Controller*  
method), 114
- `stop_buffer()` (*pymea-  
sure.instruments.keithley.Keithley2000*  
method), 139

- method), 131
- stop\_buffer() (pymeasure.instruments.keithley.Keithley2400 method), 139
- stop\_buffer() (pymeasure.instruments.keithley.Keithley2450 method), 145
- stop\_buffer() (pymeasure.instruments.keithley.Keithley2700 method), 148
- stop\_buffer() (pymeasure.instruments.keithley.Keithley6221 method), 152
- stop\_buffer() (pymeasure.instruments.keithley.Keithley6517B method), 156
- stop\_frequency (pymeasure.instruments.advantest.advantestR3767CG.AdvantestR3767CG (class in pymeasure.instruments.agilent.Agilent8257D property), 66
- stop\_frequency (pymeasure.instruments.agilent.Agilent8257D property), 68
- stop\_frequency (pymeasure.instruments.agilent.Agilent8722ES property), 70
- stop\_frequency (pymeasure.instruments.agilent.AgilentE4408B property), 70
- stop\_power (pymeasure.instruments.agilent.Agilent8257D property), 69
- stop\_ramp() (pymeasure.instruments.danfysik.Danfysik8500 method), 118
- stop\_sequence() (pymeasure.instruments.danfysik.Danfysik8500 method), 118
- stop\_step\_sweep() (pymeasure.instruments.agilent.Agilent8257D method), 69
- StoppableQThread (class in pymeasure.display.thread), 54
- StringInput (class in pymeasure.display.inputs), 51
- strip\_chart\_dat1 (pymeasure.instruments.srs.SR860 property), 212
- strip\_chart\_dat2 (pymeasure.instruments.srs.SR860 property), 212
- strip\_chart\_dat3 (pymeasure.instruments.srs.SR860 property), 213
- strip\_chart\_dat4 (pymeasure.instruments.srs.SR860 property), 213
- subsystem\_info (pymeasure.instruments.rohdeschwarz.sfm.SFM property), 195
- supply\_current (pymeasure.instruments.ami.AMI430 property), 106
- sweep\_auto\_abort() (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 96
- sweep\_status (pymeasure.instruments.oxfordinstruments.ITC503 property), 181
- sweep\_table (pymeasure.instruments.oxfordinstruments.ITC503 property), 181
- sweep\_time (pymeasure.instruments.agilent.Agilent8722ES property), 70
- sweep\_time (pymeasure.instruments.agilent.AgilentE4408B property), 70
- sweep\_time (pymeasure.instruments.rohdeschwarz.fsl.FSL property), 200
- sweep\_timing() (pymeasure.instruments.agilent.agilentB1500.AgilentB1500 method), 95
- SweepMode (class in pymeasure.instruments.agilent.agilentB1500), 102
- sync\_sequence() (pymeasure.instruments.danfysik.Danfysik8500 method), 118
- synchronous\_sweep\_source() (pymeasure.instruments.agilent.agilentB1500.SMU method), 99
- system\_number (pymeasure.instruments.rohdeschwarz.sfm.SFM property), 195
- system\_setup (pymeasure.instruments.keysight.KeysightDSOX1102G property), 161
- ## T
- TabWidget (class in pymeasure.display.widgets), 56
- target\_current (pymeasure.instruments.ami.AMI430 property), 106
- target\_field (pymeasure.instruments.ami.AMI430 property), 106
- TDS2000 (class in pymeasure.instruments.tektronix), 213
- TEDSetTemperature (pymeasure.instruments.thorlabs.ThorlabsPro8000 property), 214
- TEDStatus (pymeasure.instruments.thorlabs.ThorlabsPro8000 property), 214
- TelnetAdapter (class in pymeasure.adapters), 39
- temperature (pymeasure.instruments.agilent.Agilent34450A property), 76
- temperature (pymeasure.instruments.fluke.Fluke7341 property), 120
- temperature (pymeasure.instruments.keithley.Keithley2000 property), 131
- temperature\_1 (pymeasure.instruments.oxfordinstruments.ITC503 property), 181

`temperature_2` (`pymea-  
sure.instruments.oxfordinstruments.ITC503  
property`), 181  
`temperature_3` (`pymea-  
sure.instruments.oxfordinstruments.ITC503  
property`), 181  
`temperature_A` (`pymea-  
sure.instruments.lakeshore.LakeShore331  
property`), 166  
`temperature_B` (`pymea-  
sure.instruments.lakeshore.LakeShore331  
property`), 166  
`temperature_digits` (`pymea-  
sure.instruments.keithley.Keithley2000 prop-  
erty`), 131  
`temperature_error` (`pymea-  
sure.instruments.oxfordinstruments.ITC503  
property`), 182  
`temperature_nplc` (`pymea-  
sure.instruments.keithley.Keithley2000 prop-  
erty`), 131  
`temperature_reference` (`pymea-  
sure.instruments.keithley.Keithley2000 prop-  
erty`), 131  
`temperature_setpoint` (`pymea-  
sure.instruments.oxfordinstruments.ITC503  
property`), 182  
`text_enabled` (`pymea-  
sure.instruments.keithley.Keithley2700 prop-  
erty`), 148  
`theta` (`pymea-  
sure.instruments.srs.SR830 property`), 208  
`theta` (`pymea-  
sure.instruments.srs.SR860 property`), 213  
`ThorlabsPM100USB` (class in `pymea-  
sure.instruments.thorlabs`), 214  
`ThorlabsPro8000` (class in `pymea-  
sure.instruments.thorlabs`), 214  
`TIME` (`pymea-  
sure.instruments.agilent.agilentB1500.ADCMod-  
attribute`), 102  
`time` (`pymea-  
sure.instruments.Mock property`), 63  
`time` (`pymea-  
sure.instruments.rohdeschwarz.sfm.SFM  
property`), 195  
`time_constant` (`pymea-  
sure.instruments.ametek.Ametek7270 prop-  
erty`), 105  
`time_constant` (`pymea-  
sure.instruments.signalrecovery.DSP7265  
property`), 203  
`time_constant` (`pymea-  
sure.instruments.srs.SR510  
property`), 204  
`time_constant` (`pymea-  
sure.instruments.srs.SR830  
property`), 208  
`time_constant` (`pymea-  
sure.instruments.srs.SR860  
property`), 213  
`time_stamp` (`pymea-  
sure.instruments.agilent.agilentB1500.AgilentB1500  
property`), 95  
`timebase` (`pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
property`), 161  
`timebase` (`pymea-  
sure.instruments.srs.SR860 property`),  
213  
`timebase_mode` (`pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
property`), 162  
`timebase_offset` (`pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
property`), 162  
`timebase_range` (`pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
property`), 162  
`timebase_scale` (`pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
property`), 162  
`timebase_setup()` (`pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
method`), 162  
`to_dict()` (`pymea-  
sure.instruments.agilent.agilentB1500.QueryLearn  
static method`), 100  
`trace()` (`pymea-  
sure.instruments.agilent.AgilentE4408B  
method`), 70  
`trace_1` (`pymea-  
sure.instruments.advantest.advantestR3767CG.AdvantestR  
property`), 66  
`trace_df()` (`pymea-  
sure.instruments.agilent.AgilentE4408B  
method`), 70  
`trace_marker` (`pymea-  
sure.instruments.anritsu.AnritsuMS9710C  
property`), 111  
`trace_marker_center` (`pymea-  
sure.instruments.anritsu.AnritsuMS9710C  
property`), 112  
`trace_mode` (`pymea-  
sure.instruments.rohdeschwarz.fsl.FSL  
property`), 200  
`tracking` (`pymea-  
sure.instruments.ni.virtualbench.VirtualBench.PowerSup  
property`), 178  
`triad()` (`pymea-  
sure.instruments.keithley.Keithley2400  
method`), 139  
`triad()` (`pymea-  
sure.instruments.keithley.Keithley2450  
method`), 145  
`triad()` (`pymea-  
sure.instruments.keithley.Keithley2700  
method`), 148  
`triad()` (`pymea-  
sure.instruments.keithley.Keithley6221  
method`), 152  
`trigger` (`pymea-  
sure.instruments.hp.HP3478A prop-  
erty`), 126  
`trigger()` (`pymea-  
sure.instruments.agilent.Agilent33220A  
method`), 84  
`trigger()` (`pymea-  
sure.instruments.agilent.Agilent33500  
method`), 87  
`trigger()` (`pymea-  
sure.instruments.keithley.Keithley2400  
method`), 139

- `trigger()` (*pymeasure.instruments.keithley.Keithley2450* method), 145
- `trigger()` (*pymeasure.instruments.keithley.Keithley6221* method), 152
- `trigger()` (*pymeasure.instruments.keithley.Keithley6517B* method), 156
- `trigger_count` (*pymeasure.instruments.keithley.Keithley2000* property), 131
- `trigger_count` (*pymeasure.instruments.keithley.Keithley2400* property), 139
- `trigger_delay` (*pymeasure.instruments.keithley.Keithley2000* property), 131
- `trigger_delay` (*pymeasure.instruments.keithley.Keithley2400* property), 139
- `trigger_immediately()` (*pymeasure.instruments.keithley.Keithley2400* method), 139
- `trigger_immediately()` (*pymeasure.instruments.keithley.Keithley6221* method), 152
- `trigger_immediately()` (*pymeasure.instruments.keithley.Keithley6517B* method), 156
- `trigger_on_bus()` (*pymeasure.instruments.keithley.Keithley2400* method), 139
- `trigger_on_bus()` (*pymeasure.instruments.keithley.Keithley6221* method), 152
- `trigger_on_bus()` (*pymeasure.instruments.keithley.Keithley6517B* method), 156
- `trigger_on_external()` (*pymeasure.instruments.keithley.Keithley2400* method), 139
- `trigger_on_external()` (*pymeasure.instruments.keithley.Keithley6221* method), 152
- `trigger_ramp_to_level()` (*pymeasure.instruments.yokogawa.YokogawaGS200* method), 217
- `trigger_source` (*pymeasure.instruments.agilent.Agilent33220A* property), 84
- `trigger_source` (*pymeasure.instruments.agilent.Agilent33500* property), 87
- `trigger_source` (*pymeasure.instruments.agilent.AgilentE4980* property), 72
- `trigger_state` (*pymeasure.instruments.agilent.Agilent33220A* property), 84
- `tristate_lines()` (*pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalInputOutput* method), 170
- `TV_country` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 186
- `TV_standard` (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 187
- ## U
- `unblank_front()` (*pymeasure.instruments.srs.SR570* method), 205
- `unique_filename()` (in module *pymeasure.experiment.results*), 48
- `unit` (*pymeasure.instruments.fluke.Fluke7341* property), 120
- `unit` (*pymeasure.instruments.lakeshore.LakeShore425* property), 167
- `units` (*pymeasure.instruments.fwbell.FWBell5080* property), 121
- `units` (*pymeasure.instruments.newport.esp300.Axis* property), 168
- `UnknownProcedure` (class in *pymeasure.experiment.procedure*), 44
- `update()` (*pymeasure.display.curves.Crosshairs* method), 50
- `update_data()` (*pymeasure.display.curves.ResultsCurve* method), 50
- `update_estimates()` (*pymeasure.display.widgets.EstimatorWidget* method), 55
- `update_line()` (*pymeasure.experiment.experiment.Experiment* method), 42
- `update_parameter()` (*pymeasure.display.inputs.Input* method), 51
- `update_pcolor()` (*pymeasure.experiment.experiment.Experiment* method), 42
- `update_plot()` (*pymeasure.experiment.experiment.Experiment* method), 42
- `update_status()` (*pymeasure.experiment.workers.Worker* method), 47
- `use_absolute_position()` (*pymeasure.instruments.parker.ParkerGV6* method), 183
- `use_external_source` (*pymeasure.instruments.rohdeschwarz.sfm.Sound\_Channel* property), 197



- [use\\_front\\_terminals\(\)](#) (*pymeasure.instruments.keithley.Keithley2400* method), 139  
[use\\_front\\_terminals\(\)](#) (*pymeasure.instruments.keithley.Keithley2450* method), 145  
[use\\_rear\\_terminals\(\)](#) (*pymeasure.instruments.keithley.Keithley2400* method), 139  
[use\\_rear\\_terminals\(\)](#) (*pymeasure.instruments.keithley.Keithley2450* method), 145  
[use\\_relative\\_position\(\)](#) (*pymeasure.instruments.parker.ParkerGV6* method), 184
- ## V
- [validate\\_auto\\_range\\_terminal\(\)](#) (*pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter* method), 172  
[validate\\_channel\(\)](#) (*pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope* method), 177  
[validate\\_channel\(\)](#) (*pymeasure.instruments.ni.virtualbench.VirtualBench.PowerSupply* method), 178  
[validate\\_dmm\\_function\(\)](#) (*pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter* method), 172  
[validate\\_lines\(\)](#) (*pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalInputOutput* method), 170  
[validate\\_range\(\)](#) (*pymeasure.instruments.ni.virtualbench.VirtualBench.DigitalMultimeter* static method), 172  
[validate\\_trigger\\_instance\(\)](#) (*pymeasure.instruments.ni.virtualbench.VirtualBench.MixedSignalOscilloscope* static method), 178  
[values\(\)](#) (*pymeasure.adapters.Adapter* method), 31  
[values\(\)](#) (*pymeasure.adapters.FakeAdapter* method), 32  
[values\(\)](#) (*pymeasure.adapters.PrologixAdapter* method), 36  
[values\(\)](#) (*pymeasure.adapters.SerialAdapter* method), 34  
[values\(\)](#) (*pymeasure.adapters.TelnetAdapter* method), 40  
[values\(\)](#) (*pymeasure.adapters.VISAAdapter* method), 37  
[values\(\)](#) (*pymeasure.adapters.VXIIAdapter* method), 39  
[values\(\)](#) (*pymeasure.instruments.anaheimautomation.DPSeriesCloningAverage* method), 109  
[values\(\)](#) (*pymeasure.instruments.fwbell.FWBell5080* method), 121  
[values\(\)](#) (*pymeasure.instruments.instrument.Instrument* method), 63  
[values\(\)](#) (*pymeasure.instruments.keysight.KeysightDSOX1102G* method), 162  
[values\(\)](#) (*pymeasure.instruments.lakeshore.LakeShoreUSBAdapter* method), 164  
[values\(\)](#) (*pymeasure.instruments.rohdeschwarz.sfm.Sound\_Channel* method), 197  
[VAR1](#) (class in *pymeasure.instruments.agilent.agilent4156*), 80  
[VAR2](#) (class in *pymeasure.instruments.agilent.agilent4156*), 81  
[VARD](#) (class in *pymeasure.instruments.agilent.agilent4156*), 81  
[VARX](#) (class in *pymeasure.instruments.agilent.agilent4156*), 81  
[VectorParameter](#) (class in *pymeasure.experiment.parameters*), 46  
[version\(\)](#) (*pymeasure.instruments.attocube.anc300.ANC300Controller* property), 114  
[video\\_bandwidth](#) (*pymeasure.instruments.rohdeschwarz.fsl.FSL* property), 200  
[VirtualBench](#) (class in *pymeasure.instruments.ni.virtualbench*), 169  
[VirtualBench.DigitalInputOutput](#) (class in *pymeasure.instruments.ni.virtualbench*), 169  
[VirtualBench.DigitalMultimeter](#) (class in *pymeasure.instruments.ni.virtualbench*), 171  
[VirtualBench.FunctionGenerator](#) (class in *pymeasure.instruments.ni.virtualbench*), 172  
[VirtualBench.MixedSignalOscilloscope](#) (class in *pymeasure.instruments.ni.virtualbench*), 174  
[VirtualBench.PowerSupply](#) (class in *pymeasure.instruments.ni.virtualbench*), 178  
[VirtualBench.Direct](#) (class in *pymeasure.instruments.ni.virtualbench*), 180  
[VISAAdapter](#) (class in *pymeasure.adapters*), 36  
[vision\\_average\\_enabled](#) (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 195  
[vision\\_balance](#) (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 195  
[vision\\_carrier\\_enabled](#) (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 195  
[vision\\_carrier\\_frequency](#) (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 195  
[vision\\_clamping\\_average](#) (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 195  
[vision\\_clamping\\_enabled](#) (*pymeasure.instruments.rohdeschwarz.sfm.SFM* property), 195

`sure.instruments.rohdeschwarz.sfm.SFM` (property), 195  
`vision_clamping_mode` (`pymea-  
sure.instruments.rohdeschwarz.sfm.SFM`  
property), 195  
`vision_precorrection_enabled` (`pymea-  
sure.instruments.rohdeschwarz.sfm.SFM`  
property), 196  
`vision_residual_carrier_level` (`pymea-  
sure.instruments.rohdeschwarz.sfm.SFM`  
property), 196  
`vision_sideband_filter_enabled` (`pymea-  
sure.instruments.rohdeschwarz.sfm.SFM`  
property), 196  
`vision_videosignal_enabled` (`pymea-  
sure.instruments.rohdeschwarz.sfm.SFM`  
property), 196  
`VMU` (class in `pymea-  
sure.instruments.agilent.agilent4156`), 82  
`voltage` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
property), 76  
`VOLTAGE` (`pymea-  
sure.instruments.agilent.agilentB1500.MeasOpModes` attribute), 102  
`voltage` (`pymea-  
sure.instruments.ametek.Ametek7270`  
property), 105  
`voltage` (`pymea-  
sure.instruments.attocube.anc300.Axis`  
property), 115  
`voltage` (`pymea-  
sure.instruments.bkprecision.BKPrecision9130B`  
property), 116  
`voltage` (`pymea-  
sure.instruments.deltaelektronika.SM7045B`  
property), 120  
`voltage` (`pymea-  
sure.instruments.keithley.Keithley2000`  
property), 131  
`voltage` (`pymea-  
sure.instruments.keithley.Keithley2260B`  
property), 133  
`voltage` (`pymea-  
sure.instruments.keithley.Keithley2400`  
property), 139  
`voltage` (`pymea-  
sure.instruments.keithley.Keithley2450`  
property), 145  
`voltage` (`pymea-  
sure.instruments.keithley.Keithley6517B`  
property), 156  
`voltage` (`pymea-  
sure.instruments.keysight.KeysightN5767A`  
property), 163  
`voltage` (`pymea-  
sure.instruments.Mock` property), 63  
`voltage` (`pymea-  
sure.instruments.signalrecovery.DSP7265`  
property), 203  
`voltage_1` (`pymea-  
sure.instruments.razorbill.razorbillRP100`  
property), 186  
`voltage_2` (`pymea-  
sure.instruments.razorbill.razorbillRP100`  
property), 186  
`voltage_ac` (`pymea-  
sure.instruments.agilent.Agilent34410A`  
property), 72  
`voltage_ac` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
property), 76  
`voltage_ac` (`pymea-  
sure.instruments.hp.HP34401A`  
property), 123  
`voltage_ac_auto_range` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
property), 76  
`voltage_ac_bandwidth` (`pymea-  
sure.instruments.keithley.Keithley2000` prop-  
erty), 131  
`voltage_ac_digits` (`pymea-  
sure.instruments.keithley.Keithley2000` prop-  
erty), 131  
`voltage_ac_nplc` (`pymea-  
sure.instruments.keithley.Keithley2000` prop-  
erty), 132  
`voltage_ac_range` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
property), 76  
`voltage_ac_range` (`pymea-  
sure.instruments.keithley.Keithley2000` prop-  
erty), 132  
`voltage_ac_reference` (`pymea-  
sure.instruments.keithley.Keithley2000` prop-  
erty), 132  
`voltage_ac_resolution` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
property), 76  
`voltage_auto_range` (`pymea-  
sure.instruments.agilent.Agilent34450A`  
property), 76  
`voltage_dc` (`pymea-  
sure.instruments.agilent.Agilent34410A`  
property), 72  
`voltage_dc` (`pymea-  
sure.instruments.hp.HP34401A`  
property), 123  
`voltage_digits` (`pymea-  
sure.instruments.keithley.Keithley2000` prop-  
erty), 132  
`voltage_filter_count` (`pymea-  
sure.instruments.keithley.Keithley2450` prop-  
erty), 145  
`voltage_filter_type` (`pymea-  
sure.instruments.keithley.Keithley2450` prop-  
erty), 145  
`voltage_high` (`pymea-  
sure.instruments.agilent.Agilent33220A`  
property), 84  
`voltage_high` (`pymea-  
sure.instruments.agilent.Agilent33500` prop-  
erty), 88  
`voltage_limit` (`pymea-  
sure.instruments.ami.AMI430`  
property), 106  
`voltage_limit` (`pymea-  
sure.instruments.yokogawa.YokogawaGS200`  
property), 217  
`voltage_low` (`pymea-  
sure.instruments.agilent.Agilent33220A`

property), 84

voltage\_low (pymeasure.instruments.agilent.Agilent33500 property), 88

voltage\_name (pymeasure.instruments.agilent.agilent4156.SMU property), 80

voltage\_name (pymeasure.instruments.agilent.agilent4156.VMU property), 82

voltage\_name (pymeasure.instruments.agilent.agilent4156.VSU property), 82

voltage\_nplc (pymeasure.instruments.keithley.Keithley2000 property), 132

voltage\_nplc (pymeasure.instruments.keithley.Keithley2400 property), 140

voltage\_nplc (pymeasure.instruments.keithley.Keithley2450 property), 145

voltage\_nplc (pymeasure.instruments.keithley.Keithley6517B property), 156

voltage\_output\_off\_state (pymeasure.instruments.keithley.Keithley2450 property), 145

voltage\_range (pymeasure.instruments.agilent.Agilent34450A property), 76

voltage\_range (pymeasure.instruments.keithley.Keithley2000 property), 132

voltage\_range (pymeasure.instruments.keithley.Keithley2400 property), 140

voltage\_range (pymeasure.instruments.keithley.Keithley2450 property), 145

voltage\_range (pymeasure.instruments.keithley.Keithley6517B property), 156

voltage\_range (pymeasure.instruments.keysight.KeysightN5767A property), 163

voltage\_reference (pymeasure.instruments.keithley.Keithley2000 property), 132

voltage\_resolution (pymeasure.instruments.agilent.Agilent34450A property), 76

voltage\_setpoint (pymeasure.instruments.keithley.Keithley2260B property), 133

VSU (class in pymeasure.instruments.agilent.agilent4156), 82

VXI11Adapter (class in pymeasure.adapters), 38

## W

wait() (pymeasure.instruments.anritsu.AnritsuMS9710C method), 112

wait\_for\_buffer() (pymeasure.instruments.keithley.Keithley2000 method), 132

wait\_for\_buffer() (pymeasure.instruments.keithley.Keithley2400 method), 140

wait\_for\_buffer() (pymeasure.instruments.keithley.Keithley2450 method), 145

wait\_for\_buffer() (pymeasure.instruments.keithley.Keithley2700 method), 148

wait\_for\_buffer() (pymeasure.instruments.keithley.Keithley6221 method), 152

wait\_for\_buffer() (pymeasure.instruments.keithley.Keithley6517B method), 156

wait\_for\_buffer() (pymeasure.instruments.signalrecovery.DSP7265 method), 203

wait\_for\_buffer() (pymeasure.instruments.srs.SR830 method), 208

wait\_for\_completion() (pymeasure.instruments.anaheimautomation.DPSeriesMotorController method), 109

wait\_for\_current() (pymeasure.instruments.danfysik.Danfysik8500 method), 118

wait\_for\_data() (pymeasure.experiment.experiment.Experiment method), 42

wait\_for\_holding() (pymeasure.instruments.ami.AMI430 method), 106

wait\_for\_ready() (pymeasure.instruments.danfysik.Danfysik8500 method), 118

wait\_for\_srq() (pymeasure.adapters.PrologixAdapter method), 36

wait\_for\_srq() (pymeasure.adapters.VISAAdapter method), 37

wait\_for\_stop() (pymeasure.instruments.newport.esp300.Axis method), 168

wait\_for\_sweep() (pymeasure.instruments.anritsu.AnritsuMS9710C method), 112

- `wait_for_temperature()` (*pymea-  
sure.instruments.lakeshore.LakeShore331  
method*), 166
- `wait_for_temperature()` (*pymea-  
sure.instruments.oxfordinstruments.ITC503  
method*), 182
- `wait_for_trigger()` (*pymea-  
sure.instruments.agilent.Agilent33220A  
method*), 84
- `wait_for_trigger()` (*pymea-  
sure.instruments.agilent.Agilent33500 method*),  
88
- `wait_time()` (*pymea-  
sure.instruments.agilent.agilentB1500.AgilentB1500  
method*), 95
- `WaitTimeType` (class in *pymea-  
sure.instruments.agilent.agilentB1500*), 103
- `wave` (*pymea-  
sure.instruments.Mock property*), 63
- `waveform_abort()` (*pymea-  
sure.instruments.keithley.Keithley6221  
method*), 152
- `waveform_amplitude` (*pymea-  
sure.instruments.keithley.Keithley6221 prop-  
erty*), 152
- `waveform_arm()` (*pymea-  
sure.instruments.keithley.Keithley6221  
method*), 153
- `waveform_data` (*pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
property*), 162
- `waveform_duration_cycles` (*pymea-  
sure.instruments.keithley.Keithley6221 prop-  
erty*), 153
- `waveform_duration_set_infinity()` (*pymea-  
sure.instruments.keithley.Keithley6221  
method*), 153
- `waveform_duration_time` (*pymea-  
sure.instruments.keithley.Keithley6221 prop-  
erty*), 153
- `waveform_dutycycle` (*pymea-  
sure.instruments.keithley.Keithley6221 prop-  
erty*), 153
- `waveform_format` (*pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
property*), 162
- `waveform_frequency` (*pymea-  
sure.instruments.keithley.Keithley6221 prop-  
erty*), 153
- `waveform_function` (*pymea-  
sure.instruments.keithley.Keithley6221 prop-  
erty*), 153
- `waveform_offset` (*pymea-  
sure.instruments.keithley.Keithley6221 prop-  
erty*), 153
- `waveform_points` (*pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
property*), 162
- `waveform_points_mode` (*pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
property*), 162
- `waveform_preamble` (*pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
property*), 162
- `waveform_ranging` (*pymea-  
sure.instruments.keithley.Keithley6221 prop-  
erty*), 153
- `waveform_source` (*pymea-  
sure.instruments.keysight.KeysightDSOX1102G  
property*), 162
- `waveform_start()` (*pymea-  
sure.instruments.keithley.Keithley6221  
method*), 153
- `waveform_use_phasemarker` (*pymea-  
sure.instruments.keithley.Keithley6221 prop-  
erty*), 153
- `wavelength` (*pymea-  
sure.instruments.thorlabs.ThorlabsPM100USB  
property*), 214
- `wavelength_center` (*pymea-  
sure.instruments.anritsu.AnritsuMS9710C  
property*), 112
- `wavelength_marker_value` (*pymea-  
sure.instruments.anritsu.AnritsuMS9710C  
property*), 112
- `wavelength_max` (*pymea-  
sure.instruments.thorlabs.ThorlabsPM100USB  
property*), 214
- `wavelength_min` (*pymea-  
sure.instruments.thorlabs.ThorlabsPM100USB  
property*), 214
- `wavelength_span` (*pymea-  
sure.instruments.anritsu.AnritsuMS9710C  
property*), 112
- `wavelength_start` (*pymea-  
sure.instruments.anritsu.AnritsuMS9710C  
property*), 112
- `wavelength_stop` (*pymea-  
sure.instruments.anritsu.AnritsuMS9710C  
property*), 112
- `wavelength_value_in` (*pymea-  
sure.instruments.anritsu.AnritsuMS9710C  
property*), 112
- `wavelengths` (*pymea-  
sure.instruments.anritsu.AnritsuMS9710C  
property*), 112
- `wires` (*pymea-  
sure.instruments.keithley.Keithley2400  
property*), 140
- `wires` (*pymea-  
sure.instruments.keithley.Keithley2450  
property*), 146
- `Worker` (class in *pymea-  
sure.experiment.workers*), 46
- `write()` (*pymea-  
sure.adapters.Adapter method*), 32

[write\(\)](#) ([pymeasure.adapters.FakeAdapter](#) method), 33  
[write\(\)](#) ([pymeasure.adapters.PrologixAdapter](#) method), 36  
[write\(\)](#) ([pymeasure.adapters.SerialAdapter](#) method), 34  
[write\(\)](#) ([pymeasure.adapters.TelnetAdapter](#) method), 40  
[write\(\)](#) ([pymeasure.adapters.VISAAadapter](#) method), 38  
[write\(\)](#) ([pymeasure.adapters.VXI11Adapter](#) method), 39  
[write\(\)](#) ([pymeasure.instruments.agilent.agilentB1500.SMU](#) method), 97  
[write\(\)](#) ([pymeasure.instruments.anaheimautomation.DPSeriesMotorController](#) method), 109  
[write\(\)](#) ([pymeasure.instruments.attocube.adapters.AttocubeConsoleAdapter](#) method), 113  
[write\(\)](#) ([pymeasure.instruments.danfysik.DanfysikAdapter](#) method), 116  
[write\(\)](#) ([pymeasure.instruments.Instrument](#) method), 63  
[write\(\)](#) ([pymeasure.instruments.keysight.KeysightDSOX1102G](#) method), 163  
[write\(\)](#) ([pymeasure.instruments.lakeshore.LakeShoreUSBAdapter](#) method), 165  
[write\(\)](#) ([pymeasure.instruments.ni.virtualbench.VirtualBenchDigitalInputOutputPointer](#) method), 170  
[write\(\)](#) ([pymeasure.instruments.parker.ParkerGV6](#) method), 184  
[write\\_binary\\_values\(\)](#) ([pymeasure.adapters.PrologixAdapter](#) method), 36  
[write\\_binary\\_values\(\)](#) ([pymeasure.adapters.SerialAdapter](#) method), 34  
[write\\_binary\\_values\(\)](#) ([pymeasure.adapters.VISAAadapter](#) method), 38  
[write\\_binary\\_values\(\)](#) ([pymeasure.instruments.lakeshore.LakeShoreUSBAdapter](#) method), 165  
[write\\_raw\(\)](#) ([pymeasure.adapters.VXI11Adapter](#) method), 39  
[writeA0\(\)](#) (in module [pymeasure.instruments.comedi](#)), 65

## X

[x](#) ([pymeasure.instruments.ametek.Ametek7270](#) property), 105  
[x](#) ([pymeasure.instruments.signalrecovery.DSP7265](#) property), 203  
[x](#) ([pymeasure.instruments.srs.SR830](#) property), 208  
[x](#) ([pymeasure.instruments.srs.SR860](#) property), 213  
[x1](#) ([pymeasure.instruments.ametek.Ametek7270](#) property), 105  
[x2](#) ([pymeasure.instruments.ametek.Ametek7270](#) property), 105  
[xpointer](#) ([pymeasure.instruments.oxfordinstruments.ITC503](#) property), 182

[xy](#) ([pymeasure.instruments.ametek.Ametek7270](#) property), 105  
[xy](#) ([pymeasure.instruments.signalrecovery.DSP7265](#) property), 203  
[xy](#) ([pymeasure.instruments.srs.SR830](#) property), 208

## Y

[y](#) ([pymeasure.instruments.ametek.Ametek7270](#) property), 105  
[y](#) ([pymeasure.instruments.signalrecovery.DSP7265](#) property), 203  
[y](#) ([pymeasure.instruments.srs.SR830](#) property), 208  
[y](#) ([pymeasure.instruments.srs.SR860](#) property), 213  
[y1](#) ([pymeasure.instruments.ametek.Ametek7270](#) property), 105  
[y2](#) ([pymeasure.instruments.ametek.Ametek7270](#) property), 105  
[Yokogawa7651](#) (class in [pymeasure.instruments.yokogawa](#)), 215  
[YokogawaGS200](#) (class in [pymeasure.instruments.yokogawa](#)), 216  
[ypointer](#) ([pymeasure.instruments.oxfordinstruments.ITC503](#) property), 182

## Z

[zero\(\)](#) ([pymeasure.instruments.ami.AMI430](#) method), 107  
[zero\(\)](#) ([pymeasure.instruments.newport.esp300.Axis](#) method), 168  
[zero\\_probe\(\)](#) ([pymeasure.instruments.lakeshore.LakeShore425](#) method), 167